

# **Versatile Polymethacrylate Standards**

## **Agilent PLgel 10 $\mu$ m MIXED-B Columns**

### **Technical Overview**

#### **Introduction**

PLgel 10  $\mu$ m MIXED-B columns are designed for high MW polymer analysis and demanding eluent conditions. The PLgel 10  $\mu$ m MIXED-B spans a wide range of molecular weights, up to 10 million, with a linear calibration curve.



**Agilent Technologies**

Methyl ethyl ketone (MEK) is a good solvent for PMMA and can be used with PLgel columns. Although not ideal for use with a UV detector, MEK has a lower refractive index than THF, a property which may be used to some advantage in optimizing RI signal response (see Figures 1 and 2).

PLgel 10  $\mu$ m MIXED-B is particularly useful for molecular weight distributions where slightly higher than average MWs are encountered. The 10  $\mu$ m particle size provides good resolution with relatively low pressures for enhanced lifetimes in demanding conditions.

#### Conditions

Columns: 2 x PLgel 10  $\mu$ m MIXED-B,  
300 x 7.5 mm (p/n PL1110-6100)  
Eluent: Methyl ethyl ketone  
Flow Rate: 1.0 mL/min  
Detection: RI




Figure 1. PMMA calibration curve




Figure 2. Separation of Agilent PMMA standards

[www.agilent.com/chem](http://www.agilent.com/chem)

This information is subject to change without notice.

© Agilent Technologies, Inc. 2010

Published in UK, September 23, 2010

SI-01777



**Agilent Technologies**