
## User's Manual

Nucleic Acid Isolation Device

# QuickGene-Mini80

### Extraction Flow

#### Preparation

Prepare lysate using  
the Isolation Kits.

#### Extraction with QuickGene-Mini80 (Pages 16 to 25)

Applying lysate into the Cartridges → Wash → Elution of DNA/ RNA

# Introduction

Thank you for purchasing the Nucleic Acid Isolation Device QuickGene-Mini80 (hereafter, called "QuickGene-Mini80"). DNA or RNA can be easily extracted by using QuickGene-Mini80.

This guide contains important information regarding the safe use of QuickGene-Mini80. Read this manual carefully before using the device.

Keep this manual after reading it, for reference whenever necessary.

## About Infectious Wastes

When handling or disposing of infectious materials, follow the laboratory guideline or the law regarding infectious waste to perform proper incineration, fusion, sterilization, and/or disinfection.

When you use a third party to dispose of it, outsource this work an operator licensed to handle medical waste subject to special control, and give them the manifest of the medical waste at the same time.

Universal biological hazard symbol

### ■ Exclusion clauses

- The contents of this manual may be updated without notice.
- The FUJIFILM Corporation assumes no liability for infringement of third party copyrights or other rights arising from the use of any information in this manual.
- The FUJIFILM Corporation assumes no liability for device failure or damage due to installation, relocation, remodeling, maintenance, and repairs performed by parties other than the FUJIFILM Corporation or agencies specified by the FUJIFILM Corporation.
- The FUJIFILM Corporation assumes no liability for failure or damage due to use of third-party products.
- The FUJIFILM Corporation assumes no liability for device failure or damage due to remodeling, maintenance, repairs, etc. using repair parts other than genuine parts specified by the FUJIFILM Corporation.
- The FUJIFILM Corporation assumes no liability for device failure or damage due to failure to observe the precautions and operation procedures described in this manual.
- The FUJIFILM Corporation assumes no liability for device failure or damage due to using the device in an environment that does not meet the operational requirements (power, installation environment, etc.) described in this manual.
- The FUJIFILM Corporation assumes no liability for device failure or damage due to natural disasters such as fire, earthquakes, floods, or lightning.
- No part of the contents of this manual may be reproduced without permission from the FUJIFILM Corporation.

### ■ Trademarks

Company names and product names are trademarks or registered trademarks of their respective owners.

---

# Notations

This manual contains important information regarding the safe use of QuickGene-Mini80.

Read this information carefully before using QuickGene-Mini80, and start operation.

The meanings of safety precaution marks are as follows:

**WARNING:** "WARNING" indicates a dangerous condition that may lead to death or serious injury.

**BIOHAZARD:** "BIOHAZARD" indicates exposure to a biological agent or condition that may be infectious or harmful.

**CAUTION:** "CAUTION" indicates the possibility of danger that may lead to minor or moderate personal injury.  
"CAUTION" may also indicate an accident of property damage.

---

## Important:

"Important" shows the important notes for usage, as well as prohibited actions.

---

## Note:

"Note" indicates the notes, procedures that should be obeyed, and supplementary information for use.

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---

# Safety Information

This chapter describes safety precautions.

For your safety and that of others, follow the guidelines provided in the following pages concerning the use of QuickGene-Mini80.

## About Mark

Check all descriptions with the  mark attached.

## About Instrument

### **WARNING:** • Ignoring the following notations may lead to fire or electric shock.

- In countries other than US and Canada, do not use the power cable attached for the US and Canadian models.
- In countries other than US and Canada, use a power cable that meets your country's standard, or contact your local distributor.
- Use the AC adapter provided.
- Do not use the power cable attached for the US and Canadian models for other products.
- Do not use QuickGene-Mini80 with voltage other than the voltage specified on the device. Do not overload the power outlet with too many devices.
- Do not use QuickGene-Mini80 with a damaged power plug or a loose socket.
- If there is dust on the prongs of the power plug or on the plug socket, remove it with a dry cloth.
- When you disconnect the plug from the outlet, be sure to hold the power plug itself. Do not pull the power cable.
- For maintenance, disconnect the power plug from the outlet.
- Do not touch the power plug when you hear the crash of thunder.

### **• Do not pour any liquid on QuickGene-Mini80. Do not place any objects containing liquid on QuickGene-Mini80. Doing so may cause a device failure, fire, or electric shock.**

### **• In the event the device overheats, starts to smoke or smells strange, immediately unplug the power cable.**

### **CAUTION:** • Never attempt to remodel QuickGene-Mini80 without the manufacturer's permission. Doing so may cause fire or electric shock.

### **• Do not place or drop objects on QuickGene-Mini80. Also refrain from bumping or knocking it, as doing so may cause a failure or malfunction of QuickGene-Mini80.**

### **• If any liquid materials are left inside the device, wipe it up a soft paper tissue, etc.**

**Otherwise, QuickGene-Mini80 may be damaged.**

### **• When replacing the pressure heads, be sure to wear gloves.**

### **• If ethanol is spilled inside the device, immediately unplug the power cable.**

### **• Repairs to QuickGene-Mini80 should only be performed by such agencies as are specifically authorized by the FUJIFILM Corporation.**

### **• Only original the FUJIFILM Corporation replacement parts should be used.**

### **• If the equipment is used in a manner not specified by the manufacturer, the protection provided by the equipment may be impaired.**

## About isolation kits

Referred to as "kit/kits" hereafter in this section.

**CAUTION:**

- When handling any of the kits, refer to the respective handbook.
- Reagents in each kit should be handled by observing the safety information and precautions regarding the kit.
- Extraction should be performed in an appropriate laboratory or workplace.

isolation kits	
For DNA	DNA whole blood kit
	DNA tissue kit
	Plasmid kit
For RNA	RNA tissue kit
	RNA cultured cell kit
	RNA cultured cell HC kit
	RNA blood cell kit

**Note:**

The kits are not supplied with QuickGene-Mini80. Select the desired kit(s) and order it (them) separately.

## About Reagents Not Included in the Kits

**CAUTION:**

- When handling any reagent not included in the kits, observe the safety information and precautions of the reagent.
- Extraction should be performed in an appropriate laboratory or workplace. Some reagents should be handled in a draft chamber.

## About Cleaning/Washing Solution

**CAUTION:**

- For handling any cleaning/washing solution for the QuickGene-Mini80 and accessories, observe the safety information and precautions regarding the solution.
- Keep ethanol away from fire. It is highly flammable.

**BIOHAZARD:**

- Always wear appropriate gloves, a mask, and safety goggles when cleaning/washing.

## About Samples

**BIOHAZARD:**

- Always wear appropriate gloves, a mask, and safety goggles etc. when handling any infectious samples.

---

# Product Use Limitations

- **The FUJIFILM Corporation assumes no liability for results regardless of the application.**
- **It is the user's responsibility to validate the performance of the Isolation Kits and QuickGene-Mini80.**
- **QuickGene-Mini80 is a device for research and study only.**

When you use QuickGene-Mini80, read this manual thoroughly and follow the safety management regulations of your facilities (laboratory).

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## 1 Setting Up

1

This chapter explains how to set up QuickGene-Mini80.

## 2 Operation

2

This chapter explains the extraction procedures.

## 3 Maintenance

3

This chapter explains how to maintain QuickGene-Mini80.

## 4 Troubleshooting

4

This chapter explains the procedures that should be taken when you suspect device failures.

## 5 Main Specifications

5

This chapter explains the basic specifications of QuickGene-Mini80.

## Appendix A

A

This appendix describes after-sales services and points to note for transportation.

---

# Contents

---

## 1 Setting Up

---

1.1	Operating Environment/Use Conditions .....	1
1.2	Checking QuickGene-Mini80 and Accessories .....	3
1.3	Remove the Packing Materials.....	4
1.4	Overall View .....	5
1.5	Basic Functions.....	8
1.6	Basic Operation .....	9
1.7	Precautions on Transportation .....	10
1.8	Precautions on Storage .....	11

## 2 Operation

---

2.1	Turning the Power ON and OFF .....	13
2.2	Preparation.....	15
2.3	Extraction.....	16

## 3 Maintenance

---

3.1	Cleaning QuickGene-Mini80 .....	27
3.2	Cleaning the Parts .....	28
3.3	How to Clean/Replace the Pressure Heads .....	28

## 4 Troubleshooting

---

## 5 Main Specifications

---

## Appendix A

---

A.1	Replacement/Optional Parts.....	35
A.2	Product Guarantees .....	35
A.3	About Services .....	36
A.4	Notes for Transporting/Shipping QuickGene-Mini80.....	36

# 1 Setting Up

This chapter explains how to set up QuickGene-Mini80.

## 1.1 Operating Environment/Use Conditions

This section explains the operating environment and use conditions of QuickGene-Mini80.

### ■Operating environment


- Warning:**
- Do not use QuickGene-Mini80 to a location where it can be splashed with water. It may cause a device failure, fire, or electric shock.
  - When relocating QuickGene-Mini80, disconnect the plug from the outlet first. If the power cable is damaged, this may cause a device failure, fire, injury, or electric shock.

- Caution:**
- Do not use QuickGene-Mini80 in an unstable place such as a slanted surface or a place subject to vibrations. It may cause injury or device failure.
  - Do not use QuickGene-Mini80 in direct sunlight or close to a heating device. It may shorten the life of QuickGene-Mini80, or cause a trouble.

For correct and safe use of QuickGene-Mini80, use it in a location that meets the following conditions:

- The space shown in the following figure must be allocated.

- A location where power can be provided (This is a Class II device. Be sure to use the attached AC adapter. The  mark is the mark of the Class II device.)
- A location where the temperature is 15 to 30°C, and humidity is 30 to 80%RH (non condensing)
- A location that is flat and stable, with no vibration
- A location away from direct sunlight (Block the sunlight by closing curtains or blinds as necessary.)
- A location which is well-ventilated and not dusty
- A location where the temperature does not go up and down suddenly (Warming a cold room suddenly or moving QuickGene-Mini80 from a room with low temperature to a warm room may cause condensation inside the device, resulting in abnormal extraction.)
- A location where the temperature and humidity are kept within the specified range (far from water taps, water heaters, humidifiers, air-conditioners and heaters)

## 1 Setting Up

- A location far from objects which generate strong magnetic fields (motors, transformers, TV, audio speakers, magnets, etc.) (Bringing QuickGene-Mini80 close to any type of magnetic field may cause a malfunction.)

### ■Use conditions


1

Items		Specifications
Temperature (°C)	During operation	15 - 30
	During downtime	0 - 55
Humidity (RH)	During operation	30 - 80
	During downtime	10 - 80
Max. wet bulb temperature (°C)	During operation	29 (non condensing)
	During downtime	29 (non condensing)
Temperature gradient (°C/hr)		12 or less (non condensing)
Humidity gradient (RH/day)		30 or less (non condensing)
Altitude, operating (m)		1600 or less

## 1.2 Checking QuickGene-Mini80 and Accessories

Check that the following are included in the package.

Contact your local representative if you notice any missing part(s).

### Caution:

**Note: Use an appropriate power cable by following the descriptions below.**

**US/Canada: Use the power cable included in the package.**

**Other countries: The power cable is not included in the package. Prepare the power cable that meets your country's standard.**

**The specification of the power cable for QuickGene-Mini80 is as follows.**

#### 100-120V area

**Power Voltage: 100 to 120 V**

**Requirements for the plug/connector: AC 125 V 3 A**

**Requirements for the cable: SJT 2 x 18 AWG 60°C**

#### 200-240V area

**Power Voltage: 200 to 240 V**

**Requirements for the plug/connector: AC 250 V 3 A**

**Requirements for the cable: CENELEC 2 x 0.75 mm<sup>2</sup> 70°C**

### 1.3 Remove the Packing Materials

After unpacking QuickGene-Mini80 from the box, proceed as follows:


- (1) Remove the tape that holds the rotary switch and cushion, and then remove the cushion.

**Important: Keep the shipping box and the removed cushion. They will be needed again when transporting QuickGene-Mini80**

---

**Note: The cushion is for protecting pressure heads from being damaged.**


---

- (2) Place QuickGene-Mini80 upside down and then make sure that the pressure heads are fully inserted.

**Note: Place QuickGene-Mini80 on a soft cloth, etc. to prevent damage to the body.**


---


## 1.4 Overall View

This section explains the name and function of each part in QuickGene-Mini80.

### ■Front and side view


### ■Rear view

## 1 Setting Up


### ■Cartridge holder

1


Name	Description
Release Lever	<p>Releases the Cartridges by sliding the lever.</p> <p>All Cartridges are released by sliding the lever to the right-hand end.</p> <p>To release the No.1 to No.4 Cartridges only, stop the lever at the center position.</p> <p><b>Note: Before setting Cartridges, make sure that the release lever is set in the left-hand position.</b></p>

### ■Tube holder


### ■Label

Name	Description
ON/Off Mark	<ul style="list-style-type: none"> <li>ON Mark: Rotating the switch to the ON Mark starts pressure application.</li> <li>OFF Mark: Rotating the switch to the OFF Mark finishes pressure application.</li> </ul>
WASH Label	Indicates that tube holder and cartridge holder are set in the Wash Position (W) correctly. When the tube holder and cartridge holder are set in the Elution Position (E), this label is hidden below the tube holder.
Waste Fluid Scale	<p>The solution level changes according to the number of wash processes.</p> <p>L: Pressure application of lysate completed            W1: First wash process completed            W2: Second wash process completed            W3 → E: Third wash process completed            After wash processes are completed, move the cartridge holder to the Elution Position (E) to perform elution process.</p>

### ■Cover

This cover is for storage purposes.

**Caution:** •Do not operate the QuickGene-Mini80 with the cover on.


•Always put the cover on the QuickGene-Mini80, after cleaning the outside of the body (refer to P. 27).

•Do not use the cover inside out.

**BIOHAZARD:** Wear appropriate gloves when handling the cover.

## 1.5 Basic Functions

The following describes the main components of QuickGene-Mini80.

Component	Function
Main Unit	<p>Starts/stops putting pressurized air into Cartridges by operating the rotary unit. Pressure application automatically stops in about one minute.</p> <p>Operating state is indicated by the lamp (LED). When the power is turned on, the lamp (LED) lights. During pressure application, the lamp (LED) blinks.</p> <p>The main unit contains a pump for pressurizing air, a sensor for detecting start/stop of pressure application, and a circuit for controlling the duration of pressure application.</p>
Pressure Unit	Puts pressurized air into the Cartridges through the pressure heads. When the rotary switch is turned on, the pressure heads move down to put pressurized air into Cartridges.
Rack Unit	<p>The rack unit is comprised of the cartridge holder and tube holder. Lysate/Wash Buffer/Elution Buffer put in a Cartridge passes through the filter in the Cartridge by pressurized air and filtered solution drops into a Waste Tube or Collection Tube.</p> <p><b>Note:</b> Pressure application of lysate and wash, and that of elution operation are switched by moving the position of the cartridge holder manually.</p>

### ■Lamp (LED) indication information

The following describes the lamp (LED) indications.

Lamp (LED) Indication	Description
ON	Indicates that the power is ON and pressure application can be performed.
Blink	Indicates that pressure application is in progress.

**Note:** When pressure application is completed, the lamp (LED) indication returns to the ON state.

## 1.6 Basic Operation

The following describes the working flow of QuickGene-Mini80:


- (1) Set the Waste Tubes and Collection Tubes into the tube holder and then insert the cartridge holder into the Wash Position (W) of the tube holder.
- (2) Set the Cartridges into the cartridge holder.
- (3) Apply lysate to the Cartridges and then set the cartridge holder and tube holder in QuickGene-Mini80.
- (4) Rotate the rotary switch toward the front side to apply air pressure to the Cartridges.
- (5) Make sure that lysate is discharged completely and then return the rotary switch to the original position to terminate pressure application.
- (6) Pull out the cartridge holder and tube holder.
- (7) Apply Wash Buffer to the Cartridges and then set the cartridge holder and tube holder in QuickGene-Mini80.
- (8) Rotate the rotary switch toward the front side to apply air pressure to the Cartridges.
- (9) Make sure that Wash Buffer is discharged completely and then return the rotary switch to the original position to terminate pressure application.
- (10) Repeat Steps 6 to 9 depending on the protocol of the relevant Isolation Kits.
- (11) Pull out the cartridge holder and tube holder and then insert the cartridge holder into the Elution Position (E) of the tube holder.
- (12) Apply Elution Buffer to the Cartridges and then set the cartridge holder and tube holder in QuickGene-Mini80.
- (13) Rotate the rotary switch toward the front side to apply air pressure to the Cartridges.
- (14) Make sure that Elution Buffer is discharged completely and then return the rotary switch to the original position to terminate pressure application.

The extracted sample is eluted into the Collection Tube and the waste is discharged into the Waste Tube.

## 1.7 Precautions on Transportation

When transporting QuickGene-Mini80, note the following points.


- Be sure to hold the handles.
- Do not hold the rotary switch.
- Do not hold the tray.
- Do not insert your hands into the pressure head section.
- Remove the tube holder before transportation.

## 1.8 Precautions on Storage

When storing QuickGene-Mini80, note the following points.

- Do not keep it with the rotary switch adjusted to the position of turning on with Cartridges set in the cartridge holder. If not, the rubber packing on the pressure heads may become deformed and proper operation cannot be guaranteed.
- Put the supplied cover on QuickGene-Mini80. When the outside of the body is heavily soiled, clean it before fitting the cover.
- Wear appropriate gloves when handling the cover.

# 2 Operation

This chapter explains how to operate QuickGene-Mini80 and shows the flow of extraction.

**BIOHAZARD:**

Always wear appropriate gloves, a mask, and safety goggles during any biohazardous operations in extraction process. Even when touching the device after any operation with a biohazard risk, wear appropriate gloves and a mask since the device may be contaminated.

**Important:** Before starting extraction, put on appropriate gloves, a mask, and safety goggles if required by the operation. In the operation from preparing samples to extraction completion, be careful not to contaminate the samples with sweat, saliva, etc.


2

## 2.1 Turning the Power ON and OFF

This section explains how to turn the power of QuickGene-Mini80 on and off.


**Note:** Since QuickGene-Mini80 is not equipped with a power switch, the power is turned on/off by connecting/disconnecting the power cable.

### ■Turning on the power

- (1) Make sure that the AC adapter and power cable are connected securely, and the plug of the AC adapter is securely connected to the QuickGene-Mini80.

**Note:** QuickGene-Mini80 is a Class II device. Be sure to use the attached AC adapter. The mark is the mark of the Class II device.


- (2) Connect the plug of the power cable to the outlet.

The QuickGene-Mini80 turns on and the lamp (LED) lights.

## 2 Operation

---

### ■Turning off the power

- (1) Disconnect the plug of the power cable from the outlet.

The lamp (LED) turns off and then the QuickGene-Mini80 turns off.

**Note:** Do not place any objects near the outlet so that the plug of the power cable can be disconnected immediately in case of emergency.

---

## 2.2 Preparation

The following preparations are required for extraction operation.

### ■Items to be prepared

- Gloves
- Mask
- Safety Goggles
- Collection Tubes (CT) (\*1, \*2)
- Waste Tubes (WT) (\*1)
- Cartridges (CA) (\*1)
- Buffer (Wash Buffer, Elution Buffer) (\*1)
- Lysate (Samples)

For the preparation of samples, refer to the handbook of each kit.

\*1: Provided with the Isolation Kits.

Tube adapters are required when using Collection Tubes provided with the Isolation Kits ("When using Collection Tubes included in the Isolation Kits:"(P.17)).

\*2: Commercially-available 1.5 ml micro tubes can be used in place of Collection Tubes.

2

### ■Preparation of lysate

Prepare lysate.


Make lysate using the specified kit for each target sample. For details, refer to the handbook of the Isolation Kits.

### 2.3 Extraction


This section explains extraction procedure including before and after extraction operations. Be sure to perform extraction in the following order.

#### ■Before Extraction Operation


2

- (1) Turn on the power by connecting the power cable and then make sure the lamp (LED) lights.


- (2) Pull out the tube holder and cartridge holder.


- (3) Remove the cartridge holder from the tube holder.

**Note:** Pull up the cartridge holder vertically to prevent it becoming caught in the tube holder.



**(4) Set Collection Tubes and Waste Tubes into the tube holder.**

- **When using Collection Tubes included in the Isolation Kits:**  
Set a tube adapter into the tube holder as shown and then insert a Collection Tube into the tube adapter.


- **When using commercially-available 1.5 ml micro tubes as collection tubes:**  
Set 1.5 ml micro tubes so that the cap is positioned as shown. Tube adapters are not used.

**(5) Insert the cartridge holder into the Wash Position (W) of the tube holder.**

## 2 Operation

2


### (6) Set Cartridges into the cartridge holder.

#### Important:


- Make sure that the release lever is positioned at the left-hand end. If any Cartridges are set without the release lever positioned at the left-hand end, the Cartridges will drop into the Waste Tubes.
- Align the fins on the Cartridge with the grooves in the cartridge holder and then push in the Cartridge.
- Set Cartridges into the cartridge holder in the order of the number on the cartridge holder.

### (7) Check the location of the Tubes and Cartridges.

**Important:** Put the tubes into the tube holder in the order of the number on the tube holder. Insert the Cartridges and Tubes starting from the left, in order. If they are not in order, Buffer and/or lysate may spill. This results in an extraction failure, and Buffer and/or lysate becoming waste.

### ■ Extraction


#### <Applying Lysate to Cartridges>

##### (1) Apply lysate to each Cartridge.

Use a pipette to transfer the lysate into the Cartridges. In this operation, take care not to touch the filter in the Cartridge, the Cartridge itself or the cartridge holder.


**Note:** If any lysate is left on the rim of a Cartridge or on the cartridge holder, wipe it off with a soft paper tissue.


- (2) Set the tube holder and cartridge holder into the QuickGene-Mini80.

**Important:** After setting it, check the following:

- The cartridge holder is placed against the stopper inside the QuickGene-Mini80.
- The WASH label can be seen.

2

- (3) Rotate the rotary switch toward the front side to lower the pressure heads.


- (4) The lamp (LED) blinks and pressurized air is put into the Cartridges.

**Note:**

- When the solution in each Cartridge is discharged completely before pressure application automatically stops, return the rotary switch to the original position to terminate pressure application (Refer to Step 6).
- Pressure application automatically stops in about one minute.

## 2 Operation

2


- (5) Make sure that the lysate does not remain in the Cartridges.

**Important:** Make sure that no liquid remains in the Cartridges. Residual solution in the Cartridges may cause contamination.


**Note:**

- If any solution remains in the Cartridges even after about one minute has elapsed, return the rotary switch to the original position (Lamp (LED): ON) and then rotate it toward the front side (Lamp (LED): Blink) to put pressurized air into the Cartridges again.
- If the rotary switch is returned swiftly to the original position, it may result in the splashing of solutions. Accordingly, rotate the rotary switch slowly.
- If any solution in a Cartridge is not filtered even after pressure application is repeated several times, remove the Cartridge by pulling it up vertically and then take appropriate measures referring to the troubleshooting procedures described in the handbook of the Isolation Kits.

- (6) Return the rotary switch to the original position to raise the pressure heads.

**Note:** Do not pull out the tube holder or cartridge holder while the pressure heads are lowered.



- (7) Pull out the tube holder and cartridge holder.

## &lt;Wash&gt;

**(8) Apply Wash Buffer to each Cartridge.**


Use a pipette to apply Wash Buffer to the Cartridges. In this operation, take care not to touch the filter in the Cartridge, the Cartridge itself or the cartridge holder.

**Note:** If any Wash Buffer is left on the rim of a Cartridge or on the cartridge holder, wipe it off with a soft paper tissue.


**(9) Set the tube holder and cartridge holder into the QuickGene-Mini80 again.**

**Important:** After setting it, check the following:

- The cartridge holder is placed against the stopper inside the QuickGene-Mini80.
- The WASH label can be seen.


**(10) Repeat the wash process the necessary number of times according to the protocols.****Note:**

- Refer to Steps (3) to (8).
- The waste fluid scale indicates the number of wash processes.  
L: Pressure application of lysate completed  
W1: First wash process completed  
W2: Second wash process completed  
W3 → E: Third wash process completed  
When the wash processes are completed, move the cartridge holder to the Elution Position (E) to perform the elution process.

## 2 Operation


2

- (11) When wash process is completed, return the rotary switch to the original position to raise the pressure heads.

**Important:** Make sure that no liquid remains in the Cartridges. Residual solution in the Cartridges may cause contamination.


**Note:** Do not pull out the tube holder and cartridge holder while the pressure heads are lowered.

- (12) Pull out the tube holder and cartridge holder.

### ■Precautions for Putting DNase into Cartridges

Apply DNase directly onto the filter in each Cartridge. Be careful not to touch the filter with the pipette tip.



### <Elution>

- (13) Remove the cartridge holder and then insert it into the Elution Position (E) of the tube holder.

Cartridges are positioned above the Collection Tubes.

**Important:** Ensure that the tube holder is set in the Elution Position (E) before applying Elution Buffer.

**Note:** When moving the cartridge holder, do not slide the release lever.

**(14) Apply Elution Buffer to each Cartridge.**


Use a pipette to apply Elution Buffer to the Cartridges. In this operation, take care not to touch the filter in the Cartridge, the Cartridge itself or the cartridge holder.

**Note:** If any Elution Buffer is left on the rim of a Cartridge or on the cartridge holder, wipe it off with a soft paper tissue.


**(15) Set the tube holder and cartridge holder into the QuickGene-Mini80 again.**

**Important:** After setting it, check the following.

- The cartridge holder is placed against the stopper inside the QuickGene-Mini80.
- The WASH label is placed under the tube holder cannot be seen.

**Note:** Make sure that the WASH label is hidden below the tube holder before starting pressure application.


**(16) Repeat Steps 3 and 5.****(17) Return the rotary switch to the original position to raise the pressure heads.**

**Important:** Make sure that no liquid remains in the Cartridges.

**Note:** Do not pull out the tube holder and cartridge holder while the pressure heads are lowered.


## 2 Operation

(18) Pull out the tube holder and cartridge holder.


2

### ■Notes on Incubation

When the extraction process is followed by incubation, leave the QuickGene-Mini80 at the room temperature with the tube holder and cartridge holder set.

### ■After Extraction Operation

#### (1) Remove the cartridge holder from the tube holder and then dispose of the Cartridges.

When sliding the release lever to the right-hand end, all Cartridges will fall down.


To release the No.1 to No.4 Cartridges only, stop the release lever at the center position. In this case, to release the remaining Cartridges, slide the release lever from the center position to the right-hand end.

When you release your finger pressure from the release lever, the release lever returns to the left-hand end.

**Note:** Dispose of the Cartridges according to appropriate laws and rules.


**Dispose of infectious  
BIOHAZARD: waste according to your  
biological safety  
manuals.**


## (2) Remove the Collection Tubes.

- When using Collection Tubes included in the Isolation Kits:

Insert the Collection Tubes into the Tube Rack (optional) and then put Caps (included in the Isolation Kits) on the Collection Tubes.



**Note:** When the Tube Rack is not used, remove the Collection Tubes after putting Caps on them. The eight Caps included in the Isolation Kits can be detached by hand.

- When using commercially-available 1.5 ml micro tubes:

Put caps on 1.5 ml micro tubes and then remove them.

## (3) Remove the Waste Tubes.

**Note:** Dispose of the Waste Tubes and waste solution according to appropriate laws and rules.

BIOHAZARD:

Dispose of infectious waste according to your biological safety manuals.

# 3 Maintenance

To use QuickGene-Mini80 properly, check whether or not each part is clean and the pressure heads are in normal condition before use. If necessary, perform cleaning and/or part replacement referring to the descriptions in this chapter.

## 3.1 Cleaning QuickGene-Mini80

This section explains how to clean QuickGene-Mini80.


**BIOHAZARD:** When cleaning QuickGene-Mini80, wear appropriate gloves, a mask, and safety goggles.

**Warning:** Before starting the maintenance work, turn off the power of QuickGene-Mini80 and disconnect the power plug from the outlet. Otherwise it may cause an electric shock.

3

Parts to clean	Cleaning methods
Body (outside)	Remove dust gently with a dry, soft cloth. If the outside of QuickGene-Mini80 is heavily soiled, or if any samples that may cause infection are adhered to the outside of the QuickGene-Mini80, wipe with a soft paper tissue, soaked with 0.5% sodium hypochlorite solution or ethanol.
Body (inside)	Remove dust gently with a dry, soft cloth. If the inside of QuickGene-Mini80 is heavily soiled, or if any samples that may cause infection are adhered to the inside of the QuickGene-Mini80, wipe with a soft paper tissue, soaked with 0.5% sodium hypochlorite solution or ethanol.  <b>Note: If the tray and its circumference are heavily soiled, remove the tray before cleaning.</b>
Pressure Heads	Wipe with a soft paper tissue, soaked with ethanol.  <b>Note: If the pressure heads are heavily soiled or any abnormal conditions (damage, deformation, hardening, etc.) are found, replace them with new ones (P.28).</b>

**Caution: Keep ethanol away from fire when using. It is highly flammable.**

### 3 Maintenance

## 3.2 Cleaning the Parts

The following table describes how to clean the parts.

Parts	Cleaning methods
Cartridge Holder	Clean the parts using one of the following methods. <ul style="list-style-type: none"><li>• Clean with neutral detergent, and then wash out with distilled water. Next, soak into 3% hydrogen peroxide solution at room temperature for 10 minutes. Lastly, wash the solution out with DEPC water, and dry the parts.</li><li>• Soak in a 0.1N sodium hydroxide/0.1% EDTA solution overnight, and then wash the solution out with DEPC water, and dry the parts.</li><li>• Clean the parts with neutral detergent, then wash out with 0.5% sodium hypochlorite solution. Next, wash the solution out with DEPC water and then dry the parts.</li></ul>
Tube Holder	
Tube Adapter	

### 3


## 3.3 How to Clean/Replace the Pressure Heads

Heavily-soiled, damaged, deformed or hardened pressure heads may cause contamination or insufficient pressure application, resulting in extraction failure. Accordingly, if any abnormal conditions are found in the pressure heads, replace them promptly.

How to clean/replace the pressure heads is described below.


**Caution:** When replacing the pressure heads, wear gloves.


**(1) Turn off the power.**

Disconnect the plug of the power cable from the outlet.


**(2) Pull out the tube holder and cartridge holder.**



- (3) Place the QuickGene-Mini80 upside down and then rotate the rotary switch to eject the pressure heads.

**Note:** Before placing the QuickGene-Mini80 upside down, spread a cloth, etc. on the table to avoid damaging the body.


- (4) While holding the QuickGene-Mini80 by one hand, clean (or replace) the pressure heads. Wipe dirt off the pressure heads with a soft paper tissue, soaked with ethanol. If the pressure heads are heavily soiled or any abnormal conditions (damage, deformation, hardening, etc.) are found, replace them with new ones.

**Note:** Set the pressure heads in the direction shown in the figure, and push them completely into the body.

- (5) Rotate the rotary switch to retract the pressure heads.

## 4 Troubleshooting

This chapter explains the procedures that should be taken when you suspect device failures.

Phenomenon	Cause	Appropriate action
The lamp (LED) does not light.	<ul style="list-style-type: none"> <li>The power plug is not connected to the outlet.</li> <li>The AC adapter is not connected to QuickGene-Mini80.</li> <li>The power cable is not connected to the AC adapter.</li> </ul>	<ul style="list-style-type: none"> <li>Connect the power plug to the outlet completely.</li> <li>Connect the AC adapter to QuickGene-Mini80 completely.</li> <li>Connect the power cable to the AC adapter completely.</li> </ul>
The solution put in Cartridges does not steadily drop into the Collection or Waste Tubes.	Cartridges are incorrectly inserted into the cartridge holder.	<ul style="list-style-type: none"> <li>Insert the Cartridges completely into the cartridge holder.</li> <li>Set the Collection and Waste Tubes corresponding to the position of the Cartridge to be inserted.</li> </ul>
	Pressure heads have become deformed or hardened.	Replace the pressure heads with new ones (P.28).
	Pressure heads are incorrectly installed.	Install the pressure heads correctly (P.28).
Cartridges cannot be released/held properly.	The cartridge holder is soiled.	Clean the cartridge holder (P.28).
Pressure application does not start even after rotating the rotary switch toward the front side.	QuickGene-Mini80 may have a mechanical or electrical problem.	Contact your local representative.
QuickGene-Mini80 emits an abnormal noise.	QuickGene-Mini80 may have a mechanical or electrical problem.	Contact your local representative.
QuickGene-Mini80 emits an abnormal smell.		Turn off the power immediately by disconnecting the power cable from the outlet and then contact your local representative.

# 5 Main Specifications

The following describes the main specifications of QuickGene-Mini80.

## ■Main specifications

Items	Specifications
Model name	QuickGene-Mini80
Dimensions	276 mm (W) x 219 mm (D) x 176 mm (H)
Weight	2.8 kg (including the main body, rack set, AC adapter and power cable*) *The power cable is provided only for US and Canada.
Samples	Max. 8 samples can be processed simultaneously
Range of temperature and relative humidity	Temperature: 15 - 30°C Humidity: 30 - 80%RH (non condensing)

## ■Power supply conditions

### -AC Adapter

Items	Specifications
Rated input voltage	AC 100 to AC 240 V
Voltage fluctuation range	± 10%
Phase	Single
Power frequency	50-60 Hz
Rated input current	0.75 to 0.4 A
Output voltage	DC 12 V
Rated output current	3 A

### -Main Unit

Items	Specifications
Rated voltage	DC 12 V ± 5%
Rated current	0.5 A

## ■Environmental conditions

Items	Specifications
Usage environment	Indoor use only
Overtoltage category definition	Transient overvoltage category II
Pollution degree	Pollution Degree 2
IP	IP x 0

EMC      FCC Rule Part 15B Class A

EN61326 Class A

EN61000-3-2

EN61000-3-3

Safety     UL61010-1: Second Edition

CSA C22.2 No.1010

IEC61010: 2001

EN61010: 2001

IEC61010-2-081: 2001

CE        LV directive (73/23/EEC)

EMC directive (89/336/EEC)

## 5 Main Specifications

---

### FCC WARNING

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

This equipment has been tested and found to comply with the limits for a Class A digital device pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference, in which case, the users will be required to correct the interference at their own expense.

### Warning

This is a Class A product of Electromagnetic Interference (EMI) standard. In a domestic environment this product may cause radio interference in which case the user may be required to take adequate measures.

This Class A digital apparatus complies with Canadian ICES-003.

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



P/N: 15557927

A

\*The rack set is comprised of one cartridge holder, one tube holder and eight tube adapters.

## A.2 Product Guarantees

- The guaranty period shall be one year from the delivery of the product.
- Any part which fails to function properly under normal use, provided that all warnings and cautions in this manual are observed, during the guarantee period will be repaired by the FUJIFILM Corporation free of charge.
- If the problem is due to the causes listed below, some charge may be applied for the repair even within the guarantee period.
  - Any failures due to improper use or that are affected by other devices except for those approved by the FUJIFILM Corporation
  - Any malfunctions or damages during transportation or due to dropping of the product by a user

---

## A.3 About Services

Before you ask for device repairs, refer to “4 Troubleshooting (P.31)” to check the problems. If the problems cannot be solved, contact your local representative, or the relative customer support.

## A.4 Notes for Transporting/Shipping QuickGene-Mini80

When transporting/shipping QuickGene-Mini80 to a new location, perform the following tasks:

**Warning:** Before carrying QuickGene-Mini80, disconnect the plug from the outlet. If the power cable is damaged, this may cause a device failure, fire, injury, or electric shock.

When transporting/shipping QuickGene-Mini80, be sure to perform the following tasks:

Remove all attached parts (tube holder and cartridge holder) from QuickGene-Mini80.

Place the cushion removed at the time of unpacking onto the pressure heads and then secure the cushion and rotary switch with tape (Refer to “1.3 Remove the Packing Materials (P.4)”).

After completing the above tasks, pack QuickGene-Mini80 in its original shipping box or some other equivalent box.

---

**Nucleic Acid Isolation Device  
QuickGene-Mini80 User's Manual**

897N100125B  
Date of issuance: June 2007  
FUJIFILM Corporation  
Printed in Japan

---

- The contents of this manual may be updated without notice.
- The FUJIFILM Corporation assumes no liability for damages to third party copyrights or other rights arising from the use of any information in this manual.
- Contents of this manual are not to be reproduced without permission from the FUJIFILM Corporation.
- Manuals with missing or wrongly collated pages will be replaced free of charge.

# **FUJIFILM**

## **FUJIFILM Corporation**

FUJIFILM Corporation, Life Science Products Division  
MIDTOWN WEST, 7-3, AKASAKA 9-CHOME,  
MINATO-KU, TOKYO 107-0052, JAPAN  
Tel: +81-3-6271-2158  
Fax: +81-3-6271-3136  
E-mail: [sginfo@fujifilm.co.jp](mailto:sginfo@fujifilm.co.jp)  
<http://lifescience.fujifilm.com>

Printed on Recycled Paper

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Alle Rechte vor behalten, Copyright © FUJIFILM Corporation 2006-2007