

FLAG[®] Western Detection Kit

INSTRUCTION MANUAL

Catalog #200470

Revision A

For In Vitro Use Only

200470-12

LIMITED PRODUCT WARRANTY

This warranty limits our liability to replacement of this product. No other warranties of any kind, express or implied, including without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by Agilent. Agilent shall have no liability for any direct, indirect, consequential, or incidental damages arising out of the use, the results of use, or the inability to use this product.

ORDERING INFORMATION AND TECHNICAL SERVICES

United States and Canada

Agilent Technologies

Stratagene Products Division

11011 North Torrey Pines Road

La Jolla, CA 92037

Telephone (858) 373-6300

Order Toll Free (800) 424-5444

Technical Services (800) 894-1304

Internet techservices@agilent.com

World Wide Web www.stratagene.com

Europe

Location	Telephone	Fax	Technical Services
Austria	0800 292 499	0800 292 496	0800 292 498
Belgium	00800 7000 7000	00800 7001 7001	00800 7400 7400
	0800 15775	0800 15740	0800 15720
France	00800 7000 7000	00800 7001 7001	00800 7400 7400
	0800 919 288	0800 919 287	0800 919 289
Germany	00800 7000 7000	00800 7001 7001	00800 7400 7400
	0800 182 8232	0800 182 8231	0800 182 8234
Netherlands	00800 7000 7000	00800 7001 7001	00800 7400 7400
	0800 023 0446	+31 (0)20 312 5700	0800 023 0448
Switzerland	00800 7000 7000	00800 7001 7001	00800 7400 7400
	0800 563 080	0800 563 082	0800 563 081
United Kingdom	00800 7000 7000	00800 7001 7001	00800 7400 7400
	0800 917 3282	0800 917 3283	0800 917 3281

All Other Countries

Please contact your local distributor. A complete list of distributors is available at www.stratagene.com.

FLAG® Western Detection Kit

CONTENTS

Materials Provided	1
Storage Conditions	1
Additional Materials Required	2
Western Blotting and Antibody Probing of Electroblots.....	2
Color Development Reaction	2
Dot Blot Assay	2
Notices to Purchaser	3
Introduction	4
Western Blotting and Antibody Probing of Electroblots	4
Color Development Reaction	6
Preparation of the Color Development Solution Containing Nitroblue Tetrazolium (NBT) and 5-Bromo-4-chloro-3-indolyl Phosphate (BCIP)	6
Color Development on the Nitrocellulose Membrane.....	6
Dot Blot Assay (Optional)	7
Troubleshooting	9
Western Blotting, Antibody Probing of Electroblots, and Color Development Reaction	9
Dot Blot Assay (Optional).....	9
Preparation of Reagents	10
Stock Solutions.....	10
Standard Reagents	10
References	11
Endnotes	11
MSDS Information	11
Quick-Reference Protocol	14

FLAG® Western Detection Kit

MATERIALS PROVIDED

Materials provided ^a	Quantity
Anti-FLAG M2 antibody ^{b-d}	200 µg
Goat anti-mouse antibody, AP conjugated, immunoaffinity purified, polyclonal ^c	100 µl
Sterile, filtered fraction V bovine serum albumin (BSA) [30% (w/v) in saline]	15 ml
Nitroblue tetrazolium (NBT) ^b [75 mg/ml in 70% N,N-dimethylformamide (DMF)]	2 ml
5-Bromo-4-chloro-3-indolyl phosphate (BCIP) ^b (50 mg/ml in 100% N,N-DMF)	1.5 ml

^a Sufficient materials are provided for ten 100-mm nitrocellulose membranes.

^b Also available separately [Anti-FLAG M2 antibody, Stratagene Catalog #200471 (200 µg), #200472 (1 mg); BCIP (5-bromo-4-chloro-3-indolyl phosphate, Stratagene Catalog #300045 (1.5 ml); NBT (nitroblue tetrazolium), Stratagene Catalog #300130 (2 ml)].

^c Please see the *Certificate of Analysis* for the concentration.

^d Immediately on receipt, store the anti-FLAG M2 antibody at -20°C and the goat anti-mouse antibody at 4°C. Avoid multiple freeze-thaw cycles. Aliquot and fast freeze the primary and secondary antibodies for long-term storage.

STORAGE CONDITIONS

BSA: -20°C

Anti-FLAG® M2 Antibody: -20°C

Note For best results, store the anti-FLAG M2 antibody at -20°C.

Goat Anti-Mouse Antibody: 4°C

NBT: 4°C

BCIP: 4°C

Revision A

© Agilent Technologies, Inc. 2009.

ADDITIONAL MATERIALS REQUIRED

Western Blotting and Antibody Probing of Electoblots

Laemmli sample buffer (1×)[§]

Sodium dodecyl sulfate–polyacrylamide gel electrophoresis (SDS–PAGE) system

Nitrocellulose membranes (available from Stratagene, Catalog #420106, #420107, and #420108)

Transfer buffer (1×)[§]

Blocking solution[§] or bovine lacto transfer technique optimizer (BLOTTO)[§]

Tris-buffered saline (TBS)[§]

Color Development Reaction

Notes *Do not prepare the color development solution containing the NBT and BCIP until just prior to use in the color development reaction [see Preparation of the Color Development Solution Containing Nitroblue Tetrazolium (NBT) and 5-Bromo-4-chloro-3-indolyl Phosphate (BCIP) subsection of Color Development Reaction].*

Color development solution[§]

Tris-buffered saline (TBS)

Stop solution[§]

Whatman[®] 3MM paper

Dot Blot Assay

Strip of nitrocellulose membrane (1 × 6 cm)

Blocking solution or BLOTTO

Tris-buffered saline (TBS)

Color development solution

Whatman[®] 3MM paper

Stop solution

[§] See *Preparation of Reagents*.

NOTICES TO PURCHASER

FLAG[®] License Agreement

The enclosed DNA expression vector and/or antibody are specifically adapted for a method of producing selected protein molecules covered by one or more of the following patents owned by Sigma-Aldrich Co.: U.S. Patent Nos. (5,011,912, 4,703,004, 4,782,137 and 4,851,341; EP Patent No. 150,126 (Austria, Belgium, Switzerland, France, United Kingdom, Italy, Netherlands and Sweden); EP Patent No. 335,899 (Belgium, Switzerland, Germany, France, United Kingdom, Italy, Luxembourg and Sweden); German Patent No. P3584260.1; Canadian Patent No. 1,307,752; and Japanese Patent Nos. 1,983,150 and 2,665,359. Your payment includes a limited license under these patents to make only the following uses of these products:

A. Vector License: You may use the enclosed vector to transform cells to produce proteins containing the amino acid sequence DYKDDDDK for research purposes provided, however, such research purposes do not include binding an unlicensed antibody to any portion of this amino acid sequence nor using such proteins for the preparation of antibodies having an affinity for any portion of this amino acid sequence.

B. Antibody License: You may only use the enclosed antibody for research purposes to perform a method of producing a protein in which the protein is expressed in a host cell and purified by use of the antibody in accordance with a claim in one of the above patents in force in a country where the use actually occurs so long as: (1) you perform such method with a DNA expression vector licensed from Sigma-Aldrich Co.; and (2) you do not bind (or allow others to bind) an unlicensed antibody to any DYKDDDDK epitope of any fusion protein that is produced by use of the method.

This license does not include any rights under any other patents. You are not licensed to use the vector and/or antibody in any manner or for any purpose not recited above. As used above, the term “unlicensed antibody” means any antibody which Sigma-Aldrich Co. has not expressly licensed pursuant to Paragraph B, above. Sigma-Aldrich Co. hereby expressly retains all rights in the above listed patents not expressly licensed hereunder.

If the terms and conditions of this License Agreement are acceptable to you, then you may open the vessel(s) containing the vector and/or antibody and, through such act of opening a vessel, will have shown your acceptance to these terms and conditions.

If the terms and conditions of this License Agreement are not acceptable to you, then please return the vessel(s) unopened to Agilent for a complete refund of your payment.

For additional licensing information or to receive a copy of any of the above patents, please contact the Sigma-Aldrich Co. licensing department at telephone number 314-771-5765.

INTRODUCTION

Epitope tagging technology is a powerful method for the analysis of protein structure and function both *in vivo* and *in vitro*, and the FLAG[®] peptide (DYKDDDDK)¹ is one of the most widely used epitope tags in the field of molecular biology, cell biology, and genetics for protein labeling.^{1,2,3,4} The Stratagene FLAG[®] western detection kit includes an anti-FLAG[®] M2 antibody that can detect recombinant fusion proteins with the FLAG epitope located at either the N- or C-terminal, as well as recombinant fusion proteins with the FLAG epitope located internally. The anti-FLAG M2 antibody can detect as little as 5 ng of FLAG epitope-tagged recombinant fusion protein. The anti-FLAG M2 antibody is highly specific, with low cross-reactivity with prokaryotic and eukaryotic proteins.^{3,4}

The FLAG western detection kit utilizes the anti-FLAG M2 antibody as the primary antibody and an alkaline phosphatase (AP)-conjugated goat anti-mouse antibody as the secondary antibody. The signal is detected by a color development solution containing nitroblue tetrazolium (NBT) and 5-bromo-4-chloro-3-indolyl phosphate (BCIP). The FLAG western detection kit can be used to monitor expression of FLAG-tagged proteins in prokaryotic and eukaryotic cells, to follow FLAG epitope-tagged protein during protein purification steps, and to identify proteins following immunoprecipitation. In addition, the FLAG western detection kit provides a rapid and sensitive way to detect proteins expressed in the ESP yeast protein expression and purification system (Stratagene Catalog #251600).

WESTERN BLOTTING^{5,6} AND ANTIBODY PROBING OF ELECTROBLOTS

Note *Unless otherwise indicated, conduct all procedures at room temperature on a shaker.*

1. Boil the antigenic protein samples in 1× Laemmli sample buffer for 2 minutes and load the samples and a molecular-weight marker into separate lanes of an appropriate percentage SDS–polyacrylamide gel for analysis. Resolve the gel by SDS–PAGE.
2. Equilibrate a nitrocellulose membrane in 1× transfer buffer by soaking the membrane for 20 minutes.

Note *Never touch the nitrocellulose membranes with bare hands; always use gloves or forceps. Also, do not allow the nitrocellulose membranes to dry during any step.*

3. Electroblot the antigenic protein samples onto the nitrocellulose membrane either at 30 V and ~200 mA or at the voltage and approximate current specified by the manufacturer of the transfer apparatus for 1 hour in 1× transfer buffer.
4. Following electrophoretic transfer of the proteins onto the nitrocellulose membrane, block the membrane by incubating the nitrocellulose membrane in blocking solution, using at least 25 ml of blocking solution for an 8- × 6-cm membrane, with gentle rocking at room temperature for a minimum of 1 hour.

Notes *Blocking may also be performed overnight at 4°C.*

BLOTTO may be substituted for the blocking solution throughout the western blotting, antibody probing, and the optional dot blot assay procedures when blocking the membrane. Do not use BLOTTO to dilute the anti-FLAG M2 antibody or the goat anti-mouse antibody.

5. Thoroughly wash the nitrocellulose membrane twice with a generous amount of TBS for a minimum of 2 minutes/wash.

Notes *If nitrocellulose membranes are combined during the wash steps, less of the wash solution can be used. Separate the nitrocellulose membranes with forceps at each change to ensure adequate exposure of each membrane.*

6. For an 8- × 6-cm membrane, add 20 µg of the anti-FLAG M2 antibody to at least 10 ml of TBS and incubate the membrane in this antibody dilution at room temperature for a minimum of 1 hour.
7. Thoroughly wash the nitrocellulose membrane twice with a generous amount of TBS for a minimum of 2 minutes/wash.
8. Dilute the goat anti-mouse antibody 1:5000 in blocking solution.
9. Gently rock the nitrocellulose membrane in the goat anti-mouse antibody dilution, using at least 10 ml for an 8- × 6-cm nitrocellulose membrane, at room temperature for a minimum of 1 hour.
10. Thoroughly wash the nitrocellulose membrane three times with a generous amount of TBS for a minimum of 5 minutes/wash and then proceed to the *Color Development Reaction*.

COLOR DEVELOPMENT REACTION

Preparation of the Color Development Solution Containing Nitroblue Tetrazolium (NBT) and 5-Bromo-4-chloro-3-indolyl Phosphate (BCIP)

Prepare the color development solution containing NBT and BCIP just prior to use as outlined in the following steps.

1. Dilute 40 μ l of NBT (75 mg/ml) in 10 ml of color development solution to a final concentration of 0.3 mg/ml for a standard 8- \times 6-cm nitrocellulose membrane.
2. Add 30 μ l of BCIP (50 mg/ml) dropwise to a final concentration of 0.15 mg/ml.

Note *If a precipitate forms, warm the color development solution containing the NBT and BCIP to 37°C and pass this solution through a 0.2- μ m filter.*

3. Protect the color development solution containing the NBT and BCIP from light and leave the solution at room temperature until ready for use in step 2 of the following section.

Color Development on the Nitrocellulose Membrane

Note *Allow the color development reaction to proceed in the dark.*

1. Remove the nitrocellulose membrane from the final TBS wash (see step 11 of *Western Blotting and Antibody Probing of Electroblobs*) and blot off any excess moisture with Whatman® 3MM paper.
2. Immerse the nitrocellulose membrane in the color development solution containing NBT and BCIP. Allow the color development reaction to proceed in the dark until the positive reactions are clearly visible (usually within minutes).

Notes *Extended times for color development can result in high backgrounds. If extended times (1–5 hours) are required for the color development reaction, store the nitrocellulose membranes in the dark and check the membranes periodically to minimize the background.*

If a visible precipitate forms during the color development reaction that is not associated with specific dots or bands, replace the color development solution containing NBT and BCIP with a freshly prepared solution and allow the color development reaction to continue.

3. Remove the color development solution containing NBT and BCIP from the nitrocellulose membrane and rinse the membrane with TBS.

4. Terminate the color development reaction by immersing the nitrocellulose membrane in stop solution.
5. Following completion of the color development reaction, air dry the nitrocellulose membrane and store the membrane protected from light.

If desired, the optional *Dot Blot Assay* may be performed to quickly and simply determine if the recombinant fusion protein has been produced.

DOT BLOT ASSAY (OPTIONAL)

The dot blot assay assesses the ability of the anti-FLAG M2 antibody dilution and the goat anti-mouse antibody dilution to detect different levels of recombinant fusion protein. The dot blot assay also verifies the reaction of the NBT and BCIP with the AP-conjugated goat anti-mouse antibody.

1. Cut out a 1- × 6-cm strip of nitrocellulose membrane. Using a soft pencil, mark the strip into six equal squares measuring ~1 cm².
2. Prepare the nitrocellulose membrane with serial protein dilutions, anti-FLAG M2 antibody, and goat anti-mouse antibody as follows:
 - a. Spot 1 μl of each of the following serial protein dilutions onto separate squares of the strip: 100 ng/μl, 10 ng/μl, 5 ng/μl, and 1 ng/μl.
 - b. Spot 1 μl of the anti-FLAG M2 antibody diluted in TBS to a concentration of 10 ng/μl onto a separate square of the strip.
 - c. Spot 1 μl of the goat anti-mouse antibody diluted 1:5000 in blocking solution onto a separate square of the strip.
3. Air dry the nitrocellulose membrane for 5 minutes.
4. Immerse the strip in blocking solution for 1 hour at room temperature to block nonspecific protein-binding sites.
5. Thoroughly wash the nitrocellulose membrane two to three times for 2 minutes each in a generous amount of TBS.
6. Prepare 10 ml of the anti-FLAG M2 antibody diluted in TBS to a concentration of 2 μg/ml and incubate the strip in the anti-FLAG M2 antibody dilution at room temperature for 1 hour.

7. Thoroughly wash the nitrocellulose membrane two to three times for 2 minutes each in a generous amount of TBS.
8. Prepare 10 ml of the goat anti-mouse antibody diluted 1:5000 in blocking solution and gently rock the strip in the goat anti-mouse antibody dilution at room temperature for a minimum of 1 hour.
9. Thoroughly wash the nitrocellulose membrane three times with a generous amount of TBS for a minimum of 5 minutes/wash.
10. Prepare the color development solution containing NBT and BCIP as outlined in the following steps:
 - a. Dilute 40 μ l of NBT (75 mg/ml) in 10 ml of color development solution to a final concentration of 0.3 mg/ml.
 - b. Add 30 μ l of BCIP (50 mg/ml) dropwise to a final concentration of 0.15 mg/ml.

Note *If a precipitate forms, warm the color development solution containing the NBT and BCIP to 37°C and pass this solution through a 0.2- μ m filter.*
 - c. Protect the color development solution containing the NBT and BCIP from light and leave the solution at room temperature until ready for use in step 12 of this section.
11. Remove the strip from the final TBS wash and blot off any excess moisture with Whatman 3MM paper.
12. Immerse the nitrocellulose membrane in color development solution containing NBT and BCIP. Allow the color development reaction to proceed in the dark until the positive reactions are clearly visible. Check the membrane every minute for a positive reaction.
13. Rinse the residual precipitated dye from the strip with a generous amount of TBS.
14. Terminate the color development reaction by immersing the nitrocellulose membrane in stop solution.
15. Air dry the nitrocellulose membrane and store the membrane protected from light.

TROUBLESHOOTING

Western Blotting, Antibody Probing of Electroblobs, and Color Development Reaction

Observations	Suggestions
The signal is not visible on the nitrocellulose membrane	Ensure that the goat anti-mouse antibody or the color development solution is fresh and/or properly prepared
	To test the activity of the working dilution of the goat anti-mouse antibody, add 1 ml of the color development solution containing NBT and BCIP with 1 ml of the goat anti-mouse antibody diluted in TBS as indicated on the <i>Certificate of Analysis</i> ; if no color precipitates within 30 minutes, repeat this test with a fresh dilution of the goat anti-mouse antibody
	If color does appear, proceed to the <i>Dot Blot Assay</i> to check the protein, the anti-FLAG M2 antibody, and the goat anti-mouse antibody
	Verify the FLAG sequence in the protein
Excessive background on the nitrocellulose membrane is present	If the nitrocellulose membrane shows excessive background, check that the blocking solution or BLOTTO is properly prepared and ensure that the nitrocellulose membrane is thoroughly washed
	Check to ensure that that the anti-FLAG M2 antibody is properly diluted in TBS
	Check to ensure that that the goat anti-mouse antibody is properly diluted in the blocking solution

Dot Blot Assay (Optional)

Note *The observations in this section are addressed by performing the optional Dot Blot Assay.*

Observations	Suggestions
The protein is not visible on the dot blot strip, while the anti-FLAG M2 antibody probe is visible	Because the protein is folded, the FLAG epitope tag, which is located in the interior of the protein molecule, is not making contact with the anti-FLAG M2 antibody probe; perform western blotting noting that the denaturing characteristics of the gel should cause the protein to unfold, making the FLAG epitope tag available for detection by the anti-FLAG M2 antibody
	Verify the FLAG sequence in the protein
The anti-FLAG M2 antibody probe is not visible on the dot blot strip, while the goat anti-mouse antibody is visible	Use a freshly prepared dilution of the anti-FLAG M2 antibody probe or replace the probe completely with a new anti-FLAG M2 antibody probe
The goat anti-mouse antibody is not visible on the dot blot strip	To test the activity of the working dilution of the goat anti-mouse antibody, add 1 ml of the color development solution containing NBT and BCIP with 1 ml of the goat anti-mouse antibody diluted in TBS as indicated on the <i>Certificate of Analysis</i> ; if no color precipitates within 30 minutes, repeat this test with a fresh dilution of the goat anti-mouse antibody

PREPARATION OF REAGENTS

Stock Solutions

Use sterile concentrated stock solutions prepared with deionized, double-distilled water (ddH₂O) or the equivalent. The stock solutions necessary for all reagent preparations in the FLAG western detection kit are as follows:

Stock solutions	Sterilization method
1 M Tris-HCl (pH 7.5)	Autoclave
1 M Tris-HCl (pH 9.5)	Autoclave
1 M Tris-HCl (pH 2.9)	Autoclave
5 M NaCl	Autoclave
1 M CaCl ₂	Autoclave
1 M MgCl ₂	Filter sterilize

Standard Reagents

Note Use freshly prepared reagents. The suggested working volume is 10 ml for each 8- × 6-cm nitrocellulose membrane.

Blocking Solution 1% (w/v) BSA in TBS	BLOTTO 5 g of nonfat dry milk 1 ml of 1% (w/v) thimerosal-H ₂ O (Sigma Catalog #T 5125) Adjust the volume to 100 ml with TBS Store at 4°C for 7–10 days
TBS 50 mM Tris-HCl (pH 7.5) 150 mM NaCl	
Color Development Solution 100 mM Tris-HCl (pH 9.5) 100 mM NaCl 5 mM MgCl ₂	Laemmli Sample Buffer (2×) 250 mM Tris-HCl (pH 6.8) 8% (w/v) SDS 40% (v/v) glycerol 0.01% (w/v) bromphenol blue dye 200 mM β-mercaptoethanol
Stop Solution 20 mM Tris-HCl (pH 2.9) 1 mM CaCl ₂	Transfer Buffer (1×) 12 mM Tris-HCl (pH 8.3) 96 mM glycine 20% (v/v) methanol

REFERENCES

1. Hopp, T., Prickett, K., Price, V., Libby, R., March, C. *et al.* (1988) *BioTechnology* 6:1204-1210.
2. Chiang, C. M. and Roeder, R. G. (1993) *Pept Res* 6(2):62-4.
3. Su, X., Prestwood, A. K. and McGraw, R. A. (1992) *Biotechniques* 13(5):756-62.
4. Smit, M. J., Timmerman, H., Alewijnse, A. E., Punin, M., van den Nieuwenhof, I. *et al.* (1995) *Biochem Biophys Res Commun* 214(3):1138-45.
5. Harlow, E. and Lane, D. (Eds.). (1988). *Antibodies: A Laboratory Manual*. Cold Spring Harbor Laboratory Press, Plainview, NY.
6. Towbin, H., Staehelin, T. and Gordon, J. (1979) *Proc Natl Acad Sci U S A* 76(9):4350-4.

ENDNOTES

FLAG[®] and Anti-FLAG[®] are registered trademarks Sigma-Aldrich Co.
Whatman[®] is a registered trademark of Whatman Ltd.

MSDS INFORMATION

The Material Safety Data Sheet (MSDS) information for Stratagene products is provided on the web at <http://www.stratagene.com/MSDS/>. Simply enter the catalog number to retrieve any associated MSDS's in a print-ready format. MSDS documents are not included with product shipments.

STRATAGENE

An Agilent Technologies Division

FLAG® Western Detection Kit

Catalog #200470

QUICK-REFERENCE PROTOCOL

Western Blotting and Antibody Probing of Electoblots

- ◆ Boil the protein samples in 1 × Laemmli sample buffer for 2 minutes, load the samples and a molecular-weight marker into separate lanes of a SDS–polyacrylamide gel, and resolve by SDS–PAGE
- ◆ Equilibrate a nitrocellulose membrane in 1 × transfer buffer for 20 minutes
- ◆ Electroblot the protein samples onto the nitrocellulose membrane either at 30 V and ~200 mA or at the voltage and approximate current specified by the manufacturer of the apparatus for 1 hour in 1 × transfer buffer
- ◆ Block the nitrocellulose membrane in blocking solution with gentle rocking at room temperature for a minimum of 1 hour
- ◆ Wash the membrane twice with TBS for a minimum of 2 minutes/wash
- ◆ Add 20 µg of the anti-FLAG M2 antibody to at least 10 ml of TBS and incubate the membrane in this antibody dilution at room temperature for a minimum of 1 hour
- ◆ Wash the membrane twice with TBS for a minimum of 2 minutes/wash
- ◆ Dilute the goat anti-mouse antibody 1:5000 in blocking solution
- ◆ Gently rock the nitrocellulose membrane in the goat anti-mouse antibody dilution at room temperature for a minimum of 1 hour
- ◆ Wash the membrane three times with TBS for a minimum of 5 minutes/wash

Color Development Reaction

- ◆ Prepare the color development solution containing NBT and BCIP
- ◆ Blot off any excess moisture and immerse the nitrocellulose membrane in color development solution containing NBT and BCIP until the positive reactions are clearly visible
- ◆ Rinse the membrane with TBS
- ◆ Immerse the nitrocellulose membrane in stop solution
- ◆ Air dry the membrane and store protected from light