

# N-to-C 5.0 and C-to-N 5.0 Protein Sequencer Methods

## Technical Note

### HP 241 Protein Sequencer (N+C)

The N-to-C 5.0 and C-to-N 5.0 sequencer methods are responsible for the switching between the N-terminal and C-terminal sequencing chemistry methods on the Hewlett-Packard column-based protein sequencer.


The N-to-C 5.0 method is run when switching from the N-terminal to the C-terminal sequencing chemistry. The column method pauses the column, while the flask method prepares a coordinated injection with the on-line HPLC to place the proper HPLC solvents in the HPLC binary pumps for the C-terminal chromatography.

The C-to-N 5.0 method is run when switching from the C-terminal to the N-terminal sequencing chemistry. The column method washes the column, while the flask method, in addition to washing the flask, prepares a coordinated injection with the on-line HPLC to place the proper HPLC solvents in the HPLC binary pumps for the N-terminal chromatography. These methods are generally run with a Zitex strip in an empty RP/SAX biphasic column, but any column configuration is applicable. No column preparation is required with this method.

The methods are:

- N-to-C 5.0 (column)
- N-to-C 5.0 (flask)
- Column wash 5.0
- Flask wash 5.0

The sequence program controls all of the column and flask methods, as well as any cycle exception methods. The N-to-C 5.0 and C-to-N 5.0 sequencer methods do not use an exception method.



**Bottle Diagram**

**Drawer 1**

S2A 1R R1 2R R2

**Drawer 2**

S2/3 R2A 3R S3 R3

**Drawer 3**

L3 PTH TH R4 S5

**HPLC**

nLCA nLCB cLCA cLCB

The Routine 5.0 sequencer methods require the bottle configuration shown. The reagents and solvents are purchased from Hewlett-Packard. These methods produce less than 1mL/cycle of liquid waste.

<b>S2A</b>	Neat ethyl acetate
<b>1R</b>	Diphenylphosphorylthiocyanate (DPP-ITC) in toluene/heptane (23:27:50)
<b>R1</b>	Phenylisothiocyanate (PITC) in heptane (3:97)
<b>2R</b>	Pyridine in ethylacetate (2:98)
<b>R2</b>	Diisopropylethylamine (DIEA) in 1-propanol/water (1:6:3)
<b>S2/3</b>	Acetonitrile/toluene (23:77)
<b>R2A</b>	Ocylamine in heptane (3:97)
<b>3R</b>	Potassium trimethylsilanolate (0.1M) in methanol/t-butanol (50:50)
<b>S3</b>	Acetonitrile/toluene (15:85)
<b>R3</b>	Neat trifluoroacetic acid (TFA)
<b>L3</b>	Acetic acid in methanol/water (1:74:25)
<b>PTH</b>	Mixture of PTH amino acids (10 pmol/100 µl) in acetonitrile with DPTU as a marker
<b>TH</b>	Mixture of TH amino acids (50 pmol/100 µl) in acetonitrile
<b>R4</b>	Trifluoroacetic acid (TFA) in water (1:3)
<b>S5</b>	Phosphate buffer (pH 2.9), 0.2% ion pairing reagent

## Method Details

The steps for the N-to-C 5.0 and C-to-N 5.0 methods are described below. The steps indicate the ranges in volume that are appropriate for the various reagents/solvent deliveries. All sequencer methods are accessed by choosing the top menu item Edit/Method in the Protein Sequencer window. The sequencer column configuration used with the N-to-C 5.0 and C-to-N 5.0 methods consists of an empty reverse-phase (RP) sample column (top) mated with an empty strong anion exchange (SAX) column (bottom).

### Method (Column): Column wash 5.0

Total time: 45.5 min

All metering steps deliver the specified reagent or solvent to waste.

**Steps 2, 3:** "Flush with S2A, S3, S2/3 flushes the delivery lines and should not wet the column

#### Steps 5, 7, 8, 9, 10, 11, 12, 27, 28:

The delivery of S2A should completely fill the column, followed by drying.

#### Steps 14, 15, 25, 33, 34, 38, 39:

The delivery of S3 should completely fill the column, followed by drying.

#### Steps 17, 18, 30, 31, 41, 42:

The delivery of S2/3 should completely fill the column, followed by drying.

#### Steps 19, 21, 23:

The volume of R3 (approximately 1/2 column volume) delivered during this step should completely wet the membrane.

Step	Description	Primary Time	Temperature
1:	Wash: Dry column DOWN	30.0	55
2:	Wash: Flush with S2A, S3, S2/3	14.0	55
3:	Wash: Flush with S2A, S3, S2/3	14.0	55
4:	Wash: Meter S2A UP	60.0	55
5:	Wash: Deliver and DRY UP (closed)	30.0	55
6:	Wash: Meter S2A DOWN	60.0	55
7:	Wash: Deliver and DRY DOWN (closed)	30.0	55
8:	Wash: Flush with S2A DOWN	14.0	55
9:	Wash: Flush with S2A UP	14.0	55
10:	Wash: Flush with S2A DOWN	14.0	55
11:	Wash: Flush with S2A UP	14.0	55
12:	Wash: Flush with S2A DOWN	14.0	55
13:	Wash: Dry Column DOWN	30.0	55
14:	Wash: Flush with S3 UP	15.0	55
15:	Wash: Flush with S3 DOWN	15.0	55
16:	Wash: Dry column UP	30.0	55
17:	Wash: Flush with S2/3 DOWN	15.0	55
18:	Wash: Flush with S2/3 UP	15.0	55
19:	Wash: Flush with R3 UP	30.0	55
20:	Wash: Dry column UP	30.0	55
21:	Wash: Flush with R3 DOWN to waste	30.0	55
22:	Wash: Dry column DOWN	30.0	55
23:	!Wash: Flush with R3 DOWN to flask	30.0	55
24:	!Wash: Flush with S3	15.0	55
25:	!Wash: Dry column DOWN to flask	60.0	55
26:	!Wash: Dry column and flask	30.0	55
27:	!Wash: Flush with S2A DOWN to flask	14.0	55
28:	!Wash: Flush with S2A UP to flask	14.0	55
29:	!Wash: Dry column and flask	30.0	55
30:	!Wash: Flush with S2/3 DOWN to flask	15.0	55
31:	!Wash: Flush with S2/3 UP to flask	15.0	55
32:	!Wash: Dry column and flask	30.0	55
33:	!Wash: Flush with S3 UP to flask	15.0	55
34:	!Wash: Flush with S3 DOWN to flask	15.0	55
35:	!Wash: Dry column and flask	30.0	55
36:	!Wash: Flush with L3 and R4	20.0	55
37:	!Wash: Flush with L3 and R4	20.0	55
38:	!Wash: Flush with S3 UP	15.0	55
39:	!Wash: Flush with S3 DOWN	15.0	55
40:	Wash: Dry column DOWN	30.0	55
41:	Wash: Flush with S2/3 DOWN	15.0	55
42:	Wash: Flush with S2/3 UP	15.0	55
43:	Wash: Dry column DOWN	30.0	55
44:	Wash: Flush with S2/3 and S2A	15.0	55
45:	Wash: Dry column UP	60.0	55

! indicates system step

---

**Method (Column): N-to-C 5.0**

Total time: 20.7 min

The column pauses during this method while the N-to-C 5.0 flask method runs

Step	Description	Primary Time	Temperature
1:	wait	20.0	30
2:	!wait system	20.0	30
3:	wait	600.0	30
4:	wait	600.0	30
! indicates system step			

**Method (Flask): Flask wash 5.0**

Total time: 67.1 min

**Step 3:** After Step 3, the volume of R4 delivered to the flask should be approximately 70  $\mu\text{l}$  (+/- 5  $\mu\text{l}$ ).

**Step 5:** By the end of Step 5, the flask will have been dry for 100-150 seconds.

**Step 6:** By the end of Step 6, 70  $\mu\text{l}$  (+/- 5  $\mu\text{l}$ ) of L3 will have been delivered to the flask.

**Step 7:** By the end of Step 7, the flask will have been dry for 75-100 seconds.

**Step 8:** The delivery of L3 may not be visible.

**Step 9:** After Step 9, the volume of L3 and S5 in the flask should be 75-80  $\mu\text{l}$ . If adjustment is needed, it should be done by adjusting the metering time of S5, leaving L3 unchanged.

Step	Description	Primary Time	Temperature
1:	Wash: Deliver R4 rinse	80.0	50
2:	Wash: Empty flask	60.0	50
3:	Convert: Deliver R4	29.0	70
4:	Convert: React	400.0	70
5:	Convert: Evaporate dry	430.0	70
6:	Convert: Flush with L3	15.0	70
7:	Convert: Evaporate dry L3	180.0	70
8:	Convert: Deliver L3 dripwise	2.0	60
9:	Convert: Deliver S5 dripwise	12.0	60
10:	Convert: Mix and equilibrate	15.0	60
11:	Convert: Solubilize	20.0	60
12:	Convert: Fill loop and inject	1.8	60
13:	Convert: Empty flask	30.0	60
14:	Convert: Flush with L3	85.0	45
15:	Convert: Mix	20.0	45
16:	Convert: Empty flask	30.0	45
17:	Convert: Deliver R4 wash	40.0	45
18:	Convert: Empty flask	30.0	45
19:	Convert: Deliver R4 wash	40.0	45
20:	Convert: Empty flask	30.0	45
21:	Convert: Dry waste line	60.0	45
22:	Convert: Dry vent line	60.0	45

**TEST Delivery protocol:**

*To check the delivery volume of L3 and S5, copy Steps 8, 9, 10 and 11 to the clipboard and run from the clipboard. At the end of Step 11, remove the flask from the heating chamber. Carefully unscrew the flask cap and measure the delivered volume using a syringe. Be careful that some of the delivered volume is not retained on the tube sides or top of the flask.*

**Step 12:** Before injection the injector loop should be filled, leaving the solvent visible in both the inlet and outlet lines of the injector loop. The time for Step 12 may have to be adjusted in order to optimize the delivery.

---

**Method (Flask): N-to-C 5.0**

Total time: 11.2 min

**Step 5:** After Step 5, the volume of S5 in the flask should be 75-80 µl.

***TEST delivery protocol:***

*To check the delivery volume of S5, copy Step 5 to the clipboard and run from the clipboard. At the end of Step 5, remove the flask from the heating chamber. Carefully unscrew the flask cap and measure the delivered volume using a syringe. Be careful that some of the delivered volume is not retained on the tube sides or top of the flask*

**Step 8:** Before the injection the injector loop should be filled, leaving the solvent visible in both the inlet and outlet lines of the injector loop. The time for Step 8 may have to be adjusted in order to optimize the delivery.

Step	Description	Primary Time	Temperature
1:	Wash: Empty flask	30.0	50
2:	Wash: Flush with L3	85.0	50
3:	Wash: Mix	20.0	50
4:	Wash: Empty flask	30.0	50
5:	Wash: Deliver S5 dripwise	12.5	50
6:	Wash: Mix and equilibrate	15.0	50
7:	Wash: Solubilize	20.0	50
8:	Wash: Fill loop and inject	2.0	50
9:	Wash: Empty flask	20.0	50
10:	Wash: Flush with L3	85.0	50
11:	Wash: Mix	20.0	50
12:	Wash: Empty flask	30.0	50
13:	Wash: Dry waste line	60.0	50
14:	Wash: Dry vent line	60.0	50

Hewlett-Packard shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance or use of this material.

Information, descriptions and specifications in this publication are subject to change without notice.

Copyright © 1998  
Hewlett-Packard Company  
All rights reserved. Reproduction and adaptation is prohibited.

Printed in the U.S.A. July 1998  
5968-0127E