

Routine 3.1 PVDF Sequencer Methods

Technical Note

**HP G1005A
N-terminal Protein
Sequencing System**

TN 96-2

The N-terminal Routine 3.1 PVDF sequencer methods implement the Edman degradation chemistry (coupling, cleavage, and conversion) on Hewlett-Packard's column-based protein sequencer. These methods are appropriate for N-terminal amino acid sequence analysis of most proteins electrophoretically transferred to PVDF membrane in the nanomole to low picomole range. The methods are:


- Routine column
- Routine converter
- Routine Cycle 1 (column)
- Routine PTH Std (converter)
- Column Prep (column)

The Routine column and Routine converter methods typically control all sequence cycles with the exception of Cycle 1 for each sample. These methods are initiated by the Routine sequence program.

The method Routine Cycle 1 is used in the Routine sequence program as the column method for cycle 1. It provides an initial sample cleanup and performs a double coupling on the N-terminal amino acid. There is no specialized cycle 1 method for the converter.

The Routine PTH Std methods deliver 10 pmol/100 µl PTH-amino acid standard solution to the on-line HPLC allowing quantitation and identification of sequenced amino acids. You may choose to run a PTH amino acid standard using Routine PTH Std when scheduling your sample.

The Column Prep method is used to wash the sequencer column prior to sample loading. This is especially important for low-level (<10 pmol) samples.



The Routine 3.1 PVDF sequencer method requires the bottle configuration shown. These reagents and solvents are purchased from Hewlett-Packard (HP reagent/solvent kit G2046A). These methods produce less than 1mL/cycle of liquid waste.

- L1 neat acetonitrile
- S3A trifluoroacetic acid (TFA) in water (1:499)
- R1 phenylisothiocyanate (PITC) in heptane (3:97)
- R2A octylamine in heptane (3:97)
- R2 diisopropylethylamine (DIEA) in 1-propanol/water (1:3:6)
- S2/3 acetonitrile/toluene (23:77)
- S3 acetonitrile/toluene (15:85)
- S2A neat ethyl acetate
- R3 neat trifluoroacetic acid (TFA)
- L2 neat methanol
- Std mixture of PTH amino acids (10 pmol/100 µl) in acetonitrile with DPTU as marker in water (1:3)
- S4 acetonitrile/water (10:90)
- R4 trifluoroacetic acid (TFA) in water (1:3)

The Steps in the Routine methods are described below. All sequencer methods are accessed by choosing the top menu item Edit/Method in the Protein Sequencer window. The steps indicate the ranges in volume that are appropriate for the various reagent/solvent deliveries. The sequencer column configuration used with Routine 3.1 PVDF method consists of a membrane-compatible empty column on top mated to a strong anion exchange (SAX) sample column on bottom.

Step 2: The delivery of R2A should completely fill the column (both halves).

Steps 5, 8: The delivery of R1 should completely fill the column (both halves).

Step 11: The volume of R2 delivered during Step 11 should completely saturate the PVDF membrane(s).

Step 12: The membrane(s) should still be moist at the end of this step following removal of excess R2 and remain moist during the coupling reaction.

Step 15: The membrane(s) should be completely dry after this step.

Steps 18, 20, 22, 24: These solvent deliveries completely fill both the membrane-compatible and the SAX column, followed by drying.

Step 25: The membrane(s) should be dry after this step.

Step 28: No more than the upper frit of the membrane-compatible column should be wet.

Step 31: The volume of S4 delivered to the converter should be 50 μl ($\pm 5 \mu\text{l}$).

Steps 33, 37, 41: The PVDF column and the SAX column should be completely filled with solvent during this step. Solvent should not enter the RV6 valve before the flow is directed to the converter.

Steps 34, 38: Membrane(s) should be saturated with solvent during these steps.

Step 44: Membrane(s) should be dry after this step.

[illegible]

Method (Column): Routine 3.1

PVDF Cycle 1

Total time: 61.2 min (3670.1 sec)

Steps 2, 4, 6: Cycle 1 begins with S3A washes to remove contaminants from the PVDF sample. During these steps, the columns should be completely filled with liquid.

Step 9: The membrane(s) should be completely dry at the end of this step.

Step 11: The delivery of R2A should completely fill the column (both halves).

Steps 14, 17, 25: The delivery of R1 should completely fill the column (both halves).

Steps 20, 28: The volume of R2 delivered during Steps 20 and 28 should completely saturate the PVDF membrane(s).

Step 21, 29: The membrane(s) should still be moist at the end of these steps following removal of excess R2 and remain moist during the coupling reaction.

Steps 23, 32: The membrane(s) should be completely dry after these steps.

Steps 35, 37, 39, 41: These solvent deliveries completely fill both the PVDF and the SAX columns, followed by drying.

Step 42: The membrane(s) should be dry after this step.

Step 45: No more than the upper frit of the membrane-compatible column should be wet.

Step 48: The volume of S4 delivered to the converter should be 50 µl (±5 µl).

Steps 50, 54, 58: The PVDF column and the SAX column should be completely filled with solvent during this step. Solvent should not enter the RV6 valve before the flow is directed to the converter.

Steps 51, 55: Membrane(s) should be saturated with solvent during these steps.

Step 61: Membrane(s) should be dry after this step.

Method (Column): Routine 3.1		PVDF Cycle 1	
Total time: 61.2 min (3670.1 sec)			
Step	Description	Start Time (min)	End Time (min)
2	S3A wash	0.00	1.00
4	S3A wash	1.00	2.00
6	S3A wash	2.00	3.00
9	Dry membrane	3.00	4.00
11	R2A delivery	4.00	5.00
14	R1 delivery	5.00	6.00
17	R1 delivery	6.00	7.00
20	R2 delivery	7.00	8.00
21	R2 delivery	8.00	9.00
23	Dry membrane	9.00	10.00
25	R1 delivery	10.00	11.00
28	R2 delivery	11.00	12.00
29	R2 delivery	12.00	13.00
32	Dry membrane	13.00	14.00
35	Solvent delivery	14.00	15.00
37	Solvent delivery	15.00	16.00
39	Solvent delivery	16.00	17.00
41	Solvent delivery	17.00	18.00
42	Dry membrane	18.00	19.00
45	Upper frit wet	19.00	20.00
48	S4 delivery	20.00	21.00
50	Solvent delivery	21.00	22.00
54	Solvent delivery	22.00	23.00
58	Solvent delivery	23.00	24.00
51	Solvent delivery	24.00	25.00
55	Solvent delivery	25.00	26.00
61	Dry membrane	26.00	27.00

PVDF

Total time: 33.1min (1988.1sec)

Step 1: During Step 1, the liquid in the converter partially evaporates but must not dry down completely (10-20 μ l remains).

Step 2: During Step 2, the volume of R4 delivered to the converter should be approximately 70 μ l ($\pm$ 5 μ l).

TEST Delivery protocol:

To check the delivery volume of R4, copy Step 2 to the clipboard and run the clipboard. At the end of Step 2, remove the conversion flask from the heating chamber. Carefully unscrew the flask and measure the delivered volume.

Step 4: By the end of Step 4, the converter will have been dry for 100-200 secs.

Step 5: By the end of Step 5, 100-120 μ l of L2 will have been delivered to the converter.

Step 6: By the end of Step 6, the converter will have been dry for 75-100 secs.

Step 7: After Step 7, the volume of S4 delivered to the converter should be 75-80 µl.

TEST Delivery protocol:

To check the delivery volume of S4, copy Steps 7, 8, and 9 to the clipboard and run from the clipboard. At the end of Step 9, remove the conversion flask from the heating chamber. Carefully unscrew the flask and measure the delivered volume.

Step 10: Before injection the injector loop should be filled, leaving solvent visible in both the inlet and outlet lines of the injector loop.

[illegible]

Method (Converter): Routine

3.1/3.5 PTH Std

Total time: 34.3 min (2060.0 sec)

This Routine PTH-standard method is appropriate with any of the sequence programs: 3.1, 3.1 PVDF, or 3.5. Ten (10) pmols of PTH-amino acid standard are delivered to the conversion flask in this method.

Step 7: Step 7 should add 100 µl (± 2 µl) of PTH-amino acid standard to the dry converter.

TEST Delivery protocol:

To check the delivery volume of PTH-standard: first remove the conversion flask from the heating chamber and allow to cool.

Confirm that the PTH delivery time in the Test method is the same as the Routine PTH method delivery time. Run from the Clipboard the “Test 3.1/3.5 PTH Delivery” method. At Step 2 of the test method, carefully unscrew the flask and measure the delivered volume.

Step 8: The solvent is completely evaporated by end of the 180 sec dry.

Step 9: After Step 9, the volume of S4 delivered to the converter should be 75-80 µl.

Step 12: Before injection the injector loop should be filled, leaving solvent visible in both the inlet and outlet lines of the injector loop.

Method (Converter): Routine		3.1/3.5 PTH Std	
Total time: 34.3 min (2060.0 sec)			
Step 7:	Step 7 should add 100 µl (± 2 µl) of PTH-amino acid standard to the dry converter.	TEST Delivery protocol:	
To check the delivery volume of PTH-standard: first remove the conversion flask from the heating chamber and allow to cool.		Confirm that the PTH delivery time in the Test method is the same as the Routine PTH method delivery time. Run from the Clipboard the “Test 3.1/3.5 PTH Delivery” method. At Step 2 of the test method, carefully unscrew the flask and measure the delivered volume.	
Step 8:		The solvent is completely evaporated by end of the 180 sec dry.	
Step 9:		After Step 9, the volume of S4 delivered to the converter should be 75-80 µl.	
Step 12:		Before injection the injector loop should be filled, leaving solvent visible in both the inlet and outlet lines of the injector loop.	

Method (Column): Column Prep
3.1 PVDF

Total time: 20.5 min (1232.3 sec)

Steps 3, 14, 20, 23, 26, 28, 30, 32:
These solvents deliveries completely fill the column, followed by drying.

Steps 6, 17: These reagents deliveries should completely fill the column (both halves).

Steps 9: The delivery of R2 in Step 9 should completely fill the lower SAX column.

SAX Column		SAX Column	
Step	Description	Step	Description
1	Initial solvent delivery	1	Initial solvent delivery
2	Initial solvent delivery	2	Initial solvent delivery
3	Initial solvent delivery	3	Initial solvent delivery
4	Initial solvent delivery	4	Initial solvent delivery
5	Initial solvent delivery	5	Initial solvent delivery
6	Initial solvent delivery	6	Initial solvent delivery
7	Initial solvent delivery	7	Initial solvent delivery
8	Initial solvent delivery	8	Initial solvent delivery
9	Initial solvent delivery	9	Initial solvent delivery
10	Initial solvent delivery	10	Initial solvent delivery
11	Initial solvent delivery	11	Initial solvent delivery
12	Initial solvent delivery	12	Initial solvent delivery
13	Initial solvent delivery	13	Initial solvent delivery
14	Initial solvent delivery	14	Initial solvent delivery
15	Initial solvent delivery	15	Initial solvent delivery
16	Initial solvent delivery	16	Initial solvent delivery
17	Initial solvent delivery	17	Initial solvent delivery
18	Initial solvent delivery	18	Initial solvent delivery
19	Initial solvent delivery	19	Initial solvent delivery
20	Initial solvent delivery	20	Initial solvent delivery
21	Initial solvent delivery	21	Initial solvent delivery
22	Initial solvent delivery	22	Initial solvent delivery
23	Initial solvent delivery	23	Initial solvent delivery
24	Initial solvent delivery	24	Initial solvent delivery
25	Initial solvent delivery	25	Initial solvent delivery
26	Initial solvent delivery	26	Initial solvent delivery
27	Initial solvent delivery	27	Initial solvent delivery
28	Initial solvent delivery	28	Initial solvent delivery
29	Initial solvent delivery	29	Initial solvent delivery
30	Initial solvent delivery	30	Initial solvent delivery
31	Initial solvent delivery	31	Initial solvent delivery
32	Initial solvent delivery	32	Initial solvent delivery
33	Initial solvent delivery	33	Initial solvent delivery
34	Initial solvent delivery	34	Initial solvent delivery
35	Initial solvent delivery	35	Initial solvent delivery
36	Initial solvent delivery	36	Initial solvent delivery
37	Initial solvent delivery	37	Initial solvent delivery
38	Initial solvent delivery	38	Initial solvent delivery
39	Initial solvent delivery	39	Initial solvent delivery
40	Initial solvent delivery	40	Initial solvent delivery
41	Initial solvent delivery	41	Initial solvent delivery
42	Initial solvent delivery	42	Initial solvent delivery
43	Initial solvent delivery	43	Initial solvent delivery
44	Initial solvent delivery	44	Initial solvent delivery
45	Initial solvent delivery	45	Initial solvent delivery
46	Initial solvent delivery	46	Initial solvent delivery
47	Initial solvent delivery	47	Initial solvent delivery
48	Initial solvent delivery	48	Initial solvent delivery
49	Initial solvent delivery	49	Initial solvent delivery
50	Initial solvent delivery	50	Initial solvent delivery
51	Initial solvent delivery	51	Initial solvent delivery
52	Initial solvent delivery	52	Initial solvent delivery
53	Initial solvent delivery	53	Initial solvent delivery
54	Initial solvent delivery	54	Initial solvent delivery
55	Initial solvent delivery	55	Initial solvent delivery
56	Initial solvent delivery	56	Initial solvent delivery
57	Initial solvent delivery	57	Initial solvent delivery
58	Initial solvent delivery	58	Initial solvent delivery
59	Initial solvent delivery	59	Initial solvent delivery
60	Initial solvent delivery	60	Initial solvent delivery
61	Initial solvent delivery	61	Initial solvent delivery
62	Initial solvent delivery	62	Initial solvent delivery
63	Initial solvent delivery	63	Initial solvent delivery
64	Initial solvent delivery	64	Initial solvent delivery
65	Initial solvent delivery	65	Initial solvent delivery
66	Initial solvent delivery	66	Initial solvent delivery
67	Initial solvent delivery	67	Initial solvent delivery
68	Initial solvent delivery	68	Initial solvent delivery
69	Initial solvent delivery	69	Initial solvent delivery
70	Initial solvent delivery	70	Initial solvent delivery
71	Initial solvent delivery	71	Initial solvent delivery
72	Initial solvent delivery	72	Initial solvent delivery
73	Initial solvent delivery	73	Initial solvent delivery
74	Initial solvent delivery	74	Initial solvent delivery
75	Initial solvent delivery	75	Initial solvent delivery
76	Initial solvent delivery	76	Initial solvent delivery
77	Initial solvent delivery	77	Initial solvent delivery
78	Initial solvent delivery	78	Initial solvent delivery
79	Initial solvent delivery	79	Initial solvent delivery
80	Initial solvent delivery	80	Initial solvent delivery
81	Initial solvent delivery	81	Initial solvent delivery
82	Initial solvent delivery	82	Initial solvent delivery
83	Initial solvent delivery	83	Initial solvent delivery
84	Initial solvent delivery	84	Initial solvent delivery
85	Initial solvent delivery	85	Initial solvent delivery
86	Initial solvent delivery	86	Initial solvent delivery
87	Initial solvent delivery	87	Initial solvent delivery
88	Initial solvent delivery	88	Initial solvent delivery
89	Initial solvent delivery	89	Initial solvent delivery
90	Initial solvent delivery	90	Initial solvent delivery
91	Initial solvent delivery	91	Initial solvent delivery
92	Initial solvent delivery	92	Initial solvent delivery
93	Initial solvent delivery	93	Initial solvent delivery
94	Initial solvent delivery	94	Initial solvent delivery
95	Initial solvent delivery	95	Initial solvent delivery
96	Initial solvent delivery	96	Initial solvent delivery
97	Initial solvent delivery	97	Initial solvent delivery
98	Initial solvent delivery	98	Initial solvent delivery
99	Initial solvent delivery	99	Initial solvent delivery
100	Initial solvent delivery	100	Initial solvent delivery

Hewlett-Packard shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance or use of this material.

Information, descriptions and specifications in this publication are subject to change without notice.

**Copyright © 1996
Hewlett-Packard Company
All rights reserved. Reproduction and adaption is prohibited.**

**Printed in U.S.A. July 1996
(23) 5965-3466E**