NELLCOR PURITAN BENNETT

OPERATOR'S MANUAL

N-180 Pulse Oximeter

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N-180 Pulsioximeter

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Pulsioxímetro N-180

MANUALE PER L'USO

Saturimetro N-180

To contact Nellcor Puritan Bennett's representative: In the United States, call 800.NELLCOR or 510-463-4000; outside of the United States, call your local Nellcor Puritan Bennett representative.

Caution: Federal law (U.S.) restricts this device to sale by or on the order of a physician.

NELLCOR PURITAN BENNETT

OPERATOR'S MANUAL

N-180 Pulse Oximeter

To contact Nellcor Puritan Bennett's representative: In the United States, call 800.NELLCOR or 510-463-4000; outside of the United States, call your local Nellcor Puritan Bennett representative.

Caution: Federal law (U.S.) restricts this device to sale by or on the order of a physician.

Corporate Headquarters

Nelcor Puritan Bennett Inc.
4280 Hacienda Drive
Pleasanton, California 94588
USA
Tel. 510 463-4000
1-800-NELCOR

European Office

Nelcor Puritan Bennett Europe BV
Eembakenvatering 1
5231 DD 's-Hertogenbosch
The Netherlands
Tel. +31.73.6485200

Asia/Pacific Office

Nelcor Puritan Bennett HK Ltd.
Room 1602 Evergo House
88 Gloucester Road
Wanchai
Hong Kong
Tel. +852.2529.0333

Regional/Local Offices

Nelcor Puritan Bennett UK Ltd.
10 Tatham Business Centre
London Road
Bicester, Oxfordshire OX6 0JX
United Kingdom
Tel. +44.1869.322700

Nelcor Puritan Bennett Belgium NV/SA
Interleuvenlaan 62/8, Zone 2
B-3001 Heverlee
Belgium
Tel. +32.16.400467

Nelcor Puritan Bennett France Sarl
21 rue Albert Calmette
78355 Jouy-en-Josas Cedex
France
Tel. +33.1.34.63.06.00

Nelcor Puritan Bennett Germany GmbH
Black-&-Decker-Strasse 28
65510 Idstein
Germany
Tel. +49.6126.5330

Nelcor Puritan Bennett Italia Srl
Via dei Tulipani, 3
20090 Pieve Emanuele (MI)
Italy
Tel. +39.2.90786404

To obtain information about a warranty, if any, for this product, contact Nelcor Puritan Bennett Technical Services or your local Nelcor Puritan Bennett representative.

Durasensor, Dura-Y, Oxihand, OxiClig, Oxisensor and the Nelcor Puritan Bennett knob configuration are trademarks of Nelcor Puritan Bennett Inc.

Covered by one or more of the following U.S. Patents and foreign equivalents: 4,621,643; 4,653,496; 4,700,708; 4,770,179; 4,863,254.

CONTENTS

Figures
Tables

Safety Information	1
Warnings	1
General	1
Alarms	2
Electrical	2
Sensors	2
Measurements	3
Symbols	5
Quick Guide to Operation	7
Basic Operation	7
Alarm Functions	7
Check Alarm Limits	7
Adjust Alarm Limits	8
Adjust Alarm Volume	8
Silence Alarm Temporarily	8
Adjust Alarm Silence Period	8
Disable Audio Alarm	8
Features	9
Overview	9
Features	9
Principles of Operation	11
Oximetry Overview	11
Automatic Calibration	11
Functional versus Fractional Saturation	12
Measured versus Calculated Saturation	12
Setup	15
Unpacking and Inspection	15
List of Components	15
Testing	15
Components	15
Readout	15
Controls	16
Front Panel (100–120 VAC)	16
Front Panel (220–240 V~)	16
Rear Panel	18
Rear-Panel Switches	19

Contents

Guide to Operation	21
Basic Operation	21
Pulse Oximetry Subsystem Features	22
Pulse Tone	22
Alarm Functions	23
Operating Modes	25
Battery Operation	26
Serial Data Output	27
Nelcor Puritan Bennett Sensors	29
Selecting a Nelcor Puritan Bennett Sensor	29
Cleaning and Reuse	30
Performance Considerations	30
Connecting to Other Instruments	31
RS-232 Communication Protocol	31
Serial Data Connector Pin Assignments	31
Communications Formats	32
Setting Baud Rate	34
Maintenance	35
Service	35
Cleaning	35
Battery Testing	35
Determining Software Version	36
Technical Assistance	36
Returning the N-180	36
Troubleshooting	37
Status Messages	37
Troubleshooting Guide	38
Specifications	43
Performance	43
Environmental Requirements	44
Electrical Characteristics	44
Physical Characteristics	45

FIGURES

1	Oxyhemoglobin Dissociation Curve	13
2a	N-180 Front Panel, 100–120 VAC Version	16
2b	N-180 Front Panel, 220–240 V~ Version	16
3	N-180 Rear Panel	18
4	Rear-Panel Switches	19

TABLES

1	Default Alarm Settings	24
2	Selected <i>Nelcor Puritan Bennett</i> Sensors	29
3	Serial Communications Connector Pin Assignments	32
4	Output Format Switch Settings	33
5	Baud Rate Switch Settings.....	34

SAFETY INFORMATION

Warnings

WARNINGS

General

DANGER! Explosion hazard. Do not use in the presence of flammable anesthetics.

The N-180 is to be operated by qualified personnel only. Before use, carefully read this manual, accessory directions for use, all precautionary information, and specifications. The user must check that the equipment functions safely and see that it is in proper working condition before using it.

The N-180 is intended only as an adjunct in patient assessment. It must be used in conjunction with clinical signs and symptoms.

Do not use *Nellcor Puritan Bennett* pulse oximeters during magnetic resonance image (MRI) scanning. Adverse reactions include: potential burns to patients as a result of contact with attachments heated by the MRI RF (radio frequency) pulse; potential degradation of the MR image; and, potential reduced accuracy of oximeter measurements. Always remove oximetry devices and attachments from the MRI imaging environment before scanning a patient.

Alarms

Do not silence the audible alarm or decrease its volume if patient safety could be compromised.

Check the audible alarm silence duration before temporarily silencing the audible alarm.

Each time the monitor is used, check alarm limits to ensure that they are appropriate for the patient being monitored.

Safety Information

Electrical

Electric shock hazard. Cover to be removed only by qualified service personnel. There are no user-serviceable parts inside.

Note: Do not connect to an electrical outlet controlled by a wall switch because power to the monitor could be inadvertently turned off.

Sensors

Before use, carefully read the sensor directions for use.

Use only *Nellcor Puritan Bennett* oxygen transducers (sensors). Use of other oxygen transducers may cause improper oximeter performance.

Tissue damage can be caused by incorrect application or use of a sensor (for example, applying the sensor too tightly or applying supplemental tape). Inspect the sensor site routinely to ensure skin integrity and correct positioning and adhesion of the sensor. If skin integrity changes, move the sensor to another site.

Inspect the sensor and cable for fraying, cracking, breakage, or other damage. If defects are noted, do not use the sensor or sensor cable. Do not immerse sensor completely in water, solvents, or cleaning solutions (because the connector is not waterproof). Do not sterilize by irradiation, steam, or ethylene oxide.

Measurements

Loss of pulse signal can occur for the following reasons:

- The sensor is too tight.
- There is excessive illumination, such as from sunlight or a surgical or bilirubin lamp.
- The sensor is placed on an extremity with a blood pressure cuff, arterial catheter, or intravascular line.
- The patient is in shock, has hypotension, severe vasoconstriction or anemia, hypothermia, arterial occlusion proximal to the sensor, or cardiac arrest.

Inaccurate measurements may be caused by:

- Incorrect application or use of a sensor.
- Significant levels of dysfunctional hemoglobins, such as carboxyhemoglobin or methemoglobin.
- Significant levels of indocyanine green, methylene blue, or other intravascular dyes.
- Exposure to excessive illumination, such as surgical lamps (especially ones with a xenon light source), bilirubin lamps, fluorescent lights, infrared heating lamps, or direct sunlight. Exposure to excessive illumination can be corrected by covering the sensor with a dark or opaque material.
- Excessive patient movement.
- Venous pulsations.
- High-frequency electrosurgical interference and defibrillators.
- Placement of the sensor on an extremity that has a blood pressure cuff, arterial catheter, or intravascular line.

SYMBOLS


Attention: refer to manual

Type BF equipment; patient electrically isolated

Fuse replacement symbol (110-120 VAC)

Fuse replacement symbol (220-240 V~)

Equipotential ground

Audible alarm off

Power on/off switch

ON/STDBY switch

SpO₂%

Percentage OXYGEN SATURATION

/min

PULSE RATE, measured in beats per minute (bpm)

PULSE SEARCH indicator

Symbols

BATT IN USE

LOW BATT indicator

HIGH SAT

LOW SAT

HIGH RATE

LOW RATE

NEC

Neonatal

QUICK GUIDE TO OPERATION

Basic Operation
Alarm Functions

BASIC OPERATION

1. Select an appropriate *Nellcor Puritan Bennett* sensor and apply it to the patient, following the sensor directions for use (see section *Nellcor Puritan Bennett Sensors*). Connect the sensor to the N-180. If needed, use a *Nellcor Puritan Bennett* EC-4 or EC-3 sensor extension cable.
2. Turn on the monitor: set the N-180 rear-panel mains power switch to the ON position (the switch indicator lights up), and set the front-panel ON/STDBY switch to the ON position (see section *Guide to Operations*).

Alternatively, operate the N-180 on internal battery.

3. Check alarm limits. If necessary, adjust them to suit the patient's needs (see section *Guide to Operations*).
4. Adjust the HIGH SAT, LOW SAT, HIGH RATE, and LOW RATE alarm values by pressing the corresponding pushbuttons and adjusting the control knob (within 3 seconds).

Adjust the alarm silence period by pressing and holding the alarm disable button while adjusting the control knob.

ALARM FUNCTIONS (See Section *Guide to Operations*)

Check Alarm Limits

Press the appropriate alarm button (HIGH SAT, LOW SAT, HIGH RATE, or LOW RATE).

Quick Guide to Operation

Adjust Alarm Limits

Press the appropriate alarm button (HIGH SAT, LOW SAT, HIGH RATE, or LOW RATE), then turn the control knob until the desired setting appears.

Adjust Alarm Volume

Simultaneously press the LOW SAT and HIGH SAT buttons, turning on the alarm tone. Turn control knob until the desired setting appears in the OXYGEN SATURATION numeric display. Turning the knob varies the audible alarm tone volume.

Silence Alarm Temporarily

Press the Audible Alarm Off button. The Audible Alarm Off indicator lights steadily during the alarm-silence period.

Adjust Alarm Silence Period

Press and hold the Audible Alarm Off button and turn the control knob until the desired setting (30–120 seconds) appears in the OXYGEN SATURATION numeric display.

Disable Audio Alarm

Press and hold the Audible Alarm Off button while turning the control knob clockwise until “OFF” appears in the OXYGEN SATURATION display and the Audible Alarm Off indicator flashes.

Note: This condition will remain in effect until you re-enable the audio alarm.

Note: On 220–240 voltage N-180 monitors, when the audible alarm is disabled, a single tone sounds every 3 minutes.

To re-enable the audio alarm, press the Audible Alarm Off button again; re-enabling the audio alarm will stop the indicator flashing.

Turning the N-180 off and then back on automatically resets the alarms to default settings.

FEATURES

Overview
Features

OVERVIEW

The *Nellcor Puritan Bennett* N-180 pulse oximeter measures functional oxygen saturation of arterial hemoglobin (SpO₂), and pulse rate. The N-180 monitors SpO₂ and pulse rate continuously and noninvasively, with measurements updated at each pulse beat. An internal AC power supply provides isolated power for operating the monitor and charging its internal batteries. In addition, the N-180 provides a digital output for external devices.

FEATURES

The N-180 is ready for immediate use upon power-up, without the need for operator calibration or configuration. It offers:

- Automatic self-test and error messages.
- Automatic oximetry calibration.
- Visible and audible oximetry displays.
- An early warning system that provides an audible indicator for both SpO₂ and pulse rate; a tone signals each pulse and its pitch varies with changes in SpO₂.
- Operator-configurable audible oximetry alarms, with default alarm limits preset for adults or neonates.
- Visual alarms with pulse rate and oxygen saturation displays that change from green to red during an alarm condition.
- Battery operation for emergency or portable use.

The N-180 provides the operator with the capability to tailor the system for specific clinical applications. Capabilities include:

- Audible alarms that can be silenced temporarily or disabled; the alarm tone has adjustable volume.
-

Features

- Three oximetry operating modes that change the measurement averaging time to suit different clinical applications.

PRINCIPLES OF OPERATION

Oximetry Overview

OXIMETRY OVERVIEW

Pulse oximetry is based on two principles: that oxyhemoglobin and deoxyhemoglobin differ in their absorption of red and infrared light (i.e., spectrophotometry), and that the volume of arterial blood in tissue (and hence, light absorption by that blood) changes during the pulse (i.e., plethysmography). A pulse oximeter determines SpO₂ by passing red and infrared light into an arteriolar bed and measuring changes in light absorption during the pulsatile cycle. Red and infrared low-power light-emitting diodes (LEDs) in the oximetry sensor serve as light sources; a photodiode serves as the photo detector.

Because oxyhemoglobin and deoxyhemoglobin differ in light absorption, the amount of red and infrared light absorbed by blood is related to hemoglobin oxygen saturation. To identify the oxygen saturation of *arterial* hemoglobin, the monitor uses the pulsatile nature of arterial flow. During systole, a new pulse of arterial blood enters the vascular bed, and blood volume and light absorption increase. During diastole, blood volume and light absorption reach their lowest point. The monitor bases its SpO₂ measurements on the difference between maximum and minimum absorption (i.e., measurements at systole and diastole). By doing so, it focuses on light absorption by pulsatile arterial blood, eliminating the effects of nonpulsatile absorbers such as tissue, bone, and venous blood.

Automatic Calibration

Because light absorption by hemoglobin is wavelength dependent and because the mean wavelength of LEDs varies, an oximeter must know the mean wavelength of the sensor's red LED to accurately measure SpO₂. During manufacturing, the mean wavelength of the red LED is encoded in a resistor in the sensor.

Principles of Operation

During monitoring, the instrument's software reads this resistor and selects coefficients that are appropriate for the wavelength of that sensor's red LED; these coefficients are then used to determine SpO₂. This resistor is read when the monitor is turned on, periodically thereafter, and each time a new sensor is connected.

Additionally, to compensate for differences in tissue thickness, the intensity of the sensor's LEDs are adjusted automatically.

Functional versus Fractional Saturation

This monitor measures functional saturation — oxygenated hemoglobin expressed as a percentage of the hemoglobin that can transport oxygen. It does not detect significant amounts of dysfunctional hemoglobin, such as carboxyhemoglobin or methemoglobin. In contrast, laboratory hemoximeters report fractional saturation — oxygenated hemoglobin expressed as a percentage of all measured hemoglobin, including measured dysfunctional hemoglobins. To compare functional saturation measurements to those from an instrument that measures fractional saturation, fractional measurements must be converted as follows:

$$\text{functional saturation} = \frac{\text{fractional saturation}}{100 - (\% \text{ carboxyhemoglobin} + \% \text{ methemoglobin})} \times 100$$

Measured versus Calculated Saturation

When saturation is calculated from a blood gas partial pressure of oxygen (PO₂), the calculated value may differ from the SpO₂ measurement of a pulse oximeter. This usually occurs because the calculated saturation was not appropriately corrected for the effects of variables that shift the relationship between PO₂ and saturation (Figure 1): pH, temperature, the partial pressure of carbon dioxide (PCO₂), 2,3-DPG, and fetal hemoglobin.


Figure 1: Oxyhemoglobin Dissociation Curve

SETUP

Unpacking and Inspection
Testing
Components

UNPACKING AND INSPECTION

Notify the carrier immediately if the N-180 shipping carton is damaged. Carefully unpack the instrument and its accessories. Confirm that the items in the list of components packed with the unit are included.

Inspect each component. If any component is missing or damaged, contact Nellcor Puritan Bennett's Technical Services Department or Nellcor Puritan Bennett's representative.

List of Components

Oximeter unit, hospital-grade power cord, and *Nellcor Puritan Bennett* EC-8 8-ft. sensor extension cable.

TESTING

Verify all functions as described in section *Guide to Operations*. If a difficulty occurs, refer to the *Troubleshooting* section. If that does not resolve the difficulty, contact qualified service personnel or Nellcor Puritan Bennett's representative.

COMPONENTS

Readout

Two 3-digit green and red displays for oxygen saturation and pulse rate. Sixteen-segment display for pulse amplitude indicator. Annunciators for LOW BATT, PULSE SEARCH, HIGH SAT, LOW SAT, HIGH RATE, and LOW RATE alarms, Audible Alarm Off, BATT IN USE, and NEONATAL alarm limits.

Setup

Controls

Control knob to adjust volume and set alarm limits, and five buttons to select alarm limits and disable audio alarm. Rear-panel switches for adult/neonatal alarm settings, RS-232 format, and baud rate.

Front Panel (100–120 VAC)


Figure 2a: N-180 Front Panel, 100–120 VAC Version

Front Panel (220–240 V~)


Figure 2b: N-180 Front Panel, 220–240 V~ Version

1. **ON/STDBY switch**
2. **BATT IN USE indicator**
3. **LOW BATT indicator:** Flashes when 5 or fewer minutes of battery power remain
4. **OXYGEN SATURATION numeric display**
5. **Control knob:** Changes instrument functions and sets values
6. **HIGH SAT, LOW SAT, HIGH RATE, LOW RATE buttons:** Display alarm settings and allow the changing of alarm functions
7. **Audible Alarm Off button:** Changes alarm functions and other instrument functions
8. **Audible Alarm Off indicator:** Lights up steadily when the audio alarm has been temporarily silenced; flashes when the audio alarm has been disabled
9. **HIGH SAT, LOW SAT, HIGH RATE, LOW RATE indicators:** Flash during an alarm state
10. **PULSE RATE numeric display**
11. **Pulse amplitude indicator:** Vertical column of light bars that qualitatively indicates pulse amplitude
12. **NEONATAL indicator:** Lights when neonatal alarm limits are in effect
13. **PULSE SEARCH indicator:** Flashes when the N-180 is attempting to locate the patient's pulse
14. **Sensor or sensor extension cable connector and lock**

Note: The sensor extension cable is electrically isolated.

Setup

Rear Panel


Figure 3: N-180 Rear Panel

- 15. Mains power on/off switch: Indicator lights up when connected to mains and power is on for 220–240 V~ units
- 16. Rear panel switches (Figure 4)
- 17. SERIAL COMM connector
- 18. Equipotential ground
- 19. AC mains inlet
- 20. Fuse compartment

Rear-Panel Switches


Figure 4: Rear-Panel Switches

Eight switches for digital output, adult/neonate alarm limits.

Switch	Function
1	Default alarm settings for adult (up) or neonate (down)
2	RS-232 format
3, 4, 5, 6	not used
7, 8	Baud rate select

GUIDE TO OPERATION

Basic Operation Pulse Oximetry Subsystem Features

BASIC OPERATION

1. Plug the power cord into the N-180 power inlet and then to either a 100 to 120 VAC or 220 to 240 V~ outlet as appropriate. Use only an outlet that has a grounding connection and the original hospital-grade plug and cord or an equivalent hospital-grade plug and cord. If in doubt about the integrity of the grounding of the mains supply connection, operate the unit from its internal battery.
2. Set the rear-panel mains power switch to the ON position. The switch's indicator lights up if mains power is on.
3. Apply an appropriate *Nelcor Puritan Bennett* sensor in accordance with its directions for use. Plug the sensor into the N-180 sensor connector on the front panel. Close the sensor lock over the plug. If needed, use a *Nelcor Puritan Bennett* EC-4 or EC-8 sensor extension cable.
4. Turn on the N-180 using the ON/STDBY switch. After a few seconds, following successful completion of the automatic self-test (signified by a two-tone audio signal), the monitor displays measurements, and the pulse amplitude bar graph indicator rises and falls with each pulse beat. A "beep" signals each pulse beat, and its pitch increases and decreases to reflect changes in SpO₂.

Note: When the N-180 is turned off and back on, all operator-adjustable features return to their default status.

Guide to Operation

5. Check the alarm limits each time you use the N-180. Here is how to see a 3-second display of each alarm limit value:
 - Press HIGH SAT button and observe the applicable alarm setting in the oxygen saturation LED display.
 - Press LOW SAT button and observe the applicable alarm setting in the oxygen saturation LED display.
 - Press HIGH RATE button and observe the applicable alarm setting in the pulse rate LED display.
 - Press LOW RATE button and observe the applicable alarm setting in the pulse rate LED display.
6. To adjust the alarm limits to meet a specific patient's needs, press the appropriate button (HIGH SAT, LOW SAT, HIGH RATE, or LOW RATE), and then turn the control knob until the display shows the desired value. Oxygen saturation alarm limits may be set for any value from 20% to 100%, and the pulse rate alarm limits may be set for any value from 20 to 250 beats per minute (bpm). The upper limit must be higher than the lower limit.

PULSE OXIMETRY SUBSYSTEM FEATURES

This portion of the N-180 measures and displays SpO₂ and pulse rate, and provides data for alarms and for the digital output.

Pulse Tone

When a *Nellcor Puritan Bennett* sensor is connected to the N-180 and to a patient, each pulse beat is signaled by a "beep" that varies in pitch to reflect changes in oxygen saturation, rising as saturation increases and falling as it decreases.

Changing Pulse Tone Volume

Turn the control knob to adjust pulse tone volume. When the N-180 is turned off and back on again, the pulse tone returns to its default volume.

Alarm Functions

WARNING: When the Audible Alarm Off button has been pressed and the Audible Alarm Off indicator is flashing, no audio alarm sounds in the event of an adverse patient condition. The Audible Alarm Off button must not be used in situations in which patient safety could be compromised.

WARNING: In normal operation, the Audible Alarm Off button silences the audio alarm temporarily. While it is possible to disable the audio alarm by following the instructions in this section, this must not be done if patient safety could be compromised.

Audio and visual alarms can be set for high and low oxygen saturation, high and low pulse rate, and loss of pulse. The audio alarm volume is adjustable. The alarm tone is interrupted briefly for detected pulses. Audio alarms can be disabled for a 60-second period by pressing the Audible Alarm Off button; the disable period can be changed to 30–120 seconds, or disabled indefinitely by setting it to OFF. A visual and audio warning will alert the user when the audible alarm is disabled.

Alarm States

If the oxygen saturation level or pulse rate moves beyond alarm limits, the corresponding alarm indicator flashes and the indicator color changes from green to red. The appropriate numeric display flashes, and an audio alarm sounds (unless it has been disabled).

Note: The alarm sounds immediately if the sensor or sensor extension cables become disconnected from the monitor.

If the alarms activate because the pulse signal is lost, the OXYGEN SATURATION and PULSE RATE numeric displays flash "0," the PULSE SEARCH indicator flashes, and an audio alarm sounds (unless it has been disabled).

Guide to Operation

Default Alarm Limits

Default alarm limits are in effect when the N-180 is turned on. Table 1 shows two sets of default alarm limits, one for adults and one for neonates.

Table 1: Default Alarm Settings

Alarm Limit	Adult Setting	Neonatal Setting	Adjustable Range
High oxygen saturation	100%	95%	20-100
Low oxygen saturation	85%	80%	20-100
High pulse rate	140 bpm	200 bpm	20-250
Low pulse rate	55 bpm	100 bpm	20-250

Notes: When the operator turns the N-180 off and back on again, alarm limits return to default values.

Changing from Adult to Neonatal Alarm Limits

To change the adult/neonatal setting, set rear-panel switch 1 up for adult default limits; down for neonatal default limits.

Configuring Audio Alarm

The audio alarm function can be altered in several ways: it can be silenced temporarily, it can be disabled, and its volume can be adjusted. The visible alarm cannot be turned off.

Adjusting Audio Alarm Volume

To adjust audio alarm volume, simultaneously press the HIGH SAT and LOW SAT buttons, and then turn the control knob clockwise to increase the volume or counterclockwise to decrease it. A reference numeric display in OXYGEN SATURATION correlates with the relative volume.

Silencing the Audio Alarm Temporarily

To silence the audio alarm for 60 seconds, press the Audible Alarm Off button once. The Audible Alarm Off indicator lights steadily to show that the audio alarm has been silenced temporarily. After 60 seconds, the alarm sounds again if the alarm state continues.

Note: The alarm resumes at the end of alarm silence period only if a pulse was detected at the time the alarm was turned off.

Adjusting the Temporary Silence Period

To change the period during which the audio alarm is silenced temporarily, hold the Audible Alarm Off button down, and turn the control knob until the desired period appears in the OXYGEN SATURATION display. Release the button. This period can be set for any value between 30 and 120 seconds.

Note: When the operator turns the N-180 off and back on again, the audio alarm returns to the default temporary silence period.

Disabling the Audio Alarm

To disable the audio alarm, hold the Audible Alarm Off button down, and turn the control knob clockwise until "OFF" appears on the OXYGEN SATURATION display. The audio alarm can be reactivated by pressing the Audible Alarm Off button again.

Note: For 220–240 voltage N-180 monitors, when the audible alarm is disabled, a single tone sounds every 3 minutes and the indicator next to the Audible Alarm Off button flashes to warn the operator that the alarm will not sound.

Operating Modes

The N-180 has three operating modes that enable it to make accurate measurements despite differing levels of patient activity. In all three modes, the N-180 updates its measurements with every pulse beat. Data from the most recent beat replaces data from the earliest beat, and new averages are determined and displayed.

Mode 1, the default operating mode, uses a 5- to 7-second averaging time and is useful in situations in which the patient is relatively inactive.

Guide to Operation

Mode 2 uses a 2- to 3-second averaging time and therefore is more affected by patient motion. It is useful for special applications that require a fast response time, such as sleep studies.

Mode 3 uses a 10- to 15-second averaging time and consequently is least affected by patient motion. In this mode, pulse rate is not displayed and there is no pulse tone.

Setting the Operating Mode

When the N-180 is turned on, it is set for Mode 1. If another operating mode is better suited to the specific clinical setting, simultaneously press the HIGH RATE and LOW RATE buttons, and turn the control knob to the right until the desired number appears in the PULSE RATE display.

When the N-180 is turned off and back on again, it returns to Mode 1.

Battery Operation

If AC power is lost or disconnected, the N-180 operates on its internal battery. In this case, the BATT IN USE indicator lights. When the battery level is too low to power the instrument reliably, the LOW BATT indicator lights, and after a short period, an internal switch turns off the N-180 automatically.

Note: The monitor does not audibly alarm when the battery is low.

The battery recharges when the instrument is plugged into AC power and the mains switch is on. A minimum of 14 hours is required to recharge the battery fully.

Serial Data Output

Serial data output (Table 3) from the N-180 can be in two forms: **conversation** mode, where the N-180 responds to inputs from an external ASCII device, or **computer** mode, where the N-180 provides a periodic data string of saturation and pulse rate data, alarm conditions, and monitor status.

To select the desired mode, refer to Table 4.

NELLCOR PURITAN BENNETT SENSORS

Selecting a *Nellcor Puritan Bennett* Sensor
Cleaning and Reuse
Performance Considerations

WARNING: Use only *Nellcor Puritan Bennett* oxygen transducers. Other oxygen transducers may cause improper oximeter performance.

WARNING: Before use, carefully read the sensor directions for use.

SELECTING A NELLCOR PURITAN BENNETT SENSOR

When selecting a sensor, consider the patient's weight and activity, the adequacy of perfusion, the available sensor sites, the need for sterility, and the anticipated duration of monitoring. For more information, refer to Table 2 and your local *Nellcor Puritan Bennett* representative.

Table 2: Selected *Nellcor Puritan Bennett* Sensors

Oxygen Transducer	Model	Patient Weight
<i>Oxisensor</i> ® II or <i>Oxisensor</i> ® (sterile, single use)	N-25	<3 kg, >40 kg
	I-20	3 to 20 kg
	D-20	10 to 50 kg
	D-25(L)	>30 kg
	R-15	>50 kg
<i>OxiCliq</i> ® (sterile, single use)	A	>30 kg
<i>Dura-Y</i> ® (nonsterile, reusable)	D-YS	>1 kg
<i>Oxiband</i> ® (nonsterile, reusable)	OXI-A/N	<3 kg, >40 kg
	OXI-P/I	3 to 40 kg
<i>Durgesensor</i> ® (nonsterile, reusable)	DS-	>40 kg
	100A	
<i>Nellcor Puritan Bennett</i> ® reflectance sensor (nonsterile, limited reuse)	RS-10	>40 kg

Nellcor Puritan Bennett Sensors

CLEANING AND REUSE

Caution: Do not sterilize any *Nellcor Puritan Bennett* sensor or adhesive by irradiation, steam, or ethylene oxide, and do not immerse in water or cleaning solutions.

Clean the surface of reusable sensors by wiping with an agent such as 70% isopropyl alcohol. When using adhesive sensors, use a new adhesive sensor, wrap, or bandage for each patient.

PERFORMANCE CONSIDERATIONS

Plug sensor into front panel or sensor extension cable only.

Always select an appropriate sensor, apply it as directed, and observe all warnings and cautions.

If ambient light affects performance, ensure that the sensor is properly applied, and cover the sensor site with opaque material. Failure to do so may result in inaccurate measurements. Light sources that can affect performance include surgical lights (especially those with a xenon light source), bilirubin lamps, fluorescent lights, infrared heating lamps, and direct sunlight.

If patient movement presents a problem, verify that the sensor is properly and securely applied; move the sensor to a less active site; use an adhesive sensor that tolerates some patient motion; or use a new sensor with fresh adhesive backing.

If poor perfusion affects performance, consider using an *Ox-sensor R-15* adult nasal sensor or a *Nellcor Puritan Bennett RS-10* reflectance sensor. The R-15 obtains measurements from an artery supplied by the internal carotid, the nasal septal anterior ethmoid artery; therefore, it may obtain measurements during low peripheral perfusion. The RS-10 obtains measurements from the forehead or temple, areas that may be spared when peripheral perfusion is poor.

CONNECTING TO OTHER INSTRUMENTS

RS-232 Communication Protocol

WARNING: To ensure patient electrical isolation, connect only to other equipment with circuits that are electrically isolated.

WARNING: When connecting this monitor to any instrument, verify proper operation before clinical use. Refer to the other device's manual for full instructions. Accessory equipment connected to the monitor's analog or digital outputs must be certified according to the respective standards, i.e., IEC 950 for data-processing equipment or IEC 601.1 for electromedical equipment and must comply with the requirements for medical systems, IEC 601.1.1. The accessory equipment must be powered from separate wall outlets to avoid possible summation of earth leakage (ground) currents. If in doubt, contact Nellcor Puritan Bennett's Technical Services Department or your local Nellcor Puritan Bennett representative.

RS-232 COMMUNICATION PROTOCOL

This section contains information to allow the user to provide serial communications between the N-180 and external digital devices. For greater detail refer to the N-180 service manual, Appendix A.

Serial Data Connector Pin Assignments

The serial communications connector is a 9-pin, D-type, subminiature connector.

Table 3 lists pin assignments on the serial communications connector and the corresponding pin assignments for a 25-pin RS-232 connector.

Connecting to Other Instruments

Table 3: Serial Communications Connector Pin Assignments

Connector Pin	Signal	Direction	Output Device Pin	Output Device Signal
1	not used	none	none	none
2	Rx data	←	2	Tx data
3	Tx data	→	3	Rx data
4	DTR	→	6	DSR
5	Signal Ground	↔	7	Signal Ground
6	DSR	←	4	RTS
7	RTS	→	5	CTS
8	CTS	←	20	DTR
9	not used	none	none	none

Note: Use only shielded serial input/output cables for serial communications port connections.

Communications Formats

The RS-232 format switches are used to set the communication port format, as shown in Table 4. The serial communication format of the N-180 is eight data bits, no parity bit, and one stop bit.

Computer Mode

In this mode the N-180 transmits a data string every 10 seconds, or when a parameter has changed, whichever is sooner. Set switch 2 up.

The string transmitted from the N-180 has the format shown below.

<STX>RnnnSnnnP000LnnnHnnnOnnnAnnnMnnnT000000
Qnnn<CR><LF><checksum><ETX>

<u>Character Transmitted</u>	<u>String Received</u>
R (pulse rate)	Rnnn
S (oxygen saturation)	Sann
P (unused message)	P000
L (low rate alarm limit)	Lnnn
H (high rate alarm limit)	Hnnn
O (low saturation alarm limit)	Onnn
A (alarm status in ASCII-coded decimal; logic "1" = alarm condition)	Bit 0 = high rate Bit 1 = low rate Bit 2 = high sat Bit 3 = low sat
M monitor status (in ASCII-coded decimal);	Bit 0 = pulse search status; 1 = locked, 0 = search Bit 1 = sensor status; 1 = attached, 0 = off Bit 2 = audio alarm status; 1 = enabled, 0 = disabled Bit 3 = always 1 Bits 4 through 7 not used
T (unused message)	T000000
Q (high saturation alarm limit)	Qnnn

Table 4: Output Format Switch Settings

Format	Switch 2	
	Description	Position
Conversation mode	Request for parameter	Down
Computer mode	ASCII string containing single identifier characters plus values	Up

Connecting to Other Instruments

Setting Baud Rate

Identify the baud rate required by the device that is to be connected, and set the N-180 baud rate using the rear-panel baud rate switches (switch 7 and 8). Refer to Table 5.

Table 5: Baud Rate Switch Settings

Baud Rate	Switch Position	
	Switch 7	Switch 8
2400	Down	Down
4800	Up	Down
9600	Down	Up
19,200	Up	Up

MAINTENANCE

Service
Cleaning
Battery Testing
Determining Software Version
Technical Assistance
Returning the N-180

SERVICE

WARNING: For protection against fire, replace fuses only with those of the same type and rating.

The N-180 requires no routine service other than that which is mandated by the operator's institution. The *Troubleshooting* section discusses potential difficulties, their possible causes, and suggestions for resolving them. Complete service instructions, including preventive maintenance checks, are contained in the *Nellcor Puritan Bennett N-180 service manual*.

CLEANING

Caution: Do not immerse the N-180 in liquid or use caustic or abrasive cleaners.

To clean the N-180 surfaces, dampen a cloth with a commercial, nonabrasive cleaner and wipe the top, bottom, and front surfaces lightly. Do not spray or pour any liquid directly on the N-180 or its accessories. Do not allow any liquid to come in contact with the power connector, fuse holder, or switches. Do not allow any liquid to penetrate the connectors or openings in the chassis.

BATTERY TESTING

The battery should be tested by qualified service personnel every 2 years using the procedures described in the N-180 service manual.

Maintenance

DETERMINING SOFTWARE VERSION

To determine the N-180 main software version, simultaneously press the HIGH SAT and Audible Alarm Off buttons. Turn the control knob until a "0" appears in the OXYGEN SATURATION display. Observe the number displayed in the PULSE RATE display. Turn the control knob until a "1" appears in the OXYGEN SATURATION display. Observe the number displayed in the PULSE RATE display. The combination of the numbers observed in the PULSE RATE display represents the N-180 main software version.

For example, software version 1.10 would be indicated as follows:

OXYGEN SATURATION display	0	1
PULSE RATE display	1	10

TECHNICAL ASSISTANCE

For technical information and assistance, contact Nellcor Puritan Bennett's Technical Services Department or your local Nellcor Puritan Bennett representative. To order parts or a service manual, contact Nellcor Puritan Bennett's Customer Service Department or your local Nellcor Puritan Bennett representative. The service manual includes circuit diagrams and parts lists required by qualified personnel when servicing the N-180.

RETURNING THE N-180

If it is necessary to return the monitor to Nellcor Puritan Bennett, call Nellcor Puritan Bennett's Technical Services Department for shipping instructions.

To repack the N-180, disconnect the sensor. It is not necessary to return sensors. Pack the N-180 in the original shipping carton. In case the original carton is not available, use a suitable box with an appropriate amount of packing material.

TROUBLESHOOTING

Status Messages Troubleshooting Guide

STATUS MESSAGES

This section lists status messages, along with the actions that the operator should take. If the recommended actions do not cause the message to disappear, contact qualified service personnel or your local Nellcor Puritan Bennett representative.

The service manual, which is for use by qualified service personnel, provides additional troubleshooting information.

The error messages appear in the N-180 OXYGEN SATURATION and PULSE RATE displays.

Err 1

Defective data memory. Contact qualified service personnel.

Err 2

Defective program memory. Contact qualified service personnel.

Err 3

Defective display segment or indicator, or possibly a circuit malfunction. Contact qualified service personnel.

WARNING: Continue to use the N-180 only in an urgent situation with an "Err 3" message and only if the defective segment (or segments) has been identified. If a defective segment cannot be identified, do not use the N-180.
--

Note: The N-180 can be operated if any button is pressed while "Err 3" is showing. However, if any numeric display segment or indicator is missing, the display or warning indicators may be incorrect.

TROUBLESHOOTING GUIDE

This section discusses potential difficulties and suggestions for resolving them. If the difficulty persists, contact qualified service personnel or your local Nellcor Puritan Bennett representative.

The service manual, which is for use by qualified service personnel, provides additional troubleshooting information.

This section is divided into two troubleshooting categories: Items 1 and 2 describe general system problems. Items 3 through 8 describe general oximetry subsystem problems.

1. N-180 does not turn on.

- Check AC connections. Check that monitor is connected properly to the AC supply.
- Check battery operation. If battery is discharged, connect N-180 to an appropriate AC outlet and turn on the MAINS switch. The instrument may be operated on AC with a discharged battery. The N-180 requires a minimum of 14 hours to recharge the battery completely.
- Check AC fuses.
- Check AC MAINS switch. Make sure it is on.

2. N-180 operates on AC power but not on battery.

- The battery may be discharged. To recharge the battery, connect the N-180 to an appropriate AC power outlet. Turn ON the MAINS power ON/OFF switch and confirm that the BATT IN USE indicator is OFF. A minimum of 14 hours is required to recharge the battery completely.
- The battery pack or battery charger circuit may be defective, or the battery fuse may be open.
- Check AC MAINS switch. Make sure it is on.

3. PULSE SEARCH indicator is on; SpO₂ and pulse rate not displayed.

- *Check the patient.* The patient may be experiencing shock, hypotension, severe vasoconstriction, severe anemia, hypothermia, arterial occlusion proximal to the sensor, or cardiac arrest.
- The sensor may be improperly applied (e.g., too tight) or it may not be plugged in.
- There may be excessive ambient light; cover the sensor site with opaque material.
- The sensor may be placed on an extremity with a blood pressure cuff, arterial catheter, or intravascular line.
- The patient's perfusion may be too low for the N-180 to detect an acceptable pulse.
- Test the N-180 on someone else; try another sensor site; or try another sensor (see section *Nellcor Puritan Bennett Sensors* for suggestions).
- The sensor may be damaged; replace it.

4. The pulse amplitude indicator tracks a pulse, but there is no oxygen saturation or pulse rate.

- *Check the patient.*
- Excessive patient motion may be making it impossible for the N-180 to find a pulse pattern. If possible, keep the patient still; check whether the sensor is applied securely and properly and replace if necessary; move the sensor to a new site; use a sensor that tolerates more motion; or set the N-180 for operating mode 3.
- The sensor may be damaged; replace it.
- The patient's perfusion may be too low to allow the N-180 to measure saturation and pulse rate.

Troubleshooting

5. SpO₂ or pulse rate change rapidly; pulse amplitude indicator is erratic.
- *Check the patient.*
 - Excessive patient motion may be making it impossible for the N-180 to find a pulse pattern. If possible, keep the patient still; check whether the sensor is applied securely and properly, and replace it if necessary; move the sensor to a new site; use a sensor that tolerates more motion; or set the N-180 for operating mode 3.
 - An electrosurgical unit (ESU) may be interfering with performance:
 - Move the N-180 and the cables as far from the ESU as possible.
 - Plug the N-180 and the ESU into different AC circuits.
 - Move the ESU ground pad as close to the surgical site as possible.
 - The sensor may be damp or may have need to be replaced with a new *Nellcor Puritan Bennett* sensor.
 - If using a sensor extension cable, remove it and connect the sensor directly to the N-180.
 - If the patient weighs less than 3 kg or more than 40 kg, apply an *Oxisensor N-25* oxygen transducer to an appropriate site. This sensor has added protection against electrosurgical interference.

6. Displayed pulse rate does not correlate with that of ECG monitor.

- Excessive patient motion may be making it impossible for the N-180 to find a pulse pattern. If possible, keep the patient still; check whether the sensor is applied securely and properly and replace if necessary; move the sensor to a new site; use a sensor that tolerates more motion; or set the N-180 for operating mode 3.
- The patient may have a pronounced dicrotic notch, which may cause the pulse rate measurement to double. Try another sensor site.
- An electrosurgical unit may be interfering. Refer to discussion under item 5.

7. Oxygen saturation measurement does not correlate with the value calculated from a blood gas determination.

- The SpO₂ calculation may not have correctly adjusted for the effects of pH, temperature, PCO₂, 2,3-DPG, or fetal hemoglobin. Check whether calculations have been corrected appropriately for relevant variables (See section *Principles of Operation* for more information.) In general, calculated saturation values are not as reliable as direct laboratory hemoximeter measurements.
- Accuracy can be affected by incorrect sensor application or use, significant levels of dysfunctional hemoglobins, intravascular dyes, bright light, excessive patient movement, venous pulsations, electrosurgical interference, and placement of a sensor on an extremity that has a blood pressure cuff, arterial catheter, or intravascular line. Observe all instructions, warnings, and cautions in this manual and in the sensor's directions for use.

Troubleshooting

8. Oxygen saturation does not correlate with laboratory hemoximeter.
 - a. Fractional measurements may not have been converted to functional measurements before the comparison was made. The N-180, as well as other two-wavelength oximeters, measures functional saturation. Multi-wavelength oximeters, such as laboratory hemoximeters, measure fractional saturation. Fractional measurements must be converted to functional measurements for comparison. Refer to the equation for this conversion in the *Principles of Operation* section.
 - b. Close correlation requires simultaneous blood sampling and pulse oximeter measurements from the same arterial supply.

SPECIFICATIONS

Performance
Environmental Requirements
Electrical Characteristics
Physical Characteristics

PERFORMANCE

Range

Saturation
0-100%

Pulse Rate
20-250 bpm (beats per minute)

Accuracy

SpO ₂ ¹		
Adults	70-100%	±2 digits ²
	0-69%	unspecified
Neonates	70-100%	±3 digits ²
	0-69%	unspecified

Pulse Rate
20-250 bpm ±3 bpm²

1. The reference for accuracy testing is an N-200 oximeter and D-25 sensors that have been validated in human blood studies against a laboratory hemoximeter. Adult accuracy specification is based on testing with D-25 sensors; neonatal accuracy specification is based on testing with N-25 sensors. All accuracy testing was conducted on health adult volunteers in induced hypoxia studies. For specifications with other *Nellcor Puritan Bennett* sensors, call Nellcor Puritan Bennett's Technical Services Department or your local Nellcor Puritan Bennett representative.
2. This variation equals one standard deviation (SD). Plus or minus one SD encompasses 68% of the population.

Specifications

ENVIRONMENTAL REQUIREMENTS

Temperature

N-180

50–104°F (10–40°C) operating

32–122°F (0–50°C) storage

Sensor

Within physiologic range 82.4–107.6°F
(28–42°C) for accurate measurement.

Humidity

Any humidity/temperature combination without condensation

Altitude/Barometric Pressure

0–10,000 ft (0–3,048 m)

760–530 mmHg

ELECTRICAL CHARACTERISTICS

Protective Class

Class I: mains-supplied unit using a protective ground

Degree of Protection

Type BF: patient electrically isolated

Voltage

100–120 VAC, 50/60 Hz

220–240 V~, 50/60 Hz

Data I/O

SIP/SOP nominal output voltage ± 12 V

Specifications

Fuses

100-120 VAC
2 x T 0.25 A, slow-blow, 250 V

220-240 V~
2 x T 0.10 A, 250 V

Power Consumption

Maximum rating: 25 VA

Battery

Type
Lead-acid battery pack, 1.9 AH

Battery Life
6 hours typical on full charge; battery age and usage may affect battery performance.

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© 1996 Nellcor Puritan Bennett Inc. Tutti i diritti riservati.