
User's Guide

Nucleic Acid Isolation System QuickGene-610L

Rev. 1.0

P3X1-E380-01EN

Introduction

Thank you for using the Nucleic Acid Isolation System QuickGene-610L (hereafter, called "QuickGene-610L"). Easy operation of this device enables extraction of DNA and RNA.
This manual contains important information regarding the safe and proper use of the QuickGene-610L functions. Read this manual carefully before use.
Keep this manual after reading, for reference whenever necessary.

About Infectious Materials

When handling or disposing of infectious materials, follow the laboratory guideline or the law regarding infectious wastes to perform the incineration, fusion, sterilization, and/or disinfection.

When you use a third party to dispose it, outsource it to a licensed operator for industrial waste subject to special control, and give them the management table (manifest) for industrial waste subject to special control at the same time.

Universal biological hazard symbol

■ Exclusion clauses

- The contents of this manual may be updated without notice.
- Fuji Photo Film Co., Ltd. assumes no liability for damages to third party copyrights or other rights arising from the use of any information in this manual.
- Fuji Photo Film Co., Ltd. assumes no liability for device failures or damages due to installation, relocation, remodeling, maintenance, and repairs performed by parties other than Fuji Photo Film Co., Ltd. or agencies specified by Fuji Photo Film Co., Ltd.
- Fuji Photo Film Co., Ltd. assumes no liability for failures or damages of our products due to third-party products other than products delivered by Fuji Photo Film Co., Ltd.
- Fuji Photo Film Co., Ltd. assumes no liability for device failures or damages due to remodeling, maintenance, repairs, etc. using repair parts other than genuine parts specified by Fuji Photo Film Co., Ltd.
- Fuji Photo Film Co., Ltd. assumes no liability for device failures or damages due to failure to observe the precautions and operation procedures described in this manual.
- Fuji Photo Film Co., Ltd. assumes no liability for device failures or damages due to using the device in an environment that does not meet the operational requirements (power, installation environment, etc.) described in this manual.
- Fuji Photo Film Co., Ltd. assumes no liability for device failures or damages due to natural disasters such as fire, earthquakes, floods, or lightning.
- No part of the contents of this manual may be reproduced without permission from Fuji Photo Film Co., Ltd.

■ Trademarks

Company names and product names are trademarks or registered trademarks of their respective owners.

All Rights Reserved, Copyright © Fuji Photo Film Co., Ltd. 2005

Notations

This manual contains important information regarding the safe use of the QuickGene-610L.

Read this information carefully before using the QuickGene-610L, and start operation.

The meanings of safety precaution marks are as follows:

"WARNING" indicates a dangerous condition that may lead to death or serious injury.

BIOHAZARD

"BIOHAZARD" indicates exposure to a biological agent or condition that may be infectious or harmful.

"CAUTION" indicates the possibility of danger that may lead to minor or moderate personal injury.

"CAUTION" may also indicate an accident of property damage.

Important:

"Important" shows the important notes for usage, as well as prohibited actions.

Note:

"Note" indicates the notes, procedures that should be obeyed, and supplementary information for use.

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Safety Information

This chapter describes safety precautions.

For your safety and that of others, follow the guidelines provided in the following pages concerning the use of the QuickGene-610L.

About Instrument

WARNING:

- **Ignoring the following notations may lead to fire or electric shock.**
 - Use a power cable that meets your country's standard, or contact your local distributor.
 - Connect to a grounded outlet to reach the power source.
 - Do not touch the power plug with wet hands
 - Do not use the QuickGene-610L with voltage other than the voltage specified on the device. Do not overload the power outlet with too many devices.
 - Do not use the QuickGene-610L with a damaged power plug or a loose socket.
 - When dust is on the prongs of the power plug or the plug socket, clean with a dry cloth.
 - When you disconnect the power plug from the outlet, be sure to hold the power plug itself. Do not pull the power cable.
 - For maintenance, disconnect the power plug from the outlet.
 - Do not touch the power plug when you hear the crash of thunder.
 - When the circuit breaker trips, it may be due to a short circuit. Do not dismantle the QuickGene-610L.
- **Do not pour any liquid on the QuickGene-610L. Do not place any objects containing liquid on the QuickGene-610L. Doing so may cause a device failure, fire, or electric shock.**
- **In the event the device overheats, starts to smoke or smells strange, immediately disconnect the power cable.**

CAUTION:

- **Never attempt to remodel the QuickGene-610L without the permission of Fuji Photo Film Co., Ltd. Doing so may cause fire or electric shock.**
- **Do not place or drop objects on the QuickGene-610L. Also refrain from bumping or knocking it, as doing so may cause a failure or malfunction of the QuickGene-610L.**
- **If any liquid materials are left inside the device, wipe with a Kimwipe® or paper tissue. Otherwise, the QuickGene-610L may be damaged.**
- **The front cover will not stay on the way to the upper limit. If you release the front cover while opening it, it may cause injury or device damage.**
- **When replacing the Pump Tube, be sure to wear gloves.**
- **In the event of an ethanol spill, immediately disconnect the power cable.**
- **Repairs to the QuickGene-610L should only be performed by such agencies as are specifically authorized by Fuji Photo Films Co., Ltd.**
- **Only original Fuji Photo Films Co., Ltd. replacement parts should be used.**
- **If the equipment is used in a manner not specified by the manufacturer, the protection provided by the equipment may be impaired.**

About Installation Conditions

-  **WARNING:**
- Do not install the QuickGene-610L to a location where it can be splashed with water. Doing so may cause a device failure, fire, or electric shock.
 - When relocating the QuickGene-610L, be sure to disconnect the power plug from the outlet first. Watch your step when you relocate the QuickGene-610L. If the power cable is damaged, this may cause a device failure, fire, injury, or electric shock.

-  **CAUTION:**
- Do not install the QuickGene-610L in an unstable place such as a slanted surface or a place subject to vibrations. Doing so may cause injury or device failure.
 - Do not install the QuickGene-610L in direct sunlight or close to a heating device. Doing so may shorten the life of QuickGene-610L, or cause a failure.
 - Two or more people are necessary to relocate the QuickGene-610L. Otherwise, you may hurt your back or hands because of the weight of QuickGene-610L (24kg).

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Product Use Limitations

- Fuji Photo Film Co., Ltd. assumes no liability for results regardless of the intended use.
- It is the user's responsibility to validate the performance of QuickGene kits and QuickGene-610L for any case of use.
- The QuickGene-610L is a device for research and study only.

When you use the QuickGene-610L, read this manual thoroughly and follow the safety management regulations of your facilities (laboratory).

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2 Basic Functions

2.1	Components	19
2.2	Working Process of QuickGene-610L	20
2.3	Safety Mechanism	20
2.4	LCD Display Information.....	21
2.5	Main Specifications.....	22

3 Operation

3.1	Quick Guide.....	23
3.2	Turn the Power ON and OFF	24
3.3	Preparation	25
3.4	Preparation of Lysate (Samples)	26
3.5	Extraction Operations	27
3.6	When the Front Cover is Opened During PROCESSING.	38
3.7	Operations to Stop Extraction Process.....	39

4 Maintenance

4.1	Before Using the QuickGene-610L.....	41
4.2	When the QuickGene-610L is left for a Week or More.....	41
4.3	Cleaning the QuickGene-610L	42
4.4	Cleaning the Parts	43
4.5	How to Replace the Pump Tube.....	43

5 How to Set Parameters

5.1	Parameters	47
5.2	How to Set/Change Parameters.....	49
5.3	Parameters in EXPERT Mode.....	51
5.4	Changing Extraction Mode Names.....	53

6 Troubleshooting

6.1	Extraction Failure.....	55
6.2	Typical Failure	55
6.3	Error Messages	56
6.4	When the Buffer is Contaminated.....	57

Appendix A

A.1	Usable Power Supply Cords.....	59
A.2	Options	60
A.3	About Services.....	61
A.4	Repair Service Period of Product	61
A.5	Customer Support.....	61
A.6	Notes for Transfer QuickGene-610L.....	62

INDEX

1 Installation

This chapter explains how to install the QuickGene-610L.

1

2 Basic Functions

This chapter explains the basic functions of the QuickGene-610L.

2

3 Operation

This chapter explains the extraction procedures.

3

4 Maintenance

This chapter describes how to maintain the QuickGene-610L.

4

5 How to Set Parameters

This chapter explains how to set parameters.

5

6 Troubleshooting

This chapter describes the procedures that should be taken when you suspect device failures, and the messages shown on the LCD.

6

Appendix A

This appendix describes after-sales services and points to note for transfer.

A

1 Installation

This chapter explains how to install the QuickGene-610L.

1.1 Installation and Use Conditions

This section explains the installation conditions and operating environment of the QuickGene-610L.

■ Installation conditions


- WARNING:**
- Do not install the QuickGene-610L to a location where it can be splashed with water. Doing so may cause a device failure, fire, or electric shock.
 - When relocating the QuickGene-610L, be sure to disconnect the power plug from the outlet first. Watch your step when you relocate the QuickGene-610L. If the power cable is damaged, this may cause a device failure, fire, injury, or electric shock.

- CAUTION:**
- Do not install the QuickGene-610L in an unstable place such as a slanted surface or a place subject to vibrations. Doing so may cause injury or device failure.
 - Do not install the QuickGene-610L in direct sunlight or close to a heating device. Doing so may shorten the life of QuickGene-610L, or cause a failure.
 - Two or more people are necessary to relocate the QuickGene-610L. Otherwise, you may hurt your back or hands because of the weight of QuickGene-610L (24kg).

For correct and safe use of QuickGene-610L, install it in a location that meets the following conditions:

- The space shown in the following figure must be allocated.

1 Installation

1

- A location where power can be provided. Use an AC100±10% to AC240V±10%, 50Hz or 60Hz outlet.
- A location where the temperature is 15 to 30°C, and humidity is 30 to 80%RH (non condensing).
- A location that will support the weight of QuickGene-610L (24kg), that is flat and stable, with no vibration.
- A location away from direct sunlight. Block the sunlight by curtains or blinds as necessary.
- A location which is well-ventilated and not dusty.
- A location where the temperature does not go up and down suddenly. Warming a cold room suddenly or moving the QuickGene-610L from a room with low temperature to a warm room may cause condensation inside the device (condensing), and may have a negative influence on the sample result.
- A location far from water taps, water heaters, humidifiers, air-conditioners and heaters. (Avoid areas of high temperature and high humidity, or low temperature and low humidity.)
- A location far from objects which generate strong magnetic fields (motors, transformers, TV, audio speakers, magnets, etc.). Bringing the QuickGene-610L close to any type of magnetic field may cause a malfunction.


■ Use conditions

The following shows the operating environment of the QuickGene-610L.

Items		Specifications
Temperature (°C)	During operation	15 ~ 30
	During downtime	0 ~ 55
Humidity (RH)	During operation	30 ~ 80
	During downtime	10 ~ 80
Max. wet bulb temperature (°C)	During operation	29 (non condensing)
	During downtime	29 (non condensing)
Temperature gradient (°C/Hr)		12 or less (non condensing)
Humidity gradient (RH/day)		30 or less (non condensing)
Altitude, operating (m)		2000 or lower
Installation (Overvoltage) Category Definition		CAT II Local-level mains (wall sockets). Equipment at this level includes appliances, portable tools, and similar products. Equipment is usually cord-connected.
Pollution Degree Definition		Pollution Degree2 Normally only dry, nonconductive pollution occurs. Occasionally a temporary conductivity that is caused by condensation must be expected. This location is a typical office/home environment. Temporary condensation occurs only when the product is out of service.
Pollution Degree		Pollution Degree 2. Note: Rated for indoor use only.

1.2 Checking the QuickGene-610L and Accessories

Check that the following are included in the package.
Contact your local representative if you notice any missing part(s).


1.3 Remove the Packaging Materials

After unpacking the QuickGene-610L from the box, do as follows:

1

(1) Remove the tape that holds the front cover.


(2) Open the front cover, and remove the following items.

- Clamps for fixing the carrier both side of the Carrier
- Remove the tape that holds the tubes.

Important: Keep the shipping box and 2 clamps. They will be necessary for transporting QuickGene-610L.


CAUTION: Open the front cover as far as it will go, using both hands. The front cover will not stay open until it is fully opened. If you release the front cover while opening it, it will drop suddenly and may cause injury or damage to the device.

1.4 Overview


This section explains the name and function of each part in the QuickGene-610L.

■ Front and side view

1


■ Rear view

1 Installation


■ Internal view

1

No.	Name	Function
1	Wash Buffer Pump	Sends the wash buffer in the Wash Buffer Bottle to the Carrier.
2	Bottle Holder	Set the Wash Buffer Bottle and Elution Buffer Conical Tube.
3	Elution Buffer Pump	Sends the elution buffer in the Elution Buffer Conical Tube to the Carrier.
4	Discharge Tray Table	Set the Discharge Tray to the Discharge Tray Table.
5	Holder Stand	Set the Holder Carriage with the Collection Tube Holder and Cartridge Holder mounted, to the Holder Stand.
6	Carrier	Moves automatically during extraction operation.
7	Pressure Head	Presses out the air and uses the pressure to flush out solutions.
8	Elution Buffer Nozzle	Puts elution buffer into each cartridge.
9	Wash Buffer Nozzle	Puts wash buffer into each cartridge.
10	Sensor	Checks that the carrier is placed in the correct position. Caution: Do not touch the sensors.
11	Cartridge Holder	Set cartridges in the Cartridge Holder.
12	Holder Carriage	Set the Waste Tube in the Holder Carriage.
13	Collection Tube Holder	Set Collection Tubes in the Collection Tube Holder.

■ Operation Panel

No.	Name	Function
1	LCD	Displays messages. For the displayed messages, refer to "2.4 LCD Display Information"(→ page 21) and "6.3 Error Messages"(→ page 56).
Lamp (LED)		
2	POWER	The green light is on when the power is on. The light is off in the following cases. •When the power is off •When an abnormality is detected in the power unit
3	READY	The green light is on when the QuickGene-610L is ready for operation. The light is off during extraction process or setup.
4	BINDING	The green light blinks during binding process or when extraction process is completed. The orange light blinks when an error is detected. The light is off in other cases.
5	WASHING	The green light blinks during washing process or when extraction process is completed. The orange light blinks when an error is detected. The light is off in other cases.
6	ELUTION	The green light blinks during collection process or when the extraction process is completed. The orange light blinks when an error is detected. The light is off in other cases.
Button.		
7	START	Starts a program. The operation cannot be started while the front cover is open, or if an error has occurred. In the setup mode, changed parameters are saved and quit setup mode.
8	STOP	Halts the program of QuickGene-610L. In the setup mode, changed parameters are not saved and quit setup mode.
9	RESET	Can be used in the following cases to abort the program. • When pressing the [STOP] button during the extraction process • When the front cover is opened during extraction process and the process is paused, and then the front cover is closed In the setup mode, changed parameters are not saved and quit setup mode.

1 Installation

1

No.	Name	Function
10	DISCHARGE	Discharges the old wash buffer and elution buffer remaining in the tubes. Discharges the remaining air from the tubes, and fills the tubes with new wash buffer and elution buffer.
11	MODE	Switch the extraction mode. Each time the button is pressed, the mode is switched as the following order. DNA WHOLE BLOOD ISOLATE A ISOLATE B ISOLATE C ISOLATE D ISOLATE E ISOLATE F Extraction mode names can be changed.
12	CR MOVE	Moves the carrier from side to side. Each time the switch is pressed, the carrier moves from side to side. Can be used while the front cover is closed.
13	▲	Increments the setting value of the item displayed in the LCD at setup mode. Can only be used at setup mode. Pressing and simultaneously starts the setup.
14	▼	Decrements the setting value of the item displayed in the LCD at setup mode. Can only be used at setup mode.

1.5 Setup Parts

This section explains how to install the parts included in the package.

Install the following parts:

- Power Cable → "■Setup the power cable"(→ page 9)
- Pump Tube → "■Setup the Pump Tube"(→ page 9)
- Tube → "■Setup the tubes"(→ page 12)
- Bottle Holder → "■Setup the Bottle Holder"(→ page 14)
- Wash Buffer Cap, Elution Buffer Cap, Wash Buffer Bottle (500ml)
→ "■Setup the Wash Buffer Cap and Elution Buffer Cap"(→ page 14)
- Discharge Tray → "■Setup the Discharge Tray"(→ page 18)
- Holder Carriage, Holder Stand → "■Setup the Holder Carriage"(→ page 16)

An empty Wash Buffer Bottle and Elution Buffer Conical Tube are needed to set a Wash Buffer Cap and Elution Buffer Cap.

Prepare a Wash Buffer Bottle and an Elution Buffer Conical Tube in order to set a Wash Buffer Cap and an Elution Buffer Cap respectively. One Wash Buffer Bottle (500ml) is included in the package.

The recommended bottles are as follows:

Type	Product name
Wash Buffer Bottle (500ml)	Nalgenunc 500ml Narrow-Mouth Bottle
Elution Buffer Conical Tube	BD Falcon 50ml conical tube

■ Setup the power cable

The following explains how to install the power cable.


WARNING:

- Ignoring the following notations may lead to fire or electric shock.
- Use a power cable that meets your country's standard, or contact your local distributor.

Important: Refer to "A.1 Usable Power Supply Cords" (→ page 59), for details about the power cables that can be used in each country.

- (1) Connect one of the power plugs of the power cable to the inlet of QuickGene-610L.

- (2) Connect the other plug of power cable to the outlet.

WARNING:


Connect to a grounded outlet to reach the power source.

■ Setup the Pump Tube

The following explains how to setup the Pump Tube.

Note: When setting the Pump Tube, wear gloves.

- (1) Turn the power on.

Press the [I] side of the switch.

When the QuickGene-610L has been initialized, the LCD displays "DNA WHOLE BLOOD" (this name can be changed). It takes about a minute.

- (2) Press the [DISCHARGE] button.

Press the [DISCHARGE] button. This operation makes the drive sockets (→ page 10) of the Wash Buffer Pump and Elution Buffer Pump rotate to the position for ease to setup Pump Tubes.

- (3) Turn off the power.


Press the [O] side of the switch.

1

(4) Open the front cover.

Hold the handle and open the front cover until the upper limit is reached.



CAUTION: The front cover will not stay on the way to the upper limit. If you release the front cover while opening it, it may cause injury or device damage.

(5) Attach the Pump Tube to the pump.


Attach the Wash Buffer Pump Tube (purple) and Elution Buffer Pump Tube (transparent) to the Wash Buffer Pump and the Elution Buffer Pump respectively.


Important: When attaching the Pump Tube, wear gloves. Do not touch the drive socket in the pump without gloves.

(a) Detach the pump cover.

Remove the two screws while holding the Pump Cover so that it does not fall.

(b) Attach the Pump Tube to the rear of the Pump Cover.

Slide the connectors at the ends of the Pump Tube into the grooves of the Pump Cover.

Check that the notches on both connectors are facing the front of the Pump Cover. Slide the Pump Tube into place inside the Pump Cover by rotating the rotor.

Adjust the rotor to its original vertical position.

"WASH" or "ELUTION" is printed on the bag of the Pump Tubes. Set each type of Pump Tube correctly.

WASH: for Wash Buffer Pump


ELUTION: for Elution Buffer Pump

Important: Check that both ends of the Pump Tube are placed in the grooves.

(c) Attach the Pump Cover.

Turn the rotor vertical position, then re-attach the Pump Cover to the retainer, making sure that the drive head and drive socket are correctly mated. Secure the Pump Cover in place with the original 2 screws.

Important: Check the Pump Cover is fixed well to the Rear. Then tighten the screws until the spring washer is squeezed.

This completes the setup of the Pump Tube. For the next step, setup the tubes.


■ Setup the tubes

The following explains how to setup the tubes.

Before setting up the tubes, setup the Pump Tubes.(→ "■Setup the Pump Tube"(→ page 9))

1

Note: When setting the tubes, wear gloves.


(1) Connect the tubes to the Pumps.

- (a) Connect the upper tube (from the instrument) to the Wash Buffer Pump, and connect the lower tube (from the instrument) to the Elution Buffer Pump, by turning as shown in the figure.
- (b) Connect the thick tube (upper: purple) that has the Wash Buffer Cap to the Wash Buffer Pump, and connect the thin tube (lower: blue) that has the Elution Buffer Cap, by turning as shown in the figure.

(2) Take off the Wash Buffer Cap and Elution Buffer Cap from the tubes connected to the pumps, by turning Tube Connector.


Note: Hold the cap, and turn the Tube Connector to disconnect the tube.


To avoid contamination, place each Buffer Cap on cleaned bottle or tube.

(3) After removing the caps, remove the tubes through the holes at the lower part of the QuickGene-610L.

Remove the thick tube connected to the Wash Buffer Pump through the lower hole, and remove the thin tube connected to the Elution Buffer Pump through the lower hole.


(4) Push in the rubber packing to the holes from inside the instrument.


(5) Hook the tubes with clamps (2 positions).

Push the tubes into the clamps attached to the inside of the instrument.
Clamps are located 2positions as shown in the figure.


(6) Put the tubes along the grooves on the left panel.
Right groove is for the thin tube, Left groove is for the thick tube.


■ Setup the Bottle Holder

The following explains how to set the Bottle Holder.

(1) Attach the Bottle Holder to the QuickGene-610L.


Set the Bottle Holder at the Roller Catch on the left panel of the instrument.

(2) Set an empty Wash Buffer Bottle and Elution Buffer Conical Tube in the Bottle Holder.

■ Setup the Wash Buffer Cap and Elution Buffer Cap


The following explains how to install the Wash Buffer Cap and Elution Buffer Cap.

(1) Set the Wash Buffer Cap in the Wash Buffer Bottle, and the Elution Buffer Conical Cap in the Elution Buffer Conical Tube.

- Set the purple Wash Buffer Cap on the Wash Buffer Bottle (Large).
- Set the blue Elution Buffer Cap on the Elution Buffer Conical Tube (Small).

Note: During buffer bottle replacement, each tubes can be placed on each clamps temporarily


To avoid contamination, do not touch the end of the tubes.


(2) Connect the Wash Buffer Cap to the Tube Connector ,thick tube and connect the Elution Buffer Cap to the Tube Connect of thin tube, by turning as shown in the figure.

- Connect the thick tube for the Wash Buffer Pump (upper, purple) and the purple Wash Buffer Cap.
- Connect the thin tube for the Elution Buffer Pump (lower, blue) and the blue Elution Buffer Cap.

Note: Hold the cap, and turn the Tube Connector to disconnect the tube.

(3) Set the Bottle Holder Cover.


■ Setup the Holder Carriage

The following explains how to install the Holder Carriage.

1


(1) Attach the Holder Stand to the QuickGene-610L.

Insert the tabs of the Holder Stand in the slits on the QuickGene-610L.

Note: After attaching the Holder Stand, check that there is no gap between the bottom of Holder Stand and the QuickGene-610L.

(2) Set the Collection Tube Holder and Cartridge Holder in the Holder Carriage.


- To set the Collection Tube Holder, put its left and right latches into the notches on the Holder Carriage.
- To set the Cartridge Holder, put the left and right tabs of the Holder Carriage into the side notches of the Cartridge Holder.

(3) Check that the Collection Tube Holder and Cartridge Holder are correctly set in the Holder Carriage.

Check the following.

- The left and right long holes at the front of the Holder Carriage can be seen.
- The left and right tabs of Holder Carriage are set in the notches on the left and right sides of the Cartridge Holder.


(4) Attach the Holder Carriage to the QuickGene-610L.

Hold the handles of the Holder Carriage and attach the Holder Carriage onto the Holder Stand of the QuickGene-610L.

Check the following in order to secure the correct setting.

- The left and right hooks of the Cartridge Holder are put into the grooves on the shaft of the QuickGene-610L.
There is a little gap between the hook and the groove.


- The bottom front of the Holder Carriage matches the table front.
- The 2 long holes on the left and right at the front of the Holder Carriage match the prongs of the Table.

Note: When the Carrier is not at the left end and you cannot attach the Holder Carriage, close the front cover and turn on the QuickGene-610L. Press the [CR MOVE] button to move the Carrier to the left end, and then set the Holder Carriage.

■ Setup the Discharge Tray

1

(1) Set the Discharge Tray to the Discharge Tray Table.

Note: When the Carrier is at the left end and you cannot set the Discharge Tray, close the front cover and turn on the QuickGene-610L. Press the [CR MOVE] button to move the Carrier to the right end.

(2) Close the front cover.

1.6 Checking QuickGene-610L Operation

This section explains how to check the operation of QuickGene-610L.

After all the parts are installed, check the operation of QuickGene-610L.

To check the operation, use sterile distilled water instead of buffers and dissolved samples, and perform the same operation as the normal extraction operation.


For the extraction operation, refer to "3 Operation" (→ page 23).

2 Basic Functions

This chapter explains the basic functions of the QuickGene-610L.

2.1 Components

The following describes the main components of the QuickGene-610L:

2

No.	Component	Function
1	Operation unit	Operate the buttons to select an extraction mode and start/stop an extraction process. The LCD and LEDs display the device status.
2	Pump unit	Sends the buffers (wash buffer, elution buffer) to the Carrier unit according to the extraction program.
3	Carrier unit	Puts pressurized air and buffers (wash buffer, elution buffer) into the cartridge and performs the extraction process (binding - washing - elution).
4	Sample unit	By moving the Holder Carriage, the waste liquid and extracted liquid are discharged into the Waste Tube and Collection Tube, respectively.
5	Buffer setting unit	The buffers (wash buffer, elution buffer) are set here. The buffers set are sent to the discharge unit through the pump unit.

2.2 Working Process of QuickGene-610L

The following describes the working flow of the QuickGene-610L:

- (1) **The extraction program is started by using the operation unit.**
- (2) **The wash buffer and elution buffer are sent from the buffer setting unit to the pump unit.**
- (3) **The wash buffer and elution buffer are sent from the pump unit to the carrier unit.**
- (4) **The pressurized air, wash buffer, and elution buffer are sent from the carrier unit into the cartridge (with lysate), and the extraction process is performed.**
- (5) **The extracted sample is eluted to the Collection Tube set in the sample unit.**

The waste is discharged to the Waste Tube.

2.3 Safety Mechanism

This section explains the safety mechanisms of the QuickGene-610L.

The QuickGene-610L provides the following safety mechanisms:

- Front cover
- Earth leakage breaker
- LCD, Lamp (LED), and buzzer

■ Front cover

When the front cover is open, operations such as extraction processes cannot be performed.

Also, if the front cover is opened during the extraction process, it will be temporarily interrupted.

However, when the buffer pump is working, the pump will continue to work until the step is completed. Operation can be restarted from this point in the process by closing the front cover and pressing the start button.

Note: During the extraction process, do not open the front cover unless it is necessary.

■ Earth leakage breaker

When there is an electric leak or current surge caused by failures of the QuickGene-610L, the earth leakage breaker will cut off the power supply.

Check the breaker condition through the Breaker Checkpoint. When the breaker is OK you can see "I". If the breaker has tripped you will see "O".

If the earth leakage breaker trips, contact your local representative.

■ LCD, Lamp (LED), and buzzer

The LCD, Lamp (LED), and buzzer display the operating status of the QuickGene-610L.

2.4 LCD Display Information

The following describes the messages displayed in the LCD during normal operations.
For the messages displayed in case of errors, refer to "6.3 Error Messages" (→ page 56).

Message	Contents
FIRMWARE: VxxLxx	Shows the Software version. This appears when you turn on the power.
WAIT A MOMENT	Initializing the QuickGene-610L or disposing of the waste.
DNA WHOLE BLOOD	"DNA WHOLE BLOOD" mode is set to operate on the QuickGene-610L. The operation has not started yet.
ISOLATE A	"ISOLATE A" mode is set to operate on the QuickGene-610L. The operation has not started yet.
ISOLATE B	"ISOLATE B" mode is set to operate on the QuickGene-610L. The operation has not started yet.
ISOLATE C	"ISOLATE C" mode is set to operate on the QuickGene-610L. The operation has not started yet.
ISOLATE D	"ISOLATE D" mode is set to operate on the QuickGene-610L. The operation has not started yet.
ISOLATE E	"ISOLATE E" mode is set to operate on the QuickGene-610L. The operation has not started yet.
ISOLATE F	"ISOLATE F" mode is set to operate on the QuickGene-610L. The operation has not started yet.
PROCESSING	Working the extraction program. Or waiting to shift to wash buffer discharges in the second round (the waiting time can be set).
COMPLETED	The extraction operation is completed. The extraction result will be displayed after showing this message for one second.
START SW → RESTART	Temporarily stops. This message is displayed in the following cases: <ul style="list-style-type: none"> • Temporarily suspended state for the second round wash buffer discharge • Temporarily suspended state during the multiple times sample processing • When opening the front cover during the extraction process to interrupt it, and then closing the cover • When pressing the [STOP] button during the extraction process
1?2?3?4?5?6?	This shows the extraction result. "?" indicates one of the following marks and meanings: <ul style="list-style-type: none"> ✓ (check): Successfully extracted - (hyphen): Extraction failure (Cartridge is clogged) _ (underscore): No cartridge, or No sample

2 Basic Functions

2.5 Main Specifications

The following describes the main specifications of the QuickGene-610L.

■ Main specifications

Items	Specifications
Model name	QuickGene-610L
Dimensions	448mm (W) x 332mm (D) x 398mm (H)
Weight	24kg
Rated input voltage	AC100±10% to AC240V±10%
Power frequency	50/60 Hz
Current consumption	1.0A to 0.5A
Phase	Single(with 3P grounding pole)
Range of temperature and relative humidity	Temperature: 15 - 30°C Humidity: 30 - 80%RH (non condensing)
Samples	Max. 6 samples can be processed simultaneously
Altitude	2000m or lower (for operation)
Noise	Sound pressure: 65db (A-weighted sound pressure level)
Operation panel display	LCD 16 characters x 1 line
Sit	Indoor use only

■ Environmental conditions

EMC	FCC Rule Part15b ClassA EN61326 ClassA EN61000-3-2 EN61000-3-3
Safety	UL61010A-1 CSA C22.2 No.1010 EN61010-1 IEC61010-1
CE	LV directive (73/23/EEC) EMC directive (89/336/EEC)

FCC WARNING

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.


This equipment has been tested and found to comply with the limits for a Class A digital device pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference, in which case, the users will be required to correct the interference at their own expense.

Warning


This is a Class A product of Electromagnetic Interference (EMI) standard. In a domestic environment this product may cause radio interference in which case the user may be required to take adequate measures.

This Class A digital apparatus complies with Canadian ICES-003.

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Set the Cartridge Holder notch to the Holder Carriage tab

Close the cover,
and lock the both
fasteners



Table prong
Long hole
The bottom front matches the table front

3 Operation

- Set cartridges, Collection Tubes, and Waste Tubes in the correct positions.

Insert the cartridges and tubes starting from the left, in order.

If they are not in order, a spill of buffer and/or dissolved sample may occur. This results in an extraction failure, and buffer and/or dissolved sample waste.


3.2 Turn the Power ON and OFF

3

This section explains how to turn the power of QuickGene-610L on and off.

■ Turning on the power

(1) Check the power cable connection and the earth leakage breaker.


- Check whether the power cable is connected to the QuickGene-610L and the outlet. If not, connect it.
- Check that the earth leakage breaker has not tripped. When the breaker is OK you can see "I". If the breaker has tripped you will see "O". If you see "O", disconnect the power cable and contact your local representative.

! WARNING:

- Connect to a grounded outlet to reach the power source.
- Do not dismantle the QuickGene-610L. Doing so may cause an electric shock.

(2) Press the [I] side of the switch.

Just after the power is turned on, the LCD displays the software version for a few seconds, then changes the display to "WAIT A MOMENT". LCD displays extraction mode name, when the initialization has been completed.

■ Turning off the power

(1) Press the [O] side of the switch.

3.3 Preparation

This section explains the preparation of the QuickGene-610L for an extraction process. The following procedures must be performed as preparation of the QuickGene-610L:

- Preparation items
- Discharge

■ Preparation items

Before a test run or extraction operation, perform the following:

- Check all the parts are installed ("1.5 Setup Parts"(\rightarrow page 8)).
- Prepare the following:
 - Gloves
 - Mask
 - Wash Buffer Bottle (500ml) *The bottle provided with the kit can also be used.
 - Elution Buffer Conical Tube
 - 1.5ml micro tube
 - Waste Tube (WTL) *provided with the kit
 - Cartridge (CAL2) *provided with the kit
 - Buffer (wash buffer, elution buffer) *provided with the kit
(For a test run, use sterile distilled water)
 - Lysate (Samples) (For a test run, use sterile distilled water)
(For the preparation of samples, refer to each kit's handbook.)

Refer to "1.5 Setup Parts"(\rightarrow page 8).

3

■ Discharge

(1) Check the QuickGene-610L has been cleaned.


For clean up procedure, refer to "■Tasks required when the device is not used for a week or more"(\rightarrow page 41) and "4.4 Cleaning the Parts"(\rightarrow page 43).

(2) Turn on the power of QuickGene-610L.

"3.2 Turn the Power ON and OFF"(\rightarrow page 24)


(3) Open the front cover.

Take the handle and open the front cover as shown in the figure.

CAUTION: The front cover will not stay on the way to the upper limit. If you release the front cover while opening it, it may cause injury or device damage.

3

Elution Buffer Conical Tube

(4) Check that the Discharge Tray is set on the Discharge Tray Table.

If it is not set, set it.


"■Setup the Discharge Tray"(→ page 18)

(5) Put buffer (elution buffer) into a new Elution Buffer Conical Tube.

The Wash Buffer Bottle provided with the kit can be used as it is.

Important: Buffer (elution buffer) can be used for a day in the tube inside the QuickGene-610L.

(6) Remove the Bottle Holder Cover, set the Wash Buffer Bottle and Elution Buffer Conical Tube (prepared in Step (5)) in the Bottle Holder, and then put caps on both.

Note: Do not touch the inner tubes of the Wash Buffer Cap and Elution Buffer Cap.
In the case of contamination, please refer to "6.4 When the Buffer is Contaminated"(→ page 57).

(7) Set the Bottle Holder Cover.

(8) Close the front cover.

(9) Press the [DISCHARGE] button on the operation panel.

3.4 Preparation of Lysate (Samples)

Prepare lysate (samples).

Make lysate (samples) using the specified kit for each sample.

For details, refer to the kit's handbook.

3.5 Extraction Operations

This section explains each operation before the start of a program and after the completion of a program.

■ Operations before the start


Important: Perform the extraction operation according to the following procedures.

(1) Open the front cover.

Take the handle and open the front cover as shown in the figure.

CAUTION: The front cover will not stay on the way to the upper limit. If you release the front cover while opening it, it may cause injury or device damage.


3


(2) Take the Holder Carriage out of the QuickGene-610L.

Grip the handles, lift the Holder Carriage up, and then take it out.

Note: •When the Carrier is not at the left end and you cannot remove the Holder Carriage, close the front cover and press the [CR MOVE] button to move the Carrier to the left end.
•Hold the handles upright and lift them straight up. Failure to do so may make work difficult. Lift the Holder Carriage straight up and move it forward to remove.


(3) Take the Collection Tube Holder and Cartridge Holder out of the Holder Carriage.


(4) Set the Waste Tubes in the Holder Carriage.

Put the Waste Tubes into the bottom of the Holder Carriage.

(5) Open the cover of the Collection Tube Holder and set the Collection Tubes.


Put the Collection Tube Cap into the groove of the Collection Tube Holder to set it.


Important: Put the Collection Tubes into the Collection Tube Holder in the ascending order of the number printed on the Collection Tube Holder from "1".

(6) Close the cover of the Collection Tube Holder.

Collection Tubes are set in the Collection Tube Holder.

(7) Set the Collection Tube Holder and Cartridge Holder in the Holder Carriage.



- To set the Collection Tube Holder, put its left and right latches into the notches on the Holder Carriage.
- To set the Cartridge Holder, put the left and right tabs of the Holder Carriage into the side notches of the Cartridge Holder.

(8) Check that the Collection Tube Holder and Cartridge Holder are correctly set in the Holder Carriage.

Check the following.


- The left and right long holes at the front of the Holder Carriage can be seen.
- The left and right tabs of Holder Carriage are set in the notches on the left and right sides of the Cartridge Holder.

(9) Release the lock of the Cartridge Holder to open the cover.

Release the left and right fasteners on both sides of the Cartridge Holder and its cover.

3

(10) Set cartridges in the Cartridge Holder.

Set the cartridges firmly so that they are stable in the Cartridge Holder.



Important: •Set cartridges in the ascending order of the number printed on the Cartridge Holder from "1".
•When setting cartridges, be careful that the edge of the cartridges does not touch the Cartridge Holder.

(11) Check the location of the tubes and cartridges.

Check that the cartridges, Collection Tubes, and Waste Tubes are not incorrectly inserted as shown in the figure.

If they are not in order, the buffer and/or lysate may spill, causing a biohazard risk.


(12) Put lysate (samples) into each cartridge.

Be careful not to touch the cartridges and the Cartridge Holder when transferring the lysate (samples) into the cartridges.

Note: If any lysate is left on the edge of the cartridge, wipe it with a Kimwipe® or paper tissue.

(13) Close the Cartridge Holder Cover and lock it.



Lock the left and right fasteners on both sides.

**(14) Attach the Holder Carriage to the QuickGene-610L.**

Hold the handles of the Holder Carriage and attach the Holder Carriage onto the Holder Stand of the QuickGene-610L.

Check the following in order to secure the correct setting.

- The left and right hooks of the Cartridge Holder are put into the grooves on the shaft of the QuickGene-610L.
(There is a little gap between the hook and the groove.)

- The bottom front of the Holder Carriage matches the table front.
- The 2 long holes on the left and right at the front of the Holder Carriage match the prongs of the Table.

Note: When the Carrier is not at the left end and you cannot attach the Holder Carriage, close the front cover and turn on the QuickGene-610L. Press the [CR MOVE] button to move the Carrier to the left end, and then set the Holder Carriage.

(15) Close the front cover.

(16) Select an extraction mode.

Press the [MODE] button several times until the target extraction mode appears.

(17) Press the [START] button.


The QuickGene-610L automatically starts the extraction program.

LCD displays "PROCESSING" during extracting. Extraction progress can be checked by each lamp (LED).

- Note:**
- Do not open the front cover during the extraction process (while PROCESSING is displayed in the LCD). If the front cover is opened, the extraction process will be halted. To continue the process, refer to "3.6 When the Front Cover is Opened During PROCESSING."(→ page 38).
 - To quit the process, refer to "3.7 Operations to Stop Extraction Process"(→ page 39).
 - In some extraction modes, operations are required during processing. For details, refer to "■Optional Operation Procedure [WASH2 program]"(→ page 33).

(18) Two beeps and all the three lamps (BINDING, WASHING, ELUTION) blinking indicates process completion.

The LCD displays the extraction results.

The numerics are cartridge numbers
Number 1 to 3: Extracted successfully.
Number 4: Extraction failed.
Number 5 to 6: There were no cartridges.

Display	Meaning
✓ (check)	Successfully extracted
- (hyphen)	Extraction failure (Cartridge is clogged)
- (underscore)	No cartridge, or No sample

■ Optional Operation Procedure [WASH2 program]



When you chose optional operation mode such as "xxx PLUS" mode (a program added WASH 2 process; see "5 How to Set Parameters"(→ page 47)), QuickGene-610L pauses after first WASH and LCD displays "START SW → RESTART". Then the following steps are continued the program.

(1) Open the front cover.

Open the front cover after LCD displays "START SW → RESTART" with 8 beeps.

CAUTION: The front cover will not stay on the way to the upper limit. If you release the front cover while opening it, it may cause injury or device damage.


3

(2) Take the Holder Carriage out of the QuickGene-610L.

Grip the handles, lift the Holder Carriage up, and then take it out.


Note: Hold the handles upright and lift them straight up. Failure to do so may make work difficult. Lift the Holder Carriage straight up and move it forward to remove.

Dispense a reagent, etc.

(3) Perform the additional operations.

For example, dispense a reagent into the cartridges.


(4) Attach the Holder Carriage to the QuickGene-610L.

Hold the handles of the Holder Carriage and attach the Holder Carriage onto the Holder Stand of the QuickGene-610L.

Check the following in order to secure the correct setting.

- The left and right hooks of the Cartridge Holder are put into the grooves on the shaft of the QuickGene-610L.
There is a little gap between the hook and the groove.

- The bottom front of the Holder Carriage matches the table front.
- The 2 long holes on the left and right at the front of the Holder Carriage match the prongs of the Table.


(5) Close the front cover.

(6) Press the [START] button.


There may be no noise for a while, depending on the extraction mode or the settings, but the LCD should show a "PROCESSING".

■ Operations after the completion

(1) Open the front cover.


CAUTION: Open the front cover as far as it will go, using both hands. The front cover will not stay open until it is fully opened. If you release the front cover while opening it, it may cause injury or damage to the device.


(2) Remove the Collection Tube Holder.

Press down the tab, and remove the Collection Tube Holder.

Note: Grip the handle of the Collection Tube Holder while pressing down the tab of the cover, and then remove the Collection Tube Holder. If the Collection Tube is tilted, the extracted liquid in it may spill out. Lift the Collection Tube Holder upward, and be careful not to tilt it to pull it out.

**(3) Open the Collection Tube Holder Cover and close the Collection Tube Caps. Then, take them out.**


Note: Check the caps are put on the tubes completely.

**(4) Take the Holder Carriage out of the QuickGene-610L.**

Grip the handles, lift the Holder Carriage up, and then take it out.

Note: •When the Carrier is not at the left end and you cannot remove the Holder Carriage, close the front cover and press the [CR MOVE] button to move the Carrier to the left end.
•Hold the handles upright and lift them straight up. Failure to do so may make work difficult. Lift the Holder Carriage straight up and move it forward to remove.

(5) Remove the Cartridge Holder from the Holder Carriage, and dispose of the cartridges.


Release the lock of the Cartridge Holder to open the cover.

One by one pull out the cartridges and dispose of them into a waste container.

Note: When taking cartridges out, be careful so that the edge of the cartridges does not touch the Cartridge Holder.
If the edge of the cartridge touches the Cartridge Holder, wash the Cartridge Holder.

3


(6) Take the Waste Tubes out.

Dispose of the Waste Tubes according to appropriate laws and rules.

Dispose of the infectious wastes according to your biological safety manuals.

(7) Check whether waste liquid remains in the Holder Stand and Discharge Tray.

If waste liquid remains in them, dispose of it.

Note: When the Carrier is not at the left end and you cannot remove the Discharge Tray, close the front cover and press the [CR MOVE] button to move the Carrier to the right end.

(8) Set the Cartridge Holder and Collection Tube Holder in the Holder Carriage, and then install them to the QuickGene-610L.

"■Setup the Holder Carriage"(→ page 16)

(9) Close the front cover.

The operation has been completed.

To continue the process, refer to "3.3 Preparation"(→ page 25).

When finishing the operation, perform the following tasks:

- Turn off the power of QuickGene-610L.
- Dispose of the remaining buffer.
- Clean the dirty parts.

Refer to "4 Maintenance"(→ page 41) for maintenance details.

Important: After QuickGene-610L operation, if it will not be used again for one week or more, put 14ml of ethanol in the Wash Buffer Bottle and 7ml of ethanol in the Elution Buffer Conical Tube, then execute the discharge operation 4 times. After the discharge operation, if ethanol remains in the tubes, discharge it.

For details, refer to "■Tasks required when the device is not used for a week or more"(→ page 41).

When the amount of types of kit increases, the contents of wash buffer and elution buffer will differ from the present kits. For details about the remaining liquid in the tubes, refer to "6.4 When the Buffer is Contaminated"(→ page 57), and discharge it.

3.6 When the Front Cover is Opened During PROCESSING.

When the front cover is opened during processing, the program is halted and the LCD displays "COVER OPEN". If this occurs while the Wash Buffer Pump or Elution Buffer Pump is running, operation of the device will not be paused until the process for one cartridge is completed.

After that, the extraction process can be resumed from the operation panel.

To resume the operation after opening the front cover during extraction operation, perform the following procedures.

To quit the program, refer to "3.7 Operations to Stop Extraction Process"(→ page 39).

Note: When the front cover is opened, do not remove the cartridges or tubes in the QuickGene-610L. Also, do not move the Pressure Head.

Extraction process restarts from the position where the process was paused. Operations such as discharge may be performed even without cartridges or tubes, which may result in liquid spill.

(1) Close the front cover.

(2) Press the [START] button when "START SW → RESTART" is displayed on the LCD.

QuickGene-610L will restart from the position where the program is paused.

Restarting operations vary depending on the timing when the front cover is opened (when the operation is paused).

The following shows the timing when the front cover is opened and the restart operations.

Process status	Operations when opening the front cover	Operations at restart
During initialization	All operations stop	The Pressure Head returns to the original position, then the operations restart from the paused position.
During binding	All operations stop	<ul style="list-style-type: none">The Pressure Head returns to the original position, then moves to the paused position. Operations restart beginning with the cartridge that was in operation.If the operation is paused while applying pressure, pressure application is performed again starting with the applicable cartridge.
During washing	<ul style="list-style-type: none">All operations stop, except for Wash Buffer Pump and Elution Buffer Pump.	<ul style="list-style-type: none">If stopped while applying pressure<ul style="list-style-type: none">The Pressure Head moves to the paused position, and only pressure application restarts.If wash buffer or elution buffer was being discharged, operations restart with the next cartridge.
During elution	<ul style="list-style-type: none">During the wash buffer or elution buffer discharge process, the operations stop after discharge for one cartridge	<ul style="list-style-type: none">If stopped when pressure application was completed or with no pressure application<ul style="list-style-type: none">The Pressure Head moves to the next cartridge of the cartridge being processed, then pressure application and wash buffer discharge processes restart.

3.7 Operations to Stop Extraction Process

To stop the extraction operation, perform the following procedures:

(1) Press the [STOP] button, or open the front cover.

The LCD displays a message as follows:

If the [STOP] button is pressed: "START SW → RESTART"

If the front cover is opened: "COVER OPEN"

(2) In the case when the front cover is opened, close the front cover.

(3) When the LCD displays "START SW → RESTART", press the [RESET] button (not the [START] button) to stop the process.

The mid-process extraction results will be displayed on the LCD.

(4) Perform the operation at extraction completion.

"■Optional Operation Procedure [WASH2 program]"(→ page 33)

4 Maintenance

This chapter describes how to maintain the QuickGene-610L.

When cleaning the QuickGene-610L, wear gloves.

4.1 Before Using the QuickGene-610L

Before using the QuickGene-610L, perform the following tasks:

- Keep clean inside the instrument and parts.
"■Tasks required when the device is not used for a week or more"(→ page 41)
- Replace the Wash Buffer Bottle and Elution Buffer Conical Tube.
- Supply/Replace buffers (wash buffer and elution buffer)
The buffers must be used within one day from when they are set in the QuickGene-610L.
- Discharge
Discharge the old wash buffer, elution buffer, and air remaining in the tubes, and refill with the new wash buffer and elution buffer.

4.2 When the QuickGene-610L is left for a Week or More

4

Even after completing the operation and turning off the power, there are buffers (wash buffer, elution buffer) remaining in the QuickGene-610L. When the QuickGene-610L will not be used for a week or more, the remaining buffers may become crystallized, and cause a clog or damage to the tubes.

Discharge the remaining buffers from the QuickGene-610L and clean the inside with ethanol.

■ Tasks required when the device is not used for a week or more

If the QuickGene-610 will not be used for a week or more after the operation, perform the following procedures:

(1) Put 14ml of ethanol in the Wash Buffer Bottle and 7ml of ethanol in the Elution Buffer Conical Tube, then execute the discharge operation.

For the discharge operation, refer to "■Discharge"(→ page 25).

(2) Repeat the discharge operation 4 times, check the tube is no ethanol remaining.

(3) Dispose of the waste in the Discharge Tray.

(4) Turn off the power.

■ Tasks required when the QuickGene-610L is used again (after one week or more)

When using the QuickGene-610L again, perform the following procedures:

(1) Put 14ml of sterile distilled water in the Wash Buffer Bottle and 7ml of sterile distilled water in the Elution Buffer Conical Tube, then execute the discharge operation.

For the discharge operation, refer to "■Discharge"(→ page 25).

4 Maintenance

- (2) Repeat the discharge operation 4 times, check the tube is no water remaining.
- (3) Dispose of the waste in the Discharge Tray.
- (4) Perform the extraction operation as usual.

4.3 Cleaning the QuickGene-610L


! WARNING: Before starting the maintenance work, turn off the power of QuickGene-610L and disconnect the power plug from the outlet. Otherwise, it may cause an electric shock.

This section explains how to clean the QuickGene-610L.

Note: When cleaning the QuickGene-610L, be sure to wear gloves.

The following describes the parts to clean and cleaning methods:

4

Parts to clean	Cleaning methods
Pressure Head	Wipe with a Kimwipe® or paper tissue soaked with ethanol.
QuickGene-610L (outside)	Remove the dust lightly with a dry, soft cloth.
QuickGene-610L (inside)	Wipe with a Kimwipe® or paper tissue soaked with ethanol or 0.5% sodium hypochlorite solution. Remove the dust lightly with a dry, soft cloth.

4.4 Cleaning the Parts

The following table describes how to clean the parts.

Parts	Cleaning methods
Holder Carriage	Clean the parts using one of the following methods.
Cartridge Holder	- Clean with neutral detergent, and then wash out with distilled water. Next, dip into 3% hydrogen peroxide solution at room temperature for 10 minutes. Lastly, wash the solution out with DEPC water, and dry the parts.
Collection Tube Holder	- Soak in a 0.1N sodium hydroxide/0.1% EDTA solution overnight, and then wash the solution out with DEPC water, and dry the parts.
Bottle Holder	- Clean with neutral detergent, and then wash out with 0.5% sodium hypochlorite solution. Lastly, wash out with distilled water, and dry the parts.
Holder Stand	- Soak in a 0.1N sodium hydroxide/0.1% EDTA solution overnight, and then wash the solution out with DEPC water, and dry the parts. - Perform autoclave sterilization (121° C/20min). Note that if autoclave sterilization (121° C/20min) is performed for parts other than metal parts and packing, the parts may be deformed.
Discharge Tray	After using the QuickGene-610L, dispose of the liquid collected in the Discharge Tray, and wash the tray with water.
Wash Buffer Pump Tube Elution Buffer Pump Tube	Wipe with a Kimwipe® or paper tissue soaked with ethanol.
Wash Buffer Cap Elution Buffer Cap	Wipe with a Kimwipe® or paper tissue soaked with ethanol.
Holder Cover	Remove the dust lightly with a dry, soft cloth.

4.5 How to Replace the Pump Tube

This section explains how to replace the Pump Tube.

If the Pump Tube is continuously used without this maintenance, the Pump Tube may break and the buffer may leak. It is recommended that the Pump Tube be replaced before the required replacement date.

Note: If the QuickGene-610L will not be used for 1 year or more, the Pump Tubes must be removed from the device.

- Wash Buffer Pump Tube

Replacing the Pump Tube is necessary approximately 3 years after the beginning of QuickGene-610L use, or before performing 10,000 extraction operations.

- Elution Buffer Pump Tube

The Pump Tube has a life of about 20,000 extraction operations, which is equal to the instrument life.

Replace the Pump Tube with the spare one included in the package. If all the tubes have been used, purchase them separately as needed. Note that 20,000 of extraction operations is the instrument life. When purchasing them, contact the Customer Support.

■ How to Replace the Pump Tube

Perform the following procedures to replace the Pump Tube.

 **CAUTION:**

When replacing the Pump Tube, wear gloves.

(1) Turn the power on.

Press the [I] side of the switch.

(2) Press the [DISCHARGE] button.

Repeat the discharge operation 4 times, check the tube is no liquid remaining.

Press the [DISCHARGE] button. This operation makes the drive sockets of the Wash Buffer Pump and Elution Buffer Pump rotate to the position for ease to setup Pump Tubes.

(3) Turn off the power.


Press the [O] side of the switch.

(4) Open the front cover.



Hold the handle and open the front cover until the upper limit is reached.

CAUTION: The front cover will not stay on the way to the upper limit. If you release the front cover while opening it, it may cause injury or device damage.


(5) Remove the tube connected to the Wash Buffer Bottle, the tube connected to the Elution Buffer Conical Tube, and the tubes connected to the instrument, from the pumps.

(6) Detach the pump cover.

Loosen and remove the screws to detach the pump cover.

(7) Take off the Pump Tube from the pump, and attach a new Pump Tube to the rear of the Pump Cover.

Slide the connectors at the ends of the Pump Tube into the grooves of the Pump Cover. Check that the notches on both connectors are facing the front of the Pump Cover. Slide the Pump Tube into place inside the Pump Cover by rotating the rotor. Adjust the rotor to its original vertical position. "WASH" or "ELUTION" is printed on the bag of the Pump Tube. Set each type of Pump Tube correctly.

WASH: for Wash Buffer Pump

ELUTION: for Elution Buffer Pump

Note: Use the specified Pump Tube.

After this, the procedure is the same as in "How to install the Pump Tube". Refer to "■Setup the Pump Tube"(→ page 9).

■ How to replace the packing

Perform the following procedures to replace the packing.


4

(1) Open the front cover.

Hold the handle and open the front cover until the upper limit is reached.


CAUTION: The front cover will not stay on the way to the upper limit. If you release the front cover while opening it, it may cause injury or device damage.

(2) Take the Holder Carriage out of the QuickGene-610L.

Grip the handles, lift the Holder Carriage up, and then take it out.


Note: •When the Carrier is not at the left end and you cannot remove the Holder Carriage, close the front cover and press the [CR MOVE] button to move the Carrier to the left end.

•Hold the handles upright and lift them straight up. Failure to do so may make work difficult. Lift the Holder Carriage straight up and move it forward to remove.

- (3) Release the lock of the Cartridge Holder to open the cover. Then raise the tabs of the packing and lift it up.


Raise the tabs of the packing and lift it up.

- (4) Gently squeeze the new packing into the eight latches on the Cartridge Holder Cover to set it.


Make the packing tabs face up (when the cover is open).

4

- (5) Close the Cartridge Holder Cover and lock it.

Fasten the left and right fasteners on both sides to lock the Cartridge Holder Cover.

- (6) Attach the Holder Carriage to the QuickGene-610L.

- (7) Close the front cover.

5 How to Set Parameters

This chapter explains how to set parameters.

"Set Parameters" means to change the values (such as buffer volume, number of discharges, etc.) for each extraction mode.

The parameters set in the QuickGene-610L are stored in 7 separate memory areas.

The name of each extraction mode can be selected from 35 names.

Therefore, even if an extraction mode name is changed, the parameter values will not be changed.

- → "5.2 How to Set/Change Parameters" (→ page 49)
- → "5.4 Changing Extraction Mode Names" (→ page 53)

Important: • The setup is used to change the parameters, according to new protocol of the kit, etc. Please check the parameters for each kit's handbook.
 • In the case of wrong value for parameters, it will cause buffer spill or sample extraction failure. It is the user's responsibility to change the parameters.
 In EXPERT mode, there are some crucial parameters. Therefore, you must not change any parameters in EXPERT mode without our guidance.

5.1 Parameters

Settings of each parameter can be changed from the setup mode.

Some parameters are changed from the normal setup mode (USER mode) and others are changed from the EXPERT mode.

The following table shows the display sequence, setting range, default value, etc. for each parameter in USER mode.

Display sequence	LCD message	Setting range	Unit	Mode name (The value displays the default)						
				DNA WHOLE BLOOD	ISOLATE A	ISOLATE B	ISOLATE C	ISOLATE D	ISOLATE E	ISOLATE F
1	BIND PEAK	30 ~ 160	kpa	120	120	120	120	120	120	120
2	WASH COUNT	0 ~ 5	time(s)	3	3	3	3	3	3	3
3	WASH PEAK	30 ~ 160	kpa	90	90	90	90	90	90	90
4	WASH VOL1	0 ~ 7500	µl	7500	7500	7500	7500	7500	7500	7500
5	WASH VOL2	0 ~ 7500	µl	6500	6500	6500	6500	6500	6500	6500
6	WASH VOL3	0 ~ 7500	µl	5500	5500	5500	5500	5500	5500	5500
7	WASH VOL4	0 ~ 7500	µl	0	0	0	0	0	0	0
8	WASH VOL5	0 ~ 7500	µl	0	0	0	0	0	0	0
9	WASH DIP TM	0 ~ 510	sec	0	0	0	0	0	0	0
10	WAS2 WAIT T	0 ~ 255	min	0	0	0	0	0	0	0
11	WAS2 COUNT	0 ~ 5	time(s)	0	0	0	0	0	0	0
12	WAS2 PEAK	30 ~ 160	kpa	90	90	90	90	90	90	90
13	WAS2 VOL1	0 ~ 7500	µl	7500	7500	7500	7500	7500	7500	7500
14	WAS2 VOL2	0 ~ 7500	µl	6500	6500	6500	6500	6500	6500	6500
15	WAS2 VOL3	0 ~ 7500	µl	5500	5500	5500	5500	5500	5500	5500
16	WAS2 VOL4	0 ~ 7500	µl	0	0	0	0	0	0	0
17	WAS2 VOL5	0 ~ 7500	µl	0	0	0	0	0	0	0
18	ELUT VOL	50 ~ 500	µl	500	500	500	500	500	500	500
19	ELUT PEAK	30 ~ 160	kpa	100	100	100	100	100	100	100
20	ELUT DIP TM	0 ~ 510	sec	0	0	0	0	0	0	0

5 How to Set Parameters

The following table shows the meaning of each parameter.

Display sequence	LCD message	Parameter name	Meaning
1	BIND PEAK: 30 - BIND PEAK: 160	Peak pressure on the sample	Sets the peak pressure during pressure application. Judgment criteria is "the pressure before pressure application + the pressure set in this item". Pressure application may be stopped before reaching the peak pressure, if air pressure drops suddenly, indicating that all liquid has been discharged. Can be specified between 30 - 160 kpa in increments of 1kpa.
2	WASH COUNT: 0 - WASH COUNT: 5	Frequency of wash buffer discharges	Set frequency of wash buffer discharges Can be specified between 0 - 5 times.
3	WASH PEAK: 30 - WASH PEAK: 160	Wash buffer peak pressure	Sets the peak pressure during pressure application. Judgment criteria is "the pressure before pressure application + the pressure set in this item". Pressure application may be stopped before reaching the peak pressure, if air pressure drops suddenly, indicating that all liquid has been discharged. Can be specified between 30 - 160 kpa in increments of 1kpa.
4	WASH VOL1: 0 - WASH VOL1: 7500	1st wash buffer discharge rate	Sets the first wash buffer discharge rate. Can be specified between 0 - 7500μl in increments of 250μl. When 0μl is specified, although the Carrier moves to the wash buffer discharge position, it does not discharge the wash buffer.
5	WASH VOL2: 0 - WASH VOL2: 7500	2nd wash buffer discharge rate	Sets the second wash buffer discharge rate. Can be specified between 0 - 7500μl in increments of 250μl. When 0μl is specified, although the Carrier moves to the wash buffer discharge position, it does not discharge the wash buffer.
6	WASH VOL3: 0 - WASH VOL3: 7500	3rd wash buffer discharge rate	Sets the third wash buffer discharge rate. Can be specified between 0 - 7500μl in increments of 250μl. When 0μl is specified, although the Carrier moves to the wash buffer discharge position, it does not discharge the wash buffer.
7	WASH VOL4: 0 - WASH VOL4: 7500	4th wash buffer discharge rate	Sets the fourth wash buffer discharge rate. Can be specified between 0 - 7500μl in increments of 250μl. When 0μl is specified, although the Carrier moves to the wash buffer discharge position, it does not discharge the wash buffer.
8	WASH VOL5: 0 - WASH VOL5: 7500	5th wash buffer discharge rate	Sets the fifth wash buffer discharge rate. Can be specified between 0 - 7500μl in increments of 250μl. When 0μl is specified, although the Carrier moves to the wash buffer discharge position, it does not discharge the wash buffer.
9	WASH DIP TM: 0 - WASH DIP TM: 510	Wash buffer immersion management timer	Sets the wash buffer immersion management time. Can be specified between 0 - 510 sec. in increments of 2 seconds.
10	WAS2 WAIT T: 0 - WAS2 WAIT T:255	Waiting time	Sets the time to wait for shifting to wash buffer discharges in the second round. Can be specified between 0 - 255 min. in increments of 1 minute.

Display sequence	LCD message	Parameter name	Meaning
11	WAS2 COUNT: 0 - WAS2 COUNT: 5	Frequency of wash buffer discharges in the 2nd round	Sets the frequency of wash buffer discharges in the second round. Can be specified between 0 - 5 times.
12	WAS2 PEAK: 30 - WAS2 PEAK: 160	Wash buffer peak pressure in the 2nd round	Sets the peak pressure during pressure application.Judgment criteria is "the pressure before pressure application + the pressure set in this item". Pressure application may be stopped before reaching the peak pressure, if air pressure drops suddenly, indicating that all liquid has been discharged. Can be specified between 30 - 160 kpa in increments of 1kpa.
13	WAS2 VOL1: 0 - WAS2 VOL1: 7500	1st Wash Buffer Discharge Rate in the 2nd round	Sets the first wash buffer discharge rate in the second round. Can be specified between 0 - 7500μl in increments of 250μl.
14	WAS2 VOL2: 0 - WAS2 VOL2: 7500	2nd WBDR in the 2nd round	Sets the second wash buffer discharge rate in the second round. Can be specified between 0 - 7500μl in increments of 250μl.
15	WAS2 VOL3: 0 - WAS2 VOL3: 7500	3rd WBDR in the 2nd round	Sets the third wash buffer discharge rate in the second round. Can be specified between 0 - 7500μl in increments of 250μl.
16	WAS2 VOL4: 0 - WAS2 VOL4: 7500	4th WBDR in the 2nd round	Sets the fourth wash buffer discharge rate in the second round. Can be specified between 0 - 7500μl in increments of 250μl.
17	WAS2 VOL5: 0 - WAS2 VOL5: 7500	5th WBDR in the 2nd round	Sets the fifth wash buffer discharge rate in the second round. Can be specified between 0 - 7500μl in increments of 250μl.
18	ELUT VOL: 50 - ELUT VOL: 500	Elution buffer discharge rate	Sets the elution buffer discharge rate. Can be specified between 50 - 500μl in increments of 50μl.
19	ELUT PEAK: 30 - ELUT PEAK: 160	Elution buffer peak pressure	Sets the peak pressure during pressure application.Judgment criteria is "the pressure before pressure application + the pressure set in this item". Pressure application may be stopped before reaching the peak pressure, if air pressure drops suddenly, indicating that all liquid has been discharged. Can be specified between 30 - 160 kpa in increments of 1kpa.
20	ELUT DIP TM: 0 - ELUT DIP TM: 510	Elution buffer immersion management timer	Sets the elution buffer immersion management time. Can be specified between 0 - 510 sec. in increments of 2 seconds.

5.2 How to Set/Change Parameters


Settings of each parameter can be changed from the setup mode.

■ How to change Parameters

The following explains how to change settings of each parameter.

(1) Select the extraction mode.

Press the [MODE] button several times until the target extraction mode appears.

(2) Press the [▲] [▼] buttons simultaneously.

When the LCD displays as shown in the figure, the setup mode starts.

(3) Press the [MODE] button several times until the target parameter appears.

To return to the previous parameter, press the [DISCHARGE] button.

(4) Change the parameter settings using the [▲] [▼] buttons.

- [▲] switch: raises the setting value.
- [▼] switch: lowers the setting value.
- Holding one of the above switches down can change the setting values consecutively.


When changing other items, repeat steps (3) and (4).

(5) Press the [START] button then save the changed parameters.

The saved parameters are written into non-volatile memory. The data will not be lost even if the power is turned off.

When the LCD displays as shown in the figure at the left, the changed settings are saved and the setup mode ends.

Important: To cancel the saving parameters, press the [STOP] or [RESET] button.
The setup mode ends without saving the changed parameters.

5.3 Parameters in EXPERT Mode

In EXPERT mode, in addition to the USER mode parameters values of special parameters can also be changed.

Note: In EXPERT mode, there are some crucial parameters. Therefore, you must not change any parameters in EXPERT mode without our guidance.

The following table shows the display sequence, setting range, default value, etc. for each parameter in EXPERT mode.

Dis pla y se qu en ce	LCD message	Setting range	Unit	Mode name (The value displays the default)						
				DNA WHOLE BLOOD	ISOLATE A	ISOLATE B	ISOLATE C	ISOLATE D	ISOLATE E	ISOLATE F
1	BIND COUNT	1 ~ 5	time(s)	1	1	1	1	1	1	1
2	BIND SPEED	1 ~ 12	-	12	12	12	12	12	12	12
3	BIND PEAK	30 ~ 160	kpa	120	120	120	120	120	120	120
4	BIND UP TM	6 ~ 50	sec	15	15	15	15	15	15	15
5	BIND RETRY	30 ~ 160	kpa	160	160	160	160	160	160	160
6	BIND LOWER	3 ~ 160	kpa	90	90	90	90	90	90	90
7	BIND DOWN T	2 ~ 510	sec	26	26	26	26	26	26	26
8	BIND R DN T	2 ~ 510	sec	100	100	100	100	100	100	100
9	BIND FALL	5 ~ 60	kpa	60	60	60	60	60	60	60
10	WASH COUNT	0 ~ 5	time(s)	3	3	3	3	3	3	3
11	WASH SPEED	1 ~ 12	-	12	12	12	12	12	12	12
12	WASH PEAK	30 ~ 160	kpa	90	90	90	90	90	90	90
13	WASH UP TM	1 ~ 50	sec	15	15	15	15	15	15	15
14	WASH RETRY	30 ~ 160	kpa	90	90	90	90	90	90	90
15	WASH LOWER	3 ~ 160	kpa	80	80	80	80	80	80	80
16	WASH DOWN T	2 ~ 510	sec	20	20	20	20	20	20	20
17	WASH R DN T	2 ~ 510	sec	100	100	100	100	100	100	100
18	WASH FALL	5 ~ 60	kpa	40	40	40	40	40	40	40
19	WASH VOL1	0 ~ 7500	µl	7500	7500	7500	7500	7500	7500	7500
20	WASH VOL2	0 ~ 7500	µl	6500	6500	6500	6500	6500	6500	6500
21	WASH VOL3	0 ~ 7500	µl	5500	5500	5500	5500	5500	5500	5500
22	WASH VOL4	0 ~ 7500	µl	0	0	0	0	0	0	0
23	WASH VOL5	0 ~ 7500	µl	0	0	0	0	0	0	0
24	WASH DIP TM	0 ~ 510	sec	0	0	0	0	0	0	0
25	WAS2 WAIT T	0 ~ 255	min	0	0	0	0	0	0	0
26	WAS2 COUNT	0 ~ 5	time(s)	0	0	0	0	0	0	0
27	WAS2 SPEED	1 ~ 12	-	12	12	12	12	12	12	12
28	WAS2 PEAK	30 ~ 160	kpa	90	90	90	90	90	90	90
29	WAS2 UP TM	1 ~ 50	sec	15	15	15	15	15	15	15
30	WAS2 RETRY	30 ~ 160	kpa	90	90	90	90	90	90	90
31	WAS2 LOWER	3 ~ 160	kpa	80	80	80	80	80	80	80
32	WAS2 DOWN T	2 ~ 510	sec	20	20	20	20	20	20	20
33	WAS2 R DN T	2 ~ 510	sec	100	100	100	100	100	100	100

5 How to Set Parameters

Display sequence	LCD message	Setting range	Unit	Mode name (The value displays the default)						
				DNA WHOLE BLOOD	ISOLATE A	ISOLATE B	ISOLATE C	ISOLATE D	ISOLATE E	ISOLATE F
34	WAS2 FALL	5 ~ 60	kpa	40	40	40	40	40	40	40
35	WAS2 VOL1	0 ~ 7500	μl	7500	7500	7500	7500	7500	7500	7500
36	WAS2 VOL2	0 ~ 7500	μl	6500	6500	6500	6500	6500	6500	6500
37	WAS2 VOL3	0 ~ 7500	μl	5500	5500	5500	5500	5500	5500	5500
38	WAS2 VOL4	0 ~ 7500	μl	0	0	0	0	0	0	0
39	WAS2 VOL5	0 ~ 7500	μl	0	0	0	0	0	0	0
40	ELUT VOL	50 ~ 500	μl	500	500	500	500	500	500	500
41	ELUT COUNT	0 ~ 5	time(s)	1	1	1	1	1	1	1
42	ELUT SPEED	1 ~ 12	-	12	12	12	12	12	12	12
43	ELUT PEAK	30 ~ 160	kpa	100	100	100	100	100	100	100
44	ELUT UP TM	1 ~ 255	sec	15	15	15	15	15	15	15
45	ELUT RETRY	30 ~ 160	kpa	100	100	100	100	100	100	100
46	ELUT LOWER	3 ~ 160	kpa	65	65	65	65	65	65	65
47	ELUT DOWN T	2 ~ 510	sec	20	20	20	20	20	20	20
48	ELUT R DN T	2 ~ 510	sec	30	30	30	30	30	30	30
49	ELUT FALL	5 ~ 60	kpa	40	40	40	40	40	40	40
50	ELUT DIP TM	0 ~ 510	sec	0	0	0	0	0	0	0
51	W DS VOL S	5000 ~ 15000	μl	5000	5000	5000	5000	5000	5000	5000
52	E DS VOL S	2200 ~ 6000	μl	2200	2200	2200	2200	2200	2200	2200
53	BUZZER	ON or OFF	-	ON	ON	ON	ON	ON	ON	ON
54	DETECT PRES	1 ~ 50	kpa	15	15	15	15	15	15	15
55	DOWN PRES	10 ~ 50	kpa	50	50	50	50	50	50	50

5

■ EXPERT mode

The following explains how to shift to the EXPERT mode.

(1) Select the extraction mode whose parameter is to be changed.

Press the [MODE] button several times until the target extraction mode appears.

(2) Press the [▲] [▼] buttons simultaneously.

Setup mode starts.

(3) Press the [MODE] or [DISCHARGE] button several times until the following message is displayed on the LCD:

"EXPERT MODE? Y → ▲"

(4) Press the [▲] button.

Change the mode to the EXPERT mode.

Settings of parameters in EXPERT mode can be changed with the same procedures as in the normal setup mode. For details, refer to "■ How to change Parameters" (→ page 49).

Note: In EXPERT mode, there are some crucial parameters. Therefore, you must not change any parameters in EXPERT mode without our guidance.

5.4 Changing Extraction Mode Names

Extraction mode names displayed on the LCD can be changed.

Renaming is performed in RENAME mode.

Important: In RENAME mode, only extraction mode names are changed without changing parameter values of the corresponding modes.

The following names can be used for extraction mode name change:

Display sequence	Extraction mode name	Display sequence	Extraction mode name
1	DNA WHOLE BLOOD (Default)	19	DNA PLANT
2	ISOLATE A (Default)	20	RNA PLANT
3	ISOLATE B (Default)	21	DNA VIRUS
4	ISOLATE C (Default)	22	RNA VIRUS
5	ISOLATE D (Default)	23	SWAB
6	ISOLATE E (Default)	24	STOOL
7	ISOLATE F (Default)	25	PCR PURIFY
8	ISOLATE G	26	GEL EXTRACTION
9	ISOLATE H	27	REACTION CLEANUP
10	RNA WHOLE BLOOD	28	DYE CLEANUP
11	RNA BLOOD PLUS	29	RNA BACTERIA
12	RNA CELL	30	DNA BACTERIA
13	RNA CELL PLUS	31	BACTERIA PLUS
14	DNA TISSUE	32	DNA YEAST
15	RNA TISSUE	33	DNA CULTURE CELL
16	RNA TISSUE PLUS	34	DNA HAIR/NAIL
17	PLASMID	35	DNA BUFFYCOAT
18	PLASMID B		

5

■ How to change extraction mode names

The following explains how to change extraction mode names.

(1) Select the extraction mode whose parameter is to be changed.

Press the [MODE] button several times until the target extraction mode appears on the LCD.

(2) Press the [▲] [▼] buttons simultaneously.

Setup mode starts.

**(3) Press the [MODE] or [DISCHARGE] button several times until the following message is displayed on the LCD:
"RENAME MODE? Y → ▲"**



(4) Press the [▲] button.

RENAME mode starts.

(5) Press the [MODE] or [DISCHARGE] button several times until the extraction mode name to be applied is displayed on the LCD.

(6) Press the [START] button when the mode name to be applied is displayed on the LCD.

After the name before being changed is displayed on the LCD, a message to ask whether to save the new name is displayed.

(7) Press the [▲] button and save the new name.

When the LCD displays as shown in the figure on the left, the changed name is saved and the mode returns to the USER mode.

Important: • When not saving the new name, press the [▼] button.
The LCD display changes to the "RENAME MODE? Y → " without saving the new name.
• To cancel the new name, press the [STOP] or [RESET] button.

The setup mode ends without saving the new name, and the current extraction mode will be displayed.

6 Troubleshooting

This chapter describes the procedures that should be taken when you suspect device failures, and error messages shown on the LCD.

6.1 Extraction Failure

If samples cannot be extracted successfully, check the parameter settings.

If an inappropriate parameter value has been set, it will cause a buffer spill or sample extraction error. Be careful when changing the values. When changing parameters, follow the instructions of the kit's handbook.

For details of parameter setting, refer to "5 How to Set Parameters"(→ page 47).

If the problems cannot be solved, contact your local representative.

6.2 Typical Failure

Check the following items before contacting your local representative.

Phenomenon	Cause	Appropriate action
The QuickGene-610L does not start after turning on the power.	Is the power plug connected to the outlet?	Connect the power plug to the outlet completely.
The circuit breaker tripped.	There may be a current surge or electric leak caused by failure of the QuickGene-610L.	Contact your local representative.
The extraction has stopped in the middle of the operation.	<ul style="list-style-type: none"> •Did you open the front cover during the operation? •It sometimes stops for time adjustment during extraction work. 	<p>During extraction work (While the operation panel displays "PROCESSING"), do not open the front cover.</p> <p>If you opened the front cover during the operation and it stopped the extraction work, refer to "3.6 When the Front Cover is Opened During PROCESSING."(→ page 38).</p>
The extraction work has stopped because the wash buffer and elution buffer will not flow through the nozzle.	Did you perform the discharge operation?	Do not forget to press the [DISCHARGE] button after putting the wash buffer and elution buffer into the bottle and conical tube.
The preprocessed sample that has been dispensed in the cartridge does not steadily drop into the Collection Tube and the Waste Tube.	Is the cartridge inserted correctly into the Cartridge Holder?	<ul style="list-style-type: none"> •Insert the cartridges completely into the Cartridge Holder. •Set the Collection Tube and the Waste Tube corresponding to the cartridge inserting point.
	The position of Collection Tube Holder table may be located in wrong place.	Contact your local representative.
The QuickGene-610L emits an abnormal noise.		Contact your local representative.
The QuickGene-610L emits an abnormal smell.	A failure may be suspected.	Turn off the power at once, and remove the power cable from the outlet. Then, contact your local representative.

6 Troubleshooting

Phenomenon	Cause	Appropriate action
Wash buffer or elution buffer does not flow through the nozzle.	Is the Pump Tube completely attached to the joint?	Connect the Pump Tube completely to the joint.
	Is there buffer in the Wash Buffer Bottle and Elution Buffer Conical Tube?	When the Wash Buffer Bottle and/or Elution Buffer Conical Tube are empty, put the buffer into them.
	Has the tube connecting the device to the pump become unfastened from the device?	Contact your local representative.

6.3 Error Messages

When an error message appears on the LCD, perform the appropriate action written in the table, and contact your local representative if necessary.

The error messages displayed on the LCD are as follows:

Message	Contents	Appropriate action
MT ER CR	Carrier motor or air pump has failed.	Contact your local representative.
MT ER HD	The nozzle head motor has failed.	The message ("CR", "HD", "TB", and "WS") may appear at the same time, depending on failure type.
MT ER TB	Table motor or lift motor has failed.	
MT ER WS	Wash Buffer Pump or Elution Buffer Pump is not normal.	
POWER ALARM	The power voltage status is not normal.	
EEPROM ER	An error has occurred in the EEPROM for the settings.	Contact your local representative.
HISTORY DATA ERR	An error has occurred while writing history information.	
CARRIER POS ERR	Carrier position is not normal.	
NOZZLE POS ERR	Nozzle position is not normal.	
TABLE POS ERR	Table position is not normal.	
LIFT POS ERR	Lift position is not normal.	
WD PUMP ERR	Wash Buffer Pump position is not normal.	Close the front cover after the necessary operation.
CD PUMP ERR	Elution Buffer Pump position is not normal.	
PRESSURE ERR	Pressure sensor is not normal.	
COVER OPEN	The front cover is open.	Press the [RESET] button.
STOP SW PUSHED	The [STOP] button is pressed, and the QuickGene-610L operation is stopped.	

6.4 When the Buffer is Contaminated

When the buffer (wash buffer or elution buffer) is contaminated, wash the tube well following the procedure.

- (1) Put ethanol in the Wash Buffer Bottle and Elution Buffer Conical Tube, and perform the discharge operation 4 times.**

For the discharge operation, refer to "■Discharge"(→ page 25).

- (2) Put sterile distilled water in the Wash Buffer Bottle and Elution Buffer Conical Tube, and repeat the discharge operation 4 times.**

Appendix A

This appendix describes after-sales services and points to note for transfer.

A.1 Usable Power Supply Cords








Model	Input	Connector	Cord	Attachment plug cap
North America <*1> <*2>	100-120 Vac	IEC60320 C-22 Rated 10A, 125Vac UL, CSA Approved 	Type SJT, No.18 AWG Min. 3-Conductors(Single phase;2-current carrying conductors & ground) UL, CSA Approved	NEMA (5-15P) parallel blade Rated 10A, 125Vac UL, CSA Approved 
	200-240 Vac	IEC60320 C-22 Rated 10A, 250Vac UL, CSA Approved 	Type SJT, No.18 AWG Min. 3-Conductors(Single phase; 2-current carrying conductors & ground) UL, CSA Approved	NEMA (6-15P) tandem blade Rated 10 A, 250 Vac UL, CSA Approved 
Europe <*2>	100-240 Vac	IEC60320 C-22 Rated 10A, 250Vac <*1> <HAR>	CENELEC OC 3X1.0 square mm<*1> <HAR>	Rated 10 A, 250 Vac <*1>
Australia	100-240 Vac	IEC60320 C-22 Rated 10A, 250Vac 	Cable: AS OD 3 X1.0 square mm e.g. 	Rated 10 A, 250 Vac 
U.K <*2>	100-240 Vac	IEC60320 C-22 Rated 10A, 250Vac 	BS OC 3 X1.00 square mm	Rated 10 A. 250 Vac 
Others <*1>	Standards compliant with ground conductor and rated 10A or better			

Note: *1. Be sure that the detachable proper Supply cord has the approval of the appropriate safety agencies of the country where the equipment will be used.

Note: *2. Cable length of above Power Supply cord shall be shorter than 4.5 m.

A

A.2 Options

Name: Holder Carriage	Name: Cartridge Holder	Name: Collection Tube Holder
		
Name: Holder Stand	Name: Packing	Name: Discharge Tray
		
Name: Pump Tube	Name: Pump Tube(Purple)	Name: Bottle Holder
		
Name: Elution Buffer Cap(Blue)	Name: Wash Buffer Cap(Purple)	Name: QG610L Carriage Set
		

A

A.3 About Services

- Before you ask for device repairs, refer to "6 Troubleshooting"(→ page 55) to check the problems.
If the problems cannot be solved, contact your local Customer Support.
- If the QuickGene-610L has failed, contact your local Customer Support.

A.4 Repair Service Period of Product

The guaranteed repair service period of this product is seven years from the end of sales. Thereafter, repair services may not be provided if repair parts become out of stock. Please note this.

A.5 Customer Support

<http://lifescience.fujifilm.com>
Fuji Photo Film Co., Ltd. LIFE SCIENCE PRODUCTS DIVISION
26-30, Nishiazabu 2-Chome, Minato-ku, TOKYO 106-8620, JAPAN
Tel:+81-3-3406-2201
Fax:+81-3-3406-2158
E-mail:sginfo@fujifilm.co.jp

Subsidiaries

< United States, Canada, Mexico>
Fujifilm Medical System U.S.A.,Inc.
419 West Avenue, Stamford, CT 06902, U.S.A.
Tel:+1-203-324-2000 ext.6112 (1-800-431-1850 ext. 6112 in the U.S.)
Fax:+1-203-351-4713
E-mail:SSG@fujimed.com
URL: <http://lifescience.fujifilm.com/>

< Europe (excl. UK and Ireland)>
Fuji Photo Film (Europe) GmbH,
Heesenstr. 31, 40549 Dusseldorf, Germany,
Tel: +49-211-5089-174
Fax: +49-211-5089-139
E-mail: lifescience@fujifilmeurope.de
URL: <http://www.fujifilm.de>

< UK, Ireland >
Fuji Photo Film (U.K) Ltd.
Unit 12 St Martins way, St Martins Business centre, Bedford, MK42 OLF, U.K
Tel:+44-1234-245291
Fax:+44-1234-245293
E-mail: lifesciences@fuji.co.uk
URL: <http://lifescience.fujifilm.com/>

< China >
Fuji Photo Film (China) Investment Co., Ltd.
31st floor, Hong Kong New World Tower, No.300 Huai Hai Zhong Road, Shanghai P.R China
Tel:+86-21-3302-4655ext.363
Fax:+86-21-6384-3322
E-mail: wgxiang@fujifilm.com.cn
URL: <http://www.fujifilm.com.cn>

A

Distributors

<Australia, New Zealand>

Berthold AUSTRALIA PTY Ltd.
40 Clements Avenue, BUNDOORA Victoria 3083, Australia
Tel:+61-3-9467-6277 (1-300-300-865 in Australia)
Fax:+61-3-9467-7493
E-mail: rafael@berthold.com.au
URL: <http://berthold.com.au>

<Korea>

Shinki Hi-Tec
GUNWHA Bldg. 7-1, Yangjae, 1-dong, Secho-gu, Seoul, 137-886 Korea
Tel:+82-2-572-1600
Fax:+82-2-572-0058
E-mail: info@skhitec.co.kr
URL: <http://www.skhitec.co.kr>

<Taiwan>

HUNG CHONG CORP.
No.38, Sec. 6, Min Chuan E Road, Taipei, Taiwan
Tel:+886-2-2791-1188
Fax:+886-2-2794-2248
E-mail: fuhsing@mail.hungchong.com.tw
URL: <http://www.FUJIFILM.COM.TW>

A.6 Notes for Transfer QuickGene-610L

When transferring the QuickGene-610L, note the following points:

⚠ WARNING: When carrying the QuickGene-610L, disconnect the power plug from the outlet first. Watch your step when you relocate the QuickGene-610L. If the power cable is damaged, this may cause a device failure, fire, injury, or electric shock.

⚠ CAUTION: Two or more people are necessary to carry the QuickGene-610L. Otherwise, you may hurt your back or hands because of the weight of QuickGene-610L (more than 20kg).

When transferring the QuickGene-610L, be sure to perform the following tasks:

- Remove all attached parts from the QuickGene-610L.
- Secure the tubes inside the QuickGene-610L with tape.
"1.3 Remove the Packaging Materials"(→ page 4)
- Secure the carrier with the 2 clamps so that it does not move.
- Secure the front cover with tape so that it will not open.
"1.3 Remove the Packaging Materials"(→ page 4)
- After completing the above tasks, pack the QuickGene-610L in its original shipping box or some other equivalent box.

**Nucleic Acid Isolation System
QuickGene-610L User's Guide**

P3X1-E380-01EN

Date of issuance: November 2005

Fuji Photo Film Co., Ltd.

Printed in Japan

- The contents of this manual may be updated without notice.
- Fuji Photo Film Co., Ltd. assumes no liability for damages to third party copyrights or other rights arising from the use of any information in this manual.
- Contents of this manual are not to be reproduced without permission from Fuji Photo Film Co., Ltd.
- Manuals with missing or wrongly collated pages will be replaced free of charge.

FUJI PHOTO FILM CO., LTD.

Printed on 100% Recycled Paper

P3X1-E380-01EN

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Tous droits réservés, Copyright © Fuji Photo Film Co., Ltd. 2005
Alle Rechte vorbehalten, Copyright © 2005 Fuji Photo Film Co., Ltd.