



Identification of Pyridylethylated Cysteine Residues in Peptide Maps Using Diode-array Detection

Biopharmaceutical

Cysteine and cystine residues are important amino acid residues because they are involved in the correct formation of the tertiary and quaternary structure of peptides and proteins. The full biological activity relies on correctly synthesized intra- and intermolecular disulfide bridges between two particular cysteine residues. As a consequence of disulfide bridge formation and hence a globular structure many proteins are hardly susceptible to proteolytic attack. Therefore, the protein chemical characterization of these proteins is often not possible without prior modification of the cysteine/cystine residues.




Figure 1
Tryptic digest of a bacterial neurotoxin protein after
reduction and alkylation with vinylpyridine



Agilent Technologies
Innovating the HP Way

HPLC Performance


Typically, reduction with a variety of reagents, such as mercapto-ethanol or dithiodithreitol is carried out followed by alkylation using iodoacetic acid, iodoacetamide or vinylpyridine. Afterwards peptides are separated by HPLC and identified by Edman sequencing or MS-MS sequencing.

Prior to sequence analysis peptide maps are routinely separated by microbore or capillary HPLC. Identification of modified cysteine residues is best performed using vinylpyridine as alkylation reagent, since it adds additional spectral characteristics to this residue. As shown in figure 1, using diode-array detection vinylpyridinilated cysteine residues can be unambiguously identified in peptides separated by reversed-phase capillary HPLC by their characteristic absorption maximum at 254 nm. They can thus be easily distinguished from tyrosine and tryptophan residues which have an absorption maximum around 280 nm, even if they are present in the same peptide. With this approach collected peptide fractions can be easily analyzed further by protein sequencing or alternatively online by MS-MS sequence analysis.

Equipment

Agilent 1100 Series

- Binary pump
- Vacuum degasser
- Standard autosampler
- Diode-array detector
micro flow cell
10-mm path length, 500 nl
- Agilent ChemStation
+ 3D software
LC Packings flow splitter,
IC-100-VAR precolumn
microflow processor*



**LC Packings,
Baarsjesweg 154
1057 Amsterdam,
The Netherlands*

*LC Packings
80 Carolina Street,
San Francisco, CA 94103, USA*

© Copyright 1999 Agilent Technologies
Released 03/99
Publication Number 5968-4739E



Agilent Technologies

Innovating the HP Way