
Operating Instructions

GREAT ESCAPE
Foot and Back massager
EP518

Panasonic

Before operating this unit, please read these instructions completely.

INDEX

English

See pages 1 - 6

INDICE

Español

Vea páginas 7 - 12

IMPORTANT SAFETY INSTRUCTIONS

When using an electrical appliance, basic precautions should always be followed, including the following

Read all instructions before using the foot massager


DANGER - To reduce the risk of electric shock:

- 1 Always unplug this appliance from the electrical outlet immediately after using and before cleaning

WARNING - To reduce the risk of burns, fire, electric shock, or injury to persons:

- 1 An appliance should never be left unattended when plugged in Unplug from outlet when not in use, and before putting on or taking off parts
- 2 Do not operate under blanket or pillow Excessive heating can occur and cause fire, electric shock, or injury to person
- 3 Close supervision is necessary when this appliance is used by invalids or disabled persons
- 4 Use this appliance only for its intended use as described in this manual
Do not use attachments not recommended by the manufacturer
- 5 Do not use on or near children
- 6 Never operate this appliance if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water Call your nearest service center for repair
- 7 Do not carry this appliance by supply cord or use cord as a handle
- 8 Keep the cord away from heated surface
- 9 Never operate the appliance with the air openings blocked Keep the air opening free of lint, hair and the like
- 10 Never drop or inset any object into any opening
- 11 Do not use outdoors
- 12 Do not operate where aerosol (spray) products are being used or where oxygen is being administered
- 13 To disconnect, turn all controls to the off position, then remove plug from outlet
- 14 Keep children away from massager
- 15 Do not stand on appliance
- 16 To reduce the risk of electric shock, this appliance has a polarized plug (One blade is wider than the other) This plug will fit in a polarized outlet only one way If the plug does not fit fully in the outlet, reverse the plug If it still does not fit, contact a qualified electrician to install the proper outlet Do not change the plug in any way

SAVE THESE INSTRUCTIONS

Parts identification

- | | |
|-------------------------------|----------------------------|
| ① Tapping areas | ⑥ Intensity control button |
| ② Handle | ⑦ Power button |
| ③ Speed indicator lamps | ⑧ Control panel |
| ④ Auto program control button | ⑨ Power cord |
| ⑤ Speed control button | |

How to use


[Pre-programmed Massage]

Five professional massage techniques are repeated in sequence

1. Press OPERATION button


2. Press AUTO-PROGRAM button

- If you do not select the pre-programmed massage, you can manually select various types of massage, from slow tapping to vibration, with the speed control
- The unit simulates the 5 key massage techniques illustrated below

[Manual Massage]


1. Make sure the switch is in the OFF position and the main body is secured before plugging into the power source.
2. Use the VIBRATING function first to massage all over the soles of the feet (approx 1 minute).
3. Use the TAPPING function to concentrate the massage in the arches of the feet (approx 1 minute).
4. Use the VIBRATING function as the final massage, centering on the areas which are stiff or tense (approx 1 minute).
5. Turn switch OFF and unplug unit.

[Note]

- To prevent excessive massaging, there is a built-in timer which after 15 minutes of operation will automatically stop the motion of tapping area.


Speed Display

Examples of use


1. (Soles of Feet)

Place the massager securely on the floor with the control panel toward you, sit in a chair, and place your feet on top of the massager. Start by massaging all over the soles of the feet.


2. (Calf)

Place the massager on a bed or sofa, or other flat surface and lie down with your calves on the massager.

3. (Waist)

Place the massager on a sofa or chair with the control panel down, facing the buttons toward you. Sit down with your lower back against the massager.

Massage Effects on the Soles of Feet

There are two major types of massagers, those designed to relax the nerves, and those designed to relax the muscles. A typical massager, designed for nerves, will rub lightly and is generally equipped with a rolling mechanism. Massaging to stimulate muscles and improve circulation requires a more aggressive design using tapping and vibration.

	Massage Illustration	Effect on body
Tapping		Stimulates muscles and promotes circulation of blood and lymph fluids
Vibrating		Stimulates the nerves to relax muscles

Cautions - To avoid possible injury:

- Do not use for those who are pregnant or menstruating or persons with an illness
- Do not use for any afflicted area of your body
- Do not use on or near children
- Do not use for massaging the head, abdomen, elbows or knees
- If you feel any discomfort during use, discontinue use immediately and consult your physician
- Do not allow children to play on or around the massager
- After use, turn the power button to off
- Massaging should be limited to no more than 15 minutes. Do not massage any one point for longer than 5 minutes. Excessive massage can overstimulate the muscles and nerves and result in an adverse effect.
- To lessen the stimulation of the massager for those for whom the effect is too strong, place towels over the tapping areas
- Turn off the massager if you become sleepy during use
- Avoid using sharp or pointed objects on the massager
- Be careful not to drop lit cigarettes, ashes or matches on the massager
- Keep the massager away from direct sunlight, which can cause discoloration or fading
- Handle the massager with care. Do not stand on it, or handle it roughly
- If the massager fails to function or functions abnormally, stop use immediately and have it repaired by an authorized service center. Do not attempt to repair yourself
- Do not use to massage your feet with your shoes on
- For your safety and to prevent damage to the massager, do not poke toes or heels into the tapping area

Specifications

Power source	AC120V, 60Hz
Power consumption	20W
Speed	Approx 240-3,600 cycles/min in 5 stages, alternating tapping cycles
Weight	2.5kg

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