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[3] Ogata, Seiji; Morishige, Ken-Ichirou; Sawada, Kenjiro. Fasudil Inhibits Lysophosphatidic Acid-Induced Invasiveness of Human Ovarian Cancer Cells. International Journal of Gynecological Cancer, 2009, 19(9): 1473-1480

[4] Huaying Chen, Yang Lin, Min Han. Simultaneous quantitative analysis of fasudil and its active metabolite in human plasma by liquid chromatography electrospray tandem mass spectrometry. Journal of Pharmaceutical and Biomedical Analysis, 2010, 52(2): 242-248

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Compound	Added/(mg/L)	Found/(mg/L)	Recovery/%	RSD/%
Homopiperazine	0.5	1.373 ± 0.031	97.8	0.65
	1.0	1.896 ± 0.028	101.2	0.72
	2.0	2.850 ± 0.043	98.3	0.96

www.thermoScientific.com

©2013 Thermo Fisher Scientific Inc. All rights reserved. ISO is a trademark of the International Standards Organization. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is presented as an example of the capabilities of Thermo Fisher Scientific Inc. products. It is not intended to encourage use of these products in any manners that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

Thermo Fisher
Scientific, Austin, TX
USA is ISO 9001:2008 Certified.

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