FINO UNISEAL NANO
Versiegelungslack

REF 03036

CE 0483

DEU Gebrauchsinformation

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ENG User Information

1. Fields of application

FINO UNISEAL NANO is a nano-particle-reinforced light-curing varnish for the sealing of dentures, temporary crowns and individual trays. The advantages of FINO UNISEAL NANO are the low shrinkage during light-curing and the resulting scratch-resistant surface due to the present nano particles. This is a matrix enriched with ceramic fillers in the nano-range. FINO UNISEAL NANO seals the surface and thus protects against discolorations and the deposition of plaque.

2. Safety precautions

The product must only be used for the intended purpose described in the scope of application. FINO UNISEAL NANO contains methyl-methacrylate and is combustible!
Further safety-relevant data is given in the safety data sheet.

3. Processing instructions

Apply FINO UNISEAL NANO evenly with the brush in the lid. Take care that the resin surface is clean and dry.
Due to the slightly watery consistency FINO UNISEAL NANO is easy to apply.
The drying time of the wetted surface is approx. 30 seconds. Subsequently follows light-curing for 5 minutes in a common light-curing unit (FINO-LUX).

FINO UNISEAL NANO
Sealing varnish

REF 03036

CE 0483

ENG User Information

4. Storage

P411 Store at temperatures not exceeding 25 °C. P102 Keep out of reach of children.
 Avoid heat and sunlight!
 Close container tightly immediately after use.

5. Physical data

Boiling point:	100 °C;
Temperature of ignition:	430 °C

5.1 Composition

Contains multi-functional methyl-methacrylates, acrylates, nano-particles, photo-initiators, stabilisers

6. Delivery forms

FINO UNISEAL NANO sealing varnish
20 ml 03036

7. Guarantee

Our technical recommendations of application are based on our own experiences and tests and should only be regarded as guidelines. It rests with the skills and experience of the user to verify that the products supplied by us are suitable for the intended procedures. Our products are undergoing a continuous further development. We reserve the right of changes in construction and composition. It is understood that we guarantee the impeccable quality of our products.

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