FINOSIL PREMIUM

22/
F 22

Dubliersilikon

REF 14207/14209/14212

DEU Gebrauchsinformation

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REF 14207/14209/14212

ENG User Information

Please read carefully before use
Keep for referral purposes

1. Intended purpose

Addition-crosslinked duplicating silicone for the dental duplicating technique. Precise in detail, shrink-free and stable in storage.

2. Indication

Suitable for all types of duplication and virtually any material, e.g. plaster and metal.

2.1 Contraindication

Not known.

3. Fields of application

For processing in the dental laboratory according to the intended use.

4. Safety precautions

For safety-relevant data please see the safety data sheet.

5. Processing instructions

5.1 Preparation

Do not water model. Carefully remove contamination from the surfaces (separating agent, instant adhesives, spacer varnish etc.)

5.2 Dosing

Volume ratio component A : component B 1 : 1

5.3 Mixing

5.3.1 Manually

Place components A and B in mixing cup at a ratio of 1 : 1. Stir well with a spatula until a homogenous mixture (colour) has been achieved.

5.3.2 Automatically

Rational and clean dosage, mixing and batching in a duplicating and mixing device (e.g. FINODOSIL duplicating and mixing device, item no. 14290).

5.4 Duplicating of model

5.4.1 Manually

Pour duplicating silicone in a thin stream into the flask or duplicating form.

5.4.2 Automatically

Strength of stream is defined by the duplicating and mixing device, distance to flask or duplicating form approx. 20 cm. Optimal results are achieved through pressure duplication in a pressure pot. Allow duplicating silicone to cure (vulcanisation) until complete recovery is achieved for approx. 35 min. and devest model.

5.5 Pour

The negative mould is highly elastic, permanently dimensionally stable and can be cast several times.

Spray the negative mould with silicone wetting agent (e.g. FINO SPANN-EX silicone wetting agent, item no. 14225) for wetting the surface and preventing the formation of bubbles.

All dental plasters, modelling resins and investment materials are suitable for casting. For pressure duplication, fabricate model at same pressure as negative mould.

All dental plasters, modelling resins and investment materials are suitable for casting. For pressure duplication, fabricate model at same pressure as negative mould.

6. Storage

P404 Store in a closed container.

P411 Store at temperatures not exceeding 25 °C.

7. Physical data

FINOSIL PREMIUM 22 duplicating silicone normal-curing:

Mixing time manually approx. 1 min.; mixing time under vacuum 30-45 s; processing time 6 min.; demoulding after 30 min.; final hardness 22 Shore A; mixing ratio A : B = 1 : 1; A = transpa-white, B = blue.

FINOSIL PREMIUM F 22 duplicating silicone fast-curing:

Mixing time manually approx. 1 min.; mixing time under vacuum 30-45 s; processing time 2 min.; demoulding after 10 min.; final hardness 22 Shore A; mixing ratio A : B = 1 : 1; A = white, B = pink.

FINOSIL PREMIUM 22/F 22 duplicating silicone:

Viscosity 3000 mPa s; density 1075 g/cm³; tensile strength 4,0 N/mm²; elongation at break 500 %; Elmendorf tear 4,5 N/mm²; linear measure change 0,2 %.

7.1 Composition

Vinyl-polysiloxane-based addition curing duplicating silicone.

8. Delivery forms

FINOSIL PREMIUM 22 Duplicating
Silicone

930 ml each 14207

FINOSIL PREMIUM F 22 Duplicating
Silicone

930 ml each 14209
2 x 5172 ml 14212

9. Guarantee

Our technical recommendations of application are based on our own experiences and tests and should only be regarded as guidelines. It rests with the skills and experience of the user to verify that the products supplied by us are suitable for the intended procedures. Our products are undergoing a continuous further development. We reserve the right of changes in construction and composition. It is understood that we guarantee the impeccable quality of our products.

10. Disposal

P501 Dispose of contents/container in accordance with local/regional/national/ international regulations.

10.1 EWC waste code product

070216 WASTES FROM ORGANIC CHEMICAL PROCESSES; wastes from the MFSU of plastics, synthetic rubber and man-made fibres; waste containing hazardous silicones; hazardous waste.

11. Symbols explanation

11.1 Manufacturing

REF Item numbers

LOT Lot code

 Shelf life

 non-sterile

 Unit in the packaging

 Manufacturer

 Date of manufacture

11.2 Application

Please read carefully before use
Keep for referral purposes

Only to be used by qualified
personnel!

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