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CONTINUATION

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Data of the manual is equally had newest product as evidence, because of improvement or other change, description of manual may differ from Practical situation. Our company will reserve the right corrected at any moment, it is difficult to list the corrected place one by one.

Copyright of manual belongs to Wenzhou Sundoo Instrument Co., Ltd
Any unit and individual can't copy or plagiarize the part or full content .
Of manual without written consent of Wenzhou Sundoo Instrument Co., Ltd.

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MANUFACTURER: WENZHOU SUNDOO INSTRUMENTS CO.,LTD

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Add: Building E, No.10 High-tech Industrial Zone Wenzhou China

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Thank you for your purchasing SVM-211 Digital Biologic Microscope.
SVM-211 Digital Biologic Microscope Measuring System includes Video digital biology microscope and professional image measuring software. It is integrated traditional optic and digital video technology. The electric magnifications is 280x-2800x and the optic magnifications is 40X-400X. The image is very clear. You can capture the image very quickly and exactly by the professional measuring software. It has the image measurement, analysis, printing and so on.
The microscope is applied in Biology, Medicine, Bio-chemistry, Cosmetology, and so on.

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output connector, insert another AV terminal into the connector on the yellow MMI cable.

2. Monitor, Computer, Video Card, Measuring Software The Measuring Software can provide measuring function, that can measure width, radii, diameter, distance between two lines and angle.

Note: Above accessories should be purchased additionally.

Maintenance

1. When the microscope is not used

Please fix the protector of the lens and cover the dustproof cover. You are advised that stock the lens and base of lens into damp proof container, and place drier and prevent them to be mould.

2. When the microscope is cleaned

Don't touch the optical component with hands. Use soft cotton or brush to clean the lens. You can clean the obstinate dirt as fingerprint, grease and so on with the cotton or lens tissue, which tipped alcohol.

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dip in pure alcohol, on the bulb. Otherwise, the bulb should easy to crack or reduce the Brightness.

5. Delivery and Stock.

(a) The power must be cut off, before remove the microscope.

(b) Operation Environment

Temperature: 0℃-40℃

Relative Humidity: 85%

Operation

1. Take the microscope out of the package, it is available to be used.

2. Check the power supply

The rate input voltage is signed on the pedestal of microscope. The input voltage is according with the rate input voltage. Please connect the power supply as the voltage of power supply is same as the rate input voltage.

That will damage the microscope; the microscope is operated with the rate voltage. 2) The socket must be with earth.

3. Turn the diaphragm board to adjust the brightness.

4. Turn the focus adjusting knob to adjust the focus, turn the knob 12 to adjust focus rigor

5. Correct the screw of CCD to the CCD connector, then fasten it. You can use the Photographic system.

6. Change the light bulb

Please cut off the power before change bulb, prevent electric shock or damage the microscope. Take off the cover of the slant illumination after the bulb be cold. Pull out of the halogen bulb then change it.

Optional Accessories

1. AV/VGA, XGA-transfer card + CRT monitor

The XGA-transfer card can provide hi-definition image, display on CRT, without flash and with bright and flamboyant picture.

Connecting Step:

(a) Attach the power, Connect the CRT's signal plus with the transfer card's connector.

(b) Inert MMI signal provide with the instrument into the reciprocal connector on the transfer card, attach Ac cable with Photographic lens'

Specification

1. Parameters(view Page 11).
2. Rough adjusting range: 17mm.
3. Nice adjusting range :2mm.
4. Input voltage:220v
Transparent Illumination:Halogen bulb.

Cautions

1. Do not disassemble the microscope.
To reduce the risk of electric shock, do not remove the parts, which are not described to be removed in the manual. If the microscope is in trouble, refer servicing of quality service personnel.
2. Check the Power
 - (a)The rate input voltage is signed on the pedestal of microscope.
The input voltage is according with the rate input voltage. If the voltage of power supplied is not accord with the rate voltage, please ask the experience personnel to open the cover on the pedestal, turn The switch, which have a symbol 110V~220V. Select the voltage same as the power supply, then fix the cover on the pedestal. Or refer Service with the service personnel. That will damage the microscope; the microscope is operated with the rate voltage.
 - (b) The socket must be with earth.
3. Prevent the Fire or Scald
 - (a) During the bulb provides illumination, it is in high temperature. Pay more attention to do not touch the bulb to prevent be scalded.
 - (b) When you change the bulb, please wait for the bulb is cold.
 - (c) To prevent fire, do not take tinder close to the bulb.
4. Change the Light Bulb
 - (a) Please cut off the power before change bulb, prevent electric shock or Damage the microscope.
 - (b) The halogen bulb or fluorescent lamp is only applied with
 - (c) When you change the bulb, please wear gloves or pack the bulb with cover. Do not touch the glass part of the bulb with hand. After change the bulb, please sweep the finger mark or dirty with clean soft cloth

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Function

SVM-211 Digital Biology Microscope Measuring System includes Video Digital Microscope and professional image measuring software. It is integrated Traditional Optic and Digital Video Technology. The microscope is Inspected with not only Optica inspection but also connect with video device that provide the image for inspection of many persons in same time. The software can be used in catching the image, saving the image, Measuring sample, Printing data, and so on. The microscope is applied with biology, medicine, bio-chemistry research, and so number SONY


Hi-definition CCD Camera.

Large Stage Version.

Background Illumination.

Video Signal Output, Adapter Connected with PC Monitor.

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