APPLICATION NOTE No. 476

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[2] Yoon YJ, Kim OY & Gho YS. Extracellular vesicles as emerging intercellular communicasomes. BMB Reports vol. 47 531–539 Preprint at https://doi.org/10.5483/BMBRep.2014.47.10.164 (2014).

[3] Tetta C, Ghigo E, Silengo L, Deregibus MC & Camussi G. Extracellular vesicles as an emerging mechanism of cell-to-cell communication. Endocrine 44, 11–19 (2013).

[4] Iraci N, Leonardi T, Gessler F, Vega B & Pluchino S. Focus on extracellular vesicles: Physiological role and signalling properties of extracellular membrane vesicles. Int J Mol Sci 17, (2016).

[5] Yuana Y, Sturk A. & Nieuwland R. Extracellular vesicles in physiological and pathological conditions. Blood Rev 27, 31–39 (2013).

[6] Robbins P D & Morelli AE Regulation of immune responses by extracellular vesicles. Nat Rev Immunol 14, 195–208 (2014).

[7] Raposo G & Stoorvogel W. Extracellular vesicles: Exosomes, microvesicles, and friends. Journal of Cell Biology vol. 200 373–383 Preprint at https://doi.org/10.1083/jcb.201211138 (2013).

[8] Tamkovich SN., Tutanov OS & Laktionov PP. Exosomes: Generation, structure, transport, biological activity, and diagnostic application. Biochemistry (Moscow) Supplement Series A: Membrane and Cell Biology 2016 10:3 10, 163–173 (2016).

[9] Witwer KW et al. Updating MISEV: Evolving the minimal requirements for studies of extracellular vesicles. J Extracell Vesicles 10, (2021).

[10] Stahl PD & Raposo G. Extracellular Vesicles: Exosomes and Microvesicles, Integrators of Homeostasis. Physiology (Bethesda) 34, 169–177 (2019).

[11] Akers JC, Gonda D, Kim R, Carter BS & Chen CC. Biogenesis of extracellular vesicles (EV): exosomes, microvesicles, retrovirus-like vesicles, and apoptotic bodies. J Neurooncol 113, 1 (2013).

[12] Pan BT & Johnstone RM. Fate of the transferrin receptor during maturation of sheep reticulocytes in vitro: selective externalization of the receptor. Cell 33, 967–978 (1983).

[13] Qiu P, Zhou J, Zhang J, Dong Y & Liu Y. Exosome: The Regulator of the Immune System in Sepsis. Front Pharmacol 12, 880 (2021).

[14] Liu W. et al. Role of Exosomes in Central Nervous System Diseases. Front Mol Neurosci 12, 240 (2019).

[15] Liu J et al. The biology, function, and applications of exosomes in cancer. Acta Pharm Sin B 11, 2783–2797 (2021).

[16] Dai J et al. Exosomes: key players in cancer and potential therapeutic strategy. Signal Transduction and Targeted Therapy 2020 5:1 5, 1–10 (2020).

[17] Hoen EN, Cremer T, Gallo RC & Margolis LB Extracellular vesicles and viruses: Are they close relatives? Proc Natl Acad Sci U S A 113, 9155–9161 (2016).

[18] Lin Y, Anderson JD, Rahnama LMA, Gu Sv & Knowlton AA. Exosomes in disease and regeneration: biological functions, diagnostics, and beneficial effects. Am J Physiol Heart Circ Physiol 319, H1162–H1180 (2020).

[19] Huda MN et al. Potential Use of Exosomes as Diagnostic Biomarkers and in Targeted Drug Delivery: Progress in Clinical and Preclinical Applications. ACS Biomater Sci Eng 7, 2106–2149 (2021).

[20] Xu L, Wu LF & Deng FY. Exosome: An Emerging Source of Biomarkers for Human Diseases. Curr Mol Med 19, 387–394 (2019).

[21] Muthu S, Bapat A, Jain R, Jeyaraman N & Jeyaraman M. Exosomal therapy—a new frontier in regenerative medicine. Stem Cell Investig 8, (2021).

[22] Tejerina S, Dufey V, Hoet JF, Tacheny A & De Longueville F. The DASbox® Mini Bioreactor System as a Tool for Process Development And Stem-Cell Derived Exosome Production in Standardized Culture Conditions. (2021).

[23] Escobar Ivirico JL, Suttle A & Sha M. SciVario® twin Self-Scale-up of CHO Culture-Based Antibody Production from BioBLU ® 3c to BioBLU ® 50c Single-Use Bioreactor. www.eppendorf.group/sci-vario (2022).

[24] Escobar Ivirico JL & Sha M. Stem Cell Exosome Production on the SciVario® twin, a Flexible Controller for Your Bioprocess Needs. (2020).
[25] Ludwig N, Whiteside TL & Reichert TE. Challenges in Exosome Isolation and Analysis in Health and Disease. Int J Mol Sci 20, (2019).
[26] Sidhom K, Obi PO & Saleem AA. Review of Exosomal Isolation Methods: Is Size Exclusion Chromatography the Best Option? Int J Mol Sci 21, 1–19 (2020).
[27] Tauro BJ et al. Comparison of ultracentrifugation, density gradient separation, and immunoaffinity capture methods for isolating human colon cancer cell line LIM1863-derived exosomes. Methods 56, 293–304 (2012).
[28] Szatanek R, Baran J, Siedlar M & Baj-Krzyworzeka M. Isolation of extracellular vesicles: Determining the correct approach (review). International Journal of Molecular Medicine vol. 36 11–17 Preprint at https://doi.org/10.3892/ijmm.2015.2194 (2015).
[29] Livshts MA et al. Isolation of exosomes by differential centrifugation: Theoretical analysis of a commonly used protocol. Sci Rep 5, (2015).
[30] Chen J et al. Review on Strategies and Technologies for Exosome Isolation and Purification. Frontiers in Bioengineering and Biotechnology vol. 9 Preprint at https://doi.org/10.3389/fbioe.2021.811971 (2022).
[31] Kotmakçı M & Akbaba GE. Exosome Isolation: Is There an Optimal Method with Regard to Diagnosis or Treatment? in Novel Implications of Exosomes in Diagnosis and Treatment of Cancer and Infectious Diseases (InTech, 2017). doi:10.5772/intechopen.69407.
[32] Gurunathan S, Kang MH, Jeyaraj M, Qasim M & Kim JH. Review of the Isolation, Characterization, Biological Function, and Multifarious Therapeutic Approaches of Exosomes. Cells 8, (2019).
[33] Momen-Heravi, F et al. Current methods for the isolation of extracellular vesicles. Biological Chemistry vol. 394 1253–1262 Preprint at https://doi.org/10.1515/hsz-2013-0141 (2013).
[34] Patel GK et al. Comparative analysis of exosome isolation methods using culture supernatant for optimum yield, purity and downstream applications. Scientific Reports 2019 9:1 9, 1–10 (2019).
[35] Gupta S et al. An improvised one-step sucrose cushion ultracentrifugation method for exosome isolation from culture supernatants of mesenchymal stem cells. Stem Cell Res Ther 9, 1–11 (2018).
[36] Lane RE, Korbie D, Anderson W, Vaidyanathan R & Trau M. Analysis of exosome purification methods using a model liposome system and tunable-resistive pulse sensing. Sci Rep 5, (2015).

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*Inquire the part number for your country

www.eppendorf.cn

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Beckman Coulter® and Vi-CELL® are registered trademarks of Beckman Coulter, Inc., USA. BioCoat®, Synthemax® and Corning® are registerd trademarks of Corning Inc., USA. RoosterBio® is a registered trademark of RoostrBio Inc., USA. Zetasizer® is a registered trademark of Malvern Panalytical Ltd., UK. Thermo Fisher Scientific® is a registered trademark of Thermo Fisher Scientific Inc., USA. Eppendorf®, the Eppendorf Brand Design, BioBLU®, and ThermoMixer® are registered trademarks of Eppendorf SE, Germany. DASbox® and DASware® is a registered trademark of DASGIP Information and Process Technology GmbH, Germany.
All rights reserved, including graphics and images. Copyright © 2023 by Eppendorf SE.

Eppendorf SE reserves the right to modify its products and services at any time. This application note is subject to change without notice. Although prepared to ensure accuracy, Eppendorf SE assumes no liability for errors, or for any damages resulting from the application or use of this information. Viewing the application note alone cannot as such provide for or replace reading and respecting the current version of the operating manual.

AA476-020-02-082023