



Propaq CS®



Ambulatory Care Reference Guide

Table of Contents

Propaq CS Part Numbers, Model Numbers, Serial Numbers	3
Propaq CS Ship-with Accessories.....	9
Propaq CS Accessories Sold Separately	10
Propaq CS ECG Cables In Detail.....	17
YSI Temperature Notes.....	20
Oridion CO2 Accessories	21
Other Vendors for Accessories	22

NOTE: This reference guide is intended for the use of representatives from Welch Allyn and our ambulatory care distributors. This document in full is not intended to give to customers. However, information within this document may be shared verbally or in written form with clients and prospective customers.

Overview: Propaq CS Model Numbers and Part Numbers by Option

Propaq CS Part No.	Model 242 - No Invasive Lines	Model 244 - 1 Invasive Line	Model 246 - 2 Invasive Lines	Nellcor SpO2	Masimo SpO2	Printer	No Printer	Side stream etCO2 (SS)	Main stream etCO2 (MS) *	Both etCO2 (SS & MS) *	RS 423/ Acuity Port
9001-001288	x			x			x				x
9001-008152	x				x		x				x
9001-004720	x			x		x					x
9001-011584	x				x	x					x
9001-002146	x			x			x	x			x
9001-003004	x			x			x		x		x
9001-003862	x			x			x			x	x
9001-009010	x				x		x	x			x
9001-009868	x				x		x		x		x
9001-010726	x				x		x			x	x
9001-005578	x			x		x		x			x
9001-006436	x			x		x			x		x
9001-007294	x			x		x				x	x
9001-012442	x				x	x		x			x
9001-013300	x				x	x			x		x
9001-014158	x				x	x				x	x
9001-015016		x		x			x				x
9001-021880		x			x		x				x
9001-018448		x		x		x					x
9001-025312		x			x	x					x
9001-015874		x		x			x	x			x
9001-016732		x		x			x		x		x
9001-017590		x		x			x			x	x
9001-022738		x			x		x	x			x
9001-023596		x			x		x		x		x
9001-024454		x			x		x			x	x
9001-019306		x		x		x		x			x
9001-020164		x		x		x			x		x
9001-021022		x		x		x				x	x
9001-026170		x			x	x		x			x
9001-027028		x			x	x			x		x
9001-027886		x			x	x				x	x
9001-028744			x	x			x				x
9001-035608			x		x		x				x
9001-032176			x	x		x					x
9001-039040			x		x	x					x
9001-029602			x	x			x	x			x
9001-030460			x	x			x		x		x
9001-031318			x	x			x			x	x
9001-036446			x		x		x	x			x
9001-037324			x		x		x		x		x
9001-038182			x		x		x			x	x
9001-033034			x	x		x		x			x
9001-033892			x	x		x			x		x
9001-034750			x	x		x				x	x
9001-039898			x		x	x		x			x
9001-040756			x		x	x			x		x
9001-041614			x		x	x				x	x

**Mainstream accessories must be purchased separately.*

Propaq CS Monitor

Propaq CS monitors are shipped standard with adult, pediatric and neonatal monitoring capability. The Base Monitor includes ECG, NIBP, Temperature, Impedance Pneumography Respiration and either Nellcor Oximax® or Masimo® Motion Tolerant Pulse Oximetry. Accessories shipped with each monitor include a 3-lead ECG cable, lead wires, electrodes, Directions for Use on CD, one standard Flexi-Port adult (#11) NIBP cuff, one Flexi-Port child (#9) NIBP cuff, and one adult/pediatric NIBP hose. Included with Masimo SpO2 is one reusable pulse oximetry sensor and one extension cable. Included with Nellcor SpO2 is one reusable OxiMax sensor, one sample pack of disposable Oximax pulse oximetry sensors, and one extension cable. If you require either the 5-lead ECG cable and lead wire set or neonatal NIBP cuffs, they must be ordered separately. Monitors will accommodate (but are not shipped with) YSI temperature probes and popular invasive blood pressure transducers (5uV/V/mmHg). A three-year warranty and software license agreement are included with all Propaq monitors.

Part Number	Description
-------------	-------------

Propaq CS, Model 242, Base Monitor Configurations

9001-001288	Propaq CS, Model 242 - 8.4" Color Display - monitors Adult 3-Lead ECG, Respiration, NIBP, Nellcor OxiMax SpO2 and 2 channels of temperature, RS423 connectivity port
9001-008152	Propaq CS, Model 242 - 8.4" Color Display - monitors Adult 3-Lead ECG, Respiration, NIBP, Masimo Motion Tolerant SpO2 and 2 channels of temperature, RS 423 connectivity port

Add Printer

9001-004720	Propaq CS, Model 242 - 8.4" Color Display - monitors Adult 3-Lead ECG, Respiration, NIBP, Nellcor OxiMax SpO2, Expansion Module & Printer and 2 channels of temperature, RS423 connectivity port
9001-011584	Propaq CS, Model 242 - 8.4" Color Display - monitors Adult 3-Lead ECG, Respiration, NIBP, Masimo Motion Tolerant SpO2, Expansion Module & Printer and 2 channels of temperature, RS423 connectivity port

Add etCO2 (no Printer) - Nellcor SpO2 Versions

9001-002146	Propaq CS, Model 242 - 8.4" Color Display - monitors Adult 3-Lead ECG, Respiration, NIBP, Nellcor OxiMax SpO2, Expansion Module, 2 channels of temperature and Sidestream CO2, RS423 connectivity port
9001-003004	Propaq CS, Model 242 - 8.4" Color Display - monitors Adult 3-Lead ECG, Respiration, NIBP, Nellcor OxiMax SpO2, Expansion Module, 2 channels of temperature, and Mainstream CO2 (requires Mainstream sensor, purchased separately), RS423 connectivity port
9001-003862	Propaq CS, Model 242 - 8.4" Color Display - monitors Adult 3-Lead ECG, NIBP, Respiration, Nellcor OxiMax SpO2, Expansion Module, 2 channels of temperature and Dualstream CO2 (requires Mainstream sensor, purchased separately), RS423 connectivity port

Add etCO2 (no Printer) - Masimo SpO2 Versions

9001-009010	Propaq CS, Model 242 - 8.4" Color Display - monitors Adult 3-Lead ECG, Respiration, NIBP, Masimo Motion Tolerant SpO2, Expansion Module, 2 channels of temperature and Sidestream CO2, RS423 connectivity port
9001-009868	Propaq CS, Model 242 - 8.4" Color Display - monitors Adult 3-Lead ECG, Respiration, NIBP, Masimo Motion TolerantSpO2, Expansion Module, 2 channels of temperature, and Mainstream CO2 (requires Mainstream sensor, purchased separately), RS423 connectivity port
9001-010726	Propaq CS, Model 242 - 8.4" Color Display - monitors Adult 3-Lead ECG, NIBP, Respiration, Masimo Motion TolerantSpO2, Expansion Module, 2 channels of temperature and Dualstream CO2 (requires Mainstream sensor, purchased separately), RS423 connectivity port

Add etCO2 and Printer - Nellcor SpO2 Versions

9001-005578	Propaq CS, Model 242 - 8.4" Color Display - monitors Adult 3-Lead ECG, Respiration, NIBP, Nellcor OxiMax SpO2, Expansion Module & Printer, 2 channels of temperature and Sidestream CO2, RS423 connectivity port
9001-006436	Propaq CS, Model 242 - 8.4" Color Display - monitors Adult 3-Lead ECG, Respiration, NIBP, Nellcor OxiMax SpO2, Expansion Module & Printer, 2 channels of temperature, and Mainstream CO2 (requires Mainstream sensor, purchased separately), RS423 connectivity port

9001-007294 Propaq CS, Model 242 - 8.4" Color Display - monitors Adult 3-Lead ECG, NIBP, Respiration, Nellcor OxiMax SpO2, Expansion Module & Printer, 2 channels of temperature and Dualstream CO2 (requires Mainstream sensor, purchased separately), RS423 connectivity port

Add CO2 and Printer - Masimo SpO2 Versions

9001-012442 Propaq CS, Model 242 - 8.4" Color Display - monitors Adult 3-Lead ECG, Respiration, NIBP, Masimo Motion Tolerant SpO2, Expansion Module & Printer, 2 channels of temperature, and Sidestream CO2, RS423 connectivity port

9001-013300 Propaq CS, Model 242 - 8.4" Color Display - monitors Adult 3-Lead ECG, Respiration, NIBP, Masimo Motion Tolerant SpO2, Expansion Module & Printer, 2 channels of temperature and Mainstream CO2 (requires Mainstream sensor, purchased separately), RS423 connectivity port

9001-014158 Propaq CS, Model 242 - 8.4" Color Display - monitors Adult 3-Lead ECG, Respiration, NIBP, Masimo Motion Tolerant SpO2, Expansion Module & Printer, 2 channels of temperature and Dualstream CO2 (requires Mainstream sensor, purchased separately), RS423 connectivity port

Propaq CS, Model 244, Base Monitor Configurations

9001-015016 Propaq CS, Model 244 - 8.4" Color Display - monitors Adult 3-Lead ECG, Respiration, NIBP, 1 Invasive BP, Nellcor OxiMax SpO2 and 2 channels of temperature, RS423 connectivity port

9001-021880 Propaq CS, Model 244 - 8.4" Color Display - monitors Adult 3-Lead ECG, Respiration, NIBP, 1 Invasive BP, Masimo Motion Tolerant SpO2 and 2 channels of temperature, RS423 connectivity port

Add Printer

9001-018448 Propaq CS, Model 244 - 8.4" Color Display - monitors Adult 3-Lead ECG, Respiration, NIBP, 1 Invasive BP, Nellcor OxiMax SpO2, Expansion Module & Printer and 2 channels of temperature, RS423 connectivity port

9001-025312 Propaq CS, Model 244 - 8.4" Color Display - monitors Adult 3-Lead ECG, Respiration, NIBP, 1 Invasive BP, Masimo Motion Tolerant SpO2, Expansion Module & Printer and 2 channels of temperature, RS423 connectivity port

Add etCO2 (no Printer) - Nellcor SpO2 Versions

9001-015874 Propaq CS, Model 244 - 8.4" Color Display - monitors Adult 3-Lead ECG, Respiration, NIBP, 1 Invasive BP, Nellcor OxiMax SpO2, Expansion Module, 2 channels of temperature and Sidestream CO2, RS423 connectivity port

9001-016732 Propaq CS, Model 244 - 8.4" Color Display - monitors Adult 3-Lead ECG, Respiration, NIBP, 1 Invasive BP, Nellcor OxiMax SpO2, Expansion Module, 2 channels of temperature, and Mainstream CO2 (requires Mainstream sensor, purchased separately), RS423 connectivity port

9001-017590 Propaq CS, Model 244 - 8.4" Color Display - monitors Adult 3-Lead ECG, NIBP, 1 Invasive BP, Respiration, Nellcor OxiMax SpO2, Expansion Module, 2 channels of temperature and Dualstream CO2 (requires Mainstream sensor, purchased separately), RS423 connectivity port

Add etCO2 (no Printer) - Masimo SpO2 Versions

9001-022738 Propaq CS, Model 244 - 8.4" Color Display - monitors Adult 3-Lead ECG, Respiration, NIBP, 1 Invasive BP, Masimo Motion Tolerant SpO2, Expansion Module, 2 channels of temperature and Sidestream CO2, RS423 connectivity port

9001-023596 Propaq CS, Model 244 - 8.4" Color Display - monitors Adult 3-Lead ECG, Respiration, NIBP, 1 Invasive BP, Masimo Motion Tolerant SpO2, Expansion Module, 2 channels of temperature, and Mainstream CO2 (requires Mainstream sensor, purchased separately), RS423 connectivity port

9001-024454 Propaq CS, Model 244 - 8.4" Color Display - monitors Adult 3-Lead ECG, NIBP, 1 Invasive BP, Respiration, Masimo Motion Tolerant SpO2, Expansion Module, 2 channels of temperature and Dualstream CO2 (requires Mainstream sensor, purchased separately), RS423 connectivity port

Add CO2 and Printer - Nellcor SpO2 Versions

9001-019306	Propaq CS, Model 244 - 8.4" Color Display - monitors Adult 3-Lead ECG, Respiration, NIBP, 1 Invasive BP, Nellcor OxiMax SpO2, Expansion Module & Printer, 2 channels of temperature and Sidestream CO2, RS423 connectivity port
9001-020164	Propaq CS, Model 244 - 8.4" Color Display - monitors Adult 3-Lead ECG, Respiration, NIBP, 1 Invasive BP, Nellcor OxiMax SpO2, Expansion Module & Printer, 2 channels of temperature and Mainstream CO2 (requires Mainstream sensor, purchased separately), RS423 connectivity port
9001-021022	Propaq CS, Model 244 - 8.4" Color Display - monitors Adult 3-Lead ECG, NIBP, Respiration, 1 Invasive BP, Nellcor OxiMax SpO2, Expansion Module & Printer, 2 channels of temperature and Dualstream CO2 (requires Mainstream sensor, purchased separately), RS423 connectivity port

Add CO2 and Printer - Masimo SpO2 Versions

9001-026170	Propaq CS, Model 244 - 8.4" Color Display - monitors Adult 3-Lead ECG, Respiration, NIBP, 1 Invasive BP, Masimo Motion Tolerant SpO2, Expansion Module & Printer, 2 channels of temperature and Sidestream CO2, RS423 connectivity port
9001-027028	Propaq CS, Model 244 - 8.4" Color Display - monitors Adult 3-Lead ECG, Respiration, NIBP, 1 Invasive BP, Masimo Motion Tolerant SpO2, Expansion Module & Printer, 2 channels of temperature and Mainstream CO2 (requires Mainstream sensor, purchased separately), RS423 connectivity port
9001-027886	Propaq CS, Model 244 - 8.4" Color Display - monitors Adult 3-Lead ECG, Respiration, NIBP, 1 Invasive BP, Masimo Motion Tolerant SpO2, Expansion Module & Printer, 2 channels of temperature and Dualstream CO2 (requires Mainstream sensor, purchased separately), RS423 connectivity port

Propaq CS, Model 246, Base Monitor Configurations

9001-028744	Propaq CS, Model 246 - 8.4" Color Display - monitors Adult 3-Lead ECG, Respiration, NIBP, 2 Invasive BP, Nellcor OxiMax SpO2 and 2 channels of temperature, RS423 connectivity port
9001-035608	Propaq CS, Model 246 - 8.4" Color Display - monitors Adult 3-Lead ECG, Respiration, NIBP, 2 Invasive BP, Masimo Motion Tolerant SpO2 and 2 channels of temperature, RS423 connectivity port

Add Printer

9001-032176	Propaq CS, Model 246 - 8.4" Color Display - monitors Adult 3-Lead ECG, Respiration, NIBP, 2 Invasive BP, Nellcor OxiMax SpO2, Expansion Module & Printer and 2 channels of temperature, RS423 connectivity port
9001-039040	Propaq CS, Model 246 - 8.4" Color Display - monitors Adult 3-Lead ECG, Respiration, NIBP, 2 Invasive BP, Masimo Motion Tolerant SpO2, Expansion Module & Printer and 2 channels of temperature, RS423 connectivity port

Add etCO2 (no Printer) - Nellcor SpO2 Versions

9001-029602	Propaq CS, Model 246 - 8.4" Color Display - monitors Adult 3-Lead ECG, Respiration, NIBP, 2 Invasive BP, Nellcor OxiMax SpO2, Expansion Module, 2 channels of temperature and Sidestream CO2, RS423 connectivity port
9001-030460	Propaq CS, Model 246 - 8.4" Color Display - monitors Adult 3-Lead ECG, Respiration, NIBP, 2 Invasive BP, Nellcor OxiMax SpO2, Expansion Module, 2 channels of temperature, and Mainstream CO2 (requires Mainstream sensor, purchased separately), RS423 connectivity port
9001-031318	Propaq CS, Model 246 - 8.4" Color Display - monitors Adult 3-Lead ECG, NIBP, Respiration, Nellcor OxiMax SpO2, Expansion Module, 2 channels of temperature and Dualstream CO2 (requires Mainstream sensor, purchased separately), RS423 connectivity port

Add etCO2 (no Printer) - Masimo SpO2 Versions

9001-036446	Propaq CS, Model 246 - 8.4" Color Display - monitors Adult 3-Lead ECG, Respiration, NIBP, 2 Invasive BP, Masimo Motion Tolerant SpO2, Expansion Module, 2 channels of temperature and Sidestream CO2, RS423 connectivity port
-------------	---

9001-037324 Propaq CS, Model 246 - 8.4" Color Display - monitors Adult 3-Lead ECG, Respiration, NIBP, 2 Invasive BP, Masimo Motion TolerantSpO2, Expansion Module, 2 channels of temperature, and Mainstream CO2 (requires Mainstream sensor, purchased separately), RS423 connectivity port

9001-038182 Propaq CS, Model 246 - 8.4" Color Display - monitors Adult 3-Lead ECG, NIBP, 2 Invasive BP, Respiration, Masimo Motion TolerantSpO2, Expansion Module, 2 channels of temperature and Dualstream CO2 (requires Mainstream sensor, purchased separately), RS423 connectivity port

Add CO2 and Printer - Nellcor SpO2 Versions

9001-033034 Propaq CS, Model 246 - 8.4" Color Display - monitors Adult 3-Lead ECG, Respiration, NIBP, 2 Invasive BP, Nellcor OxiMax SpO2, Expansion Module & Printer, 2 channels of temperature and Sidestream CO2, RS423 connectivity port

9001-033892 Propaq CS, Model 246 - 8.4" Color Display - monitors Adult 3-Lead ECG, Respiration, NIBP, 2 Invasive BP, Nellcor OxiMax SpO2, Expansion Module & Printer, 2 channels of temperature and Mainstream CO2 (requires Mainstream sensor, purchased separately), RS423 connectivity port

9001-034750 Propaq CS, Model 246 - 8.4" Color Display - monitors Adult 3-Lead ECG, NIBP, Respiration, 2 Invasive BP, Nellcor OxiMax SpO2, Expansion Module & Printer, 2 channels of temperature and Dualstream CO2 (requires Mainstream sensor, purchased separately), RS423 connectivity port

Add CO2 and Printer - Masimo SpO2 Versions

9001-039898 Propaq CS, Model 246 - 8.4" Color Display - monitors Adult 3-Lead ECG, Respiration, NIBP, 2 Invasive BP, Masimo Motion Tolerant SpO2, Expansion Module & Printer, 2 channels of temperature and Sidestream CO2, RS423 connectivity port

9001-040756 Propaq CS, Model 246 - 8.4" Color Display - monitors Adult 3-Lead ECG, Respiration, NIBP, 2 Invasive BP, Masimo Motion Tolerant SpO2, Expansion Module & Printer, 2 channels of temperature and Mainstream CO2 (requires Mainstream sensor, purchased separately), RS423 connectivity port

9001-041614 Propaq CS, Model 246 - 8.4" Color Display - monitors Adult 3-Lead ECG, Respiration, NIBP, 2 Invasive BP, Masimo Motion Tolerant SpO2, Expansion Module & Printer, 2 channels of temperature and Dualstream CO2 (requires Mainstream sensor, purchased separately), RS423 connectivity port

Model Numbers

All Propaq CS models can be configured with the same optional parameters and features, except invasive blood pressure lines.

Propaq CS Model Numbers:

242: No invasive blood pressure lines

244: 1 invasive blood pressure line

246: 2 invasive blood pressure lines

The device model number is located below the handle on the back of the device, next to the serial number.

Serial Numbers: Device and Expansion Module

Propaq CS device serial number is located below the handle on the back of the device, next to the model number. Device serial numbers begin with "DA" or "GA", followed by six numbers. For example, GA108000.

If the Propaq CS device has an expansion module (required for a printer and/or CO2), there is also an expansion module serial number located on the back of the unit, near the top of the expansion module. Expansion module serial numbers are eight numbers, such as 08291015. If the customer chooses to mount a device with an expansion module on the wall, the mounting plate will cover the expansion module serial number.

Only the Propaq CS device serial number is required for service inquiries, not the expansion module serial number.

Propaq CS Ship-with Accessories

	Item Number	Item Description
<u>SPO2</u>		
Nellcor:	008-0054-01	DS-100A DURASENSOR
	660-0196-00	ASSY,8FT SPO2,DIFF. EXT CBL,NELLCOR
	680-0096-00	SENSOR PACK, SPO2, OXIMAX, A/AL
Masimo:	008-0824-00	ACCY, SPO2 CABLE, MASIMO, 8', PCO8
	008-0653-00	ADULT/PEDIATRIC SPO2 SNSOR SMPLE PACK
<u>ECG</u>		
	008-0315-00	ECG CABLE,3 LEAD,W/WIRES ATTACHED, SNAP AAMI
	008-0023-00	ECG, ELECTRODES, S/LH, PREJELLED, 1700 CLRANCE
<u>NIBP</u>		
	008-0238-00	ACCY, HOSE, ADULT, 10FT FOR NIBP CUFF
	008-1022-00	CUFF REUSE-09-1MQ, CHILD, 1 TUBE,MQ
	008-1024-00	CUFF REUSE-11-1MQ, ADULT, 1 TUBE,MQ
<u>Sample Pack</u>		
	008-1044-00	WELCH ALLYN SAMPLE PAK – Oridion Surestream Accy
<u>Paper</u>		
	008-0040-00	PAPER, THERMAL SENSITIVE,60MM W, 30M L
<u>Power Cord and Adaptor</u>		
	503-0054-02	POWER ADAPTER,120V,25W
	008-0333-00	POWER CORD,HOSP GRADE,N/S AMERICA,10A/120V
<u>Cards</u>		
	810-0137-02	CARD, EXTENDED WARRANTY, INQUIRY, PROPAQ
	810-1103-02	CARD, PROPAQ CS 3.XX QUICK REF., ENGLISH
	810-0466-03	QUICK GUIDE, CO2 WAVEFORMS, PROPAQ, ENGLISH
<u>CDs/Guides</u>		
	810-2734-00	IT, CD, DFU, PROPAQ CS 3.7X, MULTILANG
	810-0409-10	BOOKLET, PRODUCTS AND ACCESSORIES
<u>Documentation</u>		
	810-0005-01	"DEAR CUSTOMER" NOTE
	810-1088-00	TECHNICAL BULLETIN, ENCORE V.2X NIBP TESTING
	830-0487-06	DECLARATION OF CONFORMITY,PROPAQ 2XX SERIES

NOTE: *All Ship-with accessories are also available for sale individually.*

Propaq CS Accessories Sold Separately

Battery Packs

- 008-0125-00 Battery Pack (Fused) - rechargeable (double, for monitors with add-ons)
- 008-0126-00 Battery Pack (Fused) - rechargeable (single, for monitors with no add-ons)

Cables, Miscellaneous

- 008-0150-00 Cable Wraps, 1" x 6" Velcro® Cable Organizers Blue, 10/Package
- 008-0320-00 Cable, Real-Time ECG, 15', unterminated (Propaq CS & Encore only)
- 008-0634-00 Cable, Nurse Call Interface, 10', unterminated (Propaq CS only)

Carrying Cases/Shipping Containers

- 900-0820-00 Case, Protective, for Propaq CS only, with Sidestream CO2 option. Fits all monitors with Expansion Module (SpO2/CO2 and/or Printer Options)
- 900-0017-00 Case, Travel for all Propaq models
- 900-0002-00 Pouch for Accessories - with zipper, mountable on all Propaq models

CO2 Accessories

Mainstream CO2

- 008-0261-00 CO2 Airway Adapters, Multi-Use Adult/Pediatric, 3/Box (for patients with endotracheal tubes larger than 5.0mm I.D.), includes plastic storage case
- 008-0132-00 CO2 Airway Adapters, Single Patient Use, Adult/Pediatric, for Mainstream CO2 Option, box of 10
- 008-0367-00 CO2 Airway Adapters, Single Patient Use, Adult/Pediatric, for Mainstream CO2 Option, box of 24
- 008-0263-00 CO2 Airway Adapters, Single Patient Use, Low Dead-space, 10/Box (for patients with endotracheal tubes 5.0mm I.D. or smaller)
- 008-0502-01 CO2 Sensor & Airway Adapter Set (contains 1 each Mainstream CO2 Sensor with 10' Cable, Single Patient Use Airway Adapter, Storage Case, Directions for Use)

ECG Accessories, Propaq CS Monitors

All ECG cables have AAMI 6-pin connectors and are compatible with all Propaq monitors. Propaq Encore and CS monitors offer both 3- and 5-lead ECG monitoring capability.

3-Lead ECG Cables

- 008-0315-00 ECG, 3-Lead Cable, one piece shielded 10' cable with 3 attached lead wires, snap
- 008-0429-00 ECG, 3-Lead Cable, right angle, safety DIN, .060" for use with 008-0321-00, 008-0317-00 & 008-0323-00
- 008-0316-00 ECG, 3-Lead Cable, Shielded, Safety DIN .060" for use with 008-0323-00, 008-0321-00 & 008-0317-00
- 008-0321-00 ECG, 3-Lead Wire Set, for .060" safety DIN, adult pinch for use with 008-0316-00 & 008-0429-00
- 008-0317-00 ECG, 3-Lead Wire Set, for .060" safety DIN, neonatal pinch for use with 008-0316-00 & 008-0429-00
- 008-0323-00 ECG, 3-Lead Wire Set, Safety DIN .060", snap for use with 008-0316-00 & 008-0429-00
- 008-0734-00 Adapter Cable, ECG, Micropaq to Propaq, AAMI


5-Lead ECG Cables

008-0313-00	ECG, 5-Lead Cable, one piece shielded 10' cable with 5 attached lead wires, snap
008-0314-00	ECG, 5-Lead Cable, shielded, Safety DIN .060" for use with 008-0322-00 & 008-0319-00
008-0322-00	ECG, 5-Lead Wire Set, Safety DIN .060", snap for use with 008-0314-00
008-0319-00	ECG, 5-Lead Wire Set, Safety DIN, .060" adult pinch for use with 008-0314-00
008-0734-00	Adapter Cable, ECG, Micropaq to Propaq, AAMI

IBP Pressure Transducers

008-0226-01	IBP Cable for MX900 and MX860 IBP Domes for MX860, case of 50
-------------	--

008-0224-00



008-0233-00	IBP Transducer, Reuseable, MX860, case of 5
-------------	---

Mounting Systems


ADVISORY: Wall and pole mount accessories NOT intended for use in high vibration environments, such as aircraft or ambulances.

Pole Mounting


008-0205-00	Pole Mount Clamp (order monitor backplate or wraparound pan with aluminum slide separately, if not already installed)
008-0204-00	Back Plate with Aluminum Slide for monitors with Expansion Module option (for use with Pole Mount Clamp, which is purchased separately)
008-0215-00	Wraparound Pan with Aluminum Slide for monitors with SpO2 option (for use with Pole Mount Clamp, which is purchased separately)

Wall Mounting

008-0887-00	Wall Mount, Complete Kit (6" basket, tilt mount and 19" wall channel), 8" Pivot Arm, for monitors with SpO2 option, drop-on style
-------------	---




008-0888-00 Wall Mount, Complete Kit (6" basket, tilt mount and 19" wall channel), 8" Pivot Arm, for monitors with Expansion Module option, drop-on style



008-0211-01 Wall Mount, Flush Mount, for monitors with SpO2 option, drop-on style

008-0208-01 Wall Mount, Flush Mount, for monitors with Expansion Module option, drop-on style


008-0200-01 6" Wire Basket for Wall Mount



008-0203-00 19" Wall Channel


008-0890-00 Wraparound pan for monitors with SpO2 option (for use with Roll Stands and Wall Mounts)

008-0889-00 Back Plate for monitors with Expansion Module option (for use with Roll Stands and Wall Mounts)




Roll Stand, Complete Kit (handle, basket, and mount), for monitors with expansion module option, drop-on style

Roll Stand, Complete Kit (handle, basket, and mount), for monitors with expansion module option, drop-on style



Roll Stand, Complete Kit (handle, basket, and mount), for monitors with SpO2 option, drop-on style



6" Wire Basket for Roll Stand Mount

Wraparound pan for monitors with SpO2 option (for use with Roll Stands and Wall Mounts)

Back Plate for monitors with Expansion Module option (for use with Roll Stands and Wall Mounts)

Miscellaneous Mounting Accessories

Upgrade kit, drop-on style, mounting plate and adapter

Mounting Bracket for Universal Power Adapter (part number 503-0054-XX) (Requires wall mount slide, part number 008-0793-00)

NIBP Accessories

NIBP Hoses

* Interchangeable with all Propaq monitors

008-0238-00 Hose, Adult/Pediatric, 10', w/ "twist lock" cuff connector (for bladderless cuffs) *

008-0265-01 NIBP Hose, Neonatal, 8'

NIBP Disposable Cuffs - Bladderless (Hose not included)

5082-101-1 Neonate #1, disposable, single tube, box of 10

5082-102-1 Neonate #2, disposable, single tube, box of 10

5082-103-1 Neonate #3, disposable, single tube, box of 10

5082-104-1 Neonate #4, disposable, single tube, box of 10

5082-105-1 Neonate #5, disposable, single tube, box of 10

5082-241-9 Neonatal Cuff Kit, disposable, one each of sizes #1 - #5, single tube, bag of 5

NIBP Reusable Cuffs - Bladderless (Hose not included)

008-0324-00 Cuff Kit, one each of Large Adult, Small Adult and Thigh Cuff w/ "twist-lock" connector (includes WA part#'s, 5082-87-3MQ, 5082-85-3MQ, 5082-88-3MQ)

REUSE-06-1MQ Cuff, Reuse, Newbrn, 1-Tube, Twist lock connector

REUSE-07-1MQ Cuff, Reuse, Infant, 1-Tube, Twist lock connector

REUSE-08-1MQ Cuff, Reuse, Sm Child, 1-Tube, Twist lock connector

REUSE-09-1MQ Cuff, Reuse, Child, 1-Tube, Twist lock connector

REUSE-10-1MQ Cuff, Reuse, Sm Ad, 1-Tube, Twist lock connector

REUSE-11-1MQ Cuff, Reuse, Adult, 1-Tube, Twist lock connector

REUSE-11L-1MQ Cuff, Reuse, Ad Long, 1-Tube, Twist lock connector

REUSE-12-1MQ Cuff, Reuse, Lg Ad, 1-Tube, Twist lock connector

REUSE-12L-1MQ Cuff, Reuse, Lg Ad Long, 1-Tube, Twist lock connector

REUSE-13-1MQ Cuff, Reuse, Thigh, 1-Tube, Twist lock connector

SOFT-06-1MQ Cuff, Soft, Newbrn, 1-Tube, Twist lock connector, box of 20

SOFT-07-1MQ Cuff, Soft, Infant, 1-Tube, Twist lock connector, box of 20

SOFT-08-1MQ Cuff, Soft, Sm Child, 1-Tube, Twist lock connector, box of 20

SOFT-09-1MQ Cuff, Soft, Child, 1-Tube, Twist lock connector, box of 20

SOFT-10-1MQ Cuff, Soft, Sm Ad, 1-Tube, Twist lock connector, box of 20

SOFT-11-1MQ Cuff, Soft, Adult, 1-Tube, Twist lock connector, box of 20

SOFT-11L-1MQ Cuff, Soft, Ad Long, 1-Tube, Twist lock connector, box of 20

SOFT-12-1MQ Cuff, Soft, Lg Ad, 1-Tube, Twist lock connector, box of 20

SOFT-12L-1MQ Cuff, Soft, Lg Ad Long, 1-Tube, Twist lock connector, box of 20

SOFT-13-1MQ Cuff, Soft, Thigh, 1-Tube, Twist lock connector, box of 20

Power Adapters

503-0054-02 AC Adapter for Propaq CS with CO2 Option, or 25 Watt Option, North American 120 VAC, 60Hz (Power Cord not included.) (Style B Adapter - Universal)

Power Cords

All of the following hospital-grade Power Cords fit the Universal Power Adapter with Style B Connector.

NOTE: The Universal Power Adapter connects to Propaq CS Series monitors with the 25-watt power option or the CO2 Option

008-0290-00 Power Cord DC, 8 ft. Cable, for connection to fused DC power sources.

008-0333-00 Power Cord, Hospital Grade for Power Adapter for Propaq Series with CO2 Option, North American, 100 to 120 VAC

Printer Supplies

008-0040-98	Printer Paper - 1 Box, 10 rolls, standard thermal print paper
008-0040-99	Printer Paper - 100 rolls each, standard thermal print paper
008-0178-00	Printer Paper, Permanent - 1 Box, 10 rolls, premium fade-resistant thermal print paper

Pulse Oximetry Accessories

900-0021-00	Nellcor - DS-100A Durasensor (Adult), reusable
008-0742-00	Nellcor - 8' Differential Extension Cable, Pulse Oximeter Sensor
008-0654-00	Masimo - LNOP Pediatric/Slender Digit SpO2 Sensors, Disposable, box of 20
008-0655-00	Masimo - LNOP Adult SpO2 Sensors, Disposable, box of 20
008-0695-00	Masimo - Neonate/Pediatric Sensors, disposable, sensitive skin, box of 20
008-0694-00	Masimo - LNOP Neonate Sensors, Disposable, box of 20
008-0648-00	Masimo - LNOP-DCI Adult Reusable Sensor
008-0656-00	Masimo - LNOP-DCIP Pediatric Reusable Sensor
008-0692-02	Masimo - 8' Interface Cable, compatible with round side panel connector (an interface cable is required with Masimo sensors)
008-0824-00	Masimo - 8' Interface Cable, compatible with rectangle (Masimo standard) side panel connector (an interface cable is required with Masimo sensors)
LNCS-DCI	Masimo - LNCS-DCI Adult Reusable Sensor
008-1014-00	Masimo - LNCS 10' mini-d SP02 cable

Temperature Sensors & Cables

008-0234-00	Temp - Cable 10' Extension for YSI Disposable Probes, YSI #4940
008-0235-00	Temp - Cable 10' Extension for YSI Reusable Temp Probes, YSI #4010
008-0227-00	Temp - General Purpose Probe, Adult, 12 FR, YSI #4491, Disposable, box of 100 (requires 008-0452-00)
008-0419-01	Temp - Probe, Reusable, General Purpose, Adult, Autoclavable, YSI #401AC
008-0419-00	Temp - Probe, Reusable, General Purpose, Pediatric, Autoclavable, YSI #402AC
008-0230-00	Temp - Probe, Reusable, Adult, General Purpose, YSI #401
008-0229-00	Temp - Probe, Reusable, Adult, Skin Surface, YSI #409B
008-0756-00	Temp - Skin Sensor, Pediatric, Reusable, Skin Surface, YSI #427
008-0228-00	Temp - Skin Sensor, YSI #4499, Disposable, box of 100
008-0451-00	Temp - Skin Sensor, Adult, Reusable, Autoclavable, YSI #409AC
008-0236-00	Tempheart - Adhesive Pad (for use with YSI #409B, 409AC, and 427 probes), 100 Count, YSI #4012

User Documentation, Propaq CS

Kits include manuals and updates.

CO2 Waveform User Documentation

810-0466-03	Quick Guide for CO2 Waveforms, English
-------------	--

Propaq CS Monitor User Documentation

810-1103-02	Quick Reference Card, Propaq CS Software Version 3.xx, English
810-2105-01	Directions for Use Kit, Propaq CS Software Version 3.xx, English, Paper Copy
810-1441-03	Directions for Use Kit, Propaq CS Software Version 3.xx, Multi-language CD
810-1315-05	Service Manual Kit, Propaq CS, English, Paper Copy
810-1657-02	Service Manual Kit, Propaq CS, English, CD

Partners in Care Service Offerings, Propaq CS

S1-24X	Propaq CS, Comprehensive Partner Program, 1 year agreement
S1-24X-2	Propaq CS, Comprehensive Partner Program, 2 year agreement
S2-24X	Propaq CS, Biomed Partner Program, 1 year agreement
S2-24X-2	Propaq CS, Biomed Partner Program, 2 year agreement

Propaq CS ECG Cables In Detail

The Propaq CS monitor offers both 3- and 5-lead ECG monitoring capability.
All ECG cables have a 6-pin connector.

3-Lead ECG Cables

KITs:

- 008-0974-00 ECG CABLE, 3 LEAD, **NEONATAL**, **SAFETY CLIP**, AAMI
008-0329-00 ECG, 3-Lead Wire Set, Safety Clip, Neonatal, AAMI for use with 660-0023-00 (also available independently)
660-0023-00 ECG, 3-Lead Cable, Safety Clip, Neonatal, AAMI for use with 008-0329-00 (also available independently)
- 008-0976-00 ECG CABLE, 3 LEAD, **ADULT/PED**, **PINCH**, AAMI
008-0316-00 ECG, 3-Lead Cable, Shielded, Safety DIN .060", adult pinch, AAMI (also available independently)
008-0321-00 ECG, 3-Lead Wire Set, for .060" safety DIN, adult pinch, AAMI (also available independently)
- 008-0245-00 ACCY, **ADULT/PED**, CABLE/3 LEAD SET, AAMI
660-0023-00 CABLE, PATIENT, NORTH AMERICAN (also available independently)
008-0166-00 LEAD SET, HP, SHIELDED, NORTH AMERICAN - AAMI (also available independently)
008-0023-00 ECG, ELECTRODES, S/LH, PREJELLED, 1700 CLEARANCE

DIN CABLES - separate lead wires required:


- 008-0429-00 ECG CABLE, 3 LEAD, **RIGHT ANGLE**, safety DIN, .060" for use with 008-0321-00, 008-0317-00 and 008-0323-00

ECG Cable, Right Angle, 3 Lead DIN, AHA (AAMI)/IEC




REF: 008-0429-00

LOT



This cable meets the requirements of the ANSI/AAMI standard, ECG Cables and Lead Wires (ANSI/AAMI EC50-1995). This cable is reusable.

This cable complies with FDA, 21 CFR Part 606 (Medical Devices; Establishment of a Performance Standard for Electrode Lead Wires and Patient Cables).



ECG Cable



Welch Allyn Protocol, Inc.
8500 SW Creekside Place
Beaverton, OR 97008-7107 USA
www.monitoring.welchallyn.com
Phone: 503-530-7500
Phone: USA: 800-288-2500
Technical Services USA: 800-289-2501
Technical Services Fax: 503-526-4910

WelchAllyn
Welch Allyn Protocol Inc.

ECG

Propaq
Monitor


008-0316-00 ECG, 3-Lead Cable, **Safety DIN** .060" for use with 008-0323-00, 008-0321-00 & 008-0317-00

ECG Cable, 3 Lead, DIN, AHA (AAMI)



REF: 008-0316-00

LOT




This cable meets the requirements of the ANSI/AAMI standard, ECG Cables and Lead Wires (ANSI/AAMI EC53-1995). This cable is reusable.

This cable complies with FDA, 21 CFR Part 800 (Medical Devices; Establishment of a Performance Standard for Electrode Lead Wires and Patient Cables).

Welch Allyn Protocol Inc.
8500 SW Creekside Place
Beaverton, OH 97008-7107 USA
www.monitoring.welchallyn.com
Phone: 503-530-7500
Phone: USA: 800-289-2500
Technical Services USA: 800-289-2501
Technical Services Fax: 503-526-4910

WelchAllyn
Welch Allyn Protocol Inc.



DIN LEAD WIRES AVAILABLE:

008-0321-00 ECG, 3-Lead Wire Set, for .060" safety DIN, **Adult Pinch** for use with 008-0316-00 & 008-0429-00, AAMI

008-0317-00 ECG, 3-Lead Wire Set, for .060" safety DIN, **Neonatal Pinch** for use with 008-0316-00&008-0429-00, AAMI

008-0323-00 ECG, 3-Lead Wire Set, Safety DIN .060", **SNAP** for use with 008-0316-00 & 008-0429-00, AAMI

NOTE: All DIN lead wires can be used with any DIN cable! This includes the use of 3-lead DIN lead wire sets in 5-lead DIN ECG Cables.

5-Lead ECG Cables

KITs:


008-0972-00 ECG CABLE 5 LEAD, **ADULT/PED**, HP, AAMI

008-0328-00 ECG, 5-Lead Wire Set, **Snap, Adult/Ped**, HP, AAMI for use with 008-0330-00 (also available independently)

008-0330-00 ECG, 5-Lead Cable, **Snap, Adult/Ped**, HP, AAMI for use with 008-0328-00 (also available independently)


ONE PIECE CABLE – No lead wires required:

008-0313-00 ECG, 5-Lead Cable, one piece shielded 10' cable with 5 attached lead wires, **SNAP**, AAMI



DIN CABLE – separate lead wires required:

008-0314-00 ECG, 5-Lead Cable, shielded, Safety **DIN** .060" for use with 008-0322-00 & 008-0319-00, AAMI



DIN LEAD WIRES AVAILABLE:

008-0322-00 ECG, 5-Lead Wire Set, Safety DIN .060", **SNAP** for use with 008-0314-00, AAMI

008-0319-00 ECG, 5-Lead Wire Set, Safety DIN, .060", **ADULT PINCH** for use with 008-0314-00, AAMI

YSI Temperature Notes

- The Propaq CS measures all temperature options, the software/hardware is available: Esophageal, Rectal, and Skin.
- The temperature probes approved are **YSI, 400 and 700 series**.
- **Core temp** can be obtained rectally or through the esophagus.

Propaq Temperature Specifications

Connector

The T1 and T2 connectors (1/4" "phone" jacks) are compatible with both YSI 400 and YSI 700 Temp probes.

Accessories

The YSI 400 and 700 compatible probes may be used.

Display Characteristics

The Temp range is 0°C to 50°C or 32°F to 122°F.

Performance Characteristics

The Temp accuracy is from 10°C to 50°C is $\pm 0.1^{\circ}\text{C}$ exclusive of probe errors. The Temp accuracy between 10°C and 0°C is $\pm 0.2^{\circ}\text{C}$ exclusive of probe errors. In addition, the monitor provides "continuous" self-test of temperature channel accuracy.

Numeric Characteristics for Display

Numerics are recalculated at 670 msec. intervals.

**We do not sell the probes.
Customers should contact YSI directly:**

Toll Free: 800 747 5367 or (937) 427-1231 ext 770

Oridion etCO2 Accessories

etCO2 Sample Lines for Propaq CS

<u>Oridion PN:</u>	<u>Product Description:</u>
010986	Low Humidity ETCO2 -25/Bx – Intubated Sampling
010987	High Humidity ETCO2 -25/Bx – Intubated Sampling
010976	Sureline Adult, ETCO2 Nasal Cannula – 25/Bx – Non-Intubated/Nasal Sampling
010977	Sureline Pedi, ETCO2 Nasal Cannula – 25/Bx – Non-Intubated/Nasal Sampling
010979	Soft Adult O2 Delivery/ ETCO2 Nasal Cannula – 25/Bx – Non-Intubated/Nasal Sampling w/O2
010980	Soft Pedi O2 Delivery/ ETCO2 Nasal Cannula – 25/Bx – Non-Intubated/Nasal Sampling w/O2
010981	Adult Oral/Nasal O2 Sure Line w/O2 25/Bx – Non-Intubated/Oral/Nasal Sampling w/O2 delivery
010982	Pedi Oral/Nasal O2 Sure Line w/O2 25/Bx – Non-Intubated/Oral/Nasal Sampling w/O2 delivery
010177	Adult Oral/Nasal Cannula - 25/Bx – Non-Intubated/O2 Tubing/Oral/Nasal
007606	Pedi Oral/Nasal Cannula - 25/Bx – Non-Intubated/O2 Tubing/Oral/Nasal
007609	Adult CO2/O2 Cannula - 25/Bx – Non-Intubated/Nasal Sampling/O2 Tubing
007610	Pedi CO2/O2 Cannula - 25/Bx – Non-Intubated/Nasal Sampling/O2 Tubing

Miscellaneous

010993	Water Trap - 10/Bx
010994	Water Trap - 24/Bx
010991	Sample Line 7 ft – 25/Bx
010989	Adult/Pediatric Tee w/Sample Port

Reference Info

Low Humidity to be used for short term monitoring of intubated patients, emergency
High Humidity to be used for long term monitoring of intubated patients receiving heated humidified gas.

**We do not sell any of the above accessories.
Customers should contact Oridion directly:
Toll Free: 888 ORIDION or (781) 453-0500
www.oridion.com**

Other Vendors

For Additional Invasive Blood Pressure Line Accessories and Information, contact:

Smiths Medical North America

Tel: +1 800 258 5361 / +1 614 210 7300 - **Fax:** +1 614 889 2651

info.asd@smiths-medical.com

For Additional Pulse Oximetry Accessories and related information, contact:



1-949- 297-7000, or 1-800-326-4890 Option #1 for Customer Service

www.masimo.com



1-800-635-5267

<http://respiratorysolutions.covidien.com/>

For additional Mounting Option Accessories and related information, contact:



1-800-228-2555

www.gcx.com