

WelchAllyn®

Propaq LT®



Ambulatory Care Reference Guide

Table of Contents

Propaq LT Part Numbers, Model Numbers and Serial Numbers	3
Propaq LT Ship-with Accessories	6
Propaq LT Accessories Sold Separately	8
Propaq LT ECG Cables In Detail	10
Propaq LT Large Color Display.....	12
Propaq LT Mounting Options	15
Propaq LT LCDI Mounting Options	18
Propaq LT Features in Brief.....	22
Other Vendors for Accessories	23

NOTE: This reference guide is intended for the use of representatives from Welch Allyn and our ambulatory care distributors. The document in full is not intended to give to customers. However, information within this document may be shared verbally or in written form with clients and prospective customers.

Propaq LT Monitors

Propaq LT monitors are shipped complete with adult, pediatric, and neonatal monitoring capability. The Base Monitor includes ECG, NIBP, Impedance Pneumography Respiration and Nellcor Oximax® or Masimo SET® Pulse Oximetry. Accessories shipped with each monitor include a 3-lead ECG cable, electrodes, Directions for Use, and one standard Flexi-Port adult NIBP cuff. Included with Nellcor SpO2 is one reusable OxiMax sensor, and one sample pack of disposable Oximax pulse oximetry sensors. Included with Masimo SpO2 is one 10' Masimo extension cable and one LNCS DCI Reusable Sensor. A two-year warranty and software license agreement are included with all Propaq LT monitors.

Propaq LT Monitor – Standalone Devices

Nellcor SpO2

802LT0N-0E0	Propaq LT - Nellcor SpO2 (Requires charging cradle to recharge batteries. Charging cradle sold separately) - (Delivered with AAMI ECG Leads)
802LT0N-CE1	Propaq LT - Nellcor SpO2 / Charging Cradle - 120 NA (Delivered with AAMI ECG Leads)
802LT0N-UE1	Propaq LT - Nellcor SpO2 / Charging Cradle / USB port - 120 NA (Delivered with AAMI ECG Leads)

*Note: N in the model number stands for Nellcor
E in the model number stands for AAMI ECG (CE1/UE1)*

Masimo SpO2

802LT0S-0E0	Propaq LT - Masimo SpO2 - (Requires charging cradle to recharge batteries. Charging cradle sold separately) (Delivered with AAMI ECG Leads)
802LT0S-CE1	Propaq LT - Masimo SpO2 / Charging Cradle - 120 NA (Delivered with AAMI ECG Leads)
802LT0S-UE1	Propaq LT - Masimo SpO2 / Charging Cradle / USB port - 120 NA (Delivered with AAMI ECG Leads)

*Note: S in the model number stands for Masimo
E in the model number stands for AAMI ECG (CE1/UE1)*

Propaq LTR Monitor, with FlexNet Wireless Option for Acuity and Mobile Acuity Connectivity

Nellcor SpO2

802LTAN-0E0	Propaq LTR with FlexNet (802.11a) Option, Nellcor SpO2 - (Requires charging cradle to recharge batteries. Charging cradle sold separately) NA (Delivered with AAMI ECG Leads)
802LTAN-CE1	Propaq LTR with FlexNet (802.11a) Option, Nellcor SpO2 / Charging Cradle - 120 NA (Delivered with AAMI ECG Leads)
802LTAN-UE1	Propaq LTR with FlexNet (802.11a) Option, Nellcor SpO2 / Charging Cradle - USB port - 120 NA (Delivered with AAMI ECG Leads)


Masimo SpO2

802LTAS-0E0	Propaq LTR with FlexNet (802.11a) Option, (802.11A), Masimo SpO2 - (Requires charging cradle to recharge batteries. Charging cradle sold separately) (Delivered with AAMI ECG Leads)
802LTAS-CE1	Propaq LTR with FlexNet (802.11a) Option, Masimo SpO2 / Charging Cradle - 120 NA (Delivered with AAMI ECG Leads)
802LTAS-UE1	Propaq LTR with FlexNet (802.11a) Option, Masimo SpO2 / Charging Cradle / USB port - 120 NA (Delivered with AAMI ECG Leads)

Part Number Decoder:

1 = Power Plug Type A

Type A = 120V 60Hz: Flat blade attachment plug



Cradles:

0E0: 0 = No Cradle

CE1: C = With Cradle

UE1: U = USB port with Cradle

Model Numbers:

Propaq LT device model number is located on the back of the device, near the upper right corner, just above the sticker with the serial number.

Standalone Devices (no radio):

802LT0N = Propaq LT Base Model with Nellcor

802LT0S = Propaq LT Base Model with Masimo

Wireless Devices (with radio):

802LTAN = Propaq LT with 802.11a Wireless Radio and Nellcor

802LTAS = Propaq LT with 802.11a Wireless Radio and Masimo

Note: The Radio/Wireless Option is not visible on the LT device, to identify if a monitor is wireless:

1. Check the model number on back of the device; if it has an A (LTAN), it is wireless
2. Turn the monitor on; if it has a Ψ , it is wireless

Serial Numbers:

The Propaq LT **device serial number** is located on the back of the device, near the upper right corner. The device serial number begins with "**KA**" and is followed by six numbers, such as KA009768.

The Propaq LT **cradle serial number** is located on the bottom of the cradle, near the right side. The cradle serial number begins with "**KL**" or "**KV**" and is followed by six numbers, such as KL001605.

Welch Allyn's customer record should include both the device serial number and the cradle serial number.

Propaq LT Base Model

802LT0N-0E0 or 802LT0S-0E0

Propaq LT Vital Signs Monitor


Requires charging cradle to recharge battery - Charging cradle sold separately.

Propaq LT with Charging Cradle

802LT0N-CE1 or 802LT0S-CE1

Propaq LT with Charging Cradle and USB Port

802LT0N-UE1 or 802LT0S-UE1



Charging Cradle


008-0869-00

Charging Cradle

008-0873-00

Charging / Communications
Cradle with USB

NOTE: At least one Charging Cradle in every order should come with the USB Port option. The USB Port option is necessary if the user would like to custom-configure their monitor, use the Auto-Print functionality, or connect to a Large Color Display



SHIP WITH ACCESSORIES

Item Number	Item Description
-------------	------------------

SPO2

Nellcor version:

008-0054-01	DS-100A DURASENSOR
660-0196-00	ASSY,8FT SPO2,DIFF. EXTENSION CABLE,NELLCOR
680-0047-01	SENSOR PACK, SPO2, DISPOSABLE, OXIMAX, A/N

Masimo version:

008-0960-00	ACCY, SENSOR, SP02, ADULT, REUSABLE, DB-9, MASIMO
008-1010-00	ACCY, SPO2 CABLE, MASIMO DB9-LNCS, 10'

ECG

008-0880-00	ECG CABLE,3 ATTACHED LEAD,SNAP,AAMI,10'
008-0023-00	ECG,ELECTRODES,S/LH,PREJELLED,1700 CLEARANCE

NIBP

008-1024-00	CUFF REUSE-11-1MQ, ADULT, 1 TUBE, MQ
008-0864-00	HOSE NIBP ADLT/PED 6' W/TWST LOCK CNTR BLDRLESS CUFF

Power Cord and Supply

008-0913-00	POWER C PLUG A,CA,2.5A,125V,2.5M
503-0142-01	POWER SUPPLY, 15V, 24W, 100-240 VAC

USB Cord (with USB charging cradle)

008-0882-00	CABLE,USB 2.0,CHRGING/COMM CRADLE 6'
-------------	--------------------------------------

CDs/Guides

810-2711-02	KIT,CD,DFU ML,INSVC,PPQLT 1.6X
810-2775-02	KIT,CD,SW,PC UTIL 1.3X,ML,PPQLT (with USB charging cradle)
810-0409-10	BOOKLET,PRODUCTS AND ACCESSORIES
810-2212-04	INSTRUC, LANG SELECT, MULTILANG, PPQLT

Documentation/Cards/Label

810-0005-01	"DEAR CUSTOMER" NOTE
830-1067-03	DECLARATION OF CONFORMITY, PROPAQ LT
810-2816-00	ADDENDUM, PROPAQ LT 1.5X, MULTILANG
640-0647-01	LABEL, ICON DESCRIPTIONS
810-2663-00	PROPAQ LT VALUE CARD – Shipped w/Masimo Only

Others

630-0279-00	STRAP,CARRY,MOLDED,PROPAQ LT
900-0970-00	CONNECTOR PANEL PLUG,ONE,PROPAQ LT

Item Number	Item Description
-------------	------------------

LCDI ship-with cables

660-0279-00	CABLE ASSY,USB TYPE A/B,20
660-0272-00	CABLE ASSY,USB TYPE A/B, 15"
660-0269-00	CABLE ASSY,PWR,0.5M, PROPAQ LT DSPL CNVT
660-0268-00	CABLE ASSY,PWR,5.0M, PROPAQ LT DSPL CNVT


NOTE: All Ship-with accessories are also available for sale individually.

ADDITIONAL ACCESSORIES SOLD SEPARATELY

008-0290-00	Power cord DC 8' Cable for connection to fused DC power source (As in an Ambulance)
--------------------	--




008-0882-00	Cable, USB 2.0, 6' for Charging Cradle. (For connecting charging cradle to a PC for configuration or data download)
--------------------	--




503-0142-01	Power Supply, 15V, 24W, 100-240 VAC
008-0913-00	Power Cord, PLUG A , 2.5A, 125V, 2.5M
008-0394-00	Power Conditioner Isolation Transformer, Medical Grade, 120 VAC, includes Power Cord (For use w/008-0937-00)
008-0397-00	Power Conditioner Isolation Transformer, Medical Grade, 120 VAC, includes Wall Mount and Power Cord (For use w/008-0937-00)

900-0970-00	KIT, Connector covers for Propaq LT SpO ₂ , ECG, USB and Charging Ports
--------------------	--



900-0925-00	Patient Wearable Strap with Belt ; Box of 20
900-0924-00	Patient Wearable Strap with Belt ; Box of 100

Figure 40. Patient-Wearable Strap



900-0918-00	Patient Carrying Strap, Box of 5
-------------	----------------------------------



900-1008-00	Carrying Case, Propaq LT
-------------	--------------------------



ECG ACCESSORIES IN DETAIL

3-Lead ECG Cables


ONE PIECE CABLES – No lead wires required:

008-0880-00 ECG CABLE, 3 LEAD WIRES ATTACHED, **SNAP, AAMI**, 8'

DIN CABLES - separate lead wires required:


008-0892-00 ECG CABLE, 3 LEAD **DIN**, PROPAQ LT, **AAMI**, 8' (for use w/008-0323-00 lead wires)

Micropaq® and Propaq® LT ECG Cable, 3 Lead, DIN, AHA (AAMI), 8'



DIN LEAD WIRES AVAILABLE:

008-0323-00 LEAD WIRE SET, 3 LEAD,.060 SFTY **DIN SNAP AAMI**




5-Lead ECG Cables


ONE PIECE CABLES – No lead wires required:

008-0522-00 ECG CABLE, 5 LEADS WIRES ATTACHED **SNAP, AAMI**, 4'

Micropaq® and Propaq® LT Monitor ECG Cable, 5 Lead AHA (AAMI), 4'




008-0879-00 ECG CABLE, 5 ATTACHED LEAD, **SNAP, AAMI, 10'**



DIN CABLES - separate lead wires required:


008-0900-00 ECG CABLE, 5 LEAD **DIN**, MICROPAQ/PROPAQ LT, **AAMI, 2'**



008-0881-00 ECG Cable, 5-Lead Din, Micropaq, Propaq LT, 8' for use with 008-0322-00 (AAMI)

DIN LEAD WIRES AVAILABLE:

008-0322-00 LEAD WIRE SET, 5 LEAD, **DIN, SNAP** AAMI



008-0319-00 LEAD WIRE, 5L, **DIN**, ADULT **PNCH** AAMI

ECG ADAPTOR CABLE:

008-0734-00 ECG ADAPTOR CABLE, PROPAQ LT ECG CABLE TO AAMI PROPAQ CS, 6'

NOTE:

All DIN lead wires can be used with any DIN cable!

This includes the use of 3-lead DIN lead wire sets in 5-lead DIN ECG Cable

Note: Never use an ECG cable longer than 10 feet including extensions. If you use an ECG extension cable with an ECG cable longer than 4 feet, the monitor acts as though no ECG cable is connected.

LARGE COLOR DISPLAY

INTERFACE OPTION:

There are 2 different options for purchasing the Large Color Interface:

Option #1:

Customer already has a Propaq LT with a Charging Cradle/USB and wants to add the interface for a large color display. The customer will have to remove the bedrail holder from the back of cradle and mount the Interface to the cradle.

For this option you will need to order:

008-0917-00	Propaq LT Interface Box for Large Color Display
--------------------	--

This part number comes with 4 screws + 20" USB Interface Cable + Power Cable.





Figure 23. Mounting the Interface Box on the Cradle



Option #2:

Customer can order a Charging Cradle/USB with the Large Color Display Interface already mounted to it instead of a bedrail hook. Note - Propaq LT monitor will need to be ordered separately.

For this option you will need to order:

802LT0N-0E0	Propaq LT
802LT0S-0E0	<i>Requires charging cradle to recharge batteries. Charging cradle sold separately</i>
008-0931-00	Charging / Communications Cradle with USB and Large Color Display Interface

NOTES ON THE LARGE COLOR DISPLAY FLAT PANEL:

- Customer may purchase a flat panel display from us, or source their own. If they source their own display, the display must meet the specifications for which the interface software was written.

008-0937-00	17" Color Flat Panel Display for Propaq LT Interface
--------------------	--

- Customer may purchase one of two power conditioners for the flat panel display from us, or can source their own.

008-0394-00	Power Conditioner Isolation Transformer, Medical Grade
--------------------	--

008-0397-00	Power Conditioner Isolation Transformer, Medical Grade, with Wall Mount
-------------	---

Specifications for Reference Only

REVISIONS			
REV	CO NUMBER	DWN BY	CHK BY
A	N050159	JPB 17 Jun 05	DLC 20 Jun 05

DESCRIPTION:

17" Flat Panel Color Monitor

SPECIFICATION:

Size: 17 Inches (Diagonal)
 Resolution: 1,280 dots x 1,024 lines minimum
 Video Signal: Analog RGB (0.7V p-p, 75 Ω)
 Input Connections: 15 PIN D-SUB
 Power Input: 100-240VAC, 50-60Hz: removeable power cord
 Dimensions: (W x H x D) 14.5" x 16.7" x 87" (wih stand)
 Color: Gray
 Safety Marking: UL, cUL, CSA cUS, CE marked (TUV, VDE or Equiv.)
 Safety Conformance: UL 1950, CSA 22.2, No. 950, EN60950
 EMC Conformance: All models: 47CFR Part 15, Class B, EN55022 Class B, EN55024, EN61000-3-2, EN61000-3-3
 User Adjustments: Horizontal and Vertical position and size

This drawing may contain information which is the proprietary property of WELCH ALLYN.

This drawing is received in confidence and its contents may not be disclosed without the prior written consent of WELCH ALLYN.

WelchAllyn®

BEAVERTON, OR USA


DIMENSIONS ARE IN INCHES TOLERANCES ARE:	DWN BY: John P. Bender 07 Jun 05	DWG. TITLE: MONITOR, 17", FLAT PANEL	
	CHK BY: Donald L. Cox 20 Jun 05		

.XXX = \pm NA	_____	SIZE: A DRAWING NO: 008-0937-00 Rev. A	
.XX = \pm NA	_____		
ANGLES = \pm NA	_____	SCALE: NA SHEET 16 OF 25	
DO NOT SCALE DRAWING	_____		

MOUNTING OPTIONS: PROPAQ LT

IV POLE MOUNTING


008-0862-00	Propaq LT Charging Cradle IV Pole Mount
--------------------	---



MOBILE STAND MOUNTING


KIT:

008-0992-00	Propaq LT Roll Stand, Complete Kit Includes Roll Stand, Basket, Mounting Adapter and Power Adapter Bracket - NO OTHER PARTS NEED TO BE ORDERED
--------------------	--




Replacement/Individual Parts:

008-0987-00	Monitor Mounting Adapter, Propaq LT Roll Stand
--------------------	---




008-0994-00	Power Adapter Bracket, Propaq LT Roll Stand
--------------------	--




WALL MOUNTING

KITS:

008-0901-00	Wall Mount – Short Arm Kit – drop-in style, tilt/swivel. <i>This kit allows flush mounting of the monitor and power adapter to wall channel.</i> Includes wall channel, basket, mounting bracket and power adapter holder - NO OTHER PARTS NEED TO BE ORDERED
--------------------	--




008-0903-00	Wall Mount – Long Arm Kit – drop-in style, pivot block. <i>This kit allows pivot arm mounting of monitor and power adapter to wall channel.</i> Includes wall channel, basket, mounting bracket and power adapter holder - NO OTHER PARTS NEED TO BE ORDERED
--------------------	---



REPLACEMENT/INDIVIDUAL PARTS:


008-0203-00	Wall Channel, 19" (48.3 cm)
-------------	-----------------------------



008-0980-00	6" Basket with Power Adapter hook
-------------	-----------------------------------




008-0902-00	Adapter Plate, drop-in style
-------------	------------------------------



Note: this Adapter Plate allows the Propaq LT to drop on to any Propaq CS drop-on style mounting bracket!

008-0981-00	Power Supply Channel mounting bracket
-------------	---------------------------------------




MOUNTING OPTIONS: LARGE COLOR DISPLAY AND INTERFACE / PROPAQ LT

WALL MOUNTING:

KITS:


For Mounting the Flat Panel and Interface/Propaq LT together:

008-0983-00	<p>Large Color Display Interface and Flat Panel Mounting Kit.</p> <p><i>This kit allows wall mounting of monitor with LCD Interface, power supply and a Flat Panel Display.</i></p> <p>NO OTHER PARTS NEEDED – COMPLETE LCD/Interface KIT</p>
--------------------	--




For Mounting the Flat Panel Display Only

008-0982-00	<p>LCD, (Flat Panel) VESA Wall Mount Kit</p> <p><i>This is an industry-standard mounting bracket that will fit most flat-panel displays</i></p>
--------------------	---



REPLACEMENT/INDIVIDUAL PARTS:

008-0945-00	<p>Wall Mount Bracket for Large Color Display Interface – Comes with power supply holster</p>
--------------------	---



008-0203-00	<p>Wall Channel, 19" (48.3 cm) for LCD</p>
--------------------	--




008-0980-00	<p>6" Basket with Power Adapter hook</p>
--------------------	--



Monitor and Large Color Display Use Models

(A guide for ordering all necessary parts)

Use Model #1: Wall Mount for the Display, Propaq LT and Interface



Part Numbers needed for this Mounting Option

802LT0N-0E0 or 802LT0S-0E0 - Propaq LT without charging cradle

008-0931-00 - Charging cradle with USB and Large Color Display Interface attached

008-0937-00 - 17" Large Color Display (or customer can source a display that meets our specifications)


008-0983-00 - Large Color Display Interface and Flat Panel Mounting Kit

Replacement Cables if needed:

008-0947-00 - USB Cable, Type A/B, 20" for LCD Interface to Propaq LT (Gray Cable)

008-0949-00 - Power Cable, 20", for LCD Interface to Propaq LT (Black Cable)

Use Model #2: Wall Mount for the Display and Interface, Propaq LT at the Bedside



Part Numbers needed for this Mounting Option

802LT0N-UE1 or 802LT0S-UE1- Propaq LT with Charging Cradle/USB/Bedrail Hook


008-0984-00 - 15' USB + power cables and faceplate Kit

008-0917-00 - Propaq LT Interface Box for Large Color Display

008-0937-00 - 17" Large Color Display (or customer can source a display that meets our specifications)

008-0983-00 - LCD/Interface *complete mounting kit* (basket can be removed if not needed)

Use Model #3: Wall Mount for the Display, Propaq LT and Interface on a Table



Part Numbers needed for this Mounting Option

802LT0N-0E0 or 802LT0S-0E0 - Propaq LT without charging cradle

008-0931-00 - Charging cradle with USB and Large Color Display Interface attached

008-0937-00 - 17" Large Color Display

008-0982-00 - LCD, Vesa Wall Mount Kit (for Flat Panel only)

Replacement Cables if needed:

008-0947-00 - USB Cable, Type A/B, 20" for LCD Interface to Propaq LT (Gray Cable)

008-0949-00 - Power Cable, 20", for LCD Interface to Propaq LT (Black Cable)

Propaq LT Features –Brief

Size	5.4" x 7.5" x 2.0"
Weight	2lbs fully loaded
Ruggedness	75 g
Bedrail Handle	In charging cradle
Patient Modes	Adult, Ped, Neo
Battery Life	8 / 24 hrs (ECG only)
Display Type	Color LCD
Display Size	3.5"
Traces	2
Navigation	Soft Keys
Base Vitals	3/5 Lead ECG, SpO2, NIBP, RESP
Optional Parameters	None
Pulse Oximetry	Nellcor or Masimo
Noninvasive BP	Oscillometric Step Deflation, Smartcuf NIBP
Trends	Tab / Graph; 24 hours
Printer	Printing thru PC utility
Warranty	2 years

Other Vendors

For Additional Pulse Oximetry Accessories and related information, contact:



1-949- 297-7000, or 1-800-326-4890 Option #1 for Customer Service

www.masimo.com



1-800-635-5267

<http://respiratorysolutions.covidien.com/>

For additional Mounting Option Accessories and related information, contact:



1-800-228-2555

www.gcx.com