

GETINGE

GETINGE SIGNATURE SERVICE QUICK REFERENCE EQUIPMENT PARTS GUIDE



GETINGE

THE GETINGE GROUP is a leading global provider of equipment and systems that contribute to quality enhancement and cost efficiency within healthcare and life sciences. Equipment, services and technologies are supplied under the brands **ARJO** for patient hygiene, patient handling and wound care, **GETINGE** for infection control and prevention within healthcare and life science and **MAQUET** for surgical workplaces, cardiopulmonary and critical care.

To order a printed copy of the Getinge 2006 Quick Reference, use item number PGHC.

Table Of Contents

<u>Sterilizer Controls</u>	3
M/C II	3
PACS-2000 / 3000	4
<u>Sterilizers</u>	5
VINTAGE STERILIZERS	5
3500 SERIES	7
3600 SERIES	13
6900 SERIES	17
4200 (MPS) SERIES	19
100 SERIES	21
200 SERIES	23
400/500 SERIES	25
700 SERIES	27
800 SERIES	29
<u>Warming Cabinets</u>	33
5520/5525/5530	33
5540/5545/5550	34
5524/5618/5624	35
<u>20-30-36KW Steam Generators</u>	37
<u>Washers</u>	39
CASTLE 7515	39
CASTLE 7538	39
CASTLE 7550	40
2415 SONIC / RINSER DRYER	41
2460-UC SONIC CLEANER	43
2460-RD RINSER DRYER	44
GE FANUC CONTROLS	45
GE VERSAMAX CONTROLS	46
EAGLE SIGNAL 190+ CONTROLS	46
76-7700 CART WASHERS	47
78-7900 CART WASHERS	55
7642, 7652, 7662 INDEX WASHERS	59
7622 – CR22 DECONTAMINATOR	67
7626 – CR26 DECONTAMINATOR	71
4656 DISINFECTOR (DECOMAT)	75
46-X DISINFECTOR (PACS-300)	79
8666 DISINFECTOR (DECOMAT)	81
8666 DISINFECTOR (PACS-300)	87

Table Of Contents

<u>RO Water Systems</u>	89
HS-200 BOILER MATE	89
HS-450-PT	90
HS-1000	91
HS-2400	92
HS-4800	93
HS-9600	94
HS-15000	95
<u>Surgical Tables</u>	97
MAQUET 1132.0XX ALPHASTAR	97
MAQUET 1133.0XX ALPHAMAX	98
MAQUET 1130/1131/1420/1425/1532/1533	99
SHAMPAINE / ALM ST004	101
CMAX	105
<u>Surgical Lights</u>	107
ALM PRISMATIC ECL / PRC	107
ALM ANGENIEUX	108
ALM PRISMALIX (PRX)	109
ALM X'TEN & AXCEL	110
ALM ENERGIX (WPS II) CONTROLS	110
LDR (CONCEALED LIGHTING BOX)	111
21st CENTURY SERIES	112
HERAEUS 2000 SERIES	113
HERAEUS CITY SERIES	114
HERAEUS BLUE 30 & 80	114
CASTLE 2000 SERIES	115
CASTLE 2500-2900 SERIES	115
CASTLE 2300-2400 SERIES	116
GETINGE/CASTLE 500-700	116
<u>Manuals / Literature</u>	117
M/C STERILIZERS	117
PACS-2000 STERILIZERS	118
PACS-3000 STERILIZERS	119
STEAM GENERATORS/ANCILLARY	119
SONIC CLEANERS/RINSER DRYERS	120
WASHERS	120
WARMING CABINETS & AERATORS	121
OPERATING TABLES/LIGHTS	121

Sterilizer Controls

M/C CONTROLS

Power Supply Board (M/CII)	511411
Display Board	68950
Control Board	513427
Pressure Transducer M/C I	61301605830
Pressure Transducer M/C II	61301605830
RTD Temperature Probe (With Fittings)	530660

Solid-State Relays

AC Input (Yellow)	86407
AC Output (Black)	86408
DC Input (White)	518127
DC Output (Red)	56039
2.5 Amp (3600 Door)	69089

Batteries

Data Logger / Message Ctr 9V	86983
M/C I 8V	89280
M/C II 8V	511427

Circuit Breakers

2.0A White	86595
5.0A Orange	87489
7.5A White	86998
15A White	86596

Message Center

	513973
Printer Mechanism	536079
Paper Roll (8) Per Pkg.	61301600005
Ribbon Cassette	61301600006

Upgrades from MCII to PACS 2000 or 3000 are available.

Sterilizer Controls

PACS 2000 CONTROLS

Pressure Transducer	61301600636
Seal For Pressure Transducer	61301609008
RTD Temperature Probe	61301600638
CM1 Board (2A2U)	61301600314
CPU Board (2A1U)	61301600313
Input / Output Board (2A3U)	61301600263
Calibration Test Plug (Chamber, 8 pin)	61301602452
Calibration Test Plug (Jacket, 9 pin)	61301602453
Calibration Test Plug (Above Combined)	61301602980
Thermal Paper (10 Rolls)	61301601272
Take-Up Core (1) COBEX (Black)	61301603761
Take-Up Core (1) Cyber (White)	61301601611
Printer, Thermal (New, COBEX)**	61301603760
Dark Cover (Paper Compartment Door)	61301605881

*** If the printer has a toggle switch or a single black button for "Feed" & "Re-Print", order the replacement package 61301604246.*

PACS 3000 CONTROLS

On/Off Panel Switch (Green)	61301608437
Liquid Level Sensor W/Term	61301601173
RTD Temperature Probe	61301600638
RTD Load Probe	61301600231
Pressure Transducer	61301600636
Seal For Pressure Transducer	61301609008
PCB, CPU, - 1U	61301608297
PCB, CMI, - 2U	61301600314
OP 30 Assy., Display Panel	61301608168
PCB, EMI Filter, 18 VAC., - 4U	61301600312
PCB Assy., I/O I/F - 3U	61301600263
Printer, Thermal, RS-232, 24 VAC	61301603760
Thermal Paper (10 Rolls)	61301601272
Take-Up Core (1)	61301603761

Sterilizers

STRAIGHTLINE (RADIAL ARM DOORS)

Door Gaskets For Sentry & 3000 Series (Also Fits T-60)

1624 Cylindrical (Red or Grey)	44599
1616 Rectangular "	44600
2020 Rectangular "	44602
2036 Cylindrical "	44601

STRAIGHTLINE VERTIGLIDE 3100, 3300

Door Gaskets See 3500 Series (Page 5)

Door Components

The cable assembly listed below is for the door system using 2 cylindrical counter weights.

Cable	59851
Sleeve	58367
Thimble	58368

Door Switch	87318
-------------	-------

3200, 3400, ELECTRILOK

Door Gaskets

D-Shape (Blue or Red)	60864
V-Shape (Red)	517452

Door Components

Door Motor (W/brk't. Replaces 52083)	536846
Door Motor Capacitor	60788
Relay 120VAC,DPDT	68949
Cam Switch ("Open", "Closed", "Sealed")	513384

Emergency Crank Handle	26518
-------------------------------	--------------

Sterilizers

Sterilizers

3500 SERIES

Controls (See M/C II Controls Page 3)

Door Gaskets

1616 Original (Blue)	61301605017
U-Shape (Green)	513966
Oversize (Grey)	517738
2020 Original (Blue)	61301605016
U-Shape (Blue)	513967
Oversize (Grey)	517739

Door Operation Components

Pressure Switch (Door Gasket)	517492
Door Closed Switch (New Style)*	517458
Door Closed Switch (Old Style)*	511410
Door Cable (New W/Turnbuckle & Eye)	511400
Door Handle (Lexan)	511780

**The old style switch is separate from the cord. The new style switch consists of the cord and the switch as one molded assembly.*

Plumbing Components

Safety Valve (40 lbs)	61301604180
Check Valve 3/4" (ASCO)	513363
Rebuild	513369
Check Valve 3/8" (ASCO)	513361
Rebuild	513369
Strainer 3/4"	36355
Replacement Screen	513448
Steam Trap 1/2" (Jacket & Chamber)	61301602492
"Air In" Filter	61301600303
Mini-Valve (Black Knob)	11929
Stone Filter Disc	61301613275
Baffle Stud Repair Kit (Requires 3)	530625
Water Ejector	59858

Chamber Drain Screen (Pin Trap)	07181
3525 Drain Screen	87212

Sterilizers

All listed solenoid valves have green molded coils. To change from a coil with a blue casing to a green molded coil, order one of the following:

4' Cord	517360
6' Cord	517361
8' Cord	517362
10' Cord	517363

3522 VALVES

1SV Steam To Chamber	517339
Coil	517300
Rebuild	61301601739
2SV Chamber Drain (Rubber Diaphragm)	N/A
Coil	N/A
Rebuild	54190
2SV Chamber Drain (Steel Diaphragm)	541581
Coil	517303
Rebuild	61301600443
5SV Cooling Water	517338
Coil	517299
Rebuild	68164
6SV Water to Ejector	517349
Coil	517301
Rebuild	56161
5/6SV Cooling/Water to Ejector (Dual)	517356
Coil	517301
Rebuild	61301601738
8SV Steam to Jacket	517339
Coil	517300
Rebuild	61301601739
10/11SV Door Seal/Unseal (Brass)	517338
Coil	517299
Rebuild	68164
10/11SV Door Seal/Unseal (Stainless)	541582
Coil	517303
Rebuild	61301600442

Sterilizers

3523 VALVES

1SV Steam To Chamber (Steel Diaphragm)	541581
Coil	517303
Rebuild	61301600443
2SV Chamber Drain (Rubber Diaphragm)	N/A
Coil	N/A
Rebuild	54190
2SV Chamber Drain (Steel Diaphragm)	541581
Coil	517303
Rebuild	61301600443
3SV Air In	517347
Coil	517305
Rebuild	60301601740
4SV Multitherm Bypass	517339
Coil	517300
Rebuild	60301601739
5SV Cooling Water	517354
Coil	517301
Rebuild	517298
6SV Water to Ejector	517348
Coil	517301
Rebuild	56161
5/6SV Cooling/Water to Ejector (Dual)	61301601577
Coil	517301
Rebuild	61301601738
7SV Multitherm	517353
Coil	517299
Rebuild	61301601739
8SV Steam to Jacket	517339
Coil	517300
Rebuild	61301601739
10/11SV Door Seal/Unseal (Brass)	517338
Coil	517299
Rebuild	68164
10/11SV Door Seal/Unseal (Stainless)	541582
Coil	517303
Rebuild	61301600442

Sterilizers

3533 VALVES

1SV Steam To Chamber	517339
Coil	517300
Rebuild	60301601739
2SV Chamber Drain (Rubber Diaphragm)	N/A
Coil	N/A
Rebuild	54190
2SV Chamber Drain (Steel Diaphragm)	541581
Coil	517303
Rebuild	61301600443
3SV Air In	517347
Coil	517305
Rebuild	61301601740
5SV Cooling Water	517338
Coil	517299
Rebuild	68164
6SV Water to Ejector	517349
Coil	517301
Rebuild	56161
5/6SV Cooling/Water to Ejector (Dual)	517356
Coil	517301
Rebuild	61301601738
8SV Steam to Jacket	517339
Coil	517300
Rebuild	61301601739
10/11SV Door Seal/Unseal (Brass)	517338
Coil	517299
Rebuild	68164
10/11SV Door Seal/Unseal (Stainless)	541582
Coil	517303
Rebuild	61301600442

Sterilizers

3525 VALVES

1SV Steam To Chamber	517339
Coil	517300
Rebuild	61301601739
2SV Chamber Drain (Rubber Diaphragm)	N/A
Coil	N/A
Rebuild	54190
2SV Chamber Drain (Steel Diaphragm)	541581
Coil	517303
Rebuild	61301600443
2DV Drain Diaphragm	37829
Diaphragm Rebuild	517153
6SV Water to Ejector	517349
Coil	517301
Rebuild	56161
10/11SV Door Seal/Unseal (Brass)	517338
Coil	517299
Rebuild	68164
10/11SV Door Seal/Unseal (Stainless)	541582
Coil	517303
Rebuild	61301600442
12SV Pilot Air to Drain Diaphragm SS	541582
Coil	517303
Rebuild	61301600442
13SV Spray Water	517343
Coil	517301
Rebuild	54187
14SV Detergent	14029
Coil	N/A
Rebuild (plunger)	21513
16SV Vent	517339
Coil	517300
Rebuild	61301601739

Sterilizers

Sterilizers

3600 SERIES

Controls (See M/C II Controls Page 3)

Door Gasket (Red)	536651
Door Gasket (Blue)	61301602425

Door Operation Components

Hall Effect Position Switch 5V	518083
--------------------------------	--------

Power Door

Chamber Press Switch	530931
Door Motor & Cord	539384

Manual Door

Air Cylinder	541577
Water Cylinder	518146
Hook Cylinder Side (Water Only)	518386

Refer to the valves list for the solenoid and diaphragm.

Plumbing Components

Safety Valve (40 lbs)	56913
1/2" Check (ASCO)	513362
Rebuild	513369
Check Valve 3/4" (ASCO)	513363
Rebuild	513369
Check Valve 1" (ASCO)	513699
Rebuild	513700
Strainer 3/4"	36355
Replacement Screen	513448
Strainer 1"	50952
Replacement Screen	513449
Steam Trap 1/2"	61301602492
Steam Trap 3/4"	23614
Rebuild	27557
"Air In" Filter (Needs 2)	61301600303
Heat Exchanger (Brass 5 X 14)	66783
Water Ejector	60699

Chamber Drain Screen (Pin Trap)	23815
--	-------

Sterilizers

All listed solenoid valves have green molded coils. To change from a coil with a blue casing to a green molded coil, order one of the following:

4' Cord	517360
6' Cord	517361
8' Cord	517362
10' Cord	517363

3622 VALVES

1SV Steam To Chamber	61301604083
Coil	517304
Rebuild	61301604084
2SV Chamber Drain	61301604083
Coil	517304
Rebuild	61301604084
3SV Air In (Steel Diaphragm)	541581
Coil	517303
Rebuild	61301600443
5/6SV Wtr. To Drn & Ejector (Dual)	517356
Coil	517301
Rebuild	61301601738
7SV Multitherm	517339
Coil	517300
Rebuild	61301601739
8SV Steam to Jacket	517342
Coil	517299
Rebuild	54186

Manual Doors

12SV Air To Door Cylinder	517357
Coil	517303
Rebuild	68381
12SV Water To Door Cylinder	541627
Coil	517303
Rebuild	61301601837
Diaphragm Valve (Safety Shut-Off)	61301605027
Rebuild	61301603671
4 Way Valve (Air Only)	541576
4 Way Valve (Water Only)	517523

Sterilizers

3633 VALVES

1SV Steam To Chamber	61301604083
Coil	517304
Rebuild	61301604084
2SV Chamber Drain	61301604083
Coil	517304
Rebuild	61301604084
3SV Air In (Steel Diaphragm)	541581
Coil	517303
Rebuild	61301600443
4SV Multitherm Bypass	517339
Coil	517300
Rebuild	61301601739
5/6SV Wtr. To Drn & Ejector (Dual)	517356
Coil	517301
Rebuild	61301601738
7SV Multitherm	517339
Coil	517300
Rebuild	61301601739
8SV Steam to Jacket	517342
Coil	517299
Rebuild	54186

Manual Doors

12SV Air To Door Cylinder	517357
Coil	517303
Rebuild	68381
12SV Water To Door Cylinder	541627
Coil	517303
Rebuild	61301601837
Diaphragm Valve (Safety Shut-Off)	61301605027
Rebuild	61301603671
4 Way Valve (Air Only)	541576
4 Way Valve (Water Only)	517523

Sterilizers

Sterilizers

6900 SERIES

Controls (See M/C II Controls Page 3)

Door Gasket Original Round (Green)	40441
U-Shaped (Green)	513968
Oversized (Grey)	517553

Door Operation Components

"Door Closed" Internal Switch (1S/4S)	61301605325
"Door Jammed" Internal Switch (3S/6S)	61301605325
"Door Opened" Switch (2S/5S)	50982
Door Lift Screw	40458
Door Lift Nut (Brass)	40460
Timing Belt	47392
Motor Only	533596
Mounting Brackets (Need 2)	536663
Pressure Switch	517493
"U"-Bracket For New Style Pressure Sw	517494

Plumbing Components

Safety Valve (40 lbs)	56913
Check Valve 1/4" (Stainless Steel)	66838
Check Valve 3/8"	19112
Check Valve 3/4" Inline (Conbraco)	08045
Check Valve 1" Inline (Conbraco)	21815
Check Valve 1-1/4" Inline	14882
Strainer 1"	50952
Replacement Screen	513449
Steam Trap 3/4"	23614
Rebuild	27557
"Air In" Filter (Needs 2)	61301600303

Spray Arms

Outer Spinner Bearing (Short - 1/2")	41135
Inner Spinner Bearing (Long-3/4")	41136
End Plug	536958
Set Screw For End Plug	E5008

Sterilizers

6900 AMPS VALVES

1SV Steam To Chamber	517342
Coil	517299
Rebuild	54186
22SV Trap Drain (Steel Diaphragm)	541581
Coil	517303
Rebuild (Plunger)	61301600443
30M Chamber Drain (Ball Valve)	61301604372
3SV Air In	517342
Coil	517299
Rebuild	54186
6SV Water to Ejector	517343
Coil	517301
Rebuild	54187
14SV Detergent	14029
Coil	N/A
Rebuild (Plunger)	21513
15SV Steam To Sump (Steel Diaphragm)	541581
Coil	517303
Rebuild (Plunger)	61301600443
17SV Hot Water	517341
Coil	517301
Rebuild	61301601818
26SV Cold Water	517341
Coil	517301
Rebuild	61301601818
27SV Condenser Water	517340
Coil	517301
Rebuild	54198
28SV Special Water	517341
Coil	517301
Rebuild	61301601818
3-Way Air (SS)	541582
Coil	517302
Rebuild	54188

Sterilizers

4200 (MPS) SERIES

Controls (See M/C II Controls Page 3)

Door Gasket D-Shape (Blue or Red)	60866
V-Shape (Red)	517453

Door Operation Components

Door Safety Bellows & Housing (MPS)	35789
Door Motor (W/brk't. Replaces 52083)	536846
Door Motor Capacitor	60788
Relay 120VAC,DPDT	68949
Cam Switch ("Open", "Closed", "Sealed")	513384
Hook Lift Lever	27900

Plumbing Components

Safety Valve (40 lbs)	56916
Check Valve 3/4" (ASCO)	513363
Rebuild	513369
Check Valve 3/4" Inline (Conbraco)	08045
Check Valve 1" (ASCO)	513699
Rebuild	513700
Check Valve 1 1/2" Inline (Conbraco)	61301602477
Strainer 3/4"	36355
Replacement Screen	513448
Strainer 1"	50952
Replacement Screen	513449
Steam Trap 3/8"	68722
Steam Trap 3/4"	23614
Rebuild	27557
"Air In" Filter (Needs 2)	61301600303
1/4" Air Reg.	32696
Water Ejector (42" Unit)	66784
Water Ejector (76" Unit)	66762

Chamber Drain Screen (Pin Trap)	36690
--	-------

Sterilizers

All listed solenoid valves have green molded coils. To change from a coil with a blue casing to a green molded coil, order one of the following:

4' Cord	517360
6' Cord	517361
8' Cord	517362
10' Cord	517363

4200 MPS VALVES

1SV Steam To Chamber (Steel Diaphragm)	541581
Coil	517303
Rebuild	61301600443
2SV Pilot Air to Drain Diaphragm SS	541582
Coil	517303
Rebuild	61301600442
2DV Drain Diaphragm	37831
Rebuild	517154
3SV Air In (Steel Diaphragm)	541581
Coil	517303
Rebuild	61301600443
5SV Cooling Water	517343
Coil	517301
Rebuild	54187
6SV Water to Ejector (42")	517341
Coil	517301
Rebuild	61301601818
6SV Pilot Air to Water Diaphragm SS (76")	541582
Coil	517303
Rebuild	61301600442
6DV Water Diaphragm (76" Only)	37829
Rebuild	517153
8SV Steam to Jacket	517342
Coil	517299
Rebuild	54186

Sterilizers

100 SERIES

Controls (See PACS-2000 Controls Page 4)

Door Gaskets

1616 Original (Blue)	61301605017
U-Shape (Green)	513966
Oversize (Grey)	517738
2020 Original (Blue)	61301605016
U-Shape (Blue)	513967
Oversize (Grey)	517739

Door Operation Components

Gasket Pressure Switch	61301600350
Door Switch (3S,7S) Center Mount	61301601670
Door Cable (W/Turnbuckle & Eye)	511400
Door Handle	511780

Plumbing Components

Safety Valve (40 lbs)	61301604180
Check Valve 3/4" (ASCO)	513363
Rebuild	513369
Check Valve 3/8 (ASCO)	513361
Rebuild	513369
Strainer 3/4"	36355
Replacement Screen	513448
Drain Temp Switch	539126
Float Switch "Water In Drain"	61301601173
Steam Trap 1/2" (Jacket & Chamber)	61301602492
"Air In" Filter	61301600303
Mini-Valve (Black Knob)	11929
Stone Filter Disc	61301613275
Baffle Stud Repair Kit (Requires 3)	530625
Water Ejector	59858

Chamber Drain Screen (Pin Trap)	07181
--	-------

Sterilizers

100 SERIES VALVES

1SV Steam To Chamber	61301601697
Coil	61301600154
Rebuild	61301601739
2SV Chamber Drain	61301601701
Coil	61301600157
Rebuild	61301600443
3SV Air In	61301601698
Coil	61301600156
Rebuild	61301601740
4SV Multitherm Valve	61301601697
Coil	61301600154
Rebuild	61301601739
5/6SV Cooling/Water to Ejector (Dual)	61301601696
Coil	61301600153
Rebuild	61301601738
7SV Multitherm	61301601700
Coil	61301600157
Rebuild	61301601742
8SV Steam to Jacket (Black Coil)	61301609872
Coil	61301604049
Rebuild	N/A
8SV Steam to Jacket (Green Coil)	61301601697
Coil	61301600154
Rebuild	61301601739
10,11,12,13SV Door Seal/Unseal	61301601695
Coil	61301600153
Rebuild	61301601737
Black to Green Coil Cord Adapter	61301609802

Sterilizers

200 SERIES

Controls (See PACS-2000 Controls Page 4)

Door Gasket 61301600108

Door Operation Components

Door Plate Switch 61301602316

Hall Effect Position Switch 61301600639

Power Door

Chamber Press Switch 530931

Door Motor 539384

Door Motor Capacitor 518162

Door Motor Relay 61301601073

Manual Door

Air Cylinder 541577

Refer to the valves list for the solenoid and diaphragm.

Plumbing Components

Safety Valve (45 lbs) 61301601576

Check Valve 3/4" (ASCO) 513363

Rebuild 513369

Strainer 3/4" 36355

Replacement Screen 513448

Strainer 1" 50952

Replacement Screen 513449

Drain Temp Switch 539126

Float Switch "Water In Drain" 61301601173

Steam Trap 1/2" (Jacket) 61301602492

Steam Trap 3/4" (Chamber) 23614

Rebuild 27557

Steam Trap 3/4" (Chamber HC/LS Only) 61301603751

Rebuild 61301603755

"Air In" Filter (Needs 2) 61301600303

Heat Exchanger (Brass 5 X 14) 66783

Water Ejector 60699

Chamber Drain Screen Assy. (Pin Trap) 61301604624

Replacement Screen 61301600512

Sterilizers

200 SERIES VALVES

1SV Steam To Chamber	61301601701
Coil	61301600157
Rebuild	61301600443
2SV Chamber Drain	61301601702
Coil	61301600429
Rebuild	61301604084
3SV Air In	61301601701
Coil	61301600157
Rebuild	61301600443
5/6SV Cooling/Water to Ejector (Dual)	61301601696
Coil	61301600153
Rebuild	61301601738
8SV Steam to Jacket	61301603374
Coil	61301600157
Rebuild	54186

Manual Doors Only

12SV Air To Door Cylinder (Manual Door)	61301601836
Coil	61301600157
Rebuild	61301601837
Diaphragm Valve (Safety Shut-Off)	61301605027
Rebuild	61301603671
4 Way Valve	541576

Sterilizers

400/500 SERIES

Controls (See PACS-3000 Controls Page 4)

Door Gaskets

400 Series Blue Gasket With Red Stripe	61301609895
Convert From Graphite Kit	61301609371
500 Series Blue Gasket With Red Stripe	61301609896
Convert From Graphite Kit	61301609373

Door Operation Components

Door "Not Closed" Switch 3S/7S	61301608277
--------------------------------	-------------

Power Door Components

Door "Opened" Switch	61301608277
Door Motor	61301601143
Capacitor, Door Motor Run 2.5 MFD	61301608490
Resistor, PWR, 50W, 400 Ohm	61301608173
Relay, DPDT 24VAC	61301601073

Plumbing Components

Safety Valve (45 lbs)	61301608100
Check Valve 3/4" (ASCO)	513363
Rebuild	513369
Check Valve 3/8" (ASCO)	513361
Rebuild	513369
Strainer 3/4"	36355
Replacement Screen	513448
Drain Temp Switch	539126
Float Switch "Water In Drain"	61301601173
Steam Trap 1/2" (Jacket)	61301602492
Steam Trap 3/8" (Drip Leg)	68722
"Air In" Filter	61301600303
Mini-Valve (Black Knob)	11929
Water Ejector	59858

Chamber Drain Screen (Pin Trap)	61301608230
Secondary Chamber Bleed Drain Screen	61301608303

Sterilizers

400/500 SERIES VALVES

1-8-10-12SV Chmbr,Jckt, Door (Green Coil)	61301601697
Coil	61301600154
Rebuild	61301601739
1-8-10-12SV Chmbr,Jckt, Door (Black Coil)	61301609872
Coil	61301604049
Rebuild	N/A
2SV Chamber Drain	61301601701
Coil	61301600157
Rebuild	61301600443
3SV Air In (Green Coil)	61301608423
Coil	61301600156
Rebuild	61301608453
3SV Air In (Black Coil)	61301609450
Coil	61301609452
Rebuild	61301609451
5/6SV Wtr. To Drn & Ejector (Dual)	61301601696
Coil	61301600153
Rebuild	61301601738
5SV Cooling Water (Black Coil)	61301608583
Coil	61301608962
Rebuild	61301608963
6SV Water To Ejector (Black Coil)	61301609315
Coil	N/A
Rebuild	61301609505
14SV Secondary Drain (Green Coil)	61301601695
Coil	61301600153
Rebuild	61301601737
14SV Secondary Drain (Black Coil)	61301608583
Coil	61301608962
Rebuild	61301608963

Black to Green Coil Cord Adapter

61301609802

Sterilizers

700 SERIES

Controls (See PACS-3000 Controls Page 4)

Door Gasket (Blue)	61301609366
Convert From Graphite Kit	61301609375

Door Operation Components

Door Position Switch DPST 6A. ANCR	61301602316
Door Position Switch 10A NO/NC Snap	61301608795
Door Motor Control Board 90VDC	61301608595
Clutch (Torque Coupler)	61301608991
Relay, 3PDT 24VAC	61301608830
Relay, Latch ,MAG,DPDP,24VAC	61301604716

Plumbing Components

Safety Valve (45 lbs)	61301608100
Check Valve 3/4" (ASCO)	513363
Rebuild	513369
Check Valve 3/8" (ASCO)	513361
Rebuild	513369
1" Strainer	50952
Replacement Screen	513449
Drain Temp Switch	539126
Float Switch "Water In Drain"	61301601173
Steam Trap 1/2" (Jacket)	61301602492
Steam Trap 3/8" (Drip Leg)	68722
1/4" Air Reg.	32696
Heat Exchanger	61301608571
"Air In" Filter (Needs 2)	61301600303
Mini-Valve (Black Knob)	11929
Water Ejector	60699

Chamber Drain Screen (Pin Trap)	61301608230
Secondary Chamber Bleed Drain Screen	61301608303

Sterilizers

700 SERIES VALVES

1SV Steam To Chamber (Green Coil)	61301603374
Coil	61301600157
Rebuild	54186
1SV Steam To Chamber (Black Coil)	61301609872
Coil	61301604049
Rebuild	N/A
2V Chamber Drain (Pneumatic)	61301608458
Actuator Seal Kit	61301608964
Rebuild	61301608965
3V Air In (Pneumatic)	61301608458
Actuator Seal Kit	61301608964
Rebuild	61301608965
5/6SV Wtr. To Drn & Ejector (Dual)	61301601696
Coil	61301600153
Rebuild	61301601738
5SV Cooling Water (Black Coil)	61301608852
Coil	61301609378
Rebuild	61301608853
6SV Water To Ejector (Black Coil)	61301609315
Coil	N/A
Rebuild	61301609505
8SV Steam to Jacket (Green Coil)	61301603374
Coil	61301600157
Rebuild	54186
8SV Steam to Jacket (Black Coil)	61301609872
Coil	61301604049
Rebuild	N/A
10/12SV Door Seal/Unseal (Green Coil)	61301601697
Coil	61301600154
Rebuild	61301601739
10/12SV Door Seal/Unseal (Black Coil)	61301609872
Coil	61301604049
Rebuild	N/A
14SV Secondary Drain	61301608583
Coil	61301608962
Rebuild	61301608963
Air Pilot Valve (3/2 .250, 24V)	VALTE0013

Sterilizers

800 SERIES

Controls (See PACS-3000 Controls Page 4)

Door Gasket (Blue)	61301609559
Convert From Graphite Kit	61301609903

Door Operation Components

Door Gasket Pressure Switch 30-60 PSI	61301608845
"Door Open-Decel-Swing-In" Switch	61301608795
"Door Closed" Switch (3/7S)	61301602316
"Slide Out" Door Switch (25/26S)	61301608277
Door Motor Control Board 90VDC	61301608595
Clutch (Torque Coupler)	61301608991
Relay, DPDT 24VAC	61301601073
Relay, Latch ,MAG,DPDP,24VAC	61301604716

Plumbing Components

Safety Valve (45 lbs)	61301601576
Check Valve 3/4" (ASCO)	513363
Rebuild	513369
Check Valve 3/8 (ASCO)	513361
Rebuild	513369
1" Strainer	50952
Replacement Screen	513449
Drain Temp Switch	539126
Float Switch "Water In Drain"	61301601173
Steam Trap 1/2" (42" Only)	61301602492
Steam Trap 3/4" (76" Only)	61301603751
Rebuild	61301603755
1/4" Air Reg.	32696
Heat Exchanger	61301608571
"Air In" Filter (Needs 2)	61301600303
Water Ejector (42" Unit)	60699
Water Ejector (76" Unit)	66762

Chamber Drain Screen (Pin Trap)	61301609474
Secondary Chamber Bleed Drain Screen	61301608303

Sterilizers

800 SERIES 42" VALVES

1SV Steam To Chamber (Green Coil)	61301603374
Coil	61301600157
Rebuild	54186
1SV Steam To Chamber (Black Coil)	61301609872
Coil	61301604049
Rebuild	N/A
2SV Chamber Drain (Pneum)	61301609873
Actuator Seal Kit	61301609379
Rebuild	61301609893
3SV Air In (Pneum)	61301608458
Actuator Seal Kit	61301608964
Rebuild	61301608965
5SV Cooling Water	61301608852
Coil	61301609378
Rebuild	61301608853
6SV Water to Ejector	61301609315
Coil	N/A
Rebuild	61301609505
8SV Steam to Jacket (Green Coil)	61301603374
Coil	61301600157
Rebuild	54186
8SV Steam to Jacket (Black Coil)	61301609872
Coil	61301604049
Rebuild	N/A
10/12SV Door Seal/Unseal (Green Coil)	61301603374
Coil	61301600157
Rebuild	54186
10/12SV Door Seal/Unseal (Black Coil)	61301609872
Coil	61301604049
Rebuild	N/A
14SV Secondary Drain	61301608583
Coil	61301608962
Rebuild	61301608963
Air Pilot Valve (Head & Tail)	VALTE0012
Air Pilot Valve (3/2 .250, 24V)	VALTE0013
Black to Green Coil Cord Adapter	61301609802

Sterilizers

800 SERIES 76" VALVES

1SV and 8SV Solenoids to Pneum. Pkg.	61301609925
1V Steam to Chamber (Pneum)	61301609884
Actuator Seal Kit	61301608964
Rebuild	61301609904
2SV Chamber Drain (Pneum)	61301609874
Actuator Seal Kit	61301609379
Rebuild	61301609894
3SV Air In (Pneum)	61301608458
Actuator Seal Kit	61301608964
Rebuild	61301608965
5SV Cooling Water (Black Coil)	61301609872
Coil	61301604049
Rebuild	N/A
6V Water to Ejector (Pneum)	61301608458
Actuator Seal Kit	61301608964
Rebuild	61301608965
8SV Solenoid	See 1SV
8V Steam to Jacket (Pneum)	61301609884
Actuator Seal Kit	61301608964
Rebuild	61301609904
10/12SV Door Seal/Unseal (Green Coil)	61301603374
Coil	61301600157
Rebuild	54186
10/12SV Door Seal/Unseal (Black Coil)	61301609872
Coil	61301604049
Rebuild	N/A
14SV Secondary Drain (Black Coil)	61301609872
Coil	61301604049
Rebuild	N/A
Air Pilot Valve (Head & Tail)	VALTE0012
Air Pilot Valve (3/2 .250, 24V)	VALTE0013
Black to Green Coil Cord Adapter	61301609802

Sterilizers

Warming Cabinets

5520/5525/5530

Controls

Display Board (Obs.Order MP-126C)	61301603218
Heater, 120V/500W	69093
Temp Sensor (2 lead)	69200
Temp Sensor (3 Lead)	517296
Circuit Breaker, 10A	69110
Circuit Breaker, 20A	69109
Power Supply Assy.	69155
Transformer	69090
Fan, 120V	69094
Potentiometer	60768
Overheat Thermostat 212° N.O.	69107
Relay, 45A (Heater)	514107
Relay, 2.5A (Fan)	61301605584
Switch, Door	69249

Misc. Components

Door Gasket (Perimeter)	69037
Shelf / Compartment Gasket	69216
Knob	69237
Foot (5530)	14813
Foot , Screw	E5210
Catch, Magnetic (On Door)	60624
Door Handle 12"	69056
Door Handle 18"	69865

Warming Cabinets

5540/5545/5550

Controls

Display Board (Obs.Order MP-126D)	61301602894
Overheat Thermostat 212° N.O.	69107
Fuse, 3A	72441
Fuse, .25A	69425
Heater, 120V/500W	69093
Fan, 120V	69094
Circuit Breaker, 10A	69110
Circuit Breaker, 20A	69109
Power Supply Board	513013
Transformer (24V)	69090
Temp. Probe	518407
Relay, 45A (Heater)	514107
Relay, 2.5A (Fan)	61301605584
Thermostat, 212	69107
Door Switch	69249

Misc. Components

Door Gasket (Perimeter)	69037
Shelf / Compartment Gasket	69216
Catch, Magnetic (On Door)	60624
Door Handle 12"	69056
Door Handle 18"	69865

Warming Cabinets

5524/5618/5624

Temperature Probe	61301600624
Display / CPU Board	61301600382

Power Supply Board (6 pin J2 Connector)	61301601841
---	-------------

If original power supply has a 7 pin J2 connector, order the package 61301604375.

Stand-Off For Power Supply Board	61301605229
----------------------------------	-------------

Solid State Relays – Fan & Heater (4)	61301605300
---------------------------------------	-------------

Printed Circuit Board Cont/Display	61301600382
------------------------------------	-------------

Fuse, 12A, 2FU, Full size Warmer (115V)	61301600840
---	-------------

Fuse, 5A, 2FU, Half Size Warmer (115V)	61301601452
--	-------------

Fuse, 6.25A, 2FU, Full Size Warmer (230V)	61301600841
---	-------------

Fuse, 3A, 2FU, Half Size Warmer (230V)	67621
--	-------

Fuse Holder (Set Of 12)	87848
-------------------------	-------

Fan, Full & Half Size (115V) Internal	61301600439
---------------------------------------	-------------

Fan, Full & Half Size (230V) Internal	61301600440
---------------------------------------	-------------

Fan, Flat Pak (Top, Back Of Unit) 115V	61301602603
--	-------------

Fan, Flat Pak (Top, Back Of Unit) 230V	61301602604
--	-------------

Heater, Full & Half Size (115V)	61301600474
---------------------------------	-------------

Heater, Full & Half Size (230V)	61301600611
---------------------------------	-------------

6' Optional Power Cord Kit	61301605436
----------------------------	-------------

New Style Hinge Kit (To Convert From Old)	61301603505
---	-------------

Shelf Support Clip	910312600
--------------------	-----------

New Style Hinge Only	912032200
----------------------	-----------

Warming Cabinet Serial No. w/ "SW" Use The Following

5618/5624 Side Gasket (Qty.2)	61301603370
-------------------------------	-------------

5524 Side Gasket (Qty.2)	61301603368
--------------------------	-------------

5524/5618/5624 Center Gasket (Qty.2)	61301603369
--------------------------------------	-------------

Warming Cabinet Serial No. w/ "SZ" Use The Following

5524/5618/5624 Gasket (By The Foot)	R4065
-------------------------------------	-------

Warming Cabinets

20-30-36KW Steam Generators

Universal Components

Element Gasket	49943
Element Bolts(6)	49944
3-Pole Contactor 3PH	49938
4-Pole Contactor Single Phase	51012
Sight Glass With O-Rings 9-3/4"	62790
Sight Glass With O-Rings 7-7/8"	49958
Sight Glass Valves (Pair)	49957
Angled Sight Glass Valves	61250
Spark Plug Probe	49940
Probe Insulator	68308
1/4" Feed Water Solenoid Valve (CAS-20)	61245
1/2" Feed Water Solenoid Valve (CAS-30)	536586
3/8 NPT Check Valve	19112
1/2"NPT Check Valve	19114
100 PSIG Safety Valve	49949
Pressure Control (Operating)	49947
Pressure Control (Hi Limit W/man reset)	49948
Pressure Gauge	49941
Level Control Board	68307

Elements CAS-20 C,D,E (17 inch)

240 V 1ph 20KW	513944
208 V 3ph 20KW	513950
240 V 3ph 20KW	513946
380 V 3ph 20KW	513947
415 V 3ph 20KW	513943
480 V 3ph 20KW	513948
600 V 3ph 20KW	513945

Elements CAS-20A&AF (12 inch)

240 V 1ph 20KW	49935
240 V 3ph 20KW	49934
208 V 3ph 20KW	55967
380 V 3ph 20KW	68770
480 V 3ph 20KW	68718
Wire Assy From Contactor To Elements (1)	531065

20-30-36KW Steam Generators

Elements CAS-30/36

12KW

208 V 3ph 12KW	536522
240 V 1ph 12KW	536523
240 V 3ph 12KW	536524
380 V 3ph 12KW	536575
415 V 3ph 12KW	536576
480 V 3ph 12KW	536577
600 V 3ph 12KW	536578

18KW

208 V 3ph 18KW	536515
240 V 1ph 18KW	536516
240 V 3ph 18KW	536517
380 V 3ph 18KW	536518
415 V 3ph 18KW	536519
480 V 3ph 18KW	536520
600 V 3ph 18KW	536521

NOTE: Elements for all boilers ship with gasket, sealer and 6 Bolts

PROBES

1-1/2"	531032
5"	61243
5-3/4"	61241
6-3/8"	61242

PUMPS

30 & 36 KW Pump & Motor	536585
Procon Pump	58630
Motor	58629
Teel Pump	49952
CAS-30C Pump 1.6GPM	61301608522
Plastic Flo-Jet Pump	513957

Washers

CASTLE 7515

Fuse - 3A	72441
Fuse - .25A	69425
Battery - 1.5V (D)	60670
Battery Backup	530473
Det. & San. Pump Tube	511898
Viton Diaph. (Use w Tec-San)	61301600664
EPDM Diaph. (Use w Tec-Wash III)	61301600663
Post Bearing - Top Spinner	530721
Post Bearing - Bottom Spinner	530722
Spinner - Top & Bottom	541534
Spinner Cap	530720
Spray Nozzle	69657
Circulating Pump Repair Pkg.	68393
Dolly Wheels	536588
Door Gasket	69628

CASTLE 7538

Rinse Agent Replacement Tube – 21"	511898
Det. Replacement Tube – 16"	511899
Mist Eliminator Gasket	530743
Sump Drain Cyl. Assy.	78499
Replacement Stopper	72084
Flex Duct 60"	517178
Replacement Clamp	517181
Mixing Valve Gasket	531121
Nozzle - Lower Spray (Angle)	86247
Nozzle - Spray Arm (V Jet)	530657
Nozzle – Condenser (Atomizer)	78496
Spray Arm Post spinner	78355
Spray Arm (hub assy)	78387
Bearing	78356
Washer	78532
Nut	78359
End Cap	530658
Door Gasket	78392
Wiper	78637

Washers

CASTLE 7550

Controls

Door Solenoids 125VDC (4&5 Sol.)	536271
SSR's	523329
EPROM	531086
1 RTD (Wash)	531131
2 RTD (Dry)	531132
Door Switches	533174

Plumbing

Water Valve - (1,2,3 Sol)	533318
Coil	517301
Rebuild	72505
Steam Valve - Valve (6 Sol.)	517342
Coil	517299
Rebuild	54186
Cooling Water Valve (7 Sol.)	517343
Coil	517301
Rebuild	54187
Vapor Condenser Valve (8 Sol.)	541582
Coil	517302
Rebuild	54188
Sump Drain Valve - Valve (9 Sol.)	533319
Coil	517303
Rebuild	530860
Seal, Cover, Flap (In 5A Flange)	61301601037

Misc. Components

Detergent Wands (3)	533228
Single Wand	533158
Plastic Spinner Arm Nut	531126
Plastic Spinner Arm Bushing	531127
Squeeze Tubing Peristaltic Pump (21")	511898
Filter, Air, 7550	533241
Complete Door Assy.	531099
Door Spring	536506
Door Gasket	531038

Washers

2415 SONIC / RINSER DRYER

Controls

Switch, Momentary	2002030
Lense, Amber ("LID")	2024010
Lense, Amber ("FILL")	2024011
Lense, Amber ("DRAIN")	2024012
Switch, Alternating	2002019
Lense, Red ("POWER")	2024013
Lense, Amber ("DEGAS")	2024009
Fuse 1-3/4 A Fast Blow	2003053
Fuse 3A Fast Blow	2003038
Fuse 10A Fast Blow	2003045
Fuse 15A Fast Blow	2003047
Relay, 3PDT	2006005
Generator (500W)	2026000
Timer Kit (Silver Face to Black Face)	61301602255
Timer Only (Black Face)	61301602222
Contactor (2 pole)	2007003
Timer 180 sec. Delay	2005001
Heater Strip	2030016
Level Probe	49833
Isolation Transformer	541513

Lid Components

Lid Limit Switches	2002001
Gear Motor	5101005
Electric Brake	2025008
Torque Limiter Gear Motor	536029

Misc. Components

Det. Strainer	9000023
Det. Sol. Valve (1/4")	3001039
Water Sol. Valve (1/2")	3001036
Gasket Material	R3892
Injector	3023002
1/4" Needle Valve	3005021
2" Drain Valve	3002015
2" Hose (Specify Length)	R3893
5/8" Rubber Tubing (Specify Length)	3013027
1/4" Plastic Tubing (Specify Length)	R3839

Washers

Washers

2460-UC SONIC CLEANER

Controls

Digital Timer	SS20TIME0015
Controller, PLC-20 I/O	SS20PLC-0033
Light, Green	SS20LIGH0010
Light, Red	SS20LIGH0011
Circuit Breaker, 20AMP, 2 Pole	SS20CIRC0001
Sonic Generator 1200W, 40KHZ	SSKS4012
Heater 500WATT, 120V	SS20HEAT0003
Switch, Rocker "Drain"	SS10SWIT0021
Thermostat (Heater) TS1	SS20TEMP0001
Overheat Thermostat (Heater) TS2	SS20TEMP0009
Float Switch (Water Level)	SS20LEVE0005

Lid Operation Components

Lid Lift Cylinder	SS40PNEU0024
4 Way Valve (Water To Lid Cylinders)	SS30VALV0400
Valve, Flow Adj. (Lid Speed)	SS40PNEU0100
Tubing, 1/4"	SS40TUBE0019
Pressure Regulator (Water To Cylinders)	SS40REGU0008
Lid Gasket Replacement Kit	SS11116-53-15
Switch, Magnetic (Lid Closed & Footswitch)	SS20SWIT0025
Magnet	SS20MISC0022
Proximity Switch (Detergent Level)	SS20PROX0007
Switch, Front Panel Plunger	SS20SWIT0005
Kick Plate	SS13649
Latch, Front Panel Locking	SS40LATC0004
Detergent Bottle Cap	SS13689
Check Valve 1/2"	SS30VALV0022
Spring (Lid)	SS40MISC0100

Misc. Components

Valve, Hot Water 1/2"	SS30VALV0200
Lift Roller	SS12939
Check Valve 1/8"	SS30VALV0023
Detergent Pump	SS30PUMP0036

Washers

2460-RD RINSER DRYER

Controls

Digital Timer	SS20TIME0015
Controller, PLC-20 I/O	SS20PLC-0033
Light, Green	SS20LIGH0010
Light, Red	SS20LIGH0011
Circuit Breaker 30A	SS20CIRC0002
Thermocouple (Temp Probe)	SS20TEMP2008
Water Ejector (Detergent Injector)	SS30PUMP0105
Thermostat (Overtemp) TS1	SS20TEMP0023

Lid Operation Components

Lid Lift Cylinder	SS40PNEU0024
4 Way Valve (Water To Lid Cylinders)	SS30VALV0400
Valve, Flow Adj. (Lid Speed)	SS40PNEU0100
Tubing, 1/4"	SS40TUBE0019
Pressure Regulator (Water To Cylinders)	SS40REGU0008
Lid Gasket Replacement Kit	SS11116-53-15
Switch, Magnetic (Lid Closed & Footswitch)	SS20SWIT0025
Magnet	SS20MISC0022
Proximity Switch (Detergent Level)	SS20PROX0007
Switch, Front Panel Plunger	SS20SWIT0005
Kick Plate	SS13649
Latch, Front Panel Locking	SS40LATC0004
Detergent Bottle Cap	SS13689
Check Valve 1/2"	SS30VALV0022
Spring (Lid)	SS40MISC0100

Misc. Components

Valve, Hot Water 1/2"	SS30VALV0200
Lift Roller	SS12939
Check Valve 1/8"	SS30VALV0023

Washers

GE FANUC CONTROLS

Power Supply IC693PWR321	2031027
Battery	P0020693
CPU 313, 5 Slot Base	400406
CPU 331 OBSOLETE – Replaced by CPU 350	
CPU 350 IC693CPU350 *	2031029
CPU 341 IC693CPU341 *	401239
EPROM for CPU 313/331 *	400907
Base Unit, 10 Slot (includes new battery)	2031028
ASCII Module HE693ASC224x	400518
Thermocouple Module HE693THM446x	2031034
Thermocouple w/fitting	401126
Thermocouple Wire	2023030
Blank Module IC693ACC310	2031033
DC Input/AC Output Module IC693MDR390	2031031
AC Input/AC Output Module IC693MAR590	2031032
DC Input Module IC693MDL634	400756
AC Output Module IC693MDL940	400575
Touch Screen 5" B/W *	2031035
Touch Screen 6" Color *	2031036
Interface cable for Touch Screen	400190
Ribbon cable for Touch Screen	Y960479
Printer, Syntest	400868
Printer Paper	400924

**Must specify program number and washer serial number*

Washers

GE VERSAMAX CONTROLS

Power Supply (120VAC) IC200PWR102	61301605230
CPU *	P0021225
Thermocouple Module IC200ALG630	P0021189
Input Module, 24VDC, 32PT IC200MDL650	61301605120
IC200MDL330	61301605122
Carrier, I/O Box Style IC200CHS022	P0021226
Carrier, I/O Connector Style IC200CHS003	P0021227
Power Supply, 24VDC, 6.5A	61301605123

EAGLE SIGNAL 190+ CONTROLS

DAD/M	2031011
Controller, 190+	2031009
Memory Transport Module *	2031023
Battery Module	400875
Expander Board	2031010
AC Input Module	2031018
DC Input Module	2031020
AC Output Module	2031019
DC Output Module	2031021
Thermocouple Module	2031024
Thermocouple w/fitting	401126
Thermocouple Wire	2023030
Printer, Digitec	2034001
Printer Paper	2034005
Take Up Reel	OBSOLETE

**Must specify program number and washer serial number*

Washers

76-7700 / CF-4C/CS-6F CART WASHER

Controls See GE FANUC, VERSAMAX or
EAGLE SIGNAL 190+ Controls (Pages 45 & 46)

Switches

Power/Reset Switch (Allen Bradley only. Parts for GE switches no longer available)

Switch, Selector, 2 position, Green	2002047
Light Unit	2015060
Lamp, LED, 120V (Long Lasting)	2015040
Block, Contact, N.O.	2009010
Start, PB, Green, Non-Illum	2002046
Stop, PB, Red, Non-Illum	400174

Indicators

Indicator, Red	2015045
Indicator, Blue	2015050
Indicator, Yellow	2015055
Indicator, Green	2015070
Lamp, LED, 120V (Long Lasting)	2015040

Misc. Controls

Board, Liquid Level Control	2011000
Relay, Floor Tilt Control (1CR)(120V)	400181
Relay, Floor Tilt Control (2&3CR)(24VDC)	2006009
Capacitor, Floor Tilt Actuator	401071
Alarm	2025004
Rectifier	2025003
Transformer, 120V/24v, .5KVA	2010002
Transformer, 240-480/120-240, 1.5KVA	2010017
Transformer, 120/120, .5KVA	2010019

Washers

76-7700 / CF-4C/CS-6F CART WASHER Cont.

Door

Gasket, Top/Non-Hinge (Order 10')	3021015
Gasket, Hinge Side (Order 7')	400163
Latch, Brixon	7001011
Handle, Door Latch	7001010
Hinge	7002006
Cam, Hinge, Female	7002008
Cam, Hinge, Male	7002009

Switch, Proximity, 12mm (Grey Cable)	OBSOLETE
--------------------------------------	----------

Must Order 400980 & Yellow Cable below.

Switch, Proximity, (Yellow Cable)	400980
Cable, Yellow (7600)	Y960789
Cable, Yellow (7700)	Y960788
Switch, Mechanical w/Actuator	2002063
Threshold Seal	Y940230
Locking Strip, Gasket	400705
Tool, Locking Strip	9000025
Decal, "Push To Exit"	2019044
Window On Door, Lexan	Y950508
Gasket, Door Window	Y950509

Interior Light

Bulb, 25W Fluorescent	401267
Receptacle, Lamp, Porcelain	2018010
Window For Interior Light, Lexan	Y950312
Gasket, Interior Light Window	Y950311

Washers

76-7700 / CF-4C/CS-6F CART WASHER Cont.

Motor Overloads

Wash Pump 5HP, 208/240V	2008007
Wash Pump 7-1/2HP, 208/240V	2008008
Wash Pump 5HP, 480V	2008005
Wash Pump 7-1/2HP, 480V	2003057
Drying Blower, 208V	2008004
Drying Blower, 240V	2008003
Drying Blower, 480V	2008002

Contactors

Wash Pump 5HP	2007004
Wash Pump 7-1/2HP	2007012
Drying Blower	2007000
Manifold Drive	2007000

Fuses

Wash Pump Fuses (1F-3F)

20A 208V/240V	2003011
25A	2003012
15A	2003032
Wash Pump (1F-3F), 7-1/2HP, 480V – 20A	2003057
Dryer Blower (4F-6F), 208V – 7A	2003006
Dryer Blower (4F-6F), 240V – 15A	2003010
Dryer Blower (4F-6F), 480V – 6A	2003054
Manifold Motor (7F) – 7A	2003006
Exhaust Blower (8F), (Optional) – 15A	2003010
Fuse (14F) – 1A	2003036
Fuse (15F) – 15A	2003047
Fuse (16F) – 2A	2003037
Fuse (17F-23F) – 4A	400303
Fuse (24F), 7600 – 2.5A	400972
Fuse (25F-27F) – 1/2A	2003035
Fuse (28F) – 1A	2003036
Fuse (29F-30F) – 10A	2003059

Washers

76-7700 / CF-4C/CS-6F CART WASHER Cont.

Old Style Booster Heater

Booster Heater, Graham, Heliflow	3016000
Gasket Kit	9114001
Pressure Gauge	8003003
Temperature Gauge	2012006
Pressure Relief Valve	3005011
Steam Valve, 1-1/2"	3001052
Rebuild Kit	9101004
Trap, Steam, 1/2"	3006008
Trap, Steam, 3/4"	3006009
Check Valve, 1/2"	3003000
Check Valve, 3/4"	3003001

New Style Booster Heater

Booster Heater Kit	
(Replaces old ITT Flat Plate)	B010364
Heat Exchanger Only	B010180
Steam Valve, 1"	3001051
Rebuild Kit	9101009
Trap, Steam, 3/4"	3006009
Trap, Steam, 1/2", Brass Angle	STEAM0001
Check Valve, 3/4"	3003001
Vacuum Breaker, 3/4"	3005028
Water Valve, 1"	3001045
Rebuild Kit	9101048

Washers

76-7700 / CF-4C/CS-6F CART WASHER Cont.

Dryer / Exhaust

Motor, Dryer Blower	5101010
Blower, Dryer	4001000
Coil, Steam, Air Heating	3016006
Coil, Condenser (3 Coils Needed)	B020657
Motor, Damper	5101003
Hose, Blower Intake, 9"	3013033
Filter, Air Intake	9000034
Tubing, Condenser Drain, 3/4" Urebrade	3013025

Floor Lift Mechanism

Actuator, 3" Stroke (7600)	2038001
Actuator, 3" Stroke (7700)	2038002
Switch, Limit	511410
Cable, SS, 1/8" (By The Foot)	1011001
Pulley, cable	1011005
Clamp, cable	1011004
Thimble, 1/8" Cable	401343

Manifold / Drive

Motor, Manifold Drive	5101002
Hose, Manifold (By The Foot)	R3896
Clamp, Manifold Hose	3015007
Wear Sleeve, Push Pin, Manifold, Short	X740037
Wear Sleeve, Push Pin, Manifold, Long	X740038
Push Bar Assembly (7600)	Y780010-52
Push Bar Assembly (7700)	Y940182
Idler Assembly	Z830041
Bearing, Idler Sprocket	400117
Sprocket, Idler	6001008
Torque Limiter, 3/4"	6003001
Coupling, Torque Limiter	6003000
Bearing, Flange	6004009
Seal, Manifold Drive Shaft	9207000
Bearing, Drive Shaft, Nylon	6004017
Attachment, Chain w/Link	X960176
Chain, Manifold Drive, SS (By The Foot)	6002004

Washers

76-7700 / CF-4C/CS-6F CART WASHER Cont.

Plumbing

Valve, Steam, 1 1/2"	
(Sump, Booster Heater-Graham)	3001052
Kit, Rebuild	9101004
Coil	9101052
Valve, Steam, 1"	
(Drying Steam, Booster Heater-ITT)	3001051
Kit, Rebuild	9101009
Coil	9101051
Valve, Water, 1" (Sump Fill)	3001045
Kit, Rebuild	9101048
Coil	9101049
Valve, Water, 3/4"	
(Condenser Water, Non-Vented Option)	3001046
Kit, Rebuild	9101003
Coil	9101049
Valve, Water, 1/2"	
(Condenser Water, Drying System)	3001044
Kit, Rebuild	9101008
Coil	9101049
Trap, Steam, 3/4" (Drying Coil)	401112
Trap, Steam, 3/4"	
(Sump, Booster Heater, Steam Supply)	3006009
Vacuum Breaker, 1/2"	3005029
Vacuum Breaker, 3/4"	3005028
Valve, Check, 1/2"	3003000
Valve, Check, 3/4"	3003001
Tubing, Dryer Condenser Drain, 3/4"	3013025

Washers

76-7700 / CF-4C/CS-6F CART WASHER Cont.

Sump

Coil, Sump Heating, 76001, 2, 3	W950262
Coil, Sump Heating, 76004, 5	W950158
Spark Plug, Water Level Probe	2011002
Electrode, Water Level Probe	
1/2" Long	Y820134-1
3" Long	Y820134-15
Bulkhead Fitting, 90 Deg., Deterg/Rinse	400113
Dispenser, Detergent/Rinse,	
4100 Beta (7600/7700)	3020002
Multi-Time Timer	9112003
Squeeze Tube	9112006
Drain Valve	Y830020

Sump Capacities

7600-1 64" Chamber	45 gal
7600-2 72" Chamber	50 gal
7600-3 80" Chamber	55 gal
7600-4 98" Chamber	55 gal
7600-5 112" Chamber	65 gal
7700	60 gal

Washers

Washers

78-7900A CART WASHER

Controls See GE FANUC, VERSAMAX Controls
(Pages 45 & 46)

Switches

Power/Reset Switch (Allen Bradley only. Parts for GE switches no longer available)

Switch, Selector, 2 position, Green	401086
Block, Contact, N.O.	401092
Light Unit	401109
Lamp, LED, 120V (Long Lasting)	400768

Door

Gasket, Top/Non-Hinge (Order 10')	3021015
Gasket, Hinge Side (Order 7')	400163
Latch, Brixon	7001011
Handle, Door Latch	7001010
Hinge	7002006
Cam, Hinge, Female	7002008
Cam, Hinge, Male	7002009
Switch, Proximity	61301605134
Cable	61301605136
Decal, "Push To Exit"	2019044
Window On Door, Lexan	Y950508
Gasket, Door Window	Y950509

Interior Light

Bulb, 25W Fluorescent	401267
Receptacle, Lamp, Porcelain	2018010
Window For Interior Light, Lexan	Y950312
Gasket, Interior Light Window	Y950311

Washers

78-7900A CART WASHER Cont.

Motor Overloads 7800 208V*

Wash Pump (1M)	61301605168
Rinse Pump & Dryer Motor (15M-2M)	61301605164
Recirc. Pump (14M)	61301605167

Motor Overloads 7800 480V*

Wash/Recirc Pumps-Manifold (1-2-14M)	61301605164
Rinse Pump (15M)	61301605163
Dryer Motor (2M)	61301605162

Motor Overloads 7900 208V*

Wash & Recirc. Pump (1-14M)	61301605169
Rinse Pump & Dryer Motor (2-3-15M)	61301605164
Manifold (4M)	61301605165

Motor Overloads 7900 480V*

Wash Pump (1M)	61301605169
Wash Pump/Dryer (2-3-15M)	61301605164
Recirc. Pump/Manifold (4-14M)	61301605165

Motor Contactors 7800 208V*

Wash Pump (1M)	61301605159
Rinse Pump/ Dryer/Manifold (2-4-15M)	61301605154
Recirc. Pump (14M)	61301605158

Motor Contactors 7800 480V*

Wash Pump (1M)	61301605157
Recirc/Rinse Pump/Dryer/Mani (2-14-15M)	61301605154

Motor Contactors 7900 208V*

Wash Pump (1M)	61301605160
Rinse Pump & Dryer Motor (2-3-15M)	61301605154
Recirc Pump (14M)	61301605159
Manifold (4M)	61301605155

Motor Contactors 7900 480V*

Wash Pump (1M)	61301605160
Recirc/Rinse Pump/Dryer (2-3-14-15M)	61301605154
Manifold (4M)	61301605155

**Refer to manual for other supply voltages*

Washers

78-7900A CART WASHER Cont.

Dryer / Exhaust

Motor, Dryer Blower	5101010
Coil, Steam, Air Heating	D030073
Hose, 9"	3013033
Dryer Gasket	B040038

Floor Lift Mechanism

Actuator, 3" Stroke (7800)	2038001
Actuator, 3" Stroke (7900)	2038002
Proximity, Limit	61301605134
Cable, SS, 1/8" (By The Foot)	1011001
Pulley, Cable	1011005
Clamp, Cable	1011004
Thimble, 1/8" Cable	401343

Manifold / Drive

Motor, Manifold Drive	5101002
Hose, Manifold (By The Foot)	R3896
Clamp, Manifold Hose	3015007
Wear Sleeve, Push Pin, Manifold, Short	X740037
Wear Sleeve, Push Pin, Manifold, Long	X740038
Push Bar Assembly (7800 ONLY)	Y780010-52
Push Bar Assembly (7900 ONLY)	Y940182
Attachment, Chain w/Link	X960176
Chain, Manifold Drive, SS (By The Foot)	6002004

Washers

78-7900A CART WASHER Cont.

Plumbing

Valve, Steam, 1 1/2"	3001052
Kit, Rebuild	9101004
Coil	9101052
Valve, Steam, 1"	3001051
Kit, Rebuild	9101009
Coil	9101051
Valve, Water, 1"	3001045
Kit, Rebuild	9101048
Coil	9101049
Valve, Water, 1/2"	3001044
Kit, Rebuild	9101008
Coil	9101049
Trap, 3/4" TD (Watson)	P0011264
Trap, 1" TD (Watson)	P0011265
Trap, Steam, 1/2", Brass Angle	STEAM0001
Vacuum Breaker, 1/2"	3005029
Vacuum Breaker, 3/4"	3005028
Valve, Check, 1/2"	3003000
Valve, Check, 3/4"	3003001

Sump

Coil, Tank 7800	D030034
Spark Plug, Water Level	2011002
Probe 5-3/4" Long	D030374
Probe 21-1/4" Long	D030040
SS Hex Nut	P0030102
Dispenser, Detergent, Knight	P0018005
Pump, Chemical, Rinse Aid	5001024
Drain Valve	Y830020
Detergent Wand	V500480001

Sump Capacities

7800	80 gal tank
7900	100 gal tank

Washers

7642, 7652, 7662 INDEX WASHERS

Controls See GE FANUC, VERSAMAX or
EAGLE SIGNAL 190+ Controls (Pages 45 & 46)

Switches

Power/Reset Switch (Allen Bradley only. Parts for GE switches no longer available)

Switch, Selector, 2 position, Green	2002047
Light Unit	2015060
Lamp, LED, 120V (long lasting)	2015040
Block, Contact, N.O.	2009010
Start, PB, Green, Non-Illum	2002046
Stop, PB, Red, Non-Illum	400174
Switch, Load Trip (Electric Eye)	2002038
Switch, Load Trip (Micro)	2002016
Actuator	2002036
Switch, Door/Index Frame (Micro)	2002063
Arm, Switch	2002062
Switch, Index Cylinder (Proximity Sensor)	2002070
Switch, Elevator Cylinder w/Cable	400952
Switch, E-Stop, Keyed	401088
Extra Keys (2)	401508
Switch, E-Stop	401087
Switch, Door Safety, Sonic Compartment	2002065
Switch, Toggle, Bypass for Older Indexers	2002009

Indicators

Indicator, Red	2015045
Indicator, Blue	2015050
Indicator, Yellow	2015055
Indicator, Green	2015070
Lamp, LED, 120V (long lasting)	2015040

Contactors

Contactor (1M-5M, 7M) All Voltages	2007000
Contactor, Sonic Generators	2007009

Washers

7642, 7652, 7662 INDEX WASHERS Cont.

Fuses

Prewash, Wash, Rinse1, Rinse2 Pump	
Drying Blower (1F-15F) – 10A	2003009
Sonic Generators (16F) – 25A	2003012
Manifold Drive Motor (17F) – 7A	2003006
Transformer, 480V (18F-19F) – 15A	400985
110V to Control System (20F) – 8A	2003007
110V to GE FANUC (21F) – 2A	2003037
Fuse (22F) – 4A	400303
Fuse (23F-25F, 27F, 32F) – 1A	2003036
Fuse (26F, 28F-31F) – 2A	2003037

Motor Overloads

Overload Relay, 6.0-9.5A (1OL-4OL) 208V	2008005
Overload Relay, 4.0-6.0A (5OL) 208V	2008004
Overload Relay, 2.5-4.0A (1OL-4OL) 480V	2008003
Overload Relay, 1.6-2.5A (5OL) 480V	2008002

Misc. Controls

Counter	2016000
Transformer, 120VAC/24VAC, .5KVA	2010002
Transformer, 240V/480V, 1.5KVA	400970
Transformer, Sonic Generator Isolation	2010018
Transformer, PC Isolation	2010019
Gauge, Temperature, Balkamp 701-1103	8003009
Gauge, Temperature, Balkamp 701-1126	8003012

Board, Liquid Level

Single Function (Prewash, Sonic)	2011000
Dual Function (Wash, Rinse1)	2011001
Master Control Relay	400629
Socket	400188
Alarm	2025004

Washers

7642, 7652, 7662 INDEX WASHERS Cont.

Sump Capacities

Prewash	18 gal (CW)
Sonic	30 gal (HW)
Wash	18 gal (HW)
Rinse	28 gal (HW)
DI Rinse (15-30 Second Cycle)	32 GPM

Prewash Section

Spinner Manifold, Prewash	Y940347
Spark Plug, Water Level Probe	2011002
Probe Tab - 7/8 LG	Y820134-2
Probe Tab - 8 LG	Y820134-29
Drain Valve	3002015
Water Fill Solenoid	3001046
Rebuild Kit	9101003
Injector, Enzyme/Detergent	3023000
Pump, SMP1000	5001002
Pump Seal	9106002
Housing Gasket	9106001

Sonic Section

Generator, Sonic, 500 Watts	401122
Transducers	401127
Transducer Screen for Sump Bottom	W960284
Spark Plug, Water Level Probe	2011002
Probe Tab - 7/8 LG	Y820134-2
Probe Tab - 2-7/8 LG	Y820134-13
Drain Valve	3002015
Water Fill Solenoid	3001044
Rebuild Kit	9101008
Injector, Enzyme/Detergent	3023000

Washers

7642, 7652, 7662 INDEX WASHERS Cont.

Wash Section

Manifold, Complete, LH	W920323
Manifold, Complete, RH	W920491
O-Rings for Manifold	1006001
Hose, Manifold	3013031
Spark Plug, Water Level Probe	2011002
Probe Tab - 1-3/8 LG	Y820134-5
Probe Tab - 5-3/4 LG	Y820134-22
Water Fill Solenoid	3001044
Rebuild Kit	9101008
Steam Solenoid Valve	3001050
Rebuild Kit	9101036
Trap, Steam , 1/2"	3006008
Breaker, Vacuum, 3/4"	3005029
Deterg. Pump (Conc. Sensing H3100D)	3020000
90° Bulk Fitting, Det. Injection	400113
Probe, Dual Tip, Blue	3020003
Coil, Steam	Y940507
Drain Valve	3002015
Pump, SMP1000	5001002
Pump Seal	9106002
Housing Gasket	9106001

Washers

7642, 7652, 7662 INDEX WASHERS Cont.

Rinse/DI Rinse Section

Manifold, Complete, LH	W920324
Manifold, Complete, RH	W920486
O-Rings for Manifold	1006001
Hose, Manifold	3013031
Spark Plug, Water Level Probe	2011002
Probe Tab - 5-3/4 LG	Y820134-22
Coil, Steam	Y940507
Drain Valve	3002015
Water Fill Solenoid S/S	3001049
Steam Solenoid Valve	3001050
Rebuild Kit	9101036
Trap, Steam , 1/2"	3006008
Breaker, Vacuum, 1/2"	3005029
Pump, SMP1000 (Rinse 1)	5001002
Pump Seal	9106002
Housing Gasket	9106001
Pump, Final Rinse	5001006
Pump Seal	9110000
Switch, SS Float, DI Tank	2011005
Dispenser, Lube	5001007
Rebuild Kit	400914
Injection Fitting	9108010
Foot Valve	9108009

Dryer Section

Blow Off Assembly	Y860043
Hose, Blow Off Assembly	Y860043
Blower, Dryer	4001000
Motor, Blower	5101010
Motor, Damper	5101003
Steam Solenoid Valve	3001051
Rebuild Kit	9101009
Trap, Steam, 3/4"	401112
Breaker, Vacuum, 3/4"	3005028

Washers

7642, 7652, 7662 INDEX WASHERS Cont.

Pneumatic Cylinders

Index Frame, 2" Bore (Uses Switches)	8005005
Seal Kit	9103010
Index Frame, 2-1/2" Bore (Uses Switches)	8005002
Seal Kit	9103004
Index Frame, 2-1/2" Bore (Uses Sensors)	8005007
Seal Kit	401316
Valve, Flow Control, 3/8"	401081
Door Lift	8005001
Seal Kit	401315
Valve, Flow Control, 3/8"	401081
Elevator	8005003
Seal Kit	401317
Valve, Flow Control, 3/8"	401081
Filter/Regulator Assembly	401105
Regulator, Doors Down, Elevator Down	8002002

Pneumatic Valve Stacks

Norgren, Parker, 2 Valve (Old Style)	8006004
Norgren, Parker, 3 Valve (Old Style)	8006003
ASCO, 2 Valve (New Style)	Y970140
ASCO, 3 Valve 4 Coils (New Style)	Y970139
ASCO, 3 Valve 5 Coils (New Style)	B010008
Air Line, Black Vinyl, 3/8" (order by the foot)	3013022
Filter, Oil Removal	8002005
Lubricating Oil	400269

Washers

7642, 7652, 7662 INDEX WASHERS Cont.

Manifold Drive

Motor, Manifold Drive	5101004
Torque Limiter Coupling	6003000
Torque Limiter, 3/4"	6003001
Shaft, Manifold Drive	Z920258
Bearing, Manifold Driveshaft	6004000
Roller Track Weldment	Z860056

Doors

Latch, Access Doors	7001009
Door Link Assembly, Internal Doors	Y950643
Rod, All Thread, SS, 3/8X8", for Door Link	Y920377-1
Lift Shaft	6004019
Gasket, Self Gripping, Access Doors	400560
Index Frame Seals, Internal Doors	Y910441

Washers

Washers

7622 – CR22 DECONTAMINATOR

Controls (See Pages 45 & 46)

The Spirit (Grey-Pro) Controller is obsolete.

Use MP-145A GE Fanuc (Y950747) to upgrade.

For GE FANUC parts, see GE FANUC Controls.

Printer, Keltron	2034000
Paper	CHART0015
Ribbon	61301600006
RTD (Use with Spirit Controller)	2012018
Probe, Temp (Use with GE Controller)	401126
Timer, BR48A601, 0-10 min	2001008
Relay, Time Delay, 5X829C	2005001
Relay, Time Delay, Cycle Complete	2005000
Relay, Latch, 32E6CAE120	2006001
Alarm, Cycle Complete	2025004

Switches

Power/Reset Switch (Allen Bradley only. Parts for GE switches no longer available)

Switch, Selector, 2 position, Green	2002047
Light Unit	2015060
Lamp, LED, 120V (long lasting)	2015040
Block, Contact, N.O.	2009010

Wash / DI Rinse Pump Contactor	2007000
---------------------------------------	---------

Overload Relays, Wash Pump

208V	2008006
230V	2008006
460V	2008004

Overload Relays, DI Rinse Pump

208V	2008004
230V	2008004
460V	2008003

Washers

7622 – CR22 DECONTAMINATOR Cont.

Booster Heater

Gauge, Temperature	2012006
Gauge, Pressure	8003003
Valve, Pressure Relief	3005011
Control, Temperature, Sterlco	3005018

DI Holding Tank

Pump, SS Teel, 1.5HP	5001006
Valve, Float	B040061
Float, Ball – 4"	9000002
Thermostat, 180F	2012014

Door

Gasket, Self Gripping (order 7ft. per door)	3021015
Latch, Door – Kason	7001016
Hinge, Door	7002003
Handle, Access Door	7001004
Spring, Door	7002015
Tube, Door Support, Lower	Z880411
Tube, Door Support, Upper	Z880413
Switch, Door	511410
Door Interlock Assembly	W940332
Solenoid, Pull Type	2037002
Spring, Solenoid Return	400184
Window	Y950508
Gasket	Y950509
Locking Strip	400705

Washers

7622 – CR22 DECONTAMINATOR Cont.

Dryer / Exhaust

Blower/Dryer	4001001
Coil, Condenser	B020657
Contact, Drying Heater	
40A	2007007
25A	2007009
Element, Heating, 4KW, 208V, (Non-Vent)	2030034
Element, Heating, 4KW, 240/480V	2030038
Thermostat, Overheat, 300F	2012017
Motor, Damper	5101003

Manifold Drive

Motor, Manifold Drive	5101005
Contact, Manifold Drive Motor	2007003
Drive Shaft Weldment	Y900028
Torque Limiter, 5/8"	6003002
Coupling, Torque Limiter	6003000
Bearing, Drive Shaft, Old Style (Stein)	6004010
Bearing, Insert, Old Style (Stein)	6004015
Bearing, Drive Shaft	6004018
Hose, Manifold, Wash (1"ID X 42"Long)	3013031
Hose, Manifold, Rinse (3/4"ID X 42"Long)	3013032
Hose, Rinse Couple (3/4"ID X 4"Long)	Y960611
Hose, Wash Couple (1"ID X 4"Long)	Y960612
O-Ring, Manifold Tube	1006001

Washers

7622 – CR22 DECONTAMINATOR Cont.

Plumbing

Seal, Wash Pump	9104004
Rubber Ring, Wash Pump	9104003
Pump, DI, 1-1/2HP	5001006
Valve, Solenoid, 3/4" (Pre-Wash, Rinse)	3001046
Rebuild Kit	9101003
Valve, Solenoid, 1/2" (Sump Fill/Condense)	3001044
Rebuild Kit	9101008
Valve, Solenoid, 1/2" (Steam)	3001050
Rebuild Kit	9101036
Valve, Solenoid, 1" (DI Rinse)	3001043
Trap, Steam, 1/2" Inline (SS)	3006008
Trap, Steam, 1/2" Angle (Brass)	3006000
Dispenser, Detergent, Concen. Sensing	3002000
Probe	3020003
Squeeze Tube	9112006
Valve, Drain	3002015

Sump

Diverter Plate	Z960236
Electrode, Water Level (Spark Plug)	2011002
Probe, Water Level, High	Y820134
Probe, Water Level, Low	3020003
Switch, Float (GE Controls)	2011005
Bulkhead Fitting, 90 Deg., 1/4" Tubing	400288
Bulkhead Fitting, Straight, 1/4" Tubing	3020005
Coil, Steam	Y940507

Sump Capacity	18 Gallons (retained)
---------------	--------------------------

Washers

7626 – CR26 DECONTAMINATOR

Controls (See Pages 45 & 46)

The Spirit (Grey-Pro) Controller is obsolete.

Use MP-145A GE Fanuc (Y950747) to upgrade.

For GE FANUC parts, see GE FANUC Controls.

Printer, Keltron	2034000
Paper	CHART0015
Ribbon	61301600006
RTD (Use with Spirit Controller)	2012018
Probe, Temp (Use with GE Controller)	401126
Timer, BR48A601, 0-10 min	2001008
Relay, Time Delay, 5X829C	2005001
Relay, Time Delay, Cycle Complete	2005000
Relay, Latch, 32E6CAE120	2006001
Alarm, Cycle Complete	2025004

Switches

Power/Reset Switch (Allen Bradley only. Parts for GE switches no longer available)

Switch, Selector, 2 Position, Green	2002047
Light Unit	2015060
Lamp, LED, 120V (long lasting)	2015040
Block, Contact, N.O.	2009010

Wash / DI Rinse Pump Contactor	2007000
---------------------------------------	---------

Overload Relays, Wash Pump

208V	2008006
230V	2008006
460V	2008004

Overload Relays, DI Rinse Pump

208V	2008004
230V	2008004
460V	2008003

Washers

7626 – CR26 DECONTAMINATOR Cont.

Booster Heater

Gauge, Temperature	2012006
Gauge, Pressure	8003003
Valve, Pressure Relief	3005011
Control, Temperature, Sterlco	3005018

DI Holding Tank

Pump, SS Teel, 1.5HP	5001006
Valve, Float	B040061
Float, Ball – 4"	9000002
Thermostat, 180F	2012014

Door

Gasket, Self Gripping (order 7ft. per door)	3021015
Latch, Door – Kason	7001016
Hinge, Door	7002003
Spring, Door	7002015
Tube, Door Support, Lower	Z880411
Tube, Door Support, Upper	Z880412
Switch, Door	511410
Door Interlock Assembly	W940332
Solenoid, Pull Type	2037002
Spring, Solenoid Return	400184
Window	Y950508
Gasket	Y950509
Locking Strip	400705
Pulley, Door (Vertical Door Only)	6005002

Washers

7626 – CR26 DECONTAMINATOR Cont.

Dryer / Exhaust

Blower/Dryer	4001001
Blower/Dryer (Non-Vented)	400644
Blower/Exhaust (Non-Vented)	400642
Contactoer, Drying Heater	
40A	2007007
25A	2007009
Element, Heating, 4KW, 208V, Non-Vented	
Only	2030034
Element, Heating, 4KW, 240/480V	2030038
Thermostat, Overheat, 300F	2012017
Motor, Damper	5101003

Manifold Drive

Motor, Manifold Drive	5101001
Contactoer, Manifold Drive Motor	2007003
Drive Shaft Weldment	Y820056
Coupling, Drive Shaft	9204035
Coupling w/Bearings	Y810215
Push Bar	Y810133
Wear Sleeve	X820129
Spacer, Driveshaft	9202085
Chain Attachment	X960176
Torque Limiter, 3/4"	6003001
Sprocket, Torque Limiter	6001001
Sprocket, Driveshaft	6001007
Sprocket, Motor	6001005
Bearing, Drive Shaft	6004008
Hose, Manifold, Wash (1"ID X 42"Long)	3013031
Hose, Manifold, Rinse (3/4"ID X 42"Long)	3013032
Hose, Rinse Couple (3/4"ID X 4"Long)	Y960611
Hose, Wash Couple (1"ID X 4"Long)	Y960612
O-Ring, Manifold Tube	1006001
Assembly, Idler	Z830041
Bearing, Idler Sprocket	400117
Sprocket, Idler	6001008

Washers

7626 – CR26 DECONTAMINATOR Cont.

Plumbing

Seal, Wash Pump	9104004
Rubber Ring, Wash Pump	9104003
Pump, DI, 1-1/2HP	5001006
Valve, Solenoid, 3/4" (Pre-Wash, Rinse)	3001046
Rebuild Kit	9101003
Valve, Solenoid, 1/2" (Sump Fill/Condense)	3001044
Rebuild Kit	9101008
Valve, Solenoid, 1/2" (Steam)	3001050
Rebuild Kit	9101036
Valve, Solenoid, 1" (DI Rinse)	3001043
Trap, Steam, 1/2" Inline (SS)	3006008
Trap, Steam, 1/2" Angle (Brass)	3006000
Deterg. Pump (Conc. Sensing H3100D)	3020000
Probe, Dual Tip, Blue	3020003
Squeeze Tube	9112006
Valve, Drain	3002015

Sump

Diverter Plate	Z960236
Electrode, Water Level (Spark Plug)	2011002
Probe, Water Level, High	Y820134
Probe, Water Level, Low	3020003
Switch, Float (GE Controls)	2011005
Bulkhead Fitting, Straight, 1/4" Tubing	3020005
Coil, Steam	Y940507

Sump Capacity	18 Gallons (retained)
---------------	--------------------------

Washers

4656 DISINFECTOR (DECOMAT)

Controls

Circuit Board *(Use the part number printed on the circuit board to ensure a correct replacement.)*

Battery, CPU (Must Be Soldered)	V499765300
Fuse, 200mA	V499256000
Fuse, 315mA	V499223600
Fuse, 10A	V468840600
Transformer, Control Voltage	V499502700
Sensor, Temperature	V499100000
Standoffs, Ckt Board	V499082500
Plug, 9 Pin	V762320101
Pin, Male	V468555300
Pin, Female	V468555400
Black Communications Cable, 100mm	V498964100
Black Communications Cable, 500mm	V498964101
Black Communication Cable, 1500mm	V498964102
Printer, Alfa	V501865475
Cartridge, Ribbon	61301600006
Paper Roll	V479398807

Detergent System

Detergent Pump Rebuild	V499693300
Detergent Pump	V499169205
Detergent Wand, (5' Cable & Hose)	V501596800
Detergent Wand, (2.5' Cable & Hose)	V501596802
Hose, Detergent	V930830
Bottle, Detergent	V498619000

Washers

4656 DISINFECTOR (DECOMAT) Cont.

Door System

Gasket	V499224100
Gas Spring	V499179001
Hinge LH	V499176600
Hinge RH	V499176601
Repair Kit, Door Support Bracket Mount	V500654875
Locking Solenoid (Door w/o Glass)	V499168800
Locking Solenoid (Door w/Glass)	V500096500
Switch, Door Locked	V499196800
Window	V499938000
Gasket (Requires 2m per door)	V477081400

Plumbing

Valve, Single, Cold/Hot/DI	V499165402
Valve, Double, Cold	V499720300
Valve, Triple, Cold	V499722400
Valve for DW	V498988200
Valve, Single, DI (Use With Booster Pump)	V499165401
DI Boost Pump	V468918101
Rebuild Kit	V468943900
Valve, Solenoid, Steam 3/8"	V498811604
Filter, Steam Inlet	V500230000
Trap, Steam	V500367900
Pump, Circulation	V499212901
Pump, Drain (Catalog # DP025-215)	V500362601
Pump, Drain (Catalog # BE28B4-066)	V477134501
Rebuild Kit	V468943900
Valve, Check, Drain	V499105900
Plastic Hose Fitting	V498987400

Washers

4656 DISINFECTOR (DECOMAT)

Dryer

Element, Heating	V498906200
Blower	V501866600
Pre-Filter (Large)	V500576700
HEPA Filter	V498906100
Protector, Overheat	V499085700
Thermostat, Control	V481966900
Bracket, Large	V499266600
Bracket, Small	V499266500
Valve, Check	V499228000

Hoses

Hose, Inlet Water Valve to Chamber	V498625205
Hose, Drain Tempering	V498625205

Chamber

Docking Flapper Kit	V500187600
---------------------	------------

Sump Heat

Coil, Steam	V499009900
Coil, Electric	V499136500
Element Flange	V499136500
Protector, Overheat	V499085700
O-Ring, Steam or Electric Coil	V740483105

Washers

4656 DISINFECTOR (DECOMAT) Cont.

Upgrade Kits

DI Water Boost Pump	V499252576
Third Detergent Pump Kit	V499252675
Fourth Detergent Pump Kit	V499329075
Add RS-232 Port for PC Printer	V499302975
Soiled Side Printer Kit	V499326676
Steam to Electric Sump Heating	V499200277
Electric to Steam Sump Heating	V499200476
Guide Rail Kit Single Door Machine	V499295075
Guide Rail Kit Double Door Machine	V499674275

Misc. Components

Barrier Trim, 4656	Y980114
Nut, Spray Arms	V500249700
Retainer, Panel Stud	V477273900
4656 Seismic Anchor	Y960204
Adapter, 1/2" BSP Male X 1/2" NPT Male	401151
Adapter, 1/2" BSP Fem X 1/2"NPT Male	401152
Adapter, 3/4" BSP Fem X 3/4" NPT Male	401153
Adapter, 3/4 BSP Fem X 3/4 MNPT, SS	401154
Adapter, 3/4" BSP Fem X 1/2" NPT Male	V477048201
Hose, Water, SS Braid, 1/2" X 3/4"	V499765700
Washcart Wheel (4656)	V498703500
Cover, Washcart Manifold, 4656	V499218000

Washers

46-X DISINFECTOR (PACS-300)

Controls

Safety Overheat Thermostat (Sump)	V501615300
CPU/Relay Board	V500405200
Control Panel Board (Control End)	V500619100
Cable	V500890900
Control Panel Board (Remote End)	V500619200
Cable	V500890901
RTD Temp. Sensor (Dual)	V500813601

Circuit Breakers

Circulation Pump (F01)	V501442900
Sump Heat, Electric (F02)	V501443100
Dryer Heaters (F04)	V501443300

Contactors

Circulation Pump (K01)	V501443700
Sump Heat, Electric (K02)	V501443702
Dryer Heater (K04)	V501443701

Main Power Transformer

V501185200

Door Components

Door Handle	V501466500
Glass Window	V501723400
Gasket For Window (Sold By Meter)	V477081400
Door Lock Motor	V500671100
Switch, Magnetic	V499694700
Spring Cylinder	V501252100

Door Gaskets

46-2 & 46-4	V499224100
46-5	V501649900

Washers

46-X DISINFECTOR (PACS-300) Cont.

Dryer

Check Valve	V501225500
Air Hose 370mm	V500081501
Temperature Sensor	V500852000
Fan	V501165000
Fuse, 6.3A (Inline)	V502040200
Fuse Holder	V502040100
Dryer Element	V500731300
Filter 46-2	V501380400
Filter 46-4 & 46-5	V501142400
Level Switch	V499178200
Pressure Diff. Sensor	V500855100

Drain Pump	V500362601
-------------------	------------

Detergent System

Detergent Pump Rebuild	V499693300
Detergent Pump	V499169205
Detergent Wand, Standard	V500306404
Detergent Wand, (5' Cable & Hose)	V501596800
Detergent Wand, (2.5' Cable & Hose)	V501596802
Hose, Detergent	V930830
Bottle, Detergent	V498619000

Misc. Components

Nut, Spray Arms	V500249700
Washcart Wheel (4656)	V498703500

Washers

8666 DISINFECTOR (DECOMAT)

Controls

Circuit Board *(Use the part number printed on the circuit board to ensure a correct replacement.)*

Battery, CPU (Must Be Soldered)	V499765300
Fuse, 200mA	V499256000
Fuse, 315mA	V499223600
Fuse, 10A	V468840600
Transformer, Control Voltage	V499502700
Sensor, Temperature	V499100000
Standoffs, Ckt Board	V499082500
Plug, 9 Pin	V762320101
Pin, Male	V468555300
Pin, Female	V468555400
Black Communications cable, 100mm	V498964100
Black Communications cable, 500mm	V498964101
Black Communication cable, 1500mm	V498964102
Printer, Alfa	V501865475
Cartridge, Ribbon	61301600006
Paper Roll	CHART0015

Detergent System

Detergent Pump Rebuild	V499693300
Detergent Pump	V499169205
Detergent Wand, (5' Cable & Hose)	V501596800
Detergent Wand, (2.5' Cable & Hose)	V501596802
Hose, Detergent	V930830
Detergent Wand, 2 ft Long w/5 ft Cord (Deco Cont.)	V500480001
Detergent Wand, 2 ft Long w/5 ft Cord (AB Cont.)	V500480000
Bottle, Detergent	V498619000

Washers

8666 DISINFECTOR (DECOMAT) Cont.

Door System

Door Gasket	V477331300
Motor, Door Sealing, Linear (Old Style)	V498983700
Motor, Seal Gear Motor (New Style)	V500037075
Capacitor for Gear motor	V500315400
Motor, Door Up/Down	V499081300
Capacitor for Up/Down Motor	V499689900
Chain, Door Drive	V499030900
Master Link, Chain	V499031000
Cable, Door Lift, SS	V477336601
Pulley, Door Lift Cable	V477336400
Leaf Spring, Door Open (Manual Only)	V499061800
Door Guide (Round Nylon "Roller")	V498956800
Window, Small (Old Style Door)	V477331400
Packing for Small Window (Req. 2m per door)	V477081400
Window, Large Smoked Glass (New Style Door)	V499998000
Packing for Large Window	V499990200

Door Switches

Switch, Safety, Door Up/Down (Power Door Only)	V500278500
Switch, Door Position, Up, Down, Sealed, Unsealed	V471597601

Washers

8666 DISINFECTOR (DECOMAT) Cont.

Drying System

Element, Heater	V476907801
Blower	V501866600
HEPA Filter	V498906100
Pre-Filter (Small - Before Fan)*	V498929602

* Pre-Filter Assembly – (New Style (Fits between blowers and HEPA filter. Remove small pre-filters from blower housing when using this filter assembly.)

V500551975

Pre-Filter (Large)	V500560600
Thermostat, Overheat Protection	V499085700
Thermostat, Control	V481966900
Timer, Fan Delay	V498766400
Switch, Water in Drying System (Mach. after SN 81511)	V500047800
Switch, Water in Drying System (Mach. before SN 81511)	V499013300
Valve, Non-Return	V501140300

Hoses

Hose, Dryer L=370	V500081501
Hose, Dryer L=560	V500081502
Hose, Door Drain, SS	V500058801
Hose, Door Drain, CS	V500058802
Hose, Inlet Water Valve to Chamber	V498625205
Hose, DI Water Valve to Chamber	V498977100
Hose, Drain, Flexible Grey Tubing	V150196700
Hose, Condenser Spray	V498625203
Hose, Drain Tempering Water	V499207400

Washers

8666 DISINFECTOR (DECOMAT) Cont.

Plumbing

Valve, Water, Hot/Single DI (W/O Economy Tank)	V499541602
Valve, Water, Cold, Double Solenoid	V498941001
Valve, Water, Cold, Triple Solenoid	V499743400
Valve, Steam to Sump	V499535200
Rebuild Kit	ZVALS0101
Filter, Steam Inlet	V500230000
Trap, Steam	V500367900
Pump, Circulation (Mach. Before SN81511)	V477330701
Pump, Circulation (Mach. After SN 81511)	V500956201
Pump, Drain/DI Boost/Economy Tank	V468918101
Rebuild Kit	V468943900
Valve, Check, Drain	V499105900
Plastic Hose Fitting	V498987400
Valve, Drain, DI Tank	V500141300

Sump Heat

Coil, Steam (Mach. Before SN 81511)	V499009900
Coil, Steam (Mach. After SN 81511)	V500283700
Coil, Electric	V499136500
Element Flange	V500154300
Protector, Overheat (For Electric Only)	V481950101
O-Ring, Electric or Steam Coil	V501823300

Washers

8666 DISINFECTOR (DECOMAT) Cont.

Upgrade Kits

DI Water Boost Pump	V499001076
Third Detergent Pump Kit	V499120875
Fourth Detergent Pump Kit	V499426175
Add RS-232 Port for PC Printer	V499287975

The following are for machine SN below 81511

Clean Side Printer Kit	V499378775
Soiled Side Printer Kit	V499206775
Steam to Electric Sump Heating	V499333200
Electric to Steam Sump Heating	V499108875

The following are for machine SN 81511 and above

Clean Side Printer Kit (No Unloader)	V500068076
Soiled Side Printer Kit	V500068075
Steam to Electric Sump Heating	V500237100

Misc. Components

Barrier Trim, 8666	Y980115
Nut, Spray Arms	V500249700
Spray Arm, Wash Cart	V499235101
Diffuser, Vent, 8666 (Non-Vented Units)	Y970393
Stud, Panel (1/4Turn) For 8666	V477273800
Retainer, Panel Stud	V477273900
Seismic Anchor	Y960205
Adapter, 1/2" BSP Male X 1/2" NPT Male	401151
Adapter, 1/2" BSP Fem X 1/2"NPT Male	401152
Adapter, 3/4" BSP Fem X 3/4" NPT Male	401153
Adapter, 3/4 BSP Fem X 3/4 MNPT, SS	401154
Adapter, 3/4" BSP Fem X 1/2" NPT Male	V477048201
Hose, Water, SS Braid, 1/2" X 3/4"	V499765700
Large Washer	V499215200
Small Washer	V499189300
Washcart Wheel	V498646800
Washcart Wheel (Tapered Stainless/AGS)	V500035200
Bearing for Tapered Wheel	V500165200
Cover, Washcart Manifold	V498677400

Washers

Washers

8666 DISINFECTOR (PACS-300)

Controls

Safety Overheat Thermostat (Sump)	V499085700
CPU/Relay Board	V500405200
Control Panel Board (Control End)	V500618900
Cable	V500890900
Control Panel Board (Remote End)	V500619000
Cable	V500890901
RTD Temp. Sensor (Dual)	V500813601

Circuit Breakers

F01-F05

PKZM 0-4	V500976100
PKZM 0-6.3	V500845100
PKZM 0-10	V500782500
PKZM 0-16	V499573500
PKZM 0-25	V500317000
F10 (Dryer Fans)	V500266102
F11 (Controls)	V500810300

Contactors K01-K05

V500584900

Main Power Transformer

V499858400

Dryer

Check Valve	V500937700
Temperature Sensor	V500852000
Fan	V501061700
Dryer Element	V500731300
Hepa Filter	V500745800
Pressure Sensor	V500855100

Door Components

Door Gasket	V477331300
Glass Window	V499998000
Gasket For Window (Sold By Meter)	V499990200
Door Lock Motor	V500037075
Switch ("Up", "Down", "Sealed")	V500859900

Washers

8666 DISINFECTOR (PACS-300) Cont.

Drain Pump	V468918101
Rebuild	V468943900

Circulation Pump (208V-240V)	V500956201
Mechanical Seal	V498644701
O-Ring	V499040400

Detergent System

Detergent Pump Rebuild	V499693300
Detergent Pump	V499169205
Detergent Wand, (5' Cable & Hose)	V501596800
Detergent Wand, (2.5' Cable & Hose)	V501596802
Hose, Detergent	V930830
Bottle, Detergent	V498619000

Booster Tank (Rinse Water Pre-Heat)

Heater Element (208V)	V500821700
O-Ring	V501823300
Level Sensor	V500984900
Fill Solenoid Valve	V499541600
Temperature Sensor	V500852000
Overheat Protector	V499085700

Sump

Heater Element (Electric 208V)	V500821700
O-Ring	V501823300
Steam Coil	V500283700
O-Ring	V501823300
Steam Trap	V500367900
Solenoid Valve	V499535200
Strainer	V500230000
Overheat Protector	V499085700

Misc. Components

Nut, Spray Arms	V500249700
Washcart Wheel	V498646800

RO Water Systems

Universal Parts

Low Pressure Switch	61301605490
Switch Contact Block	61301605504
Control Box Transformer	61301605506
Control Box Circuit Board	61301605508
Gauge 0-300 Gly. 1/4" Pnl. Mt.	61301605512
Wrench for 10" BB Housing	61301605513
Wrench for Standard Housing	61301605454

HS-200 BOILER MATE

Kit Qty

Membrane Filter	3	61301605455*
Filter Sediment 5 micron	4	61301605449*
Filter Carbon Block 5 micron	4	61301605447*
Filter GAC/KDF	4	61301605450*
Filter Calcite	4	61301605516*
*Part of PM Kit		61301605518
Permeate Pump		61301605453
Flow Restrictor 750 ml		61301605473
5 Gal Storage tank		61301605471
Valve Check 1/2" Brass		61301605475
Valve Ball 1/2" Brass		61301605474

RO Water Systems

HS-450-PT

Kit Qty

Membrane Filter 2 X 26	1	61301605535*
Control Box Motor Relay 30 Amp	1	61301605507*
Filter Calcite 10"	4	61301605448*
O-Ring for Standard Housing	1	61301605452*
Filter Sediment 5 micron 10" BB	4	61301605460*
Filter Carbon Block 10 micron 10" BB	4	61301605456*
Filter GAC/KDF 10" BB	4	61301605458*
O-Ring for 10" BB Housing	3	61301605468*
*Part of PM Kit		61301605525
Processor Pump 125 GPH		61301605527
Processor Motor 1/2 HP		61301605496
Coupler Pump/Motor Lovejoy		61301605495
Solenoid 1/2" FPT (24VAC)		61301605488
Solenoid Plastic Body 1/2"		61301605487
24v Coil for Plastic Solenoid		61301605514
Valve Check 1/4" FPT Plastic		61301605481
Valve Check 1/4" FPT Brass		61301605475
Pressure Tank Switch 40-60psi		61301605491
Valve Needle 3/8" FPT		61301605485
Valve Check 1/2" FPT Brass		61301605476
Vessel Closure Kit 2"		61301605534
Fuse 20 Amp		61301605502
Gauge 0-100 psi Gly. 1/4" Pnl. Mt.		61301605511
Housing Filter 10" Standard		61301605451
Housing Filter 10" BB		61301605470
Blending Valve 1/2"		61301605484

RO Water Systems

HS-1000

Kit Qty

Membrane Filter 2.5 X 40	1	61301605463*
Control Box Motor Relay 30 Amp	1	61301605507*
Filter Calcite 10"	4	61301605448*
O-Ring for Standard Housing	1	61301605452*
Filter Sediment 5 micron 10" BB	4	61301605460*
Filter Carbon Block 10 micron 10" BB	4	61301605456*
Filter GAC/KDF 10" BB	4	61301605458*
O-Ring BB Housing	3	61301605468*
*Part of PM Kit		61301604729

Processor Pump 190 GPH	61301605493
Processor Motor 1/2 HP	61301605496
Solenoid 1/2" FPT (24VAC)	61301605488
Solenoid Plastic Body 1/2"	61301605487
24v Coil for Plastic Solenoid	61301605514
Valve Check 1/4" FPT Plastic	61301605481
Valve Check 1/4" FPT Brass	61301605475
Valve Needle 3/8" FPT	61301605485
Vessel Closure Kit 2.5"	61301605466
Fuse 20 Amp	61301605502
Housing Filter 10" Standard	61301605451
Housing Filter 10" BB	61301605470
Blending Valve 1/2"	61301605484

HS-1000-PT Only

Coupler Pump/Motor lovejoy	61301605495
Pressure Tank Switch 40-60psi	61301605491
Valve Check 1/2" FPT Brass	61301605476
Gauge 0-100 psi Gly. 1/4" Pnl. Mt.	61301605511

HS-1000-FT Only

Valve Check 3/4" FPT Brass	61301605477
Valve Check 1 1/4" PVC	61301605479
Horizontal Reed Float	61301605509
Tank for Repressurizer (RP-5)	61301605471
Injector for Repressurizer	61301605498
Repressurizer Pump/Motor	61301605492
Gauge 0-100psi 1/4" Btm. Mt.	61301605510
Valve Check PVC 40 1/2"	61301605480

RO Water Systems

HS-2400

Kit Qty

Membrane Filter 4X40	1	61301605464*
Control Box Motor Relay 30 Amp	1	61301605507*
Filter Calcite 10"	8	61301605448*
O-Ring for Standard Housing	2	61301605452*
Filter Sediment 5 micron 20" BB	4	61301605461*
Filter Carbon Block 10 micron 20" BB	4	61301605457*
Filter GAC/KDF 20" BB	4	61301605459*
O-Ring for 20" BB Housing	3	61301605469*
*Part of PM Kit		61301604730

Processor Pump 330 GPH	61301605494
Processor Motor 3/4 HP	61301605497
Solenoid 1/2" FPT (24VAC)	61301605488
Valve Check 1/4" FPT Plastic	61301605481
Valve Check 1/4" FPT Brass	61301605475
Valve Needle 3/8" FPT	61301605485
Valve Check 1/2" FPT Brass	61301605476
Valve Ball 1/2" (PT Only)	61301605517
Vessel Closure Kit 4"	61301605467
Fuse 20 Amp	61301605502
24v Coil for Plastic Solenoid	61301605514
Plastic Solenoid Body 1/2"	61301605487
Housing Filter 10" Standard	61301605451
Housing Filter 20" BB	61301605465
Blending Valve 1/2"	61301605484

HS-2400-PT Only

Tank Switch 40-60psi	61301605491
Gauge 0-100 psi Gly. 1/4" Pnl.	61301605511

HS-2400-FT Only

Horizontal Reed Float	61301605509
Valve Check PVC 80 11/4"	61301605479
Valve Check PVC 40 1/2"	61301605480
Valve Check PVC 80 1"	61301605478
Tank for Repressurizer (RP-12)	61301605472
Injector for Repressurizer	61301605498
Pump/Motor for Repressurizer	61301605498

RO Water Systems

HS-4800

Kit Qty

20" BB 5M Sediment filter	4	61301605461*
10" all Calcite feeder	8	61301605448*
BB "O" Ring (for housing)	1	61301605468*
Slim "O" Ring (for calcite feeder)	2	61301605452*
GAC/KDF mix bed In-Out head 1"	1	61301605536*
*Part of PM Kit		61301604731
Vertical Pump/Motor for Processor		61301605500
Brass Solenoid 3/4" 24v		61301605489
Membrane Filter 4 x 40		61301605464
Membrane Closure Kit 4"		61301605467
Valve Check PVC 1"		61301605478
TDS Monitor Panel Mt. 0-1999ppm		61301605515
Solenoid Plas 1" 24v		61301605486
Gauge 0-100 Gly. Pnl. Mt.		61301605511
Valve Labcock 1/4"		61301605482
Housing Filter 20" BB		61301605465
Gauge 0-100psi 1/4"		61301605510
Sediment 20" BB 1Micron		61301605462
Valve PVC 1/2" Panel Mt.		61301605483
Control Box Motor Relay		61301605505
Fuse 30amp		61301605503
O-Ring for 20" Housing		61301605469

HS-4800-PT Only

Tank Pressure Switch	61301605491
----------------------	-------------

HS-4800-FT Only

Valve Check Brass 1/2"	61301605476
Horizontal Reed Float	61301605509
Tank for Repressurizer (RP-12)	61301605472
Injector for Repressurizer	61301605498
Pump/Motor for Repressurizer	61301605498

RO Water Systems

HS-9600

Kit Qty

20" BB 5M Sediment filter	4	61301605461*
10" all Calcite feeder	8	61301605448*
BB "O" Ring (for housing)	1	61301605468*
Slim "O" Ring (for calcite feeder)	2	61301605452*
GAC/KDF mix bed In-Out head 1"	1	61301605536*
*Part of PM Kit		61301605520
Vertical Pump/Motor for Processor		61301605500
Brass Solenoid 3/4" 24v		61301605489
Membrane Filter 4 x 40		61301605464
Membrane Closure Kit 4"		61301605467
Valve Check PVC 1"		61301605478
TDS Monitor Panel Mt. 0-1999ppm		61301605515
Solenoid Plas 1" 24v		61301605486
Gauge 0-100 Gly. Pnl. Mt.		61301605511
Tank Pressure Switch		61301605491
O-Ring for 20" Housing		61301605469
Valve Labcock 1/4"		61301605482
Housing Filter 20" BB		61301605465
Gauge 0-100psi 1/4"		61301605510
Sediment 20" BB 1Micron		61301605462
Valve PVC 1/2" Panel Mt.		61301605483
Control Box Motor Relay		61301605505
Fuse 30amp		61301605503

HS-9600-FT Only

Horizontal Reed Float	61301605509
Valve Check Brass 1/2"	61301605476
O-Ring for 20" Housing	61310605469
Pump/Motor for Repressurizer	61301605499
Tank for Repressurizer	61301605472
Valve Check PVC 80 1 1/4"	61301605479
Solenoid 1/2" Plastic Body	61301605487
Coil for Plastic Solenoid	61301605514
Valve PVC 40 1/2"	61301605480

RO Water Systems

HS-15000

Kit Qty

20" BB 5M Sediment filter	4	61301605461*
10" all Calcite feeder	8	61301605448*
BB "O" Ring (for housing)	1	61301605468*
Slim "O" Ring (for calcite feeder)	2	61301605452*
GAC/KDF mix bed In-Out head 1"	2	61301605536*
*Part of PM Kit		61301605521
Pump/Motor for Processor		61301605501
Control Box Contactor		61301605505
Fuse 30amp		61301605503
TDS Monitor 0-1999 Pnl. Mt.		61301605515
Valve PVC 1/2" Pnl. Mt.		61301605483
Valve Labcock 1/4"		61301605482
Vessel Closure Kit 4"		61301605467
Membrane Filter 4 x 40		61301605464
Tank Pressure Switch 40-60psi		61301605491
Solenoid Plastic 1" 24v		61301605486
Solenoid Brass 3/4" 24v		61301605489
Valve Check PVC 1"		61301605478
Housing Filter 20" BB		61301605465
O-Ring for Housing 20" BB		61301605469

HS-15000-PT Only

Gauge 0-100psi Gly. Pnl. Mt.	61301605511
------------------------------	-------------

HS-15000-FT Only

Horizontal Reed Float	61301605509
Pump/Motor for Repressurizer	61301605499
Tank for Repressurizer	61301605472
Valve Check PVC 80 11/4"	61301605479
Coil for Plastic Solenoid	61301605514
Valve PVC 40 1/2"	61301605480

RO Water Systems

Surgical Tables

MAQUET 1132.0XX ALPHASTAR

Electronic Components

Battery 12V	227.0154
Charger Board	970.2083
Transformer	3113.0513
Power Cord	3113.3669

Potentiometers

Leg & Back (Short Shaft, Older)	970.2243
Leg & Back (Long Shaft, Newer)	3113.6373
Tilt	3113.0624
Longitude & Height	3113.9754

Hand Control, Cable Connected	1132.90C0
Spiral Cable With Plugs	200.4293
Hook	7421.8394
Counter Sunk Screw (For Hook)	1118.0379
Tallow-Drop Screw (For Back Cover)	1136.0499
Infrared Hand Control	1132.91C0

Mechanical Components

Base Cover (1132.01)	3113.0943
Base Cover (1132.03XX W/ Red Stripe)	3113.4563
Left Leg Cover Plate Conv. Kit	3114.3619
Right Leg Cover Plate Conv. Kit	3114.3629
Body-Bound Rivet	104.1214
Hydraulic Fluid (1ltr Bottle)	3112.6369

Surgical Tables

MAQUET 1133.0XX ALPHAMAX

Electronic Components

Battery 12V	227.0154
Power Cord	3113.3669
Transformer	3113.6253
Charger Board	970.2873

1133.02F2 Potentiometers

Long. Shift & Height	3113.6594
Trend. & Tilt	3113.7384
Back & Leg	3113.7394

1133.02F3 Potentiometers

Longitude Shift	231.0614
Leg	231.0634
Back	231.0644
Height	231.0654
Trendelenburg	3113.0624

Cable Connected Hand Control	1133.90A0
Spiral Cable With Plugs	200.4293
Hook	7421.8394
Counter Sunk Screw (For Hook)	1118.0379
Tallow-Drop Screw (For Back Cover)	1136.0499
Infrared Hand Control	1133.91A0

Column Covers

Lowest (Steel) (1 Side)	3113.7461
Lower Middle (Steel) (1 Side)	3113.7471
Upper Middle (Steel) (1 Side)	3113.7481
Upper (Steel) (1 Side)	3113.7491
Top Cover, Blue (2 Sides)	8601.3353

Mechanical Components

Body-Bound Rivet	104.1414
Twin Swivel Castor	501.3392
Hydraulic Fluid (1ltr Bottle)	3112.6369
Base cover	8601.3113

Surgical Tables

MAQUET 1130/1131/1420/1425/1532/1533

Electronic Components

Power Cord Reel	3112.0604
Transformer	3112.1293
Battery (x4)	227.0054
Charging Unit	970.1582
4AT Fuse	217.0564
Infrared Charging Unit	3110.26B9
Battery For Infrared Hand Control	227.0064
Spiral Cable for Hand Control (With Plugs)	200.4293
Clip For Hand Control	506.0434

Mechanical Components

Snap-In Plate Connector (On Table)	8600.5784
Top Snap-In Plate Connector (On Plate)	8600.5124
Bottom Snap-In Plate Connector (On Plate)	8600.5134
Column Bellows	705.0191
Bellows Fixation Strip	8600.3864

Surgical Tables

Surgical Tables

SHAMPAINE / ALM ST004

Electrical Components

Power Entry

Fuse, 4AP, Type ET, 4.0A	B063224
Power Input Receptacle	B063371
Fuse Box For Receptacle (Drawer)	B063114
Power Cord Clip	P063373
Power Cord, 14 FT.	B063089
Power Cord, 20 FT.	B06308902
Power Input Fuse (4A)	B06311802

Batteries

4900B & 5100 B/E, ST004B Only

12V, 24 Amp Hr	303368
----------------	--------

4900, 4900E & ST004E Only

12V 4 Amp Hr.	B061957
---------------	---------

Controls (General)

Micro Switch, Long Lever	B061710
--------------------------	---------

4900B/E & ST004 Controls Only

Control Board	B063082S
Display Board	B063068S
Service Switch (4900B, ST004B)	B06334601
Service Switch (4900E, ST004E)	B06334602
Foot Control Switch Assembly	B062293S

5100B/E Controls Only

Foot Control Switch Assembly	B063614S
Back Sensor	B063593
CPU Board (On Column)	B063512S
Control Board (In Base)	B063509
Leg Sensor	B063592
Tilt & Trend Sensor Assembly	B063542S

Surgical Tables

SHAMPAINE / ALM ST004 Cont.

Hand Pendants

5100B/E	B063551S
4900B/E, 4900, 4800E, 1900RC, ST004	A06200806S
4900 Pendant Hook	P056841
1800RC	A058224

Low Voltage Supply Boxes

5100B/E	B063527S
4900E, ST004E	B063245S
4900B, ST004B	B063123S
1800RC	A058223

Hydraulic Flow Controls

Flow Control Valve	B056985
Pump/Motor DC (4900B/E,5100B/E,ST004)	A06285751
Pump/Motor AC (All Others)	302469
Kepner Check Valve	B056515
Kepner O-Rings (2 Per Valve)	B104450
Solenoid Valve (Complete)	B062048
Solenoid Coil 24VDC	B062049

Lift Rebuild Kits

4900B/E,5100B/E,ST004 (15" Lift)	A063006
4800,4900,1900M,19000RC,1700RC	A058603
1800RC, 2605NLB	A058604

Hydraulic Fluid (All Tables)	A032384
------------------------------	---------

Misc. Tools

Brass O-Ring Removal Tool	Grainger 4F418
Brass Coupling (Use For Balancing)	B052950
Swivel Fitting (Requires 2)	B054332
Hydraulic Hose 26"	A06200005

Surgical Tables

SHAMPAINE / ALM ST004 Cont.

PADS

1800RC,1900RC,4800E,4900,4900B/E,5100B/E,ST004

4 Piece 2" Set (Soft Memory Foam)	B063323
3 Piece 2" Set (Soft Memory Foam)	B06332301
4 Piece 4" Set (Soft Memory Foam)	B06332302
Deluxe Head Rest 2" (Soft Memory Foam)	B063323-1
Body Section 2" (Soft Memory Foam)	B063323-2
Leg Section 2" (Soft Memory Foam)	B063323-3
4 Piece 2" Set	B052998
Leg Section 2"	B055601
Body Section 2"	B055604
Deluxe Head Rest Pad 2"	B062210

2605,2605NLB Pads

2 Piece 1-1/4" Set (Soft Memory Foam)	302909
Leg Section (Soft Memory Foam)	302910
Body Section 1-1/4" (Soft Memory Foam)	302911
Head Rest 2-1/2"	B063296
Head Rest 4"	B063296-04
Body Section 1"	B054532
Foot Section 28"x33"x5"	B061780
Leg Section 2.75"	B054531

Common Mechanical Components

T-Handle For Lock Extension Attachments	A06186301
Foot Pad	B054963
Plug For Foot Pad	P005246
Set-Screw For Foot Pad	B105314

Deluxe Head Rest Assembly	302676
----------------------------------	---------------

Surgical Tables

Surgical Tables

CMAX

Electrical Components

Battery Set	342212555
Hydraulic Pump	342211555
Sliding Motor Assembly	342231555
Sliding Power Board	342232555
Switch Kit, Leg Section	342274555
Switch, Kidney	342252555
Sensor, Back Section	342242555
Sensor, Leg Section	342263555
Sensor, Sliding	542200196
Electronic Unit (PC Board)	542200903
Up/Down Sensor	542200194
Override Hand Control	542200905
Power Supply/Battery Charger	542200902
Fuse Set	342213555
Main Hand Control	542200904
Main Hand Control Cable	542200420
Sensor, Tilt	688200019
Sensor, Trend	688200019
Foot Controls	542200912
Power Cord, Mains	689400009

Hydraulic Components

Main Valve Block	342228555
Floor Lock Valve Block	342217555
Foot Cylinder Assembly	342215555
Trend Cylinder Assembly	342223555
Tilt Cylinder Assembly	342224555
Leg Section Assembly	342261555
Back Section Assembly	342241555
Lift Cylinder Assembly	342225555
Hydraulic Screws Set	342291555
Hydraulic Oil (2L/0.53 Gallon)	5420809.98

Surgical Tables

CMAX Cont.

Base And Column Components

Base Covers	342210555
Column, Top Plastic Covers Set	342222555
Column, Stainless Steel Covers Set	342221555
Column, Base Plastic Cover Kit	342220555
Casters Kit	342145555
Pad Kit	342216555
Trend Bellow Kit	342227555
Tilt Bellow Kit	342226555
Seat Plastic Cover Kit	342230555
Back Section, Table Top	342240555
Kidney Section, Table Top	342250555
Seat Section, Table Top	342260555
Leg Support, Top Plate	342270555
Head Section Support	342280555
Head Section Assembly	542200811

Pads

2" Leg Plate Pad	1006.P7L2
3" Leg Plate Pad	1006.P7L3
2" Table Top Pad	1006.P7T2
3" Table Top Pad	1006.P7T3
2" Headrest Pad	1006.P7H2
2" Table Top Pad	1006.P7P2

Surgical Lights

ALM PRISMATIC ECL / PRC

Bulb 24V 120W	ALM84000
----------------------	----------

Bulb Holders

Single Bulb Holder	ALM84001
--------------------	----------

Double Bulb Holder	ALM84002
--------------------	----------

Triple Bulb Holder	ALM84003
--------------------	----------

Red Knob For Bulb Holder (Must Order Fiber Board Assy)	
--	--

Fiber Board Assy	ALM85897
------------------	----------

Sterilizable Handle	ALM84004
----------------------------	----------

1 Dozen Sterilizable Handles	ALM84005
------------------------------	----------

Controls

Dimmer Control (Potentiometer)	ALM85017
--------------------------------	----------

Knob	ALM85795
------	----------

Fuse Holder	ALM85101
-------------	----------

Fuse	ALM85102
------	----------

Red Indicator Light	ALM85015
---------------------	----------

Transformer 150 KVA TAPS	ALM85024
--------------------------	----------

Transformer 400 KVA TAPS	ALM85025
--------------------------	----------

Misc. Components

Male / Female Rotating Contact	ALM85898
--------------------------------	----------

Wave Washer Contact Set (Light head)	ALM85202
--------------------------------------	----------

Gray Handle Assy	ALM85145
------------------	----------

Copper Segment	ALM85165
----------------	----------

Cupola Adjust Tool (Spanner Wrench)	ALM85713
-------------------------------------	----------

Surgical Lights

ALM ANGENIEUX

Bulbs

AX4 70W 24V	ALM353092
AX4 100W 24V	ALM373431
AX4 150W 24V	ALM363508
AX8 250W 24V	ALM376554
AX10 100W 24V	ALM386762
AX14 150W 24V	ALM363508

Bulb Holder

ALM362798

Sterilizable Handle

ALM300049

(1 Dozen) Sterilizable Handles

ALM300050

Controls

Dimmer Control (Potentiometer)	ALM85017
Knob	ALM85795
Fuse Holder	ALM85101
Fuse	ALM85102
Red Indicator Light	ALM85015
Transformer 150 KVA TAPS	ALM85024
Transformer 400 KVA TAPS	ALM85025

Misc. Components

Spring Arm Adjustment cover	ALM379426
Copper Segment	ALM85165
Copper stop spring arms	ALM304867

Surgical Lights

ALM PRISMALIX (PRX)

Bulb 100W 23V	ALM50067
Package of 10 Bulbs	ALM50170

Bulb Holders

Single Bulb Holder	ALM50363
Double Bulb Holder	ALM50364
Triple Bulb Holder	ALM50365

Sterilizable Handle	ALM50258
(5) Sterilizable Handles	ALM50171
Fixed Focus Camera Handle	ALM50259
Video Zoom Camera Handle (Square)	ALM50369

Devon adapters (PRX)

*Devon soft glove adapter ALM50240 (must be ordered with each handle adapter listed below)

PRX 4000	ALM50223
PRX 6000	ALM50224
PRX 8000	ALM50225
Dummy VZ For Video Zoom	ALM50226
Flat Screen Yoke Adapter	ALM50336

Controls (WPS I)

Fuse 6.3A MINI 5x20mm	ALM50262
WPS Solid State Power Module	ALM50059
Potentiometer	ALM50188
LED (Green)	ALM50190
Control Knob (Gray)	ALM50187
CFF Power Supply Board (In Main Arm)	ALM50174

Misc. Components

Gray Plug With Bumper	ALM50237
Copper Segment 4000/6000	ALM85165
Brake Pin Repair Kit	ALM50347
Cupola Adjust Tool (Spanner Wrench)	ALM50315

Surgical Lights

ALM X'TEN & AXCEL

Bulb 100W 24V	ALM386762
----------------------	-----------

Bulb Holders

X'TEN	567804999
-------	-----------

AXCEL	567501950
-------	-----------

Sterilizable Handle	ALM50258
----------------------------	----------

(5) Sterilizable Handles	ALM50171
--------------------------	----------

Fixed Focus Camera Handle (X'TEN Only)	ALM50259
--	----------

Video Zoom Camera Handle (X'TEN Only)	ALM50369
---------------------------------------	----------

Devon adapters

X'TEN	567501253
-------	-----------

AXCEL	567501254
-------	-----------

ALM ENERGIX (WPS II) CONTROLS

Transformer 750VA	567314007
-------------------	-----------

Transformer 450VA	567314002
-------------------	-----------

Dimmer Keypad	567314001
---------------	-----------

Blank Keypad	567314011
--------------	-----------

Control Overlay Sticker	567314019
-------------------------	-----------

Regulator Board 180W	567314004
----------------------	-----------

Regulator Board 400W	567314006
----------------------	-----------

Filter Board	567314003
--------------	-----------

Control Board	567314009
---------------	-----------

Surgical Lights

LDR (CONCEALED LIGHTING BOX)

Linear Motor (Actuator)	ALM87082
Door Skin (Bottom)	ALM87078
Door Skin (Top W/ Ramp & Guides)	ALM87091
Tape Switch 4ft-6" (With Channel)	ALM87088
Tape Switch 1ft-6" (With Channel)	ALM87092
Relay, 24V	ALM87089
Circuit Breaker, W28XQ1A, .25A	ALM87090
Limit Switch For Door	ALM87097
Rocker Switch	ALM87055
Fuse, 6A (MDL 6)	ALM87068
On/Off Switch (W/ Amber Light)	ALM87069
Docking Pin For PRC	ALM87071
Docking Pin For AX or PRX	ALM87093

Surgical Lights

21st CENTURY SERIES

Bulbs

72131(Century) 75W 12V	ALM701257
72151(Centurion) 90-100W 12V	ALM701256
72181(Centura) 150W 20V	ALM701255
Surgilite II (Old) 150W 120V	ALM71800
Surgilite II A (After s/n 059601)150W 120V	ALM71809
7111(11"Refl) 100W 120V	ALM71727
7131(13"Refl) 150W 120V	ALM71727
7151(15"Refl) 150W 120V	ALM71727
7241(Spot) 250W 120V	ALM70240
7291(Spot) 300W 82V	ALM72829
Surgilite I 150W 120V	ALM86010

Bulb Holders

Centurion / Centry	ALM700088
Centura (Must Order ALM700973 & ALM71089)	

Sterilizable Handle

	ALM71282
Post For Sterilizable Handle	ALM701144
Sterilizable Handle For Surgilite II	ALM71822

Misc. Components

Diaphragm Adjust. Handle (Centura Only)	ALM701206
Side Cap Cover For Spring Arm	ALM702045
Toggle Switch	ALM71001
Switch Cover	ALM702040
Transformer For Centura	ALM701223
Transformer For Centurion, Centry	ALM701250
Cover For Transformer	ALM701159

Surgical Lights

HERAEUS 2000 SERIES

Bulb (22.8V - 50W)	567901955
Package of 8 Bulbs	567901954

Sterilizable Handles (Pkg. of 3)	5605-0756
Metal Sterilizable Handles (2007 Only)	5605-0755
Handle Repair Kit 2003/4/5 (White)	5605-0751
Handle Repair Kit 2003/4/5 (Green)	5605-3136

Controls

Dimmer Knob SLV600 PM	1400-1308
Dimmer, Wall Slide w/ Rocker Sw.	SLV600PM
Dimmer, Nova, Wall (Slide to Off)	0007-5601

Misc. Components

Adjustment Rods - 2002/3/4/5 (3)	5605-2343
Adjustment Rods - 2002/3/4/5 (4)	5605-2344
Lens 2001/2/3/4 std	5605-0734
Lens 2001/2/3/4/5i/5std	5605-2609
Lens 2004/5ixl	5605-2962
Locking Assy. London Lens	5605-0435
Plug & Socket, Radial 2001 - 5	5605-0713
Plug & Socket, Radial, London	5605-2305
Reflector 2000i/ixl	5605-0726
Reflector 2000 Std.	5605-2961
Covers, Spring Arm 2001/2/3/4/5	5605-1771

Surgical Lights

HERAEUS CITY SERIES

Bulbs

Hamburg (24V - 25W)	5601-6164
Hamburg (110V - 30W)	5601-6170
Hamburg (125V - 30W)	5601-6238
Hamburg (32V - 40W)	5601-6579
Hamburg (50V - 30W)	5601-8505
Hamburg (24V - 40W)	5601-8550
Boston, London (125V - 50W)	5601-6578
Boston, London (24V - 50W)	5601-6678
Boston, London (50V - 50W)	5601-6685
London (32V - 50W)	5601-6813
Amsterdam (22.8V - 75W)	5601-8366
London, Frankfurt (22.8V - 50W)	5601-8566
Melbourne, Oslo, Heidelberg (22.8V-40W)	5601-8769

Handles

"Universal"	5605-2596
Heidelberg	5605-2595
London, Oslo, Hamburg, Hospital	5605-2315

HERAEUS BLUE 30 & 80

Bulbs

Blue 30 (24V - 50W)	5605-3010
Blue 80 (21.5V - 130W)	5605-3026
Bulb Socket	5605-3073

Sterilizable Handles (Pkg. of 3)	5605-0756
---	-----------

On/Off Switch	5605-3016
Blue 30 Lens Cover	5605-3013
Blue 30S Lens Cover	5605-3084
Blue 80 Lens Cover	5605-3067

Surgical Lights

CASTLE 2000 SERIES

Bulb 24V 60W	516767
Bulb Socket	67501

Sterilizable Handle	516912
---------------------	--------

Controls

Dimmer, Incandescent	517022
Circuit Breaker 7.5A	69997

Misc. Components

1/4 Turn Lens Retainer	516854
1/4 Turn Snap-In Receptacle	516855
Heat Filter Cylinder	516768
Retaining Spring for Heat Filter	65064
Retaining Hook for Heat Filter	65032
Lens	516780
Reflector	516766

CASTLE 2500-2900 SERIES

Bulbs

Bulb 120V, 150W (2500-2900, Old)	49455
Bulb, 120V, 100W (2500-2900, Old)	60976
Bulb, 24V, 150W (2500-2900, New)	68078
Bulb, 24V, 100W (2500-2900, New)	68077
Bulb Socket	67501

Controls

Dimmer, Solid State	67969
Circuit Breaker 7.5A	69997

Gas Spring, 2 Pod, Vertical (Blue)	68906
Gas Spring, 2 Pod, Horizontal (Brown)	67503
Gas Spring, 5 Pod	67421

Heat Filter Cylinder	516768
Retaining Spring for Heat Filter	65064
Retaining Hook for Heat Filter	65032
Reflector, Surgery (Dichroic)	65012

Surgical Lights

CASTLE 2300-2400 SERIES

Bulbs

Bulb 150W 120V (Surgery Light)	49455
Bulb 100W 120V (Exam Light)	60976
Bulb 60W 24V (Mobile Field Light)	516767
Bulb Socket	67501

Sterilizable Handle	65196
----------------------------	-------

Misc. Components

Gas Spring For Balance Arm	65024S
----------------------------	--------

Heat Filter Cylinder	516768
----------------------	--------

Retaining Spring for Heat Filter	65064
----------------------------------	-------

Retaining Hook for Heat Filter	65032
--------------------------------	-------

Lens	65014
------	-------

Lens (Center Mount Handle)	57614
----------------------------	-------

Reflector, Exam (Aluminized)	65011
------------------------------	-------

Reflector, Surgery (Dichroic)	65012
-------------------------------	-------

Dimmer Switch 120V	67969
--------------------	-------

GETINGE/CASTLE 500-700

Bulb 150W 24V	303930
----------------------	--------

Bulb Holder 500	303932
-----------------	--------

Bulb Holder 700	303933
-----------------	--------

Sterilizable Handle	516912
----------------------------	--------

Overlay For Single Dimmer	304506
---------------------------	--------

LED PC Board (Behind Overlay)	304521
-------------------------------	--------

Dimmer PC Board (Behind LED Board)	303951
------------------------------------	--------

Circuit Breaker 8A	305050
--------------------	--------

Transformer 800VA	303306
-------------------	--------

Burnt Bulb Indicator PC Board	303950
-------------------------------	--------

Overlay For Burn Bulb Indicator	305014
---------------------------------	--------

M/C STERILIZERS

3400 Parts	30837
3422 Service	30836
3422-23 Operator	30835
3433 Operator	30845
3471 Operator	30838
3500 Parts	350065
3500 Service	350063
3522-23 Operator	350059
3525 Operator	350060
3533 Operator	350061
3571 Operator	350062
3600 Parts	350116
3622 / 3633 Service	350115
3622 / 3633 Operator	350113
3671 Operator	350119
3671 Service	350278
4200 Parts	350106
4200 Service	350105
4233 Installation	350101
4233 Operator	350102
4271 Operator	350107
6900 Operator	350088
6900 Parts	350090
6900 Service	350089
EO/CO2 Dual Load Station Parts	350564
EO/CO2 Dual Load Station Service	350427
M/C II Controls Service	350064
M/C II Message Center Operator	350110

PACS-2000 STERILIZERS

122/123/133 Installation	61301600294
122/123/133 Operator	61301600298
122/123/133 Service	61301600299
122/123/133 Parts	61301600300
100 LS/HC Installation	61301603358
100 LS Operator	61301603355
100 LS Supervisor	61301603357
100 HC Operator	61301603363
100 HC Supervisor	61301603365
100 LS/HC Service	61301603359
100 LS/HC Parts	61301603361
222 / 233 Installation	61301600862
222 / 233 Operator	61301600866
222 / 233 Service	61301600867
222 / 233 Parts	61301600868
233 LS/HC Installation	61301603648
233 LS Operator	61301603645
233 LS Supervisor	61301603647
233 HC Operator	61301603652
233 HC Supervisor	61301603654
233 LS/HC Service	61301603649
233 LS/HC Parts	61301603650

PACS-3000 STERILIZERS

400-500 Installation	61301608205
400-500 LS Operator	61301608206
400-500 LS Supervisor	61301608207
400-500 HC Operator	61301608208
400-500 HC-LS Service	61301608209
400-500 HC-LS Parts	61301608210

700 HC-LS Installation	61301608214
700 LS Operator	61301608683
700 LS Supervisor	61301608684
700 HC Operator	61301608685
700 HC-LS Service	61301608215
700 HC-LS Parts	61301608216

833 Installation	61301609130
833 HC Operator	61301609128
833 LS Operator	61301609127
833 Service	61301609133
833 Parts	61301609398

STEAM GENERATORS

CES-20	30422
CAS-20A & CAS-20AF	30893
CAS-20C	350208
CAS-20D	350510
CAS-20E	350582
CAS-30A & CAS-36A	350646
CAS-30C Service	400/500 Service
CAS-30C Parts	400/500 Parts

ANCILLARY

7510 Steam Gun Operator	350100
SS-10 Service	350255
SS-10 Parts	350256
GSS-10 User, Service & Parts	61301605235

SONIC CLEANERS/RINSER DRYERS

7946 & 7947 Service & Parts	30451
2415-RD Operator	350491
2415-RD Parts	350604
2415-RD Service	350601
2415-UW Operator	350490
2415-UW Parts	350594
2415-UW Service	350593
2460-RD User+Service+Parts	61301605109
2460-UC User+Service+Parts	61301605110

WASHERS

7515 Parts	30931
7515 Service	30930
7535 / 7936 / 7937 Operator	30436
7535 / 7936 / 7937 Service & Parts	30437
7538 Operator	350054
7538 Parts	350055
7538 Service	350056
7550 Operator	350457
7550 Parts	350375
7550 Service	350374
7620 Operator	350527
7662 Installation	350541
7662 Operator	350518
7700 Installation	61301603904
7800 Installation	61301604981
7800 Operator	61301604768
7800 Service	61301605323
7900 Installation	61301604805
7900 Operator	61301604804
7900 Service	61301605323
4656 (I) (O) (S) (P)	61301603745
8666 (I) (O) (S) (P)	61301603744

WARMING CABINETS & AERATORS

5520-25-30 Operator	30923
3441 Operator	30924
4141 Operator	30925
5520 / 5525 / 5530 Service	30926
3341 / 3441 / 4141 Parts	30927
5520 / 5525 / 5530 Parts	30927
5540-45-50 Operator	350264
5540 / 5545 / 5550 Service	350265
35-36-4241 Operator	350268
3541 / 3641 / 4241 Parts	350270
5540 / 5545 / 5550 Parts	350270
5524 / 5618 / 5624 (I) (O) (S) (P)	61301601083
5524 / 5618 / 5624 (I) (O) (S) (P) NEW	61301603163

OPERATING TABLES

1900 Operator	350466
All Maquet Manuals	ALM89300
5100B-E Operator	350519
5100B/E Parts	350552
5100B/E Service	350520
4900B/E (ST004) Operator	350388
4900B/E (ST004) Parts	350390
4900B/E (ST004) Service	350389
2605 Operator	350461

CASTLE LIGHTS

2500-2800 Parts	
2000 Installation	350178
2000 Operator	350179
2000 Parts	350181
2000 Service	350180
2300/2400 Operator	30705
2300/2400 Service	30706
500-700 Installation & Operator	303189
500-700 Parts	30975
500-700 Service	30706

GETINGE

Getinge USA, Inc.
1777 East Henrietta Road
Rochester, New York 14623-3133
U.S.A.
Phone: (800) 475-9040
Fax: (585) 272-5033
info@getingeusa.com www.getingeusa.com