

# Varian 310-MS LC/MS

The Varian 310-MS LC/MS provides maximum performance for pharmaceutical, toxicological, environmental, and other applications requiring quantitation in a heavy matrix. The API inlet is sensitive and robust and the unique MS analyzer consistently delivers superior signal-to-noise.

The 310-MS is flexible, affordable, and easy to use making it an ideal system for LC/MS.

## Specifications

- Analyzer mass range: 10 to 2000 Da
- Scan rates: 500 to 6000 Da/sec
- Resolution: unit mass over the entire range
- Linear dynamic range: up to  $10^5$ , compound-dependent based upon analytical methods
- Mass axis stability:  $\pm 0.1$  Da over 24 hours
- Manifold temperature: independent control; 30 °C to 50 °C
- Ion detector: positive or negative ions, constant 5 kV post acceleration voltage
- Triple quadrupole modes: Q1MS and/or Q3MS for Selected Ion Monitoring (SIM) or Full Scan, Precursor Scan, Product Scan, Neutral Loss Scan, Selected Reaction Monitoring (SRM)
- Turbomolecular pump: dual stage, 310/400 L/sec, air cooled
- High efficiency electron multiplier: 240 ms positive to negative switching time
- Adjustable dwell times: <5 ms

## Atmosphere Pressure Ionization (API) Interface

- SelecTemp™ allows temperature programming of the drying gas throughout the analysis
- SelecFlow™ allows electronic flow control of the drying gas throughout the analysis
- Drying gas temperature: 400 °C maximum at 8 L/min
- Drying gas flow rate: up to 10 L/min
- Spray chamber temperature: up to 60 °C
- Spray needle: off-axis from the capillary axis
- Quadrupole ion guide: 6° off-axis for the capillary axis
- Capillary voltage: up to 300 V

## Electrospray Interface (ESI)

- LC flow range: 1 to 1000  $\mu$ L/min
- Needle voltage: up to 6 kV
- Nebulizing gas flow: up to 2 L/min
- Spray plate voltage: up to 800 V
- Spray needle: wide range of adjustments for x-y positioning; independent adjustment of the inner liquid capillary needle to the outlet of the nebulizing gas

## Mass filters (Q1 and Q3)

- Quadrupole with pre- and post-filters (guides)
- Lensless design between the quadrupoles

## Collision cell (Q2)

- 180° curved cell with pre- and post-filter region:
  - In single and triple quadrupole modes, the curvature of Q2 minimizes noise by positioning the multiplier off-axis from the source.
  - In MS/MS modes, the long, 185 mm path ensures high dissociation efficiency.

## Vacuum System

- Foreline pump: one Varian HS 652 universal voltage pump
- Pump includes:
  - oil mist eliminator
  - oil return kit

## Utilities and Environment

- Power requirements for the MS: 100-240 Vac, 50/60 Hz  $\pm 3$  Hz, 1.0 kW
- Power requirements for the foreline pump: 100-240 Vac, 50/60 Hz, 2.8 kW
- Venting for the API spray chamber: up to 12 L/min
- Venting for the foreline pumps: up to 2 L/min
- Nitrogen (+ESI) nebulizing and drying gas: up to 12 L/min, regulated at 80 psi
- Air (-ESI) nebulizing gas: up to 2 L/min, regulated at 80 psi
- Humidity: 20% to 80% relative humidity (without condensation)
- Temperature: 16 °C to 30 °C

NOTICE: This document contains references to Varian. Please note that Varian, Inc. is now part of Agilent Technologies. For more information, go to [www.agilent.com/chem](http://www.agilent.com/chem).



**Agilent Technologies**

# Varian 310-MS LC/MS

## User Upgradeable Options

The APCI option is a separate needle/chamber that uses the same hinge/clip mounting system as the ESI chamber. Changing from ESI to APCI takes only a few minutes to disconnect and connect the appropriate electrical, gas and LC lines and does not require venting the instrument.

## APCI Specifications

- LC flow range: 100 to 2000  $\mu$ L/min
- Corona current: -50 to +20  $\mu$ A
- Nebulizing gas flow: up to 2 L/min
- Auxiliary gas flow: 4 L/min typical
- Auxiliary gas temperature: 550 °C maximum; optimum is compound, mobile phase and flow dependant
- Spray plate voltage: up to 800 V
- Capillary voltage: up to 300 V
- Drying gas temperature: 400 °C maximum at 8 L/min
- Drying gas flow rate: up to 10 L/min
- Spray chamber temperature: up to 60 °C
- Spray needle: off-axis from the capillary axis
- Quadrupole ion guide: 6° off-axis for the capillary axis

## Chromatography Options

The MS workstation fully controls a wide range of chromatography products including:

- AutoSamplers: ProStar™ 400, 410, 420, 430, and HTS PAL
- Solvent Delivery Modules: 212-LC; ProStar 210, 220, 230, and 240; PrepStar™ 218, SD-1, and SD-2
- Column Valve Module: ProStar 500
- Detectors: ProStar 325, 340 UV-Visible Absorbance, 335 Photodiode Array, 363 Fluorescence and PL-ELS 2100 ELSD
- Fraction Collectors: ProStar 701 and 704

## 310-MS LC/MS Dimensions

- 50 cm (w) x 38 cm (h) x 71 cm (d) / 9.5 in. x 15 in. x 28 in., 108 kg/237 lb

## Performance Specifications

Mode	Test Standard	S/N*
+ESI SIM	5 pg Reserpine	$\geq 10:1$
+ESI MS/MS	5 pg Reserpine	$\geq 20:1$

\*Based on RMS noise, 0.7 PWHM resolution



GC • LC • MS • GPC/SEC • AA • ICP • ICP-MS • UV-Vis-NIR • FT-IR • Fluorescence • Dissolution • NMR • MRI • FTMS • Consumables • Data Systems

Varian, Inc. reserves the right to revise these specifications without notification.  
SelectTemp, SelectFlow, ProStar, PrepStar, Varian and the Varian logo are trademarks or registered trademarks of Varian, Inc. in the U.S. and other countries. © 2008 Varian, Inc.

Varian, Inc.  
www.varianinc.com  
North America: 800.926.3000, 925.939.2400  
Europe The Netherlands: 31.118.67.1000  
Asia Pacific Australia: 613.9560.7133  
Latin America Brazil: 55.11.3845.0444  
Other sales offices and dealers throughout the world—  
check our Web site

SI-1121 03/08 Printed in U.S.A.