

The Agilent 6460 Triple Quadrupole LC/MS achieves sensitivity and resolution specifications after autotune without manual user intervention.

Parameter	Measure	Specification
Electrospray ESI MS/MS sensitivity	1 pg of reserpine injected on column, with Agilent Jet Stream technology, quantifying on the transition $m/z$ 609 to 195.	Signal-to-noise ratio >1000:1
Mass range		m/z 5-3,000
Polarity switching		Switch from positive ion mode to negative ion mode in 100 ms
Mass resolution		Three settings (FWHM) Unit: 0.7 u; Wide: 1.2 u; Widest: 2.5 u
Mass accuracy		0.1 Da from 5–1,000 <i>m/z</i> 0.01% from 1,000–2,000 <i>m/z</i> 0.02% from 2,000–3,000 <i>m/z</i>
Mass stability		<0.1 u in 24 h
Dynamic range		>6.0 × 10 <sup>6</sup>
Scan rate		5,200 u/s
Scan modes		MS scan, MS/MS product ion scan, MRM, MS/MS neutral loss/gain scan and precursor ion scan, SIM
Minimum MRM dwell time		1 ms
MRM transitions		500 per time segment, >10,000 MRMs possible in a method
Collision cell		High-pressure, hexapole with linear acceleration
Cross-talk		None detectable

## **General System Specifications**

Parameter	Measure	
Single point of control	Single-point data system method capability with full control of Agilent 1100 and 1200 Series HPLC systems and 6460A Triple Quadrupole LC/MS System	
Time programming	<ul> <li>Polarity change in time segment</li> <li>Scan and SIM or MRM (plus other modes of data collection)</li> <li>Solvent divert through calibrant delivery system valve</li> </ul>	
Wide range of ionization sources	<ul> <li>Electrospray (ESI)</li> <li>Nanospray with HPLC-Chip Cube MS interface</li> <li>Multimode source (simultaneous ESI and APCI)</li> </ul>	
Autotune	Automated optimization of ion optics and mass axis calibration using Agilent Jet Stream Technology	
Solvent declustering	Countercurrent gas	
Detector	High-energy conversion dynode and high-gain electron multiplier horn	
Vacuum system	Two turbomolecular pumps with one mechanical pump	

## **Ordering Information**

Description Product Number

## **LCMS Triple Quad HS Bundle**

Includes G6460 Triple Quadrupole LC/MS System, workstation data system (PC) and printer, MassHunter Workstation software with 2 data analysis software licenses. Includes 24 hours installation and familiarization training, and image recovery service.

G6460AA

Agilent shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

Information, descriptions, and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc. 2009

Printed in the USA March 3, 2009 5989-9905EN

