



SPECIFICATIONS

708-DS Dissolution Apparatus

The 708-DS Dissolution Apparatus lets you choose exactly the options you need to support simple, manual dissolution testing, as well as integrated automation for ultimate productivity and throughput. Regardless of the options you choose, you can be sure you will have a dissolution apparatus that meets your expectations for productivity, robustness and ease of use. Be confident knowing that as your needs change, you can add options to your existing apparatus.



708-DS Specifications

Condition		Performance may vary depending on environmental conditions (temperature, humidity, altitude, etc.)	
Evaporation		Less than 1% evaporative loss under specific conditions	
Ambient temperature		5 to 40 °C	
Humidity (non-condensing)		Not more than 80% RH	
Requirements	Voltage	Current 708-DS (operating)	Current Wide Input Range Heater (operating)
	90-250 V, 50-60 Hz	2.5-1.0 A	13.0-4.5 A
	115 V, 60 Hz	2.0 A	10.0 A
	230 V, 50 Hz	1.0 A	5.0 A
Water bath	Temperature Range	Probe Accuracy	
	Ambient +5 to 55 °C	±0.1 °C	
Spindle	Speed Range	Speed Accuracy	Speed Selection
	10-250 RPM	±1% over 25 RPM ±2% 10-25 RPM	Via touch screen
Sampling		Manifold option available for automated sampling	
Display		LCD with integral touch screen	
Spindle shaft material		Stainless steel	
Drive unit lift		Motorized drive	
Optional features		Manual drive unit lift, Dosage Delivery Module (DDM), Autosampling, AutoTemp, resident sampling cannulas, handheld temperature probe, printer	
Dimensions		Width: 67.31 cm (26.5 in.) Operating Height: 67.95 cm (39 in.) Clearance Height: 99.06 cm Depth: 76.2 cm (30 in.)	
Weight		54.4 kg (120 lb) machine dry with vessels and paddles	

To learn more about
Dissolution Testing, visit us at
www.agilent.com/lifesciences/dissolution

This information is subject to change without notice.
© Agilent Technologies, Inc. 2011
Printed in U.S.A., October 17, 2011
5990-6752EN

The Measure of Confidence



Agilent Technologies