

NOTICE: This document contains references to Varian.
Please note that Varian, Inc. is now part of Agilent
Technologies. For more information, go to
www.agilent.com/chem.



Agilent Technologies

APPLICATION NOTE

Beta Blockers

Application 2002 LC

Column : Polaris 3 μ C18-A, 50 x 2.0mm
Part Nr.: A2001050X020
Eluent : A: H₂O + 0.1% NH₄OH, pH 10.0 (25 mM)
B: MeCN A:B 65 : 35
Flow : 0.6 ml/min
Detection : UV @ 215nm
Temp : 40°C

Peak Identification

1. pindolol
2. metoprolol
3. propranolol

