





Agilent G2614A ALS Tray and G2615A Bar Code Reader

Information Sheet

At Agilent, we are continuously striving to improve our products in order to better satisfy our customers' needs. The enclosed G2614A tray has recently been improved to provide even more reliable vial handling. In addition, Agilent recently introduced the new G2615A bar code reader designed specifically for the G2614A tray. This reader mounts securely to the G2614A tray and alleviates the need to adjust a positioning bracket.

The new G2615A bar code reader replaces the G1926A bar code reader to offer the greatest reliability with your G2614A tray. If you did not purchase a new bar code reader with this system, you can very easily add one at a later date. Your new tray has been factory calibrated based on the fixed positioning of the G2615A bar code reader for simple installation and optimal performance.



Agilent Technologies



Agilent Technologies

© Agilent Technologies, Inc.

Printed in USA, March 2003



G2614-90115