

NOTICE: This document contains references to Varian.
Please note that Varian, Inc. is now part of Agilent
Technologies. For more information, go to
www.agilent.com/chem.



Agilent Technologies

APPLICATION NOTE

Goldenseal

Application 2012 LC

Column : Polaris 5 μ C18-A, 150 x 4.6mm
Part Nr.: A2000150X046
Eluent : A: H₂O + 0.1% TFA,
B: MeCN + 0.1% TFA
Gradient : 20 - 65%B in 10 min.
Flow : 1.0 ml/min,
Temp : 35°C
Detection : UV @ 220nm,

Peak Identification

1. hydrastine,
2. berberine

