

NOTICE: This document contains references to Varian.
Please note that Varian, Inc. is now part of Agilent
Technologies. For more information, go to
www.agilent.com/chem.



APPLICATION NOTE

Sulfa drugs

Application 2006 LC

Column : Polaris 5 μ C18-A, 150 x 4.6mm
Part Nr.: A2000150X046
Eluent : 20mM K₂HPO₄, pH 9.0: MeCN 92:8
Flow : 1.5 ml/min,
Detection : UV @ 270nm,

Peak Identification

1. sulfanilic acid
2. sulfamerazine
3. sulfisoxazole
4. sulfapyridine

